

**rea**<sup>®</sup>

# Azalea<sup>®</sup> Tall



No limits on comfort for  
taller users



Yes, you can.

**Rea Azalea Tall** is specially designed to meet the needs of tall users who require a Tilt in Space wheelchair with a longer seat support. Adapted from the **Rea Azalea**, the **Rea Azalea Tall** boasts all the advantages of a reliable, tilting wheelchair and offers a unique weight-shifting mechanism. The **Rea Azalea Tall**: at the top of its class for comfort and stability. The improved Flex 3 backrest adds individual adaptability. 5 frame colours and three upholsteries are available.

## Real compatibility

The **Rea Azalea Tall** has all the features of the **Rea Azalea** such as Tilt in Space weight-shifting for a really stable tilt. All accessories such as rear wheels for self-propelling or smaller wheels for transport and a different seat system also fit the **Rea Azalea Tall**.

## Easy maneuverability

The adjustable push bar has been designed with the carer in mind, it offers a good ergonomic grip and can be adjusted to suit the needs of a carer.



### Longer seat plate

An extra 50 mm of seat depth is designed to support longer legs.



### Chassis length

An extra 50 mm in length enables greater forward stability.



### Seat height

An extra 50 mm of height allows for longer leg-rests.

For other Rea Azalea accessories, please consult the Rea Azalea brochure or visit [www.invacarearea.com](http://www.invacarearea.com) to see videos of them

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width ▼	Seat depth ▼	Seat height ▼	Backrest height ▼	Armrest height ▼	Length of footrests ▼	Seat angle ▼	Backrest angle ▼
<b>Azalea Tall</b>	390-450/ 440-500/ 490-550 mm (390-490/ 440-540/ 490-590 mm)	480 - 550 mm	500 mm	Backrest plate: 560 - 790 mm Tension adjustable backrest: 540 - 715 mm	240 - 360 mm	330 - 627 mm	-1° to 23° (+1° to 25°)	0° - 30°
	Driving unit width ▼	Total height ▼	Total length inc. Footrests ▼	Total weight ▼	Max. User weight ▼	Transport weight ▼	Crash tested ▼	Maximum slope for use of parking brake ▼
<b>Azalea Tall</b>	SW + 220 mm (SW + 230 mm)	1000 - 1300 mm	950 - 1085 mm	34 kg (SW 440 mm)	Max 135 kg	18,7 kg <sup>1</sup>	ISO 7176-19	7°

## Frame colour



## Upholstery



1. Transport weight increases with 4,2 Kg when Power tilt & recline option is used



Tel: +44 61 487 70 70  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
 Neuhofweg 51  
 4147 Aesch BL  
 Switzerland

Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare International. The **Rea Azalea** is designed with DSS (Dual Stability system). It provides the user with enhanced stability and safety when being tilted or the backrest is reclined.

