

# Softform Premier MaxiGlide



A high risk mattress designed to help prevent the build up of shear and friction



Yes, you can.


The *Invacare **Softform Premier MaxiGlide*** is the latest Invacare mattress to incorporate the patented “Glide” mechanism. The **MaxiGlide** provides the user with exceptional pressure reduction and high levels of comfort and support, while uniquely conforming to the profile of an electric ward bed, thereby significantly reducing the risk of shear and friction forces.



### The unique Glide mechanism

The unique Patented Glide system allows the mattress to fully conform simultaneously to the Patient and bed frame, optimising the properties of a profiling bed. This significantly reduces shear and friction forces, particularly in the heel and sacral areas, which would normally be experienced by a patient on most static foam mattresses.

### Features and Options




#### Protective flap to conceal the zip

Reduces the risk of possible fluid ingress.



#### Two-way stretch, PU water proof cover

Provides exceptional recovery and durability.



#### Pressure Ulcer prevention

As well as offering the patented Glide mechanism, the mattress also provides an excellent pressure reducing surface, designed to maximise the body contact area and prevent the build up of pressure ulcers.

### Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.



247 kg

80°

**Softform  
Premier  
MaxiGlide**

Manufactured to European Standards:  
Interior Foam: BS 5852 Crib 5,  
Cover: BS 7175 Crib 5

### Colour

Baltic Blue



**Tel: +44 61 487 70 70**  
hqueurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



**Yes, you can.®**