

AQUATEC®

Sorrento



Shower and Bathtub Collection



The outstanding modular and modern *Aquatec Sorrento* range is a complete Shower and Bathtub collection, the Shower Chair, Stools as well as the Bathtub Swivel Seat and the Bathtub Seat include great function with great style and design. The non-slipping soft touch surfaces and the soft backrest provide a secure and very comfortable seat for the user.

The wide adjustment range enables an easy fit to every user. The all-corrosion free, modular shower chair range is the best choice for those who need assistance in the bathroom for their personal hygiene.

Adjustable soft backrest

The soft and adjustable backrest of the *Aquatec Sorrento* provides a very comfortable feeling for the user. With the wide adjustment range every user will feel pleased due to an easy personal adaptation.

Non-slipping Soft Touch Surface

The *Aquatec Sorrento* is the first shower chair with the innovative TPE Soft Touch Surface on the seat and the armrests. Not only the great look but also the comfortable feeling will persuade the user.



Sorrento as Shower Chair



Sorrento as Shower Stool



Sorrento as Bathtub Swivel Seat



Sorrento as Bathtub Seat

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	<i>Sorrento Shower Chair</i>	<i>Sorrento Shower Stool w/ armrests</i>	<i>Sorrento Shower Stool</i>	<i>Sorrento Shower Chair w/o recess</i>	<i>Sorrento Shower Stool w/ armrests, w/o recess</i>	<i>Sorrento Shower Stool, w/o recess</i>	<i>Sorrento Bathtub Swivel Seat</i>	<i>Sorrento Bathtub Seat</i>
Shower Chair	1532633			1532730				
Shower Stool w/armrests		1532660			1532729			
Shower Stool			1532661			1532728		
Bathtub Swivel Seat							1532693	
Bathtub Seat								1532715

Measurements								
Seat width (mm)	400	400	400	400	400	400	400	400
Seat depth (mm)	400	400	400	400	400	400	400	400
Seat height (mm)	460 - 610	460 - 610	460 - 610	460 - 610	460 - 610	460 - 610		
Total width (mm)	520	520	520	520	520	520	710 - 760	710 - 760
Total depth (mm)	430	430	430	430	430	430	650	560
Total height (mm)	850 - 990	655-805	460-610	850 - 990	655-805	460-610	520	670 - 725
Width between armrests (mm)	450	450		450	450		450	450
Armrest length (mm)	280	280		280	280		280	280
Armrest height (mm)	195	195		195	195		195	195
Backrest width (mm)	500			500			500	500
Backrest height (mm)	380			380			380	380
Backrest cover height (mm)	280			280			280	280
Width hygiene recess (mm)	140	140	140				140	140
Depth hygiene recess (mm)	170	170	170				170	170
High load capacity	135 kg	130 kg	130 kg					
Weight	5,5 kg	4 kg	3,5 kg	5,5 kg	4 kg	3,5 kg	6 kg	4,5 kg
Colour								



Tel: +41 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland



Yes, you can.®