

Storm^{®4} Max



Much more than just
a larger power chair!

LINX[®]



The *Invacare Storm⁴ Max* is designed specifically for the needs of clients with a larger bodyshape. Its longer chassis supports the shifted centre of bodymass for improved weight distribution and swift mobility. Extra strong arm- and legrest support the clients needs to push themselves back to re-position or do some type of pressure relief.

The *G-Trac™* option is also available, realizing the ultimate driving performance and traction control. The stylish *Storm⁴ Max* completes the *Storm⁴* family and conserves the trouble-free servicing.

Ease of use

Just as the entire *Invacare Storm⁴* family is characterised by its ease of use, so is the *Max*. Batteries and electronics are easily accessible, making servicing simple. The depth of the seat can be simply adjusted in 3 steps which are indicated on the seat frame. An adapter to further reduce the seat depth is offered as an option.

Tilt and recline

For individual comfort and positioning the *Storm⁴ Max* offers powered tilt and recline options. The backrest recliner has a range of 0° - 25° with a high pivot point in order to reduce shearing forces.

Features and options



LiNX controls

The new *Invacare LiNX* control system has been specially designed with ease of use, simplicity and reliability in mind. It uses the latest technology to allow for accurate control and a predictable driving experience for the user.



Armrest

The strong and stylish armrest can be mounted in several positions as well as adjusted in width.



Footboard and Calf Pads

Flip-up footboard supports various positions of the feet. It's strong enough to allow the clients to push themselves back to reposition or do some type of pressure relief.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.


	Seat width	Seat depth	Seat height	Backrest height	Length of footrests	Seat angle	Backrest angle	Total length inc. footrests
<i>Storm⁴ Max</i>	490 - 610 mm ¹	560, 595, 630 mm ²	480, 510 mm	570 mm ^{3,4} , 630 mm ⁴	360 - 480 mm	0°- 15° ⁵ , 0°- 9° ⁶	0 - 25°	1390 mm
	Driving unit width	Turning diameter	Max. climbable obstacle height	Max. safe slope	Motor power	Driving range	Speed	Max. user weight
<i>Storm⁴ Max</i>	630 mm ⁷	2080 mm	60 mm	6°(10,5 %)	350 W	up to 31 km ⁸	6 km/h	220 kg

1. Additional 100 mm adjustment by each armrest
2. Adapter to reduce seat depth by 130 mm available
3. For small seat frame (490 mm) only
4. Additional 50 mm of height adjustment


5. Powered seat tilt
6. Man. settable seat tilt
7. Width of wheelbase; Overall dimension: min 670 - max. 855 mm to 975 mm including armrest width adjustment
8. According to ISO 7176-4:2008



Standard colours



Colour options



Tel: +41 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®