

Torill Side Rails



Flexibility, safety and independence

HOMELY DESIGN | AVAILABLE IN 2 LENGTHS | HEIGHT ADJUSTABLE WITH
3 DIFFERENT HEIGHTS | COMPATIBLE WITH THE INVACARE SB® 755



Available in 2 lengths

Toril Long - 980 mm

- Both rails, when placed in their highest position, offer full protection and reassurance.
- Can be used to support transfers when one rail is placed in its middle position.



Torill Short - 760 mm

- Opening in the middle is ideal for mobilisation and support.
- Reduced risk of shear as patients do not need to move further down the bed when mobilising.



Height adjustable - 3 heights

Lowest position

- Offers full access for daily nursing care.
- Unrestricted access when making and changing bed linen.
- Rounded design ensures optimal comfort when patient is sat on the edge of the bed.



Middle position

- Ideal height when side rail is used as a support to get in and out of bed.
- Provides additional protection and reassurance when the patient is sleeping or resting.



Highest position

- Full protection for safety and comfort.
- Can be used with mattress heights between 120-160 mm.



Frame colour and wood work

Aluminium Grey

Beech

Light Oak

White Beech



For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

QR
Code

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

