

Invacare® | **H263/H273 Lima**

NEW



The foldable shower commode that comes with everything

The *Invacare Lima* is the perfect choice if you are looking for a space-saving and fully equipped shower commode. All accessories are already included and the easy-to-handle folding mechanism provides a perfect storage solution for institutional and home care demands.

The rust resistant aluminum framed commodes are ideal in the shower for users who need assistance for their personal hygiene.



Yes, you can.®



Features and Options

Space-saving foldable design

See how easy a shower commode can be folded. The *Invacare H263 Lima* and *Invacare H273 Lima* are ready for storage or transport in just a few seconds. This saves space in the bath- or repository room.

Variable seating positions

The padded 4-position seat of the *Invacare Lima* provides maximum comfort for its user and can be used over every toilet. It can be easily adapted and perfect adjustment for individual patient needs is possible.



Wheel options

The choice of 5" castors or 24" self propelling wheels provides a custom fit solution.



All corrosion free

The aluminum frame and stainless steel hardware makes it completely rust resistant.



Soft backrest

The fabric backrest shows maximum comfort to the patient.

Technical data

H263	58 cm	47 cm	57 cm	104,5 cm
H273	58 cm	47 cm	68 cm	111,5 cm
H263	100 cm	21 cm	14 kg	135 kg
H273	100 cm	30 cm	18,5 kg	135 kg

Colour



White



Blue

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Lima - EU - 10/2010



Yes, you can.®

