

Invacare® | *Topaz*



Strength and Comfort

The *Invacare Topaz* wheelchair is designed specifically to suit the unique needs of bariatric patients. The *Topaz* combines the strength and durability of a triple chrome-plated, carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster bearing blocks and dual heavy-duty crossbraces, this chair provides not only ease in maneuverability but stability that will be trusted for years to come.



Yes, you can.®



Features and Options

Secure and easy to use

The dual axle positions allow quick and simple seat-to-floor height variations for the 24" rear wheels: 44,5 and 49,5 cm. The heavy-duty triple inner liners keep seat and back from stretching.


Strong and durable arm rests

The heavy-duty armrests are designed to accommodate weight capacity and provide proper positioning.



Heavy-duty front riggings

It provides proper positioning for feet with footrest ①, one piece footboard ③ and elevating legrest ②.



Heavy-duty armrests

- ① Full-Length Fixed Height
- ② Full-Length Adjustable Height
- ③ Desk-Length Fixed Height
- ④ Desk-Length Adjustable Height



Reinforced frame

It include gusseted sides, caster journals and dual heavy-duty crossbraces.

Technical data

Invacare Topaz 51 / 56 / 61 / 66 / 71 / 76 cm	45,5 / 51 cm	39 / 44,5 / 49,5 cm	40,5 / 43 / 45,5 cm	23 cm	Seat width + 24 cm

Invacare Topaz 83 - 93 cm	125 cm	37 cm	40 kg	318 kg	Bariatric use

Invacare International Sarl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Topaz - EU - 08/2010

Frame colour



Black

Upholstery



Black vinyl



Yes, you can.®