



Cetus[®]

PREMIUM PERFORMANCE | ADVANCED COMFORT |
COMPLETE INDEPENDENCE



Yes, you can.™

The *Invacare Cetus* scooter has been specifically designed to offer a durable, premium product for those who desire complete independence to travel further afield. With high performance to support those with an active, independent lifestyle, the *Cetus* can travel up to 44 km. With advanced comfort, adjustability and front and rear suspension, the *Cetus* will ensure a smooth ride across almost all terrains.

Complete independence

Longer distances

The *Cetus* is the ideal scooter for those who live life to the full. With a maximum range of up to 44 km, the *Cetus* is able to support users on those longer journeys. With speeds up to 15 kmp/h and a high-low speed switch on the control panel to effortlessly switch between pavement and road speeds, the *Cetus* is ideal for all environments.



Suitable for all users

With a 22" wide captain swivel seat, the *Cetus* is built to support all shapes and sizes. The *Cetus* can be configured easily with a built-in sliding seat mechanism, back angle adjustment and height adjustment. With a higher maximum user weight of 226 kg, the *Cetus* is truly a diverse product for a range of users.



All terrain mobility scooter

Perfect for those rougher terrains, the *Cetus* is equipped with 15" tyres and front and back suspension to provide a smoother ride for those looking for a more adventurous journey.

5 FEATURES















Available accessories





- ▶ Front basket
- ▶ USB charger adaptor
- ▶ Oxygen cylinder holder
- ▶ Mobile phone holder
- ▶ Rollator holder
- ▶ Protective storage cover
- ▶ Cane holder
- ▶ Front and rear bumpers
- ▶ Electric seat raiser

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

						
	Backrest width ▼	Seat depth ▼	Seat to floor height* ▼	Backrest height** ▼	Driving unit width ▼	Total length ▼
Cetus	559 mm	510 mm	670 mm	500 mm	720 mm	1600 mm

						
	Empty weight*** ▼	Max. user weight ▼	Battery capacity ▼	Turning diameter ▼	Max. climbable obstacle height ▼	Ground clearance ▼
Cetus	178 kg	226 kg	2 x 12 V/100Ah	1860 mm	90 mm	110 mm

					
	Max. safe slope ▼	Max. driving distance ▼	Max. driving speed ▼	Motor capacity ▼	
Cetus	9°	44 km	15 kmp/h	800W	*Seat post adjusted at lowest hole ** Without headrest height ***Including batteries

Colour of shrouds ▼



Silver

Upholstery ▼



Black

Tel: +44 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®

