



Cetus[®]

PREMIUM PERFORMANCE | ADVANCED COMFORT |
COMPLETE INDEPENDENCE



Yes, you can.

The *Invacare Cetus* scooter has been specifically designed to offer a durable, premium product for those who desire complete independence to travel further afield. With high performance to support those with an active, independent lifestyle, the *Cetus* can travel up to 44 km. With advanced comfort, adjustability and front and rear suspension, the *Cetus* will ensure a smooth ride across almost all terrains.

Complete independence

Longer distances

The *Cetus* is the ideal scooter for those who live life to the full. With a maximum range of up to 44 km, the *Cetus* is able to support users on those longer journeys. With speeds up to 15 kmp/h and a high-low speed switch on the control panel to effortlessly switch between pavement and road speeds, the *Cetus* is ideal for all environments.



Suitable for all users

With a 22" wide captain swivel seat, the *Cetus* is built to support all shapes and sizes. The *Cetus* can be configured easily with a built-in sliding seat mechanism, back angle adjustment and height adjustment. With a higher maximum user weight of 226 kg, the *Cetus* is truly a diverse product for a range of users.



All terrain mobility scooter

Perfect for those rougher terrains, the *Cetus* is equipped with 15" tyres and front and back suspension to provide a smoother ride for those looking for a more adventurous journey.

5 FEATURES















Available accessories





- ▶ Front basket
- ▶ USB charger adaptor
- ▶ Oxygen cylinder holder
- ▶ Mobile phone holder
- ▶ Rollator holder
- ▶ Protective storage cover
- ▶ Cane holder
- ▶ Front and rear bumpers
- ▶ Electric seat raiser

Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

| |  |  |  |  |  |  |
|--------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Cetus | Backrest width ▼ 559 mm | Seat depth ▼ 510 mm | Seat to floor height* ▼ 670 mm | Backrest height** ▼ 500 mm | Driving unit width ▼ 720 mm | Total length ▼ 1600 mm |

| |  |  |  |  |  |  |
|--------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Cetus | Empty weight*** ▼ 178 kg | Max. user weight ▼ 226 kg | Battery capacity ▼ 2 x 12 V/100Ah | Turning diameter ▼ 1860 mm | Max. climbable obstacle height ▼ 90 mm | Ground clearance ▼ 110 mm |

| |  |  |  |  | |
|--------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
| Cetus | Max. safe slope ▼ 9° | Max. driving distance ▼ 44 km | Max. driving speed ▼ 15 kmp/h | Motor capacity ▼ 800W | *Seat post adjusted at lowest hole ** Without headrest height ***Including batteries |

Colour of shrouds ▼



Silver

Upholstery ▼



Black

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®

