



DOLOMITE

**Melody**



Makes you wanna dance...



*Dolomite* **Melody** is a neat, stylish rollator with features that far outstrip the competitors. Low weight, even though the frame is made of steel. Simple adjustment thanks to clear, safe functions. Easy to handle as it can be folded and stands in the folded position. Robust wireless braking system. Everything with a high standard

**Melody** turns life into a dance. – Are you going to join in?



#### Wireless, integrated brakes

No brake cable to get caught up. Components in stainless steel for long life. Patented braking mechanism with tried and tested function.

#### Stands in folded position

Standing when folded makes the rollator easier to park and store on a small surface. Push the seat to get the rollator to operating position.

Easier to walk and reach things

The seat is secure and comfortable. It can be folded up to gain more space for walking and to make it easier opening doors or reaching things for example kitchen tops.



#### Infinite adjustment

Precise height and angle of the push handles can be set. In this way, the rollator can be customized for different arm lengths and hand positions.

The unique adjustment possibilities with three height models and a wide range of accessories makes *Dolomite Melody* into one of the most flexible rollators on the market.

### Features and options



#### Tray

For transportation of a cup of coffee or similar. Can be moved to a table. Art no. 70311



#### Cane holder

Makes the cane easily accessible. Art no. 12606



#### Backrest

A comfortable and safe support to lean against. Mounted on the handles and therefore follows the height adjustment.

Art no. 12628 **Melody M-L**

Art no. 12629 **Melody S**



## Features and options



### Kerb-climber

"Pedal" that simplifies passing kerbs or similar. Mounted on the right side.  
Art no. 12563



### Wheel guard

Protects clothing from being dirty from wheel splash.  
Art no. 12544 (pair)



### Slow-down brake

A safety feature preventing the walker from rolling away. Perfect for gait training and one handed solutions.  
Art no. 1452396-L (left)  
1452396-R (right)



### Basket

Positioned high on the rollator so that the user will not lose balance when handling the basket.  
Art no. 12540: **Melody M-L**  
Max load capacity 10 kg.  
Art no. 12542: **Melody S**  
Max load capacity 7 kg.



### Hanging bag

Provides extra room when the basket is quickly filled. Can be combined with oxygen cylinder holder.  
Art no. 12523



### Drip-feed holder

Holds drip feed on place. Height adjustable. No tools needed.  
Art no. 12547 (Not available for **Melody S**)



### Screw instead of knob

Replace height adjustment knob with the screw when adjustment is not desired.  
Art no. 13655 (pair)



### Tall push handles

Adjust the maximum height of push handles by 40 mm. Fits **Melody S** only.  
Art no. 13548-S (pair)



### Adapt handle set

Can be adapted to completely optional hand position. Can only be ordered as complete product.  
Art no. 12010-44-21 ADAPT  
**Melody L**  
12012-44-21 ADAPT  
**Melody M**  
Not available for **Melody S**.



### HP-handle set

Adaptation for one-handed/stroke patients. Can only be ordered as complete product.  
Art no. 12010-44-21 HPL (left):  
**Melody L**  
12010-44-21 HPR (right):  
**Melody L**  
12012-44-21 HPL (left):  
**Melody M**  
12012-44-21 HPR (right):  
**Melody M**  
Not available for **Melody S**.



### One-hand brake

Makes it possible to activate both brakes with one hand. Can only be ordered as complete product.  
Art no. 12010-44-21L (left)  
**Melody L**  
12010-44-21R (right)  
**Melody L**  
12012-44-21L (left)  
**Melody M**  
12012-44-21R (right)  
**Melody M**  
Not available for **Melody S**.



### Oxygen cylinder holder

The oxygen cylinder is centrally positioned on the rollator and thereby provides good maneuverability despite its weight. Fits most cylinders.  
Art no. 12683

## Three heights (push handles and seat)



### MELODY S

Art no. 12014-44-21



### MELODY M

Art no. 12012-44-21



### MELODY L

Art no. 12010-44-21

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please refer to your local Invacare website.

	Width between outside wheels	Depth	Width between wheels	Overall height	Seat height	Width between handles	Wheel diameter
<b>Melody S</b>	545 mm	645 mm	485 mm	660-810 mm	485 mm	370 mm	200 mm
<b>Melody M</b>	630 mm	650 mm	550 mm	710-910 mm	530 mm	420 mm	200 mm
<b>Melody L</b>	630 mm	680 mm	550 mm	790-990 mm	610 mm	420 mm	200 mm

	Load capacity	Max. user weight	Width folded
<b>Melody S</b>	7.3 kg	135 kg	200 mm
<b>Melody M</b>	7.7 kg	135 kg	260 mm
<b>Melody L</b>	7.9 kg	135 kg	320 mm

## Frame colour

RAL 1035



Please notice colors may vary slightly from those displayed above.

QR  
Code

Tel: +41 61 487 70 70  
hqeurope@invacare.com  
www.invacare.eu.com

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Melody - EU - 10/2018