# Orion<sup>®</sup> Pro

Comfort, safety and reliability with added performance



Now available with Zeta<sup>™</sup> technology







### **Even greater performance**

Ideal for longer trips into town or afternoons driving around the park, the new ergonomically designed Orion® PRO is a powerful upgrade to our standard Orion® METRO model.

With enhanced battery performance, superior suspension and large 12" wheels, the new Orion\* PRO is ideal for users who need that extra performance.

# Yes, you can:

## Advanced suspension, smoother drive

> The advanced suspension system on the Orion® PRO has been designed to bring better performance and give a smoother ride over a variety of surfaces and terrains. In addition, the Orion® PRO features a powerful motor and 12" pneumatic tyres to ensure that whether you are in the city or the country, you get where you need to go in comfort and style.

The Orion® PRO has been built using an array of quality components, which have enabled us to make it both reliable and easy to service. The electronics and motor have additional protection to guard against water and corrosion, safeguarding the vital workings of the scooter against any unnecessary malfunctions. The top cover can also be easily removed to check the motor, batteries, electronics and cables.

# Performance with security and safety

> The Orion® PRO has been subjected to German independent TÜV laboratory testing, ensuring that users can be confident that they are driving a scooter that adheres to the latest European Medical Device Standards.

In addition, an on-board electronic speed reduction system automatically reduces speed around bends for confident and comfortable manoeuvres. The two-step disengaging lever prevents the scooter from free-wheeling if it is knocked out of drive mode.

# Get connected with new Zeta™ technology

> The Orion® PRO is now available with optional Zeta™ technology, making it one of our first connected scooters. With innovation at the core, this connectivity enables remote battery checks, activity tracking and servicing notifications that can be managed via our InvaHome app.







#### **Features and options**







with four wheels



with three wheels



Finite tiller movement

Users can easily adjust the Our ergonomically tiller with a lever to suit their needs.



Ergonomic steering

designed steering system makes for easier, more responsive control.



Hand brake

Ensures immediate braking if required.



Brake light

Clear, bright brake lighting warns people that the scooter is slowing down.



Disengaging lever

Two step disengaging lever prevents the scooter from accidentally free-wheeling.



Auto-speed reduction

Automatic speed reduction reduces speed around bends for added safety.



#### Customise your scooter

Choose from 10 shroud colours and over 25 accessories for customising your scooter to your individual style. To find out more, visit our dedicated scooter website at: www.invacarescooter.com.



#### Feel confident with Zeta™ Technology

This option is a first in scooter connectivity and is packed with exciting and useful features to help you feel more confident.

#### **Technical data**

#### For more information about this product, including the user manual, please visit our website: www.invacare.eu

#### Backrest width ▼ Seat depth 🔻 Backrest height ▼ Comfort and Deluxe: Comfort: 500 mm Min. 440 mm 3 wheels 510 mm 470 mm Deluxe: 490 mm Max. 510 mm Premium: 460 mm Premium: 630 mm Comfort and Deluxe: Comfort: 500 mm Min. 440 mm 4 wheels 510 mm 470 mm Deluxe: 490 mm Max. 510 mm Premium: 460 mm Premium: 630 mm Driving unit width 🔻 Total length 🔻 The heaviest part ▼ 3 wheels 650 mm 1300 mm 126 kg 62 kg 4 wheels 1320 mm 136 kg 62 kg 650 mm Max. User Max. Climbable weight 🔻 capacity V diameter 🔻 obstacle height 🔻 2 x 12 V/75 Ah (C20) 127 kg 3 wheels 2500 mm 100 mm AGM 2 x 12 V/75 Ah (C20) 4 wheels 160 kg 2750 mm 100 mm AGM W Max. Driving Max. Driving Ground Max Safe Motor clearance \ capacity 🔻 distance 🔻 speed ▼ slope ▼ 10 km/h: 240 W / 600 W <sup>1</sup> 10 / 12.8 km/h: 550 W / 1600 W <sup>1</sup> 15 km/h: 550 W / 1600 W <sup>1</sup> 3 wheels 100 mm 10 ° 52 km 15 km/h 6 km/h 240 W / 500 W <sup>1</sup> 10 km/h: 240 W / 600 W <sup>1</sup> 12 km/h: 550 W / 1500 W <sup>1</sup> 6 / 10 / 12 / 12.8 / 10 ° 4 wheels 100 mm 52 km 15 km/h 12.8 km/h: 550 W / 1600 W <sup>1</sup> 15 km/h: 550 W / 1600 W

Shroud colours







Silver









Crystal Copper



Upholstery



Black Vinyl

1. Peak performance

**Tel: +41 61 487 70 70** hqeurope@invacare.com www.invacare.eu.com

Invacare International GmbH Benkenstrasse 260 4108 Witterswil Switzerland

