

ScanBeta[®] NG



A safe, configurable bed that changes with your child



Invacare **ScanBeta_{NG}** paediatric beds are available in two sizes and are perfect for children aged from 3 to 12 years old. The **ScanBeta_{NG}** conforms to the Nordic Requirement Specifications¹ and comes complete with a lockable handset to avoid unintentional use.



Safety Built-in

The *Invacare ScanBeta_{NG}* in a standard configuration has been rigorously tested to **Nordic Requirement Specifications¹**. Designed with safety as a priority, the bed meets all the current regulations and safety criteria.

Most Popular Configurations

The 4-section profiling mattress platforms for the two standard sizes have been constructed according to anthropometric measurements in order to maximize a child's comfort.

The **ScanBeta_{NG}** is currently available in mattress platform sizes of:

- W 700 x L 1600mm
- W 800 x L 1800mm

Ergonomic, Safe & Comfortable

The height is adjustable from 300 - 900 mm, enabling easy access for children to transfer in and out of bed. For care givers, functions for safety and comfort such as height adjustment and lockable handset are included.

Features and options



ScanBeta_{NG} can be configured with a selection of 800mm high sides*

The **ScanBeta_{NG}** equipped with 4 braking castors and Beech veneer bed ends as standard is available in a variety of configurations:

- Folding bar panels on both sides of the bed (both mattress platform sizes)
- Folding perspex panels on both sides of the bed (both mattress platform sizes)
- Folding bar panels on one side with fixed bars situated on one side (700 x 1600mm mattress platform only)

* For more details on configurations, please contact your local Invacare Business Development Manager.



The optimized lengths of the sections are based on anthropometrical measurements of children. The 2 sizes of the **ScanBeta_{NG}** standard models have been developed specifically for children aged between 3 and 12 years old.



A slim and ergonomic handset comes with an individual lockable function for increased safety.



Lifting Pole

The handle can be adjusted both in height and depth (optional).



Upholstery

20mm deep upholstery for the bed ends and all versions of sides/side rails (optional).



Safety

Folding panels have a safe locking mechanism which can only be opened by releasing two latches.



Invacare **Dacapo Mattresses****

A range of mattresses that have been specifically designed to fit the **ScanBeta_{NG}** are available.

Configuration of the *Invacare ScanBeta_{NG}* mattress support

4-section profiling

Motor on back rest and thigh section. Manual operated leg section.



Available models of the *Invacare ScanBeta_{NG}*

Product Code	Mattress Platform Size (W x L)	Sides
1575474_S1	700 x 1600 mm	2 x folding panels
1639219_S2	800 x 1800 mm	2 x folding panels
1639218_S2	700 x 1600 mm	1 x folding panel, 1 x fixed bar
1639118_S2***	800 x 1800 mm	2 x folding perspex panels
1639137_S2***	700 x 1600 mm	2 x folding perspex panels

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk.

	Width outside Width inside	Length outside Length inside	Height adjustment	Mattress support dimensions	Backrest angle	Thigh angle	Legrest lift
700 x 1600 mm	780 mm 720 mm	1680 mm 1620 mm	300 - 900 mm	670 - 150 - 220 - 520 mm	0 - 73°	0 - 31°	0 - 24°
800 x 1800 mm	880 mm 820 mm	1880 mm 1820 mm	300 - 900 mm	750 - 190 - 270 - 550 mm	0 - 73°	0 - 27°	0 - 26°

Max. user weight



2 years warranty

All sizes

SWL 100 kg
70 kg

1. Nordic Requirement Specifications - Adjustable beds for disabled children. Version 19.06.01.
<https://tinyurl.com/y8a5ujcv>

** Please note, the Invacare Dacapo Mattresses are tested to the standards EN597-1 and EN597-2.

*** Please note these configurations are subject to extended lead times of up to 6 weeks.

Frame colour and wood work

Pearl white
(metal)

Beech (wood)



Please note colors may vary slightly from those displayed above.



Visit www.safepatienthandling.co.uk

Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print.

Invacare reserves the right to alter product specifications without prior consultation. ScanBeta NG - UK - 06/2018 - 1629925