



Invacare® | **XLT**

Get Perfect Support - Remain Active!



The new *Invacare XLT* is the ideal active, rigid chair for day to day use. Designed to be perfectly adjusted to fit the body the *XLT* offers its users a wide range of various levels of support, enabling and even improving independence and mobility.

Top Configurability made Simple!

The *Invacare XLT* has been designed with the individual in mind. Measured, constructed and fitted to suit the precise requirements of each user, the *XLT* is the ideal chair to help regain, retain and improve both mobility and independence. Complex chair-specific solutions and spare parts are a thing of the past with the new *XLT*. Benefit from the extensive and readily available range of Invacare options which can be fitted to the *XLT* without any hassle.

Stability, Rigidity and Security for an Active Life

Made entirely from titanium the *XLT* gives its user stability and rigidity without compromising on weight. In combination with the wide range of options, features and configuration possibilities offered with the *XLT* the user will benefit from a chair which offers the maximum security and support needed when going about busy day to day life.



Proven Quality

We have made no compromises when it comes to the *XLT*. Designed and built at the Rea factory in Sweden, users will benefit from the strong tradition of *Invacare* quality fused to create one ultimately unique chair: The *Invacare XLT*.



Features and options



Armrests

Wide range of armrests for best support.



Anti tipper device

The Anti-tip device can easily be folded away by the foot. The device is also angled outwards, to give space for the feet for the attendant.



Standard wheel

with optional Marathon Plus tyres.



Swing Away

for easy transfer and maneuverability as well as comfort.



Front Castor Attachments

for stability and performance.



Titanium Foottube

to improve comfort as well as lightness.



Tyres

① solid, ② pneumatic profiled or puncture proof, ③ pneumatic non marking, ④ low profile, ⑤ pneumatic profiled or puncture proof.



Frontwheels

① 75 x 24 hard (3"),
② 100 x 27 medium hard (4")
③ 125 x 27 medium hard (5"),

④ 140 x 37 soft (5,5")
⑤ 150 x 27 medium hard (6"),
⑥ 150 x 30 soft (6"),

⑦ 150 x 32 pneumatic (6"),
⑧ 180 x 40 pneumatic (7"),
⑨ 180 x 40 soft (7")


Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Seat width	Seat depth	Seat height	Backrest height	Length of footrests	Seat angle	Backrest angle
<i>Invacare</i> XLT Active	355 - 480 mm	400 / 450 / 500	400 - 520 mm	200 - 490 mm	380 - 505 mm	0 - 14°	± 10°
<i>Invacare</i> XLT Dynamic	355 - 500 mm	400 / 450 / 500	400 - 500 mm	200 - 490 mm	355 - 485 mm	0 - 14°	± 10°
<i>Invacare</i> XLT Swing	330 - 505 mm	400 / 450 / 500	390 - 520 mm	200 - 490 mm	280 - 510 mm	0 - 14°	± 10°
	in increments of 15 mm	on 40 mm adj.					

	Driving unit width	Total height	Total length inc. Footrests	Total weight	Max. User weight	Mass of the heaviest part	Crash-tested
<i>Invacare</i> XLT Active	SW + 200 mm	610 - 1050 mm	830 - 1030 mm	9.8 kg	135 kg	6 kg	ISO 7176-19
<i>Invacare</i> XLT Dynamic	SW + 200 mm	610 - 1050 mm	730 - 880 mm	9.8 kg	135 kg	6 kg	
<i>Invacare</i> XLT Swing	SW + 200 mm	610 - 1050 mm	820 - 1230 mm	11.3 kg	135 kg	6 kg	


Frame colours



Please notice colors may vary slightly from those displayed above.

* Available as special orders

Upholstery



Scan here to access the website.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
 Benkenstrasse 260
 4108 Witterswil
 Switzerland

