

Alber | **viaplus**



Pushing comfort chairs the easy way

Comfort wheelchairs offer users great comfort but sometimes pushing these wheelchairs can be a strenuous job for the attendant. **Alber viaplus**, the push and brake aid for Invacare comfort wheelchairs, offers attendants the assistance they need when pushing this type of wheelchair. What's more, the automatic braking function eliminates downhill worries and enables the safe negotiation of slopes and ramps.





Fits perfectly to Rea™ Azalea® and Clematis®

Unlike many other pushing aids, the **viaplus** does not impede the adjustability of a tilt in space wheelchair (with a seat height of 45 cm). All of the wheelchair's functions remain accessible. That's why the **viaplus** can simply be left on the wheelchair after a walk without restricting the use of its functions.

Easy to operate

Controlling the **viaplus** is as easy as operating a single, handy rotary switch. This switches the drive on or off, sets the required speed or activates the freewheel mode. **viaplus** is therefore exceptionally easy to operate for older users.

Personal freedom

viaplus not only benefits the attendant: It also gives the wheelchair user the opportunity for a fuller social life because frequent excursions with the wheelchair are no longer a problem now that handling the wheelchair is so much easier for the attendant. **viaplus** is also used in residential homes. With it, pushing wheelchairs becomes an effortless task

and lightens the load for the staff. The handy controls of the **Alber viaplus** fits onto the handle of the wheelchair. The wheelchair is pushed forwards by light thumb pressure on the operating lever, making fatigue-free operation possible even on longer journeys. The chair is put into reverse by pulling the operating lever upwards.

Features and Options



Drive unit

With puncture proof solid rubber tyres, automatic spring system for constant traction, pedal for lowering/lifting the wheels.



Control unit

Central, handy rotary switch for simple control of all functions.



Battery pack

Powerful battery pack with charging socket.

Features and Options



Battery charger

Automatic charger with indicator light for charging status.



Bracket set for wheelchair

Available for Rea™ **Clematis**® and **Azalea**®.



Battery box


Solid battery box for carrying the battery pack.



Versatility

No limitation of seat/backrest adjustability (seat height = 45 cm).

Technical data

							
viaplus	max. 5.5 km/h	16 km	24V/12Ah	9 kg	10%	20 kg (drive unit: 10,7 kg)	Azalea* : 135 kg, Clematis* : 125 kg

Ulrich Alber GmbH
 Vor dem Weißen Stein 21
 72461 Albstadt · Germany
 Tel.: +49 7432 2006 0
 Fax: +49 7432 2006 299
 Email: info@alber.de
 www.alber.de

© 2010 Invacare® International Sarl
 All rights reserved.

All information quoted is believed to be correct at time of print.
 Invacare® reserves the right to alter product specifications without prior consultation.

Alber viaplus- EU - 10/2010



Yes, you can.®