


Invacare® Dress Toileting Slings



- en** Addendum **da** Tilføjelse **de** Addendum **es** Apéndice **fi** Lisäys
- fr** Addendum **it** Addendum **nl** Aanvulling **no** Tillegg **pt** Adenda
- sv** Tillägg

1 Figures



3.2 Options / Ekstraudstyr / Optionen / Opciones / Lisävarusteet / Options / Opzioni / Opties / Ekstrautstyr / Opções / Tillval

- Extension loops / Forlængerløkker / Verlängerungsschlaufen / Anillas de extensión / Jatkesilmukat / Boucles d'extension / Occhielli di estensione / Verlenglussen / Forlængelsesløkker / Presilhas de extensão / Förlängningsöglor
- Head support/ Hovedstøtte / Kopfstütze / Reposacabezas / Päätuki / Appuie-tête / Poggiatesta / Hoofdsteun / Hodestøtte / Apoio de cabeça / Huvudstöd
- Loop marker/ Loop-markør / Schlaufenmarkierung / Marcador de anillas / Silmukkamerkki / Marqueur de boucle / Indicatore di occhielli / Lusmarkering / Løkkemarkør / Marcador de presilha / Ögelmarkör

7 Technical Data / Tekniske Data / Technische Daten / Datos Técnicos / Tekniset Tiedot / Caractéristiques Techniques / Dati Tecnici / Technische Specificaties / Tekniske Data / Características Técnicas / Teknisk Data

Dimensions / Mål / Abmessung / Dimensiones / Mitat / Dimensioni / Afmetingen / Dimensjoner / Dimensões / Mått

Variant / Ausführung / Variante / Versio	[mm]	XS	S	M	L	XL	
Dress Toileting Low	A	840	870	950	1020	1100	
	B	620	710	800	900	1080	
	C	160	200	230	230	270	
	D	380	430	495	600	760	
Dress Toileting High	A	840	870	950	1020	1100	
	B	620	710	800	900	1080	
	C	160	200	230	230	270	
	D	380	430	495	600	760	
	E	400	400	465	465	465	

60147940-A 2024-09-25 1/1



Invacare UK Operations Limited
Unit 4, Pencoed Technology Park,
Pencoed
Bridgend CF35 5AQ
UK



Invacare Portugal, Lda
Rua Estrada Velha 949
4465-784 Leça do Balio
Portugal



Yes, you can.®