

Medley Ergo

Providing clients comfort and support



An advanced range of care beds

for increased comfort and improved clinical outcomes



Yes, you can.®

Designed to care for the client and caregiver

Our Medley Ergo range of beds are an easy to handle solution that help deliver excellent client care and provide maximum comfort.

Designed to increase independence and improve clinical outcomes for varying needs, with a choice between our Medley Ergo or Medley Ergo Low and the availability of Medley Ergo Select for a more homely design with bed end covers.

Built-in features and functionalities ensure the client as well as the caregiver can operate the bed easily. Optimised dimensions of the mattress platform help to reduce shear and friction and in combination with one of our range of pressure redistributing mattresses, can provide optimal pressure care and micro climate management.



1



Easy Handling

An optional **transportation fitting** is available so that the bed can easily be **collapsed, moved and re-installed** in a variety of locations.

2



Comfort


In all resting and sitting positions with **minimised risk of slipping downwards** during profiling with optimised mattress base dimensions.

3



Support

Easy repositioning of the client with hand controls that help to **increase independence** and **support the care giver**. Available with or without lockable functions.



› **Mattress support**
Available in easy to clean steel.

› **Side rails**
Choose from 3 options that vary in length, height and finish.
Full length - Lisa wood finish in beech or Aria in steel
3/4 length - Scala, collapsible in steel

› **Select design**
A more homely design can be achieved with retro-fitted bed end covers for our Medley Ergo Select option.

› **Hand control**
With or without lockable functions.



› **Double height positions for bed ends**
Adjustable height ranges for
Medley Ergo: 33 - 73 cm or 40 - 80 cm
Medley Ergo Low: 21 - 61 cm or 28 - 68 cm

Medley Ergo Low


- › **Positioning the bed at a lower level** increases security for clients and provides peace of mind when going to sleep. Caregivers can raise the bed to a working height when required.
- › Clients can get **in and out of bed easier** increasing independence.
- › Adjustable heights from **21 - 61 cm** and **28 - 68 cm**.

Options to choose from:

- Light and easily adjustable lifting Pole
- 15 cm Extension kit for increasing mattress support length by 15 cm
- Transport fitting for easy storage and transport
- Longer side rail to fit extended bed


Helping to deliver optimal client care

A combination of advanced functionalities, sleek features and uncompromised quality for enhanced client care.



Shear and Friction Reduction

Optimised dimensions of the mattress base help reduce shear and friction and increase comfort, according to expert recommendations for anthropometrical measurements and pressure mapping.



Easy to store and transport

Tool-less assembly / disassembly in combination with our optional transport fitting enables easy transportation and storage when the bed is collapsed and placed on the dedicated transporter.



Safety and quality assured

Safety standard IEC 60601-2-52 compliant, our trusted Medley Ergo range can be relied upon for its high safety standards and optimal quality, helping to deliver reliable long-term support for clients at home or in long-term care environments.



Technical data ▼


For more information about this product, including the user manual, please visit our website: www.invacare.eu




Width outside
Width inside ▼



Length outside
Length inside ▼




Height adjustment ▼




Mattress support
dimensions ▼


	Width outside Width inside ▼	Length outside Length inside ▼	Height adjustment ▼	Mattress support dimensions ▼
Medley Ergo	1020 mm 900 mm	2240 mm 2090 mm	400 - 800 mm or 330 - 730 mm	800 - 230 - 300 - 670 mm
Medley Ergo Low	1020 mm 900 mm	2240 mm 2090 mm	280 - 680 mm or 210 - 610 mm	800 - 230 - 300 - 670 mm




Backrest angle ▼



Thigh angle ▼



Legrest lift ▼



Tilt option¹ ▼

	Backrest angle ▼	Thigh angle ▼	Legrest lift ▼	Tilt option ¹ ▼
Medley Ergo	0 - 70°	0 - 24°	0 - 18°	11°
Medley Ergo Low	0 - 70°	0 - 24°	0 - 18°	11°



Total weight ▼



Max. user weight ▼



Heaviest part of the
product ▼

	Total weight ▼	Max. user weight ▼	Heaviest part of the product ▼
Medley Ergo	64 kg	SWL 180 kg 145 kg	17 kg
Medley Ergo Low	64 kg	SWL 180 kg 145 kg	17 kg

Complies to the safety standard, IEC 60601-2-52

1. Please notice that the tilt function is only practised for users under medical supervision.

NB: The bed must not be used by clients under 12 years of age, or by clients with body size equivalent to an average 12 year old or smaller.

Silver

Beech

Please note colours may vary slightly from those displayed above.



Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

ESG ENVIRONMENTAL SOCIAL & GOVERNANCE



Yes, you can.®