

NordBed™ Accessories



To fit all needs!



Yes, you can.

The Invacare NordBed range features a wide range of accessories which gives you the possibility to customise the bed both in terms of design and taste. This allows perfect fittings to personal environment with supporting accessories which tailors the bed to individual client needs.

In the catalogue you will also find compatibility references for beds, side rails and bed ends.

Bed ends

Modern bed end range, homely / non-institutional atmosphere, fitting individual taste and preferences. Available in a variety of wood finishes and some colours. All bed ends are easy removable without the use of tools.




Una

Can be used in combination with full length side rails, split side rails or side panels



Una / Una Low

Can be used in combination with side panels



Nora

Can be used in combination with full length side rails, split side rails or side panels



Astrid

Removable middle plate for easy access. Can be used in combination with full length side rails, split side rails or side panels



Nanna

Can be used in combination with side panels and the Ane side rails



Runa

Can be used in combination with Ane side rail range and a side panels.

Bed end range	Beech	Oak	Umber	Beech/pale blue	Pale blue	White
Una bed end	✓	✓	✓	✓	✓	✓
Nora bed end	✓	✓	✓		✓	✓
Astrid bed end	✓	✓	✓		✓	✓
Nanna bed end	✓	✓	✓		✓	✓
Runa bed end	✓	✓	✓		✓	✓

Side panels

Our range of side panels promote the homely aesthetic atmosphere



Wooden side panel

Compatible with Una, Astrid and Nora bed ends



Wooden side panel

Compatible with Runa bed end and Ane side rails



Wooden side panel

Compatible with Nanna bed ends and Ane side rails

Side rails

Side rails can be used for multiple reasons: for safety to ensure the client is not falling out of bed, for repositioning purposes in the bed, as a support for transfer in and out of bed etc.

The Invacare range of side rails fit any of these needs without jeopardizing the healthy working condition for the care staff



Freya

Full length side rails in wood 200 cm, 210 cm and 220 cm



Embla

Full length aluminium side rails 200 cm



Embla extendable

Full length side rails in aluminium. Stepless adjustable length 200 - 220 cm



Torill Short

Adjustable in 3 heights



Torill Long

Adjustable in 3 heights



Ane

3/4 length steel side rails

Side rails



Ane

3/4 length steel side rails for high mattresses



Freya side rail height extender

The Freya side rails can be extended with + 6 cm or + 10 cm for safe use of high mattresses



Embla side rail height extender

The Embla side rails can be extended with + 6 cm or + 10 cm and extended in length 200 - 220 cm

Side rail range	Una	Nora	Astrid	Nanna	Runa
Freya	✓	✓	✓		
Embla	✓	✓	✓		
Torill	✓*	✓*	✓*		
Ane				✓	✓

* Bed ends without guides for full length side rail

Side rail covers

Our side rail cover range serves several purposes: the padded versions are used to prevent clients hurting themselves on the sides of the bed, the netting covers prevent the client putting their arms or legs out through the openings in the side rails and yet still allow clear visibility.

All covers can remain on the side rails when the side rails are operated.

All covers can be easily removed and washed at 90 degrees.



Padded side rail covers

To be used with Freya & Embla side rails



Padded side rail covers

To be used with Freya & Embla side rails with height extender



Net cover

To be used with Freya and Embla side rails



Padded side rail cover

To be used with Ane side rails

Repositioning aids

Our range of stable support handles and mobilization accessories promote client independence, by supporting the clients to reposition themselves in the bed and to get safely in-and out of bed.



Support Handle

H40 x W30 cm



Support Handle

H40 x W90 cm



Swivel support handle

The swivel support handles make it possible to swing the legs out over the handles before the grip is folded out. Can be used alone or in pairs



Stand-up support

Makes getting in and out of bed easier



Sit-up support

Supports the client to get to true sitting position and can support getting in and out of bed



Lifting pole

Catching handle can be adjusted in both height and depth

Other accessories



Torill side rail extension kit



Cover tubes

To be mounted on the bed frame to cover the tubes when the bed is extended




Handset holder



Distance bumper

Prevents walls from being damaged when the bed is moved. Can be placed on the base frame at the end of the bed or at the side of the bed



Extension kit +10 cm

Available on the NordBed Essential

Electrical accessories

A range of handsets are available for the NordBed range to accommodate different needs



Handset for Optimo and Ultra

With anti-trendelenburg / with and without trendelenburg



Handset for Optimo and Ultra

With anti-trendelenburg / with and without trendelenburg and under bed light



LCD Handset

With RememberMe™ functionality (Optimo and Ultra)



Handset for Essential

The standard handset



Handset for Essential

Available with anti-trendelenburg and optional trendelenburg



LED lamp

Attractive, modern design, very low power consumption (4W). Can be mounted on the lifting pole bracket on all Invacare beds. Comes in aluminium.

Optimo & Ultra																															
Mattress height cm	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30									
Freya (200 cm, 210 cm and 220 cm)																															
Embla/Embla Extended																															
Freya/Embla side rail with height extender + 6 cm																															
Freya/Embla side rail with height extender + 10 cm																															
Ane																															
Ane High																															
Torill Long/Torill Short																															

Essential																															
Mattress height cm	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30									
Freya (200 cm, 210 cm and 220 cm)																															
Embla/Embla Extended																															
Freya/Embla side rail with height extender + 6 cm																															
Freya/Embla side rail with height extender + 10 cm																															
Ane																															
Ane High																															
Torill Long/Torill Short																															

Tel: +41 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®



© 2023 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. EU - 04/23