



Yes, you can.®

EUROPEAN PRESCRIPTION FORM

ACTION³ JUNIOR-117-161212 Sweden

Invacare® Action³ JUNIOR



Local Homologation Number :
ISO , MDA , TÜV Logos
Crossbar Folding Mechanism, Seat Angle Possibility : 1,5°, 5°, 8,5° (Delivered Always 5°).
Rear Wheel Position : Standard, Mediuml, Passive (Delivered Always In Standard or Passive).
80°hanger delivered with 0° abduction (can be easily adjusted to 30°)
Heel Straps and Tool Kit. Maximum Weight : 80 Kg.
Price In Local Currency & Euros (If Required) of The Standard Product.

INVACARE France OPERATIONS SAS
Route de Saint Roch, 37230 Fondettes, France
Tel: 33 (0)247 626 466 Fax: 33 (0)247 626 488

Codification of Local Prescription Form:

Action3 Junior = CNB
Action3 Junior Evolutive = CNC

Customer Name :
Dealer :
Account No :
Address :
Tel No :
Delivery address :
Order No :
Order Qty :
Order Date :

OR DEALER STAMP



UPDATE : In Green: New modification ; In Blue: Previous Modification.

☐ Please make your selection with a cross in the box : ◆ Standard feature ○ Free of Charge (NC9999) Non Conformance with Option 9999.

➔		Customer weight =	Kg	
Code		Description		
MODEL				
CNA0000	☐	ACTION ³ JUNIOR	(NC: 0140, 0141, 0142, 0143, 0150, 0151, 0152, 0153, 0160, 0161, 0162, 0170, 0171, 0335, 0402, 0515, 0565)	
CNA0001	◆	ACTION ³ JUNIOR EVOLUTIVE	(NC: 0180, 0190, 0200, 0240, 0250, 0251, 0252, 0330, 0340, 1470)	
SPECIAL				
N	◆	NO SPECIAL	(NC: 0991)	
CNASPEC	☐	SPECIAL		
BASE MODEL FOR PRICE LIST				
CNA0025	◆	ACTION ³ JUNIOR		
CNA0030	☐	ACTION ³ JUNIOR FIXED BACKREST (ONLY FOR FRANCE)		
CNA0035	☐	ACTION ³ JUNIOR FOLDING BACKREST (ONLY FOR FRANCE)		
CNA0040	☐	ACTION ³ JUNIOR RECLINING BACKREST (ONLY FOR FRANCE)		
PROPULSION TYPE				
CNA0100	◆	SELF PROPULSION	(NC : 1100, 1110, 1170)	
CNA0110	☐	TRANSIT PROPULSION (SEAT ANGLE 1,5° ONLY) WITH RH TIPPING AID	(NC : 0280, 1060, 1070, 1075, 1080, 1090, 1180, 1190, 1195, 1275, 1375, 1390, 1391, 1392, 1430, 1440)	
SEAT				
STANDARD SEAT VERSION				
SEAT WIDTH				
CNA0180	☐	305 MM SEAT WIDTH	(NC: 0250, 0252, 0520, 0570, 0890, 1190, 1195)	
CNA0190	☐	330 MM SEAT WIDTH	(NC: 0240, 0251, 0530, 0580, 0900)	
CNA0200	☐	355 MM SEAT WIDTH	(NC: 0240, 0251, 0530, 0580, 0900)	
SEAT DEPTH				
CNA0240	☐	300 MM SEAT DEPTH		
CNA0250	☐	350 MM SEAT DEPTH		
CNA0251	☐	ADJUSTABLE SEAT DEPTH 250 TO 300 MM	(NC: 0320, 0340)	
CNA0252	☐	ADJUSTABLE SEAT DEPTH 300 TO 350 MM	(NC: 0320, 0340)	
EVOLUTIVE SEAT VERSION				
SEAT WIDTH				
CNA0140	☐	205 MM SEAT WIDTH INCLUDING STABILISER KIT (0881)	(NC: 0170, 0171, 0480, 0511, 0513, 0520, 0530, 0561, 0563, 0570, 0580, 0881, 0890, 0900, 1375, 1470)	
CNA0141	☐	230 MM SEAT WIDTH INCLUDING STABILISER KIT (0881)	(NC: 0170, 0171, 0480, 0511, 0513, 0520, 0530, 0561, 0563, 0570, 0580, 0881, 0890, 0900, 1375, 1470)	
CNA0142	☐	255 MM SEAT WIDTH INCLUDING STABILISER KIT (0881)	(NC: 0170, 0171, 0480, 0511, 0513, 0520, 0530, 0561, 0563, 0570, 0580, 0881, 0890, 0900, 1375, 1470)	
CNA0143	☐	280 MM SEAT WIDTH INCLUDING STABILISER KIT (0881)	(NC: 0170, 0171, 0480, 0511, 0513, 0520, 0530, 0561, 0563, 0570, 0580, 0881, 0890, 0900, 1375, 1470)	
CNA0150	◆	305 MM SEAT WIDTH	(NC: 0160, 0161, 0162, 0520, 0570, 0890)	
CNA0151	☐	330 MM SEAT WIDTH	(NC: 0160, 0161, 0162, 0530, 0580, 0900)	
CNA0152	☐	355 MM SEAT WIDTH	(NC: 0160, 0161, 0162, 0530, 0580, 0900)	
CNA0153	☐	380 MM SEAT WIDTH	(NC: 0160, 0161, 0162, 0530, 0580, 0900)	
SEAT DEPTH				
CNA0160	☐	ADJUSTABLE SEAT DEPTH 200 TO 250 MM DELIVERED AT 230 MM	(NC: 1060, 1070, 1075)	
CNA0161	◆	ADJUSTABLE SEAT DEPTH 250 TO 300 MM DELIVERED AT 280 MM		
CNA0162	☐	ADJUSTABLE SEAT DEPTH 300 TO 350 MM DELIVERED AT 330 MM		
CNA0170	☐	ADJUSTABLE SEAT DEPTH 275 TO 325 MM DELIVERED AT 300 MM		
CNA0171	☐	ADJUSTABLE SEAT DEPTH 325 TO 375 MM DELIVERED AT 350 MM		
OTHER SEAT INFORMATIONS				
SEAT HEIGHT				
CNA0260	◆	460 MM SEAT HEIGHT		
CNA0270	☐	435 MM SEAT HEIGHT		
CNA0280	☐	410 MM SEAT HEIGHT	(NC: 1375, 1430, 1440)	
SEAT UPHOLSTERY TYPE				
CNA0320	☐	WITHOUT SEAT UPHOLSTERY		
CNA0330	☐	NYLON WITH VELCRO STRIPS SEAT UPHOLSTERY		
CNA0335	◆	NYLON WITH VELCRO STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE		
CNA0340	☐	PADDED NYLON SEAT UPHOLSTERY		
BACKREST				
BACKREST TYPE				
CNA0370	☐	FIXED 0° BACKREST	3 HEIGHTS: 370/400/430 MM - DELIVERED AT 370 MM (SW 205 TO 280 MM) AT 400 MM (SW 305 TO 380 MM)	
CNA0380	◆	FIXED 7° BACKREST	3 HEIGHTS: 370/400/430 MM - DELIVERED AT 370 MM (SW 205 TO 280 MM) AT 400 MM (SW 305 TO 380 MM)	
CNA0381	◆	FIXED 7° BACKREST 510 MM	3 HEIGHTS: 370/400/430 MM - DELIVERED AT 370 MM (SW 205 TO 280 MM) AT 400 MM (SW 305 TO 380 MM)	
CNA0390	☐	RECLINING (RATCHET) 0/30° BACKREST	5 HEIGHTS: 350/370/390/410/430 MM - DELIVERED AT 370 MM (SW 205 TO 280 MM) AT 390 MM	(NC: 1190, 1195)
CNA0395	☐	RECLINING (GAS SPRING) 0/30° BACKREST	400 MM HEIGHT W/O BACKBRACE	(NC: 1190, 1195)
CNA0400	☐	FOLDING BACKREST	3 HEIGHTS: 370/400/430 MM- DELIVERED AT 370 MM (SW 205 TO 280 MM) AT 400 MM (SW 305 TO 380 MM)	
CNA0401	☐	FOLDING BACKREST 510 MM	2 HEIGHTS: 400/430 MM- DELIVERED AT 400 MM	
CNA0402	☐	FOLDING BACKREST	2 HEIGHTS: 380/500 MM- DELIVERED AT 500 MM; HEIGHT ADJUSTABLE PUSHHANDLE REAR SET (150 MM); INCLUDING BACK BRACE	(NC : 0480, 0855, 0902, 1426)
CNA0420	☐	FOLDING BACKREST	400 MM HEIGHT; HEIGHT ADJUSTABLE PUSHHANDLE (100 MM)	(NC : 0855)
CNA0425	☐	FIXED 7° BACKREST	400 MM (with 0475) 430 MM (with 0480) HEIGHT; HEIGHT ADJUSTABLE PUSHHANDLE (150 MM)	(NC : 0855)
BACKREST UPHOLSTERY TYPE				
CNA0465	☐	WITHOUT BACK UPHOLSTERY		(NC: 0855)
CNA0475	☐	WIDTH & HEIGHT ADJUSTABLE NYLON SLING TYPE		
CNA0480	◆	CONTOUR (TENSION) ADJUSTABLE NYLON TYPE		(NC: 0855)
HANGERS/FOOTPLATES				
HANGERS RH TYPE (HEIGHT ADJ. BY STEP OF 10 MM)				
CNA0500	◆	WITHOUT RH HANGER	(NC : 0600, 0605, 0610, 0611, 0880, 0882, 0883, 0884, 0885)	
CNA0505	☐	90° STRAIGHT RH REMOVABLE HANGER	UL STD 125 TO 260 MM	(NC: 0565, 0600, 0610, 1005)
CNA0510	☐	80° SHORT RH SWING IN/OUT HANGER	0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC: 0561, 0562, 0563, 0565, 0610)
CNA0511	☐	80° SHORT RH SWING IN/OUT HANGER	30° ABDUCTION POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC: 0560, 0562, 0563, 0565, 0610, 0611, 0881)
CNA0512	☐	80° LONG RH SWING IN/OUT HANGER	0° STANDARD POSITION; UL STD 300 TO 450 MM	(NC: 0560, 0561, 0563, 0565, 0610)
CNA0513	☐	80° LONG RH SWING IN/OUT HANGER	30° ABDUCTION POSITION; UL STD 300 TO 450 MM	(NC: 0560, 0561, 0562, 0565, 0610, 0611, 0881)
CNA0520	☐	ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL STD 275 TO 375 MM		(NC: 0565, 0600, 0611, 0880, 0886)
CNA0530	☐	ELEVATING LEGREST (0° TO 80°) RH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM		(NC : 0565, 0600, 0611, 0880, 0886)
HANGERS LH TYPE (HEIGHT ADJ. BY STEP OF 10 MM)				
CNA0550	◆	WITHOUT LH HANGER	(NC : 0611, 0620, 0625, 0630, 0880, 0882, 0883, 0884, 0885)	
CNA0555	☐	90° STRAIGHT LH REMOVABLE HANGER	UL STD 125 TO 260 MM	(NC: 0620, 0630, 1005)
CNA0560	☐	80° SHORT LH SWING IN/OUT HANGER	0° STANDARD POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC: 0630)
CNA0561	☐	80° SHORT LH SWING IN/OUT HANGER	30° ABDUCTION POSITION; UL STD 220 TO 320 MM BLACK COLOUR ONLY	(NC: 0611, 0630, 0881)
CNA0562	☐	80° LONG LH SWING IN/OUT HANGER	0° STANDARD POSITION; UL STD 300 TO 450 MM	(NC: 0630)
CNA0563	☐	80° LONG LH SWING IN/OUT HANGER	30° ABDUCTION POSITION; UL STD 300 TO 450 MM	(NC: 0611, 0630, 0881)
CNA0570	☐	ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (INCLUDING CALFPAD); UL ST 275 TO 375 MM		(NC : 0611, 0620, 0880, 0886)
CNA0580	☐	ELEVATING LEGREST (0° TO 80°) LH SWING IN/SWING OUT HANGER (WITHOUT CALFPAD, WITH CALFSTRAP); UL STD 235 TO 295 MM		(NC : 0611, 0620, 0880, 0886)
FOOTPLATES RH TYPE				
CNA0600	☐	FLIP UP RH FOOTPLATE	(NC: 0882, 0883, 0884, 0885, 0886)	
CNA0605	◆	FLIP UP, ANGLE & DEPTH ADJUSTABLE RH FOOTPLATE		
CNA0610	☐	FLIP UP AND ANGLE ADJUSTABLE RH FOOTPLATE	(NC: 0882, 0883, 0884, 0885, 0886)	
CNA0611	☐	ONE PIECE FOOTBOARD	(NC : 0620, 0625, 0630, 0882, 0883, 0884, 0885)	
FOOTPLATES LH TYPE				
CNA0620	☐	FLIP UP LH FOOTPLATE	(NC: 0882, 0883, 0884, 0885, 0886)	
CNA0625	◆	FLIP UP, ANGLE & DEPTH ADJUSTABLE LH FOOTPLATE		
CNA0630	☐	FLIP UP AND ANGLE ADJUSTABLE LH FOOTPLATE	(NC: 0882, 0883, 0884, 0885, 0886)	



Important Information:
Stabiliser kit (see seating option type; ref. 0881) is a standard feature for a better front and lateral stability for these seat widths.



September 2016 (1/2)

ACTION³ JUNIOR-117-161212 Sweden

UPDATE : In Green: New modification ; In Blue: Previous Modification.
Please make your selection with a cross in the box : Standard feature Free of Charge (NC9999) Non Conformance with Option 9999.

		Customer weight =	Kg		
		Code	Description		
ARMRESTS					
ARMRESTS RH FORM					
CNA0640	◆	REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED	DESK, STD ARMPAD (60 MM)	ADJ. 180 TO 280 MM	(NC : 0750, 0745, 0746, 0755)
CNA0645	□	REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED	JUNIOR, STD ARMPAD (60 MM)	ADJ. 150 TO 250 MM	(NC : 0740, 0746, 0750, 0755)
CNA0646	□	REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED	NARROW ARMPAD JUNIOR (50 MM)	ADJ. 150 TO 250 MM	(NC : 0740, 0745, 0750, 0755, 1470)
CNA0650	□	MUDGUARD RH		HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC: 0740, 0745, 0746, 0755, 1201, 1470, 1471, 1472)
CNA0655	□	SKIRTGUARD RH	JUNIOR		(NC: 0740, 0745, 0746, 0750, 1470, 1471, 1472)
ARMRESTS LH FORM					
CNA0740	◆	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED	DESK, STD ARMPAD (60 MM)	ADJ. 180 TO 280 MM	
CNA0745	□	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED	JUNIOR, STD ARMPAD (60 MM)	ADJ. 150 TO 250 MM	
CNA0746	□	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED	NARROW ARMPAD JUNIOR (50 MM)	ADJ. 150 TO 250 MM	(NC: 1470)
CNA0750	□	MUDGUARD LH		HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC: 1201, 1470, 1471, 1472)
CNA0755	□	SKIRTGUARD LH	JUNIOR		(NC: 1470, 1471, 1472)
SEATING OPTIONS TYPE					
CNA0820	◆	NO SEATING OPTION			(NC : 0825, 0851, 0855, 0860, 0865, 0870, 0880, 0881, 0882, 0883, 0884, 0885, 0886, 0890, 0900, 0902)
CNA0851	□	SEAT CUSHION MSS (ESSENTIAL JUNIOR)			
CNA0855	□	BACKREST EXTENSION (230 MM)			(NC 0860, 0865)
CNA0860	□	HEADREST STANDARD MULTI - ADJUSTABLE	INCLUDING WIDTH ADJUSTABLE BACKBRACE		(NC: 0865, 0902)
CNA0865	□	HEADREST ROUND MULTI - ADJUSTABLE	INCLUDING WIDTH ADJUSTABLE BACKBRACE		(NC: 0902)
CNA0870	◆	ADJUSTABLE HEELS STRAP			
CNA0880	□	ADJUSTABLE CALF STRAP			
CNA0881	□	FRONT STABILISER KIT	+ 100 MM OUTWARD		
CNA0882	□	ANKLE GRIPS BRACKET, PAIR			(NC: 0883, 0884, 0885)
CNA0883	□	ANKLE GRIPS "XS", PAIR	INCLUDING BRACKET	STRAP'S CIRCUMFERENCE 155 TO 200 M	(NC: 0884, 0885)
CNA0884	□	ANKLE GRIPS "S", PAIR	INCLUDING BRACKET	STRAP'S CIRCUMFERENCE 160 TO 215 M	(NC: 0885)
CNA0885	□	ANKLE GRIPS "M", PAIR	INCLUDING BRACKET	STRAP'S CIRCUMFERENCE 175 TO 230 MM	
CNA0886	□	FOOTPLATE ADJUSTMENT OPTIMIZER KIT	(FRONT OR REAR SET ADJUSTMENT)		
CNA0890	□	ADDITIONAL ELEVATING LEGREST (0° TO 80°) SWING IN/OUT HANGER PAIR (INCL. ADJ FOOTPLATES); UL STD 275 TO 375 MM			(NC: 0900)
CNA0900	□	ADDITIONAL ELEVATING LEGREST (0° TO 80°) SWING IN/OUT HANGER PAIR (INCL. ADJ FOOTPLATES, WITHOUT CALFPADS, WITH CALFSTRAPS); UL STD 235 TO 295 MM			
CNA0902	□	WIDTH ADJUSTABLE BACK BRACE			
COLOURS					
CNA0994	■	VIBRATING ORANGE	FOR UK ONLY; UNTIL END OF APRIL		
CNA0952	◆	PEARL WHITE			
CNA0953	◆	MATT BLACK			
CNA0954	□	MATT RED			
CNA0955	□	INTENSE YELLOW			
CNA0956	□	POLAR BLUE			
CNA0957	□	CANDY PINK			
CNA0960	■	DARK RED	FOR UK ONLY; UNTIL END OF APRIL		
CNA0970	■	SAPPHIRE BLUE	FOR UK ONLY; UNTIL END OF APRIL		
CNA0980	■	ANTHRACITE BLACK	FOR UK ONLY; UNTIL END OF APRIL		
CNA0990	■	BUBBLE GUM PINK	FOR UK ONLY; UNTIL END OF APRIL		
CNA0991	□	SPECIAL COLOUR (ON DEMAND ONLY)			
FRONT WHEELS					
CNA1000	◆	125 X 25 MM (5" x 1")	SOLID TYRE		
CNA1001	□	125 X 25 MM (5" x 1") FLASHING LIGHTS	SOLID TYRE		
CNA1005	□	125 X 45 MM (5"x 1" 3/4)	SOLID TYRE		
REAR WHEELS					
REAR WHEEL TYPE					
CNA1060	□	560 X 35 MM (22" x 1.3/8")	PNEUMATIC TYRE	SPOKED RIM	(NC: 1170)
CNA1070	◆	555 X 35 MM (22" x 1.3/8")	SOLID TYRE	SPOKED RIM	(NC: 1170)
CNA1075	□	550 X 25 MM (22" x 1")	BLACK SOLID TYRE	SPOKED RIM	(NC: 1170)
CNA1090	□	515 X 35 MM (20" x 1.3/8")	SOLID TYRE	SPOKED RIM	(NC: 1170, 1280)
CNA1100	□	305 X 50 MM (12" x 2")	PNEUMATIC TYRE	PLASTIC RIM	(NC : 1140, 1141, 1180, 1250, 1251, 1252, 1253, 1254, 1255, 1256, 1380, 1390)
CNA1110	□	305 X 45 MM (12" x 2")	SOLID TYRE	PLASTIC RIM	(NC : 1140, 1141, 1180, 1250, 1251, 1252, 1253, 1254, 1255, 1256, 1380, 1390)
HANDRIM TYPE					
CNA1140	◆	ALUMINIUM PRE-ANODIZED HANDRIM			
CNA1141	□	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC			
REAR WHEELS AXLE TYPE					
CNA1170	□	FIXED AXLE			
CNA1180	◆	QUICK RELEASE AXLE			
REAR WHEELS POSITION					
CNA1190	◆	STANDARD REAR WHEELS POSITION	25MM FRONTWARD/BACKREST CANE AXLE		
CNA1195	□	MEDIUM REAR WHEELS POSITION	ON LINE/BACKREST CANE AXLE		
CNA1200	□	PASSIVE REAR WHEELS POSITION	25MM BACKWARD/BACKREST CANE AXLE		
CNA1201	□	DUAL AMPUTEE REAR WHEELS POSITIONS	TRANSIT=75MMBACKWARD/BACKREST CANE AXLE, ADJUSTABLE TO 50 MM BY THE DEALER & SELF PROP=60 MM BACKWARD/BACKREST CANE AXLE, ADJUSTABLE TO 85 MM		(NC: 1275, 1280, 1375)
PUMP					
CNA1220	◆	WITHOUT PUMP			
CNA1230	□	PUMP AND SUPPORT			
CHASSIS OPTIONS TYPE					
CNA1240	◆	NO CHASSIS OPTION			(NC : 1260, 1261, 1262, 1253, 1254, 1255, 1256)
CNA1260	■	SPOKE GUARDS SMILEY PAIR	FOR UK ONLY; UNTIL END OF APRIL		
CNA1261	■	SPOKE GUARDS SPIDER PAIR	FOR UK ONLY; UNTIL END OF APRIL		
CNA1262	■	SPOKE GUARDS FLOWER PAIR	FOR UK ONLY; UNTIL END OF APRIL		
CNA1253	□	SPOKE GUARDS CLEAR PAIR			
CNA1254	□	SPOKE GUARDS PIXEL; PAIR			
CNA1255	□	SPOKE GUARDS SNOWFLAKES; PAIR			
CNA1256	□	SPOKE GUARDS COLOR ME!; PAIR	INCLUDING A SET OF PENCILS		
BRAKES TYPE					
CNA1270	◆	PUSH TO LOCK BRAKES			
CNA1275	□	PULL TO LOCK BRAKES			
CNA1280	□	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)			(NC : 1426)
SAFETY OPTIONS TYPE					
CNA1340	□	NO SAFETY OPTION			(NC : 1345, 1360, 1361, 1370, 1375, 1380, 1390)
CNA1345	□	TRANSPORTATION KIT INCLUDING POSTURE BELT			(NC : 1360)
CNA1360	□	POSTURE BELT (CRASH TESTED);			
CNA1361	□	PEDIATRIC BELT (NOT CRASH TESTED);			
CNA1370	◆	REAR ANTI-TIPPERS PAIR WITH TIPPING AID			(NC : 1430, 1440)
CNA1375	□	TRANSFER KIT (100 MM (4") WHEELS)	ANTI-BARRIER		
CNA1380	□	REFLECTORS KIT (FOR SPOKE WHEELS ONLY)			
CNA1390	□	BRAKE EXTENSION PAIR			(NC: 1391, 1392)
ATTENDANT OPTIONS TYPE					
CNA1420	◆	NO ATTENDANT OPTION			(NC : 1426, 1430, 1440)
CNA1426	□	STROLLER PUSH HANDLES PAIR			
CNA1430	□	RH TIPPING AID	MULTI-ANGLE ADJUSTABLE		
CNA1440	□	LH TIPPING AID			
ACCESSORIES TYPE					
CNA1460	◆	NO ACCESSORY			(NC : 1461, 1470, 1471, 1472)
CNA1461	◆	USER TOOL KIT			
CNA1470	□	ACRYLIC TRAY			
CNA1471	□	RH SWING AWAY ACRYLIC TRAY			
CNA1472	□	LH SWING AWAY ACRYLIC TRAY			
FREIGHT & PRICE					
FREIGHT PRICE					
TOTAL PRICE					