



AVIVA FX40

MPS Maxx Stander



Our new power stander for a different perspective

A revolutionary new assistive design with AVIVA FX40 base



Yes, you can.®

A new perspective

Embrace a new perspective with the AVIVA FX40 MPS Maxx power stander and enjoy your surroundings from a whole new level. Interact with friends, laugh and swap stories. Experience the benefits of standing.

Meet. Greet. Swap Stories.



Personalised fit for you

Our MPS Maxx is a clinically advanced seating solution designed for optimal comfort and support whether you are seated or standing. Our multi-award winning Matrix seating combines with a wide range of our latest options and accessories to ensure a personalised fit for you.

Matrix LIBRA seat cushions

The award-winning Matrix Libra seat cushion has a contoured design to maximise the surface area for increased protection. Made with high resilient foam and a dual layer Flo-tech sac for enhanced shear management.

* Product specifications may vary from those displayed.



Matrix ELAN® headrests

Our highly adjustable Matrix Elan headrest is designed with dual-density moulded foam that provides immersion, comfort, protection and helps prevent bottoming out.

Lateral trunk supports

A choice of adjustable lateral supports for additional postural support and comfort.

Backrest options

A range of backrest options including our Matrix series are available for maximum postural support, pressure relief and exceptional microclimate control. Featuring built in contoured support for increased lateral stability.

Bodypoint®

For optimal function, mobility and stability when seated or standing we offer a choice of chest and positioning belts. This helps ensure both pelvis and upper body are secure for more control and improved functionality.

Hip guides and lateral hip/knee supports

A range of supports in various sizes and attachments for added positioning and ease of transfer.

Effortless positioning

Our MPS Maxx Stander is a powered multi-positioning solution for increased independence and accessibility. Effortlessly reposition for pressure redistribution, functional reach or to socialise with friends. The benefits of standing are endless.



45°

seat tilt with centre of gravity shift for pressure management.

7"

seat elevation for ease of access, transfers and increased independence.

178°

open angle standing function.

170°

of power recline with Extended Shear Reduction (ESR).

Active anti-tippers

connect to the base for stability when climbing or descending obstacles.



* Product specifications may vary from those displayed.

Exceptional drive performance

Take time to discover the new **AVIVA FX40** base. A new and exciting addition to our popular AVIVA range. This new front wheel drive power base maintains the true DNA of the AVIVA series. Enjoy a premium drive experience.

Design. Technology. Performance.

AVIVA FX40



>>> LiNX G-Trac Gyro

Travelling long distances or at higher speeds has been made even easier with LiNX G-Trac Gyro technology as standard. A gyro module within the chassis ensures the **AVIVA FX40** follows the control input, providing reassurance and confidence.



4Sure® Suspension System

A new front wheel drive design that improves control and traction across different terrains and uneven surfaces. The suspension isolates the seating system from vibrations and shocks to provide a high level of comfort.



Drive Performance

Articulating front anti-tippers, true shock absorbers and elastomer dampeners all work together to maximise both comfort and safety. The seat plane stays as level as possible, even when tackling kerbs and other obstacles.



4-bar Linkage

The remarkable drive performance is further enhanced by a patented 4-bar linkage with floating pivot points creating a smooth ride, superior traction, and powerful climbing even at low speeds.



LINX[®] Controls

LINX[®] Technology is an insight inspired control system designed to provide an excellent driving experience. Featuring personalised programming to ensure you feel confident, safe and in control.

Wireless maintenance updates and programming can be carried out by your prescriber remotely and in real time.

Easy-to-use integrated function

The M290 is an easy-to-use, interactive touch screen seating control unit for caregiver control. Mounted at the rear of the chair for easy caregiver access.



3.5" Touchscreen technology

The REM400 remote has an impressive 3.5" colour touchscreen with joystick, that requires minimal force to operate. Also included on the large screen is a speedometer and an odometer.



Connect to multiple devices with Built-in Bluetooth[®]

The REM400 remote connects to multiple devices simultaneously. It features an integrated mouse mover and switch control for speakers, phones and laptops.

Display-only version available

The REM500 is an optional display-only version of the REM400, retaining the 3.5" interactive touch screen, and allows specialty controls to be used to operate the power wheelchair.

* Product specifications may vary from those displayed.

Technical data

For more information about this product, including the user manual, please visit our website: www.invacare.eu.com

Seat width ▼	Seat depth ▼	Seat height ▼	Backrest height ▼	Armrest height ▼	Length of footrests ▼	Seat angle ▼
405 - 510 mm	405 - 510 mm	470 or 495 mm	305 - 710 mm	241 - 330 mm	203 - 430 mm	0 - 45°
Backrest angle ▼	Total height ▼	Total length inc. footrests ▼	Total length without footrests ▼	Driving unit width ▼	Max. user weight ▼	Battery capacity ▼
170° ESR	1115 - 1445 mm	1250 - 1305 mm	1070 mm	604 mm	113 kg	60 Ah
Turning diameter ▼	Max. climbable obstacle height ▼	Max. safe slope ▼	Motor capacity ▼	Speed ▼	Driving range ▼	Total weight ▼
1360 mm	83 mm	9°	2 x 340 W	10 or 12 km/h	Up to 27 km	190 - 210 kg

Colours

Select from a broad range of coloured frame shrouds and wheel inserts.

Frame shroud colours



Rim insert colours



Tel: +41 61 487 70 70
 hqeuropa@invacare.com
 www.invacare.eu.com

Invacare International GmbH
 Neuhofweg 51
 4147 Aesch BL
 Switzerland



Yes, you can.®