



**Yes, you can:**

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

July 2026  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX20 Modulite

### Standard Features

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm - 510mm (16" - 20"), stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £6,500**



Please make your selection with a cross in the box     Free of charge     standard feature     POA special request

Customer weight = **KG**

Code	Description	
<b>BASE MODEL*</b>		
JBC0000	AVIVA RX20 Modulite	£6,500
JBZ1300	Special AVIVA RX20 Modulite (Extended Lead Time). Please provide details in comments section.	POA
<b>CHASSIS WIDTH*</b>		
JBC2810	Narrow (560 mm/22")	◆
<b>SEATING TYPE*</b>		
JBC2100	Modulite Seat Unit	◆
<b>SEAT FRAME*</b>		
JBC2111	Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	(NC 2156) ◆
JBC2112	Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.	(NC 2156) ○
JBX2102	Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	(NC 2122, 2123) £110
<b>SEAT WIDTH*</b>		
JBC2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	(NC 2174, 2175, 2176, 2148) ○
JBC2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	○
JBC2123	480 mm Seat Plate	(NC 2704, 2706, 2148) ○
JBC2124	530 mm Seat Plate	(NC 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2106, 2148) ○
JBC2104	550 mm Seat Plate	(NC 2132, 2134, 2142, 2143, 2144, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2183, 2192, 2193, 2194, 2198, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 1445, 0556, 0589, 0526, 2701, 2702, 2703, 2704, 2705, 2706, 1404, 1405) ○
<b>BODY FIT KIT</b>		
JBC2106	Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)	(NC 2192, 2194) ○
<b>SEAT DEPTH*</b>		
JBC2132	410 mm	(NC 2197, 2198, 2701, 2702, 2703, 2704, 2705, 2706) ○
JBC2134	460 mm	○
JBC2136	510 mm	○
<b>SEAT CUSHION*</b>		
JBC2142	Matrx Flo-tech Lite Cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) £120
JBC2143	Matrx Flo-tech Contour Visco; (2-way stretch material)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) £169
JBC2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)	£410
JBC2144	Matrx Vi Cushion (Mesh fabric)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983) £400
JBC2148	Standard Cushion (90mm high) - For 550mm Seat Plate Only	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) ○
JBC2140	Without seat cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) ○
<b>1st POSITIONING ACCESORIES FOR MATRX LIBRA</b>		
JBC1980	Tapered Wedges (Pair)	£30
JBC1981	Abductor	£23
JBC1982	Lateral Wedge	£20
JBC1983	Anterior Seat Wedges	£27

2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
JBC1980		Tapered Wedges (Pair)	£30
JBC1981		Abductor	£23
JBC1982		Lateral Wedge	£20
JBC1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
JBC1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACKREST SUPPORT*			
JBC2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	(NC 2173, 2052, 2053) ○
JBC2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	○
PUSH HANDLES FOR BACKPOSTS			
JBC2108		Backposts without Push Handles	(NC 2219, 2220, 2227, 2228, 2245, 2291, 2292) ○
BACKREST UPHOLSTERY / CUSHION*			
JBC2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 2230, 2247, 2052, 2053) ◆
JBC2179		Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).	(NC 2192, 2193, 2194, 2199, 2155, 2230, 2247) £155
JBC2172		Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2128) £250
JBC2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2128) £400
JBC2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	£353
JBC2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£353
JBC2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£494
JBC2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£558
JBC2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£546
JBC2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) £475
JBC2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) £575
JBC2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2230, 2247, 2151, 2152, 2155) ○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)			
JBC2380		405 mm (16 inch) MAC BH	(NC 1445) ○
JBC2381		455 mm (18 inch) MAC BH	(NC 1445) ○
JBC2382		505 mm (20 inch) MAC BH	(NC 1445) ○
JBC2383		560 mm (22 inch) MAC BH	(NC 1445) ○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)			
JBC2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445) £100
JBC2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445) £100
JBC2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445) £100

<b>MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)</b>				
JBC2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445)	£100
JBC2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445)	£100
JBC2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445)	£100
<b>MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)</b>				
JBC2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
JBC2394		Backrest matches seat width	(NC 1445)	○
JBC2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○
<b>UPHOLSTERY COLOUR*</b>				
JBC2151		Black Nylon Fabric TR18		◆
JBC2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26		○
JBC2155		Black Spacer Fabric TR36	(NC 2171)	○
JBC2156		Black Nylon (For 550mm SW Only)		○
JBX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2173, 2174, 2175, 2176, 2177, 2178)	○
<b>ARMRESTS*</b>				
JBC2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
JBC2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm		£40
JBC2198		Padded Armrest with padded side guard (Water Fall Type)	(NC 2195, 2196, 2197, 2199)	£120
JBC2183		Inrigged Armrests. Reduces SW by 25 mm each side		£30
JBC2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	(NC 1445, 2179)	£65
JBC2192		Following Armrest. Height adjustable from 230 - 300 mm	(NC 2179, 2106)	£185
JBC2194		Following Armrest. Height adjustable from 300 - 360 mm	(NC 2179, 2106)	£185
<b>ARMPADS*</b>				
JBC2195		Short Armpad (325 mm)		◆
JBC2196		Long Armpad (375 mm)	(NC 2198)	○
JBC2197		Long Ergonomic Armpad (397 mm)	(NC 2198, 2230, 2247)	○
JBC2199		Soft Armpad (REA Azalea Wide Pad)		£95
JBX3000		Please select with Option 2198 - Padded Armrest with padded side guard (Water Fall Type)	(NC 2181, 2183, 2193, 2192)	○
<b>HIP SUPPORTS LH</b>				
JBC2184		Hip Support LH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2186		Hip Support LH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2188		Hip Support LH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
<b>HIP SUPPORTS RH</b>				
JBC2185		Hip Support RH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2187		Hip Support RH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBC2189		Hip Support RH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>				
JBC2164		Mounting hardware for hip supports (6,1" / 155 mm high)		○
JBC2165		Mounting hardware for hip supports (8,1" / 206 mm high)		○

<b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b>				
JBC1350		Neckrest depth adjustable	(NC 2247)	£80
JBC1432		Headrest sideways adjustable	(NC 2247)	£225
JBC1427		Headrest with cheek support	(NC 2247)	£225
JBC2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	(NC 2247, 2170)	£325
JBC2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	(NC 2247, 2170)	£325
JBC2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware	(NC 2247, 2170)	£440
JBC1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171, 2174, 2175, 2176, 2177, 2178)	○
JBC1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		○
<b>LATERAL SUPPORTS</b>				
JBC1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£165
<b>MANUAL / POWERED SEATING OPTIONS*</b>				
JBC1091		Manually settable Seat Tilt 0 - 15°		○
JBC1092		Powered "fixed pivot" Tilt 20°	(NC 2201)	£400
JBC1093		Powered Tilt (0° - 30°) with CoG shift	(NC 2201)	£1,100
JBC2052		Powered Tilt (0° - 50°) with CoG shift	(NC 2201, 2171)	£1,450
JBC1094		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected	(NC 2201, 2220)	£1,900
JBC2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG. 6 point Tie downs must be selected	(NC 2201, 2220, 1910, 1912, 2171)	£2,500
<b>SEAT TO GROUND HEIGHT*</b>				
JBC2857		435 mm		○
JBC2858		460 mm		○
JBC2859		485 mm		○
<b>BACKREST RECLINER*</b>				
JBC0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°		◆
JBC0977		Powered Recliner 30°		£500
<b>LEG REST TYPE*</b>				
JBC0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 1244, 1416, 0501)	£370
JBC0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416, 0501)	£1,620
JBC0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 2365, 0559, 0516, 0517, 0501)	◆
JBC0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0559, 0516, 0517, 0554, 0555, 0501)	£120
JBC0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£235
JBC0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£660
JBC0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517)	○

FOOTPLATE TYPE*			
JBC2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£105
JBC0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBC0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)	£135
JBC0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBC0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBC0504		Angle Adjustable Footplate pair (NC 0526)	£75
JBC0513		Angle and depth adjustable Footplate pair	£120
JBC0501		Without footplates for option 0501 without Leg Rests (NC 0503, 0533, 0531)	○
AMPUTEE SUPPORT			
JBC0554		Amputee Support RH	£195
JBC0555		Amputee Support LH	£195
LINX CONTROLS*			
JBC2201		DLX-REM110 (Drive Only) (NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908)	◆
JBC2203		DLX-REM211 (Drive & Actuators) (NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 1908)	○
JBC2204		DLX-REM216 (Drive, Actuators & Lights) (NC 2232, 2233, 2230, 2247, 2257, 2240, 2270)	£200
JBC2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets) (NC 2220)	£425
JBC2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick). (NC 2220, 2270, 1709, 1710, 1207, 1208)	£765
REM400 TOGGLE OPTION			
JBC2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL) (NC 2203, 2204, 2206)	£315
CONTROL BRACKETRY*			
JBC1207		Standard Bracket RH (NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	○
JBC1208		Standard Bracket LH (NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	○
JBC1707		Swing Away RH Bracket (NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	£135
JBC1708		Swing Away LH Bracket (NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	£135
JBC1709		Matrx Resolve Swing in/out Bracket RH (NC 2283, 2287, 2207, 2265, 2275, 1404)	£245
JBC1710		Matrx Resolve Swing in/out Bracket LH (NC 2284, 2288, 2208, 2266, 2276, 1405)	£245
JBC2213		Centre Mount Swing Away to RH Nucleus Bracket (NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405)	£495
JBC2214		Centre Mount Swing Away to LH Nucleus Bracket (NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405)	£495
JBC2209		REM500 Swing Away RH Multi-Adjustable Bracket (NC 2203, 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404)	○
JBC2210		REM500 Swing Away LH Multi-Adjustable Bracket (NC 2203, 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405)	○
JBC2227		Bracket RH mounted for Attendant (NC 2220, 2128, 2108, 2219, 1221, 1312)	£65
JBC2228		Bracket LH mounted for Attendant (NC 2220, 2128, 2108, 2219, 1222, 1363)	£65
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
DUAL CONTROLS			
JBC2220		LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228, 9259)	£395
JBC2219		LiNX DLX-ACU200 (Compact Dual Control) (NC 2108, 1221, 1222, 2227, 2228, 9259)	£395
JBC2245		LiNX DLX-IDC (Intuitive Dual Control) (NC 2108, 2259, 2260, 9259)	£745
UNZ9259		mo-vis Scoot Control (Only available via Special) (NC 2108, 2259, 2260, 0000, 2259, 2260, 1221, 1222)	£2,315

<b>DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)</b>				
JBC2259		Dual Control Standard Bracket RH	(NC 2291, 2245, 1312)	○
JBC2260		Dual Control Standard Bracket LH	(NC 2292, 2245, 1363)	○
JBC1221		IDC Dual Control Detachable Bracket RH	(NC 2291, 2220, 2219, 2227, 1312)	○
JBC1222		IDC Dual Control Detachable Bracket LH	(NC 2292, 2220, 2219, 2228, 1363)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
<b>COMPACT REMOTES</b>				
JBC2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£825
JBC2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0000, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£3,190
JBC2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215, 2220)	£4,100
JBC2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215, 2220)	£2,285
<b>COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)</b>				
JBC2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBC2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBC2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBC2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
JBC2285		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JBC2286		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JBC2283		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBC2284		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBC2287		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBC2288		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2289		Individual mounting option per request (Special). Requires additional information/specification.		POA
<b>CHIN CONTROLS</b>				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,015

<b>CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)</b>			
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)	£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)	£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)	£2,280
<b>HEAD CONTROLS</b>			
JBC2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £4,220
JBC2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220) £4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0000, 2220) £3,315
<b>FOOT CONTROL</b>			
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220) £2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220) £2,740
<b>FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)</b>			
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.	POA
<b>TRAY CONTROL</b>			
JBC2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128) £6,535
<b>TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)</b>			
JBC2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415) £495
JBC2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339) £495
JBZ2269		Individual mounting option per request (Special). Requires additional information/specification.	POA
<b>WIRELESS ACCESSORIES FOR ASL CONTROLS</b>			
JBC2249		Wireless Triple Switch Receiver (ASL557-3)	£340
JBC2250		Wireless Mouse Emulator (ASL558)	£810
<b>LINX 10-WAY SWITCH</b>			
JBC2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
<b>LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)</b>			
JBC2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404) ○
JBC2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405) ○
JBC2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 2108) £65
JBC2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 2108) £65
<b>LINX USB CHARGER PORT</b>			
JBC2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £160
<b>LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)</b>			
JBC2281		RH Mount	○
JBC2282		LH Mount	○
<b>LINX INPUT MODULE</b>			
JBC2273		DLX-IN200 Module (Jack port for external power button)	£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	(NC 2201, 2203, 2204, 2220) £785

<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
JBC2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
JBC2300		To be enclosed only	○
JBZ2301		To be connected and programmed (Special)	POA
<b>BUDDY BUTTONS</b>			
JBC1685		Red Buddy Button	£90
JBC1686		Yellow Buddy Button	£90
JBC1687		Green Buddy Button	£90
JBC1690		Blue Buddy Button	£90
<b>EGG SWITCH</b>			
JBC2254		Egg Switch, black	£90
<b>PIKO BUTTONS</b>			
JBC1711		Red Piko Button	£160
JBC1712		Yellow Piko Button	£160
JBC1713		Green Piko Button	£160
JBC1714		Orange Piko Button	£160
JBC1715		Blue Piko Button	£160
JBC1716		Black Piko Button	£160
<b>BODYPOINT JOYSTICK KNOBS</b>			
JBC2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JBC2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JBC2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JBC2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
<b>5:1 SWITCH ADAPTOR</b>			
JBC1681		5:1 Switch Adapter ( 9-pin sub D to 5 jack ports)	£120
<b>SHROUD COLOUR*</b>			
JBC2871		Mid range - Black/Anthracite Shroud	○
JBC2872		Mid range - Silver Shroud	○
JBC2873		Mid range - Bengal Red Shroud	○
JBC2874		Mid range - Mercury Blue Shroud	○
JBC2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300
<b>RIM COLOUR*</b>			
JBC2882		Mid range - Black Rim Inserts	○
JBC2883		Mid range - Silver Rim Inserts	○
JBC2884		Mid range - Bengal Red Rim Inserts	○
JBC2885		Mid range - Mercury Blue Rim Inserts	○
JBC2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
<b>FORKS*</b>			
JBC2631		Single sided fork	£250
JBC2632		Double sided fork Without castor rim insert	◆

<b>CASTORS*</b>			
JBC0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBC0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JBC0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JBC0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	◆
<b>CASTOR LOCKS</b>			
JBC1490		Castor Lock	£195
<b>DRIVE WHEELS*</b>			
JBC0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBC0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
<b>C.T.C SUSPENSION*</b>			
JBC2870		CTC Suspension Control - Traction - Comfort	◆
<b>SUSPENSION SETTING*</b>			
JBC1910		Soft (NC 1094)	◆
JBC1912		Medium (NC 1094)	○
JBC1911		Firm	○
JBC1915		Very firm	○
<b>MOTOR TYPE*</b>			
JBC1294		Indoor setting, reduced speed approx. 3 km/h	○
JBC1282		6 km/h	◆
JBC1908		10 km/h Package (NC 2203, 1437, 1254) Always combined with LED Light system with Indicators, 30 Lux headlights	£515
<b>BATTERY TYPE*</b>			
JBC1259		60 Ah Gel	◆
<b>G TRAC</b>			
JBC1215		G-Trac Module	£400
<b>CHARGER*</b>			
JBC1343		Opti Charge Charger 8A (IP44)	◆
JBC1313		Opti Charge Charger 10A (IP44)	£155
<b>ANTI-TIPPERS*</b>			
JBC1205		Anti-tipper	◆
<b>LIGHTS</b>			
JBC1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors (NC 2203, 1908)	£450
JBC1254		Without Light system	○
<b>KERB RIDER</b>			
JBC1244		Kerb Rider with park position feature (NC 0556, 0589, 1294)	£215
<b>TRANSPORTATION KIT*</b>			
JBC1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JBC1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBC1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory. (NC 1091, 1092, 1093)	○

<b>EASY ACCESS FRONT TIE DOWNS</b>			
JBC2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
<b>FIXATION BELT*</b>			
JBC1329		Lap belt with Auto Lock Standard mount on backrest	○
JBC2113		Lap belt with Auto Lock Flexible mount in side rails	◆
JBC2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBC2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBC2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBC2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBC2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBC2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
<b>REAR VIEW MIRROR</b>			
JBC1415		Rear view mirror RH	(NC 2217, 2285, 2209, 2213, 2265, 2275, 2283, 2287) £45
JBC1339		Rear view mirror LH	(NC 2218, 2286, 2210, 2214, 2266, 2276, 2284, 2288) £45
<b>BAG HOOKS</b>			
JBC2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2172, 2173, 2227, 2228, 2233, 2240, 2257, 2291, 2292) £80
<b>CANE HOLDER</b>			
JBC1312		Cane holder RH	(NC 2227, 2259, 1221) £45
JBC1363		Cane holder LH	(NC 2228, 2260, 1222) £45
<b>WHEEL FENDERS</b>			
JBC1416		Fender for Castors	£80
JBC1417		Fender for Drive wheels	£160
<b>TRAY</b>			
JBC1404		Tray, Swing Away RH	(NC 1707, 2209, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217) £200
JBC1405		Tray, Swing Away LH	(NC 1708, 2210, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218) £200

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

JBC0000	AVIVA RX20 Modulite
JBC2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
JBC2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
JBC2123	480 mm Seat Plate
JBC2124	530 mm Seat Plate
JBC2106	Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)
JBC2132	410 mm
JBC2134	460 mm
JBC2142	Matrx Flo-tech Lite Cushion
JBC2143	Matrx Flo-tech Contour Visco; (2-way stretch material)
JBC2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)
JBC2144	Matrx Vi Cushion (Mesh fabric)
JBC2148	Standard Cushion (90mm high) - For 550mm Seat Plate Only
JBC1980	Tapered Wedges (Pair)
JBC1981	Abductor
JBC1982	Lateral Wedge
JBC1983	Anterior Seat Wedges
JBC1980	Tapered Wedges (Pair)
JBC1981	Abductor
JBC1982	Lateral Wedge
JBC1983	Anterior Seat Wedges
JBC1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
JBC2108	Backposts without Push Handles
JBC2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
JBC2179	Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).
JBC2172	Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBC2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBC2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBC2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
JBC2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
JBC2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBC2380	405 mm (16 inch) MAC BH
JBC2381	455 mm (18 inch) MAC BH
JBC2382	505 mm (20 inch) MAC BH

JBC2383	560 mm (22 inch) MAC BH
JBC2385	90 mm (3.5 inch) MAC Contour Depth LH
JBC2386	150 mm (6 inch) MAC Contour Depth LH
JBC2387	230 mm (9 inch) MAC Contour Depth LH
JBC2389	90 mm (3.5 inch) MAC Contour Depth RH
JBC2390	150 mm (6 inch) MAC Contour Depth RH
JBC2391	230 mm (9 inch) MAC Contour Depth RH
JBC2393	Backrest -25 mm (-1 inch) to seat width
JBC2394	Backrest matches seat width
JBC2395	Backrest +25 mm (+1 inch) to seat width
JBC2151	Black Nylon Fabric TR18
JBC2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
JBC2155	Black Spacer Fabric TR36
JBC2156	Black Nylon (For 550mm SW Only)
JBC2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
JBC2198	Padded Armrest with padded side guard (Water Fall Type)
JBC2183	Inrigged Armrests. Reduces SW by 25 mm each side
JBC2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
JBC2192	Following Armrest. Height adjustable from 230 - 300 mm
JBC2194	Following Armrest. Height adjustable from 300 - 360 mm
JBC2195	Short Armpad (325 mm)
JBC2196	Long Armpad (375 mm)
JBC2197	Long Ergonomic Armpad (397 mm)
JBC2199	Soft Armpad (REA Azalea Wide Pad)
JBC2184	Hip Support LH, pad size 140 x 100 mm
JBC2186	Hip Support LH, pad size 230 x 100 mm
JBC2188	Hip Support LH, pad size 330 x 100 mm
JBC2185	Hip Support RH, pad size 140 x 100 mm
JBC2187	Hip Support RH, pad size 230 x 100 mm
JBC2189	Hip Support RH, pad size 330 x 100 mm
JBC1350	Neckrest depth adjustable
JBC1432	Headrest sideways adjustable
JBC1427	Headrest with cheek support
JBC2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
JBC2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
JBC2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
JBC1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBC1445	Lateral Support with extended arm, multi-adjustable pad
JBC1091	Manually settable Seat Tilt 0 - 15°
JBC1092	Powered "fixed pivot" Tilt 20°
JBC1093	Powered Tilt (0° - 30°) with CoG shift
JBC2052	Powered Tilt (0° - 50°) with CoG shift
JBC1094	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected

JBC2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG. 6 point Tie downs must be selected
JBC0977	Powered Recliner 30°
JBC0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
JBC0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
JBC0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
JBC0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
JBC0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBC0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBC0501	Without Leg Rests
JBC2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JBC0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
JBC0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
JBC0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBC0512	Fixed Footplate pair
JBC0504	Angle Adjustable Footplate pair
JBC0513	Angle and depth adjustable Footplate pair
JBC0501	Without footplates for option 0501 without Leg Rests
JBC0554	Amputee Support RH
JBC0555	Amputee Support LH
JBC2201	DLX-REM110 (Drive Only)
JBC2203	DLX-REM211 (Drive & Actuators)
JBC2204	DLX-REM216 (Drive, Actuators & Lights)
JBC2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JBC2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JBC2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JBC1207	Standard Bracket RH
JBC1208	Standard Bracket LH
JBC1707	Swing Away RH Bracket
JBC1708	Swing Away LH Bracket
JBC1709	Matrx Resolve Swing in/out Bracket RH
JBC1710	Matrx Resolve Swing in/out Bracket LH
JBC2213	Centre Mount Swing Away to RH Nucleus Bracket
JBC2214	Centre Mount Swing Away to LH Nucleus Bracket
JBC2209	REM500 Swing Away RH Multi-Adjustable Bracket
JBC2210	REM500 Swing Away LH Multi-Adjustable Bracket
JBC2227	Bracket RH mounted for Attendant
JBC2228	Bracket LH mounted for Attendant
JBC2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
JBC2219	LiNX DLX-ACU200 (Compact Dual Control)
JBC2245	LiNX DLX-IDC (Intuitive Dual Control)
UNZ9259	mo-vis Scoot Control (Only available via Special)

JBC2259	Dual Control Standard Bracket RH
JBC2260	Dual Control Standard Bracket LH
JBC1221	IDC Dual Control Detachable Bracket RH
JBC1222	IDC Dual Control Detachable Bracket LH
JBC2231	LiNX DLX-CR400 (Compact Remote)
JBC2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JBC2233	ASL MEC (Proportional MEC Joystick ASL130)
JBC2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
JBC2207	RH Mount
JBC2208	LH Mount
JBC2265	Centre Mount Swing Away to RH (Nucleus)
JBC2266	Centre Mount Swing Away to LH (Nucleus)
JBZ2215	Individual mounting option per request (Special). Requires additional information/specification.
JBC2285	Tray Mount RH (ASL603, 100mm x 150mm)
JBC2286	Tray Mount LH (ASL603, 100mm x 150mm)
JBC2283	Centre Mount Swing Away to RH (Nucleus)
JBC2284	Centre Mount Swing Away to LH (Nucleus)
JBC2287	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JBC2288	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JBZ2289	Individual mounting option per request (Special). Requires additional information/specification.
JBC2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
JBC2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
JBC2257	Four-Switch Proximity Array (ASL106)
JBC2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBC2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBC2217	Bracket RH, mounted for user
JBC2218	Bracket LH, mounted for user
JBC2291	Bracket RH, mounted for attendant
JBC2292	Bracket LH, mounted for attendant
JBC2253	USB Charger Port DLX-USB02-A
JBC2273	DLX-IN200 Module (Jack port for external power button)
JBC1686	Yellow Buddy Button
JBC1687	Green Buddy Button
JBC1690	Blue Buddy Button
JBC2254	Egg Switch, black
JBC1711	Red Piko Button
JBC1712	Yellow Piko Button
JBC1713	Green Piko Button
JBC1714	Orange Piko Button
JBC1715	Blue Piko Button
JBC1716	Black Piko Button
JBC1910	Soft
JBC1912	Medium
JBC1294	Indoor setting, reduced speed approx. 3 km/h

JBC1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBC1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBC1254	Without Light system
JBC1244	Kerb Rider with park position feature
JBC2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBC2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBC2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBC2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBC2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBC2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBC1415	Rear view mirror RH
JBC1339	Rear view mirror LH
JBC2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
JBC1312	Cane holder RH
JBC1363	Cane holder LH
JBC1416	Fender for Castors
JBC1404	Tray, Swing Away RH
JBC1405	Tray, Swing Away LH



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):