



**Yes, you can:**

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

**AVIVA RX20 Standard**

**Standard Features**

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Modulite Standard Seat
- One piece seat plate (380mm/15" or 430mm/17" SW), additional 50mm/2" arm rest width adjustment
- Seat depths from 410mm - 510mm (16" - 20")
- Tension adjustable backrest
- Standard foam cushion (50mm thick)
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0°- 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM211 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction – Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



**Standard Base price: €4,826**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

| Code                  | Description   |        |
|-----------------------|---|--------|
| <b>BASE MODEL*</b>    |   |        |
| JBB0000               | AVIVA RX20 Standard   | €4,826 |
| JBZ1300               | Special AVIVA RX20 Standard (Extended Lead Time). Please provide details in comments section. | POA    |
| <b>CHASSIS WIDTH*</b> |   |        |
| JBB2810               | Chassis Narrow (560 mm/22")   | ◆      |
| <b>SEATING TYPE*</b>  |   |        |
| JBB2105               | Modulite Standard Seat (Standard seat plate and backrest frame in fixed width)                | ◆      |

| <b>SEAT FRAME*</b>  |                          |   |      |
|---|--------------------------|---|------|
| JBB2110   | <input type="checkbox"/> | One Piece Seat Plate  | ◆    |
| <b>SEAT WIDTH*</b>  |                          |   |      |
| JBB2121   | <input type="checkbox"/> | 380 mm Seat Plate. Additional 50mm arm rest width adjustment.   | ○    |
| JBB2122   | <input type="checkbox"/> | 430 mm Seat Plate. Additional 50mm arm rest width adjustment.   | ◆    |
| <b>SEAT DEPTH*</b>  |                          |   |      |
| JBB2132   | <input type="checkbox"/> | 410 mm (NC 2197)  | ○    |
| JBB2134   | <input type="checkbox"/> | 460 mm  | ○    |
| JBB2136   | <input type="checkbox"/> | 510 mm  | ○    |
| <b>SEAT CUSHION*</b>  |                          |   |      |
| JBB2141   | <input type="checkbox"/> | Standard Foam Cushion, 50 mm height   | ◆    |
| JBB2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite  | €148 |
| JBB2143   | <input type="checkbox"/> | Matrx Flo-tech Contour Visco; (2-way stretch material)  | €207 |
| JBB2140   | <input type="checkbox"/> | Without seat cushion  | ○    |
| <b>BACKREST SUPPORT*</b>  |                          |   |      |
| JBB2161   | <input type="checkbox"/> | Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs   | ◆    |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |   |      |
| JBB2171   | <input type="checkbox"/> | Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery. | ◆    |
| JBB2170   | <input type="checkbox"/> | Without Backrest Upholstery (Back Posts/Back Frame only) (NC 2151)  | ○    |
| <b>UPHOLSTERY COLOUR*</b>   |                          |   |      |
| JBB2151   | <input type="checkbox"/> | Black Nylon Fabric TR18   | ◆    |
| JBX1002   | <input type="checkbox"/> | No Upholstery Colour required - No backrest upholstery required   | ○    |
| <b>ARMRESTS*</b>  |                          |   |      |
| JBB2191   | <input type="checkbox"/> | T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm  | ◆    |
| JBB2181   | <input type="checkbox"/> | T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm   | €56  |
| JBB2183   | <input type="checkbox"/> | Inrigged Armrests. Reduces SW by 25 mm each side  | ○    |
| <b>ARMPADS*</b>   |                          |   |      |
| JBB2195   | <input type="checkbox"/> | Short Armpad (325 mm)   | ◆    |
| JBB2196   | <input type="checkbox"/> | Long Armpad (375 mm)  | ○    |
| JBB2197   | <input type="checkbox"/> | Long Ergonomic Armpad (397 mm)  | ○    |
| JBB2199   | <input type="checkbox"/> | Soft Armpad (REA Azalea Wide Pad)   | ○    |
| <b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b> |                          |   |      |
| JBB1350   | <input type="checkbox"/> | Neckrest depth adjustable   | €103 |
| JBB1495   | <input type="checkbox"/> | Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts                    | ○    |
| <b>MANUAL / POWERED SEATING OPTIONS*</b>                                |                          |   |      |
| JBB1092   | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20°  | ◆    |
| <b>SEAT TO GROUND HEIGHT*</b>   |                          |   |      |
| JBB2857   | <input type="checkbox"/> | 435 mm  | ○    |
| JBB2858   | <input type="checkbox"/> | 460 mm  | ○    |
| JBB2859   | <input type="checkbox"/> | 485 mm  | ○    |

| <b>BACKREST RECLINER*</b>  |                          |  |      |
|--|--------------------------|--|------|
| JBB0975  | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120°   | ◆    |
| <b>SIDE MOUNTED LEG RESTS</b>  |                          |  |      |
| JBB0526  | <input type="checkbox"/> | Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm                           | ◆    |
| <b>SIDE MOUNTED FOOTPLATE TYPE</b>   |                          |  |      |
| JBB0512  | <input type="checkbox"/> | Fixed Footplate pair   | ◆    |
| JBB0513  | <input type="checkbox"/> | Angle and depth adjustable Footplate pair  | €149 |
| <b>LINX CONTROLS*</b>  |                          |  |      |
| JBB2203  | <input type="checkbox"/> | DLX-REM211 (Drive & Actuators) (NC 1437, 1908)   | ◆    |
| JBB2204  | <input type="checkbox"/> | DLX-REM216 (Drive, Actuators & Lights)   | €191 |
| <b>CONTROL BRACKETRY*</b>  |                          |  |      |
| JBB1207  | <input type="checkbox"/> | Standard Bracket RH (NC 1404)  | ○    |
| JBB1208  | <input type="checkbox"/> | Standard Bracket LH (NC 1405)  | ○    |
| JBB1707  | <input type="checkbox"/> | Swing Away RH Bracket (NC 1404)  | €167 |
| JBB1708  | <input type="checkbox"/> | Swing Away LH Bracket (NC 1405)  | €167 |
| JBB2227  | <input type="checkbox"/> | Bracket RH mounted for Attendant (NC 2220, 2219, 1312)   | €81  |
| JBB2228  | <input type="checkbox"/> | Bracket LH mounted for Attendant (NC 2220, 2219, 1363)   | €81  |
| JBZ2239  | <input type="checkbox"/> | Individual mounting option per request (Special). Requires additional information/specification. | POA  |
| <b>DUAL CONTROLS</b>   |                          |  |      |
| JBB2220  | <input type="checkbox"/> | LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 2273, 2227, 2228)                    | €497 |
| JBB2219  | <input type="checkbox"/> | LiNX DLX-ACU200 (Compact Dual Control) (NC 2227, 2228)   | €497 |
| <b>DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)</b>        |                          |  |      |
| JBB2259  | <input type="checkbox"/> | Dual Control Standard Bracket RH (NC 2291, 1312)   | ○    |
| JBB2260  | <input type="checkbox"/> | Dual Control Standard Bracket LH (NC 2292, 1363)   | ○    |
| <b>LINX 10-WAY SWITCH</b>  |                          |  |      |
| JBB2251  | <input type="checkbox"/> | 10-Way Switch DLX-SW10-A for direct Powered Seating Access                                       | €601 |
| <b>LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)</b> |                          |  |      |
| JBB2217  | <input type="checkbox"/> | Bracket RH, mounted for user (NC 1404)   | ○    |
| JBB2218  | <input type="checkbox"/> | Bracket LH, mounted for user (NC 1405)   | ○    |
| JBB2291  | <input type="checkbox"/> | Bracket RH, mounted for attendant (NC 2259)  | €81  |
| JBB2292  | <input type="checkbox"/> | Bracket LH, mounted for attendant (NC 2260)  | €81  |
| <b>LINX IN200 MODULE</b>   |                          |  |      |
| JBB2273  | <input type="checkbox"/> | DLX-IN200 Module (Jack port for external power button)   | €303 |
| <b>FURTHER LINX COMPONENTS / ACCESSORIES</b>                                       |                          |  |      |
| JBB2302  | <input type="checkbox"/> | Not Necessary (NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716)                    | ◆    |
| JBB2300  | <input type="checkbox"/> | To be enclosed only  | ○    |
| JBZ2301  | <input type="checkbox"/> | To be connected and programmed (Special)   | POA  |
| <b>BUDDY BUTTONS</b>   |                          |  |      |
| JBB1685  | <input type="checkbox"/> | Red Buddy Button   | €110 |
| JBB1686  | <input type="checkbox"/> | Yellow Buddy Button  | €110 |
| JBB1687  | <input type="checkbox"/> | Green Buddy Button   | €110 |
| JBB1690  | <input type="checkbox"/> | Blue Buddy Button  | €110 |
| <b>EGG SWITCH</b>  |                          |  |      |
| JBB2254  | <input type="checkbox"/> | Egg Switch, black  | €110 |

| PIKO BUTTONS             |  |  |      |
|--------------------------|--|--|------|
| JBB1711                  |  | Red Piko Button  | €198 |
| JBB1712                  |  | Yellow Piko Button   | €198 |
| JBB1713                  |  | Green Piko Button  | €198 |
| JBB1714                  |  | Orange Piko Button   | €198 |
| JBB1715                  |  | Blue Piko Button   | €198 |
| JBB1716                  |  | Black Piko Button  | €198 |
| BODYPOINT JOYSTICK KNOBS |  |  |      |
| JBB2721                  |  | 3 inch U-Shaped, Joystick Handle, Bodypoint PC101A   | €79  |
| JBB2722                  |  | 4 inch U-Shaped, Joystick Handle, Bodypoint PC102A   | €79  |
| JBB2723                  |  | Domed Rubber, Joystick Handle, Bodypoint PC107A  | €79  |
| JBB2724                  |  | Mushroom Grooved, Joystick Handle, Bodypoint PC110A  | €79  |
| SHROUD COLOUR*           |  |  |      |
| JBB2871                  |  | Mid range - Black/Anthracite Shroud  | ○    |
| JBB2872                  |  | Mid range - Silver Shroud  | ○    |
| JBB2873                  |  | Mid range - Bengal Red Shroud  | ○    |
| JBB2874                  |  | Mid range - Mercury Blue Shroud  | ○    |
| JBB2875                  |  | Mid range - Midnight Purple Shroud   | ○    |
| JBZ0299                  |  | Premium - Individual Shroud Colour. RAL number required<br>Extended lead lime: + 15 days.  | €268 |
| RIM COLOUR*              |  |  |      |
| JBB2882                  |  | Mid range - Black Rim Inserts  | ○    |
| JBB2883                  |  | Mid range - Silver Rim Inserts   | ○    |
| JBB2884                  |  | Mid range - Bengal Red Rim Inserts   | ○    |
| JBB2885                  |  | Mid range - Mercury Blue Rim Inserts   | ○    |
| JBB2886                  |  | Mid range - Midnight Purple Rim Inserts  | ○    |
| JBZ2887                  |  | Premium - Individual Rim Insert Colour. RAL number required<br>Extended lead lime: + 15 days. Price Included with Premium Shroud Kit | ○    |
| FORKS*                   |  |  |      |
| JBB2632                  |  | Double sided fork Without castor rim insert  | ○    |
| CASTORS*                 |  |  |      |
| JBB0828                  |  | 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre   | ◆    |
| CASTOR LOCKS             |  |  |      |
| JBB1490                  |  | Castor Lock  | €248 |
| DRIVE WHEELS*            |  |  |      |
| JBB0704                  |  | DW 14" (355 mm) x 3" (75 mm), Pneumatic  | ◆    |
| JBB0705                  |  | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert   | ○    |
| C.T.C SUSPENSION*        |  |  |      |
| JBB2870                  |  | CTC Suspension Control - Traction - Comfort  | ◆    |
| SUSPENSION SETTING*      |  |  |      |
| JBB1910                  |  | Soft   | ◆    |
| JBB1912                  |  | Medium   | ○    |
| JBB1911                  |  | Firm   | ○    |
| JBB1915                  |  | Very firm  | ○    |

| <b>MOTOR TYPE*</b>         |  |  |                               |
|----------------------------|--|--|-------------------------------|
| JBB1294                    |  | Indoor setting, reduced speed approx. 3 km/h   | ○                             |
| JBB1282                    |  | 6 km/h   | ◆                             |
| JBB1908                    |  | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights  | (NC 2203, 1437, 1254)<br>€683 |
| <b>BATTERY TYPE*</b>       |  |  |                               |
| JBB1259                    |  | 60 Ah Gel  | ◆                             |
| <b>CHARGER*</b>            |  |  |                               |
| JBB1343                    |  | Opti Charge Charger 8A (IP44)  | ◆                             |
| JBB1313                    |  | Opti Charge Charger 10A (IP44)   | €191                          |
| <b>ANTI-TIPPERS*</b>       |  |  |                               |
| JBB1205                    |  | Anti-tipper  | ◆                             |
| <b>LIGHTS</b>              |  |  |                               |
| JBB1437                    |  | LED Light system with indicators 30 Lux headlights, included with 10km/h motors  | (NC 2203, 1908)<br>€602       |
| JBB1254                    |  | Without Light system   | ○                             |
| <b>KERB RIDER</b>          |  |  |                               |
| JBB1244                    |  | Kerb Rider with park position feature  | (NC 1294)<br>€266             |
| <b>TRANSPORTATION KIT*</b> |  |  |                               |
| JBB1450                    |  | Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering. | €736                          |
| JBB1494                    |  | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.   | ◆                             |
| <b>FIXATION BELT*</b>      |  |  |                               |
| JBB1329                    |  | Lap belt with Auto Lock Standard mount on backrest   | ◆                             |
| <b>REAR VIEW MIRROR</b>    |  |  |                               |
| JBB1415                    |  | Rear view mirror RH  | (NC 2217)<br>€56              |
| JBB1339                    |  | Rear view mirror LH  | (NC 2218)<br>€56              |
| <b>CANE HOLDER</b>         |  |  |                               |
| JBB1312                    |  | Cane holder RH   | (NC 2227, 2259)<br>€53        |
| JBB1363                    |  | Cane holder LH   | (NC 2228, 2260)<br>€53        |
| <b>WHEEL FENDERS</b>       |  |  |                               |
| JBB1416                    |  | Fender for Castors   | €110                          |
| JBB1417                    |  | Fender for Drive wheels  | €213                          |
| <b>TRAY</b>                |  |  |                               |
| JBB1404                    |  | Tray, Swing Away RH  | (NC 1707, 2217)<br>€210       |
| JBB1405                    |  | Tray, Swing Away LH  | (NC 1708, 2218)<br>€210       |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**

|         |   |
|---------|---|
| JBB2151 | Black Nylon Fabric TR18   |
| JBB2197 | Long Ergonomic Armpad (397 mm)  |
| JBB2203 | DLX-REM211 (Drive & Actuators)  |
| JBB1707 | Swing Away RH Bracket   |
| JBB1708 | Swing Away LH Bracket   |
| JBB2227 | Bracket RH mounted for Attendant  |
| JBB2228 | Bracket LH mounted for Attendant  |
| JBB2220 | LiNX DLX-REM050 (Dual Control used as secondary remote)                                     |
| JBB2219 | LiNX DLX-ACU200 (Compact Dual Control)  |
| JBB2259 | Dual Control Standard Bracket RH  |
| JBB2260 | Dual Control Standard Bracket LH  |
| JBB2217 | Bracket RH, mounted for user  |
| JBB2218 | Bracket LH, mounted for user  |
| JBB2291 | Bracket RH, mounted for attendant   |
| JBB2292 | Bracket LH, mounted for attendant   |
| JBB2273 | DLX-IN200 Module (Jack port for external power button)                                      |
| JBB1686 | Yellow Buddy Button   |
| JBB1687 | Green Buddy Button  |
| JBB1690 | Blue Buddy Button   |
| JBB2254 | Egg Switch, black   |
| JBB1711 | Red Piko Button   |
| JBB1712 | Yellow Piko Button  |
| JBB1713 | Green Piko Button   |
| JBB1714 | Orange Piko Button  |
| JBB1715 | Blue Piko Button  |
| JBB1716 | Black Piko Button   |
| JBB1294 | Indoor setting, reduced speed approx. 3 km/h  |
| JBB1908 | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights |
| JBB1437 | LED Light system with indicators 30 Lux headlights, included with 10km/h motors             |
| JBB1254 | Without Light system  |
| JBB1312 | Cane holder RH  |
| JBB1363 | Cane holder LH  |
| JBB1404 | Tray, Swing Away RH   |
| JBB1405 | Tray, Swing Away LH   |



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

**Special Instructions**

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):