



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 HD Modulite HD

Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO7176 - 19
- Modulite HD Seating System
- HD telescopic adjustable seat widths from 530 to 550 mm with an additional 50mm by arms
- Seat depths from 510mm - 580mm, stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0°- 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction – Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: €7,720



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBH0002	AVIVA RX40 HD Modulite HD	€7,720
JBZ1300	Special AVIVA RX40 HD Modulite HD (Extended Lead Time). Please provide details in comments section.	POA
CHASSIS WIDTH*		
JBH2830	Chassis Wide (635 mm/25")	◆
SEATING TYPE*		
JBH2101	Invacare® Modulite HD Seat Unit (NC 2273)	<input type="radio"/>

SEAT FRAME*			
JBH2102		HD Telescopic Seat Frame and Plate with Side Rails. Width 530 to 550 mm	○
SEAT WIDTH*			
JBH2124		530 mm Seat Plate	○
JBH2104		550 mm Seat Plate (NC 0556, 0589, 1404, 1405)	○
SEAT DEPTH*			
JBH2136		510 mm	○
JBH2137		580 mm	○
SEAT CUSHION*			
JBH2148		Seat Cushion (90 mm high)	○
JBH2140		Without seat cushion	○
BACKREST SUPPORT*			
JBH2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrix Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	◆
BACKREST UPHOLSTERY / CUSHION*			
JBH2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery. (NC 2230, 2247)	◆
UPHOLSTERY COLOUR*			
JBH2156		Black Nylon Fabric TR18	◆
ARMRESTS*			
JBH2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	○
JBH2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	€50
ARMPADS*			
JBH2195		Short Armpad (325 mm)	◆
JBH2196		Long Armpad (375 mm)	○
JBH2197		Long Ergonomic Armpad (397 mm) (NC 2230, 2247)	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBH1350		Neckrest depth adjustable (NC 2230, 2247)	€103
JBH1432		Headrest sideways adjustable (NC 2230, 2247)	€285
JBH1427		Headrest with cheek support (NC 2230, 2247)	€285
JBH2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware (NC 2230, 2247)	€454
JBH2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware (NC 2230, 2247)	€454
JBH2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware (NC 2230, 2247)	€615
JBH1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
MANUAL / POWERED SEATING OPTIONS*			
JBH1091		Manually settable Seat Tilt 0 - 15° (NC 1401)	○
JBH1092		Powered "fixed pivot" Tilt 20° (NC 2201)	€507
JBH1093		Powered Tilt (0° - 30°) with CoG shift (NC 2201)	€1,348
JBH1094		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected (NC 2201, 2220)	€2,377

SEAT TO GROUND HEIGHT*			
JBH2857		435 mm	○
JBH2858		460 mm	○
JBH2859		485 mm	○
BACKREST RECLINER*			
JBH0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBH0977		Powered Recliner 30° (NC 2201)	€634
LEG REST TYPE*			
JBH0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm (NC 0517, 0512, 0504, 0513, 1244, 1416)	€468
JBH0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm (NC 2201, 0512, 0504, 0513, 1244, 1416)	€2,060
JBH0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools (NC 0559, 0516, 0517)	◆
JBH0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad (NC 2365, 0559, 0516, 0517)	€292
JBH0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad (NC 2201, 2365, 0559, 0516, 0517)	€841
JBH0501		Without Leg Rests (NC 0504, 0513, 2365, 0559, 0516, 0517)	○
FOOTPLATE TYPE*			
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	€130
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	€170
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	€339
JBH0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBH0504		Angle Adjustable Footplate pair	€40
JBH0513		Angle and depth adjustable Footplate pair	€149
JBH0501		Without footplates for option 0501 without Leg Rests (NC 0503, 0533, 0531)	○
LINX CONTROLS*			
JBH2201		DLX-REM110 (Drive Only) (NC 2837, 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908)	◆
JBH2203		DLX-REM211 (Drive & Actuators) (NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908, 1909)	○
JBH2204		DLX-REM216 (Drive, Actuators & Lights) (NC 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258)	€191
JBH2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets) (NC 2220)	€428
JBH2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick). (NC 2220, 2270, 1207, 1208, 1707, 1708, 1709, 1710)	€974
REM400 TOGGLE OPTION			
JBH2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL) (NC 2203, 2204, 2206)	€397

CONTROL BRACKETRY*				
JBH1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	○
JBH1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	○
JBH1707		Swing Away RH Bracket	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	€167
JBH1708		Swing Away LH Bracket	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	€167
JBH1709		Matrx Resolve Swing in/out Bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 1404)	€281
JBH1710		Matrx Resolve Swing in/out Bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 1405)	€281
JBH2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405)	€626
JBH2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405)	€626
JBH2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2203, 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404)	○
JBH2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2203, 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405)	○
JBH2227		Bracket RH mounted for Attendant	(NC 2220, 2128, 2219, 1221, 1312)	€81
JBH2228		Bracket LH mounted for Attendant	(NC 2220, 2128, 2219, 1222, 1363)	€81
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
JBH2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 2273, 2205, 2206, 1221, 1222, 2227, 2228)	€497
JBH2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1221, 1222, 2227, 2228)	€497
JBH2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2259, 2260)	€943
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2259, 2260, 2259, 2260, 1221, 1222)	€3,020
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
JBH2259		Dual Control Standard Bracket RH	(NC 2291, 2245, 1312)	○
JBH2260		Dual Control Standard Bracket LH	(NC 2292, 2245, 1363)	○
JBH1221		IDC Dual Control Detachable Bracket RH	(NC 2291, 2220, 2219, 2227, 1312)	○
JBH1222		IDC Dual Control Detachable Bracket LH	(NC 2292, 2220, 2219, 2228, 1363)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
JBH2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€903
JBH2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€1,129
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€2,630
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0002, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	€4,160
JBH2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215)	€4,537
JBH2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215)	€2,529

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
JBH2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBH2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBH2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	€626
JBH2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	€626
JBZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
JBH2285		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JBH2286		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JBH2283		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	€626
JBH2284		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	€626
JBH2287		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	€626
JBH2288		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	€626
JBZ2289		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2203)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX- IN500 module needed)	(NC 2201, 2203, 2203)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2203)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX- IN500 module needed)	(NC 2201, 2203, 2203)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2203)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2203)	€2,630
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		€1,255
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		€1,820
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		€3,005
HEAD CONTROLS				
JBH2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427)	€4,673
JBH2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427)	€5,472
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427, 0002)	€9,135
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2171, 2197, 1350, 1432, 1427, 0002)	€5,435
FOOT CONTROL				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2203)	€3,800
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2203)	€3,615
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification		POA
TRAY CONTROL				
JBH2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128)	€7,236

TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)			
JBH2275	<input type="checkbox"/>	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415) €626
JBH2276	<input type="checkbox"/>	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339) €626
JBZ2269	<input type="checkbox"/>	Individual mounting option per request (Special). Requires additional information/specification.	POA
WIRELESS ACCESSORIES FOR ASL CONTROLS			
JBH2249	<input type="checkbox"/>	Wireless Triple Switch Receiver (ASL557-3)	€428
JBH2250	<input type="checkbox"/>	Wireless Mouse Emulator (ASL558)	€1,027
LiNX 10-WAY SWITCH			
JBH2251	<input type="checkbox"/>	10-Way Switch DLX-SW10-A for direct Powered Seating Access	€634
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
JBH2217	<input type="checkbox"/>	Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404) ○
JBH2218	<input type="checkbox"/>	Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405) ○
JBH2291	<input type="checkbox"/>	Bracket RH, mounted for attendant	(NC 2259, 1221) €81
JBH2292	<input type="checkbox"/>	Bracket LH, mounted for attendant	(NC 2260, 1222) €81
LiNX USB CHARGER PORT			
JBH2253	<input type="checkbox"/>	USB Charger Port DLX-USB02-A	€207
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
JBH2281	<input type="checkbox"/>	RH Mount	○
JBH2282	<input type="checkbox"/>	LH Mount	○
REMOTE STOP SWITCH			
JBH2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	(NC 2203, 2204) €1,319
LiNX INPUT MODULE			
JBH2273	<input type="checkbox"/>	DLX-IN200 Module (Jack port for external power button)	€303
UNZ9249	<input type="checkbox"/>	DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	(NC 2201, 2203, 2203) €814
FURTHER LiNX COMPONENTS / ACCESSORIES			
JBH2302	<input type="checkbox"/>	Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
JBH2300	<input type="checkbox"/>	To be enclosed only	○
JBZ2301	<input type="checkbox"/>	To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JBH1685	<input type="checkbox"/>	Red Buddy Button	€110
JBH1686	<input type="checkbox"/>	Yellow Buddy Button	€110
JBH1687	<input type="checkbox"/>	Green Buddy Button	€110
JBH1690	<input type="checkbox"/>	Blue Buddy Button	€110
EGG SWITCH			
JBH2254	<input type="checkbox"/>	Egg Switch, black	€110
PIKO BUTTONS			
JBH1711	<input type="checkbox"/>	Red Piko Button	€198
JBH1712	<input type="checkbox"/>	Yellow Piko Button	€198
JBH1713	<input type="checkbox"/>	Green Piko Button	€198
JBH1714	<input type="checkbox"/>	Orange Piko Button	€198
JBH1715	<input type="checkbox"/>	Blue Piko Button	€198
JBH1716	<input type="checkbox"/>	Black Piko Button	€198

BODYPOINT JOYSTICK KNOBS			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	€79
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	€79
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	€79
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€79
5:1 SWITCH ADAPTOR			
JBH1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	€152
SHROUD COLOUR*			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	€269
RIM COLOUR*			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
FORKS*			
JBH2631		Single sided fork	€285
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
CASTORS*			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
CASTOR LOCKS			
JBH1490		Castor Lock	€248
DRIVE WHEELS*			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆

MOTOR TYPE*			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254, 2837) €683
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2201, 2203, 1437, 1254, 2837) €820
BATTERY TYPE*			
JBH1259		60 Ah Gel	◆
JBH1260		73.5 Ah Gel	(NC 1450) €116
G TRAC			
JBH1215		G-Trac Module	€406
CHARGER*			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	€191
ANTI-TIPPERS*			
JBH1205		Anti-tipper	◆
LIGHTS			
JBH1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors	(NC 2203, 1908, 1909) €602
JBH2837		LED Light system extra bright with indicators 60 Lux headlights	(NC 2203) €180
JBH1254		Without Light system	(NC 1909) ○
KERB RIDER			
JBH1244		Kerb Rider with park position feature	(NC 0556, 0589, 1294) €266
TRANSPORTATION KIT*			
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	€736
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 1091, 1092, 1093) ○
FIXATION BELT*			
JBH1329		Lap belt with Auto Lock Standard mount on backrest	○
JBH2113		Lap belt with Auto Lock Flexible mount in side rails	○
REAR VIEW MIRROR			
JBH1415		Rear view mirror RH	(NC 2217, 2285, 2209, 2213, 2265, 2275, 2283, 2287) €53
JBH1339		Rear view mirror LH	(NC 2218, 2286, 2210, 2214, 2266, 2276, 2284, 2288) €53
BAG HOOKS			
JBH2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2227, 2228, 2233, 2240, 2257, 2291, 2292) €97
CANE HOLDER			
JBH1312		Cane holder RH	(NC 2227, 2259, 1221) €53
JBH1363		Cane holder LH	(NC 2228, 2260, 1222) €53
WHEEL FENDERS			
JBH1416		Fender for Castors	€12
JBH1417		Fender for Drive wheels	€213

TRAY			
JBH1404		Tray, Swing Away RH	(NC 1707, 2209, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217) €210
JBH1405		Tray, Swing Away LH	(NC 1708, 2210, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218) €210

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBH0002	AVIVA RX40 HD Modulite HD
JBH2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
JBH2197	Long Ergonomic Armpad (397 mm)
JBH1350	Neckrest depth adjustable
JBH1432	Headrest sideways adjustable
JBH1427	Headrest with cheek support
JBH2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
JBH2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
JBH2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
JBH1091	Manually settable Seat Tilt 0 - 15°
JBH1092	Powered "fixed pivot" Tilt 20°
JBH1093	Powered Tilt (0° - 30°) with CoG shift
JBH1094	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected
JBH0977	Powered Recliner 30°
JBH0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
JBH0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
JBH0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
JBH0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBH2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JBH0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
JBH0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0512	Fixed Footplate pair
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH2201	DLX-REM110 (Drive Only)
JBH2203	DLX-REM211 (Drive & Actuators)
JBH2204	DLX-REM216 (Drive, Actuators & Lights)
JBH2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JBH2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JBH2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JBH1207	Standard Bracket RH
JBH1208	Standard Bracket LH
JBH1707	Swing Away RH Bracket
JBH1708	Swing Away LH Bracket
JBH1709	Matrx Resolve Swing in/out Bracket RH
JBH1710	Matrx Resolve Swing in/out Bracket LH
JBH2213	Centre Mount Swing Away to RH Nucleus Bracket

JBH2214	Centre Mount Swing Away to LH Nucleus Bracket
JBH2209	REM500 Swing Away RH Multi-Adjustable Bracket
JBH2210	REM500 Swing Away LH Multi-Adjustable Bracket
JBH2227	Bracket RH mounted for Attendant
JBH2228	Bracket LH mounted for Attendant
JBH2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
JBH2219	LiNX DLX-ACU200 (Compact Dual Control)
JBH2245	LiNX DLX-IDC (Intuitive Dual Control)
JBH2259	Dual Control Standard Bracket RH
JBH2260	Dual Control Standard Bracket LH
JBH1221	IDC Dual Control Detachable Bracket RH
JBH1222	IDC Dual Control Detachable Bracket LH
JBH2231	LiNX DLX-CR400 (Compact Remote)
JBH2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JBH2233	ASL MEC (Proportional MEC Joystick ASL130)
JBH2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount Swing Away to RH (Nucleus)
JBH2266	Centre Mount Swing Away to LH (Nucleus)
JBZ2215	Individual mounting option per request (Special). Requires additional information/specification.
JBH2285	Tray Mount RH (ASL603, 100mm x 150mm)
JBH2286	Tray Mount LH (ASL603, 100mm x 150mm)
JBH2283	Centre Mount Swing Away to RH (Nucleus)
JBH2284	Centre Mount Swing Away to LH (Nucleus)
JBH2287	Centre Mount, Swing Away to RH, Nucleus (Incl Tray ASL603, 100mm x 150mm)
JBH2288	Centre Mount, Swing Away to LH, Nucleus (Incl Tray ASL603, 100mm x 150mm)
JBZ2289	Individual mounting option per request (Special). Requires additional information/specification.
JBH2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
JBH2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
JBH2257	Four-Switch Proximity Array (ASL106)
JBH2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2258	Remote Stop Switch (ASL504)
JBH2273	DLX-IN200 Module (Jack port for external power button)
JBH1686	Yellow Buddy Button
JBH1687	Green Buddy Button
JBH1690	Blue Buddy Button
JBH2254	Egg Switch, black
JBH1711	Red Piko Button

JBH1712	Yellow Piko Button
JBH1713	Green Piko Button
JBH1714	Orange Piko Button
JBH1715	Blue Piko Button
JBH1716	Black Piko Button
JBH1294	Indoor setting, reduced speed approx. 3 km/h
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBH2837	LED Light system extra bright with indicators 60 Lux headlights
JBH1254	Without Light system
JBH1244	Kerb Rider with park position feature
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.
JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):