



Yes, you can:

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Power Retail Prescription Form

July 2026  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX40 HD Modulite

### Standard Features

- Max User Weight 160kg (25 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm - 510mm (16" - 20"), stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short amp pads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £7,250**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
JBH0002	AVIVA RX40 HD Modulite (NC 2193)	£7,250
JBZ1300	Special AVIVA RX40 HD Modulite (Extended Lead Time). Please provide details in comments section.	POA
<b>CHASSIS WIDTH*</b>		
JBH2830	Wide (635 mm/25")	<input checked="" type="checkbox"/>
<b>SEATING TYPE*</b>		
JBH2100	Modulite Seat Unit	<input type="radio"/>

<b>SEAT FRAME*</b>			
JBH2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	○
JBH2112		Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.	○
JBX2102		Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	(NC 2122, 2123) £110
<b>SEAT WIDTH*</b>			
JBH2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	(NC 2174, 2175, 2176, 2148) ○
JBH2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	○
JBH2123		480 mm Seat Plate	(NC 2702, 2704, 2706, 2148) ○
JBH2124		530 mm Seat Plate	(NC 2174, 2175, 2176, 2177, 2178, 2702, 2704, 2706, 2106, 2148) ○
JBH2104		550 mm Seat Plate	(NC 2132, 2134, 2142, 2144, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2192, 2193, 2194, 2198, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 1445, 0556, 0589, 2701, 2702, 2703, 2704, 2705, 2706, 1404, 1405) ○
<b>BODY FIT KIT</b>			
JBH2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)	○
<b>SEAT DEPTH*</b>			
JBH2132		410 mm	(NC 2196, 2197, 2198) ○
JBH2134		460 mm	(NC 2148) ○
JBH2136		510 mm	○
<b>SEAT CUSHION*</b>			
JBH2142		Matrx Flo-tech Lite Cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) £120
JBH2147		Matrx Libra Cushion, 90 mm thick (Extended lead time)	£410
JBH2144		Matrx Vi Cushion (Mesh fabric)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983) £400
JBH2148		Standard Cushion (90mm high) - For 550mm Seat Plate Only	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) ○
JBH2140		Without seat cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) ○
<b>1st POSITIONING ACCESORIES FOR MATRX LIBRA</b>			
JBH1980		Tapered Wedges (Pair)	£30
JBH1981		Abductor	£23
JBH1982		Lateral Wedge	£20
JBH1983		Anterior Seat Wedges	£27
<b>2nd POSITIONING ACCESORIES FOR MATRX LIBRA</b>			
JBH1980		Tapered Wedges (Pair)	£30
JBH1981		Abductor	£23
JBH1982		Lateral Wedge	£20
JBH1983		Anterior Seat Wedges	£27
<b>POSITIONING ACCESSORIES FOR MATRX VI</b>			
JBH1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
<b>BACKREST SUPPORT*</b>			
JBH2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	(NC 2172, 2173) ○
JBH2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	◆

<b>PUSH HANDLES FOR BACKPOSTS</b>			
JBH2108		Backposts without Push Handles	(NC 2219, 2227, 2228, 2245, 2291, 2292) ○
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBH2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 2155, 2230, 2247) ◆
JBH2179		Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247) £155
JBH2172		Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128) £250
JBH2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128) £400
JBH2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 2121) £353
JBH2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 2121) £353
JBH2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 2121) £494
JBH2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£558
JBH2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£546
JBH2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) £475
JBH2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) £575
JBH2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2230, 2247, 2151, 2152, 2155) ○
<b>MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)</b>			
JBH2380		405 mm (16 inch) MAC BH	(NC 1445) ○
JBH2381		455 mm (18 inch) MAC BH	(NC 1445) ○
JBH2382		505 mm (20 inch) MAC BH	(NC 1445) ○
JBH2383		560 mm (22 inch) MAC BH	(NC 1445) ○
<b>MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)</b>			
JBH2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445) £100
JBH2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445) £100
JBH2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445) £100
<b>MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)</b>			
JBH2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445) £100
JBH2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445) £100
JBH2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445) £100
<b>MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)</b>			
JBH2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445) ○
JBH2394		Backrest matches seat width	(NC 1445) ○
JBH2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445) ○
<b>UPHOLSTERY COLOUR*</b>			
JBH2151		Black Nylon Fabric TR18	◆
JBH2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171) ○
JBH2155		Black Spacer Fabric TR36	(NC 2171) ○
JBH2156		Black Nylon (For 550mm SW Only)	○
JBX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) ○

<b>ARMRESTS*</b>			
JBH2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	◆
JBH2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	£40
JBH2198		Padded Armrest with padded side guard (Water Fall Type)	(NC 2132, 2195, 2196, 2197, 2199) £120
JBH2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).	(NC 1445, 2179) £65
JBH2192		Following Armrest. Height adjustable from 230 - 300 mm	(NC 1445, 2179) £185
JBH2194		Following Armrest. Height adjustable from 300 - 360 mm	(NC 1445, 2179) £185
<b>ARMPADS*</b>			
JBH2195		Short Armpad (325 mm)	(NC 2198) ◆
JBH2196		Long Armpad (375 mm)	(NC 2132, 2198) ○
JBH2197		Long Ergonomic Armpad (397 mm)	(NC 2132, 2198, 2230, 2247) ○
JBH2199		Soft Armpad (REA Azalea Wide Pad)	(NC 2198) £95
JBX3000		Please select with Option 2198 - Padded Armrest with padded side guard (Water Fall Type)	(NC 2181, 2193, 2192) ○
<b>HIP SUPPORTS LH</b>			
JBH2184		Hip Support LH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £250
JBH2186		Hip Support LH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £250
JBH2188		Hip Support LH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £250
<b>HIP SUPPORTS RH</b>			
JBH2185		Hip Support RH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £250
JBH2187		Hip Support RH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £250
JBH2189		Hip Support RH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £250
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
JBH2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JBH2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○
<b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b>			
JBH1350		Neckrest depth adjustable	(NC 2230, 2247) £80
JBH1432		Headrest sideways adjustable	(NC 2230, 2247) £225
JBH1427		Headrest with cheek support	(NC 2230, 2247) £225
JBH2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	(NC 2230, 2247, 2170) £325
JBH2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	(NC 2230, 2247, 2170) £325
JBH2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware	(NC 2230, 2247, 2170) £440
JBH1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2170, 2171, 2174, 2175, 2176, 2177, 2178) ○
JBH1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
<b>LATERAL SUPPORTS</b>			
JBH1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194) £165

<b>MANUAL / POWERED SEATING OPTIONS*</b>			
JBH1091		Manually settable Seat Tilt 0 - 15°	(NC 1401) ○
JBH1092		Powered "fixed pivot" Tilt 20°	(NC 2201) £400
JBH1093		Powered Tilt (0° - 30°) with CoG shift	(NC 2201) £1,100
JBH1094		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected	(NC 2201) £1,900
<b>SEAT TO GROUND HEIGHT*</b>			
JBH2857		435 mm	○
JBH2858		460 mm	○
JBH2859		485 mm	○
<b>BACKREST RECLINER*</b>			
JBH0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBH0977		Powered Recliner 30°	£500
<b>LEG REST TYPE*</b>			
JBH0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 1244, 1416) £370
JBH0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416) £1,620
JBH0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0559, 0516, 0517, 0501) ◆
JBH0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 2365, 0559, 0516, 0517, 0501) £235
JBH0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0501) £660
JBH0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517) ○
<b>FOOTPLATE TYPE*</b>			
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£105
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	£135
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBH0512		Fixed Footplate pair	(NC 0533, 0531) ◆
JBH0504		Angle Adjustable Footplate pair	£75
JBH0513		Angle and depth adjustable Footplate pair	£120
JBH0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531) ○
<b>LINX CONTROLS*</b>			
JBH2201		DLX-REM110 (Drive Only)	(NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908) ◆
JBH2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 1908, 1909) ○
JBH2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2270) £200
JBH2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	£425
JBH2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 2270, 1709, 1710) £765
<b>REM400 TOGGLE OPTION</b>			
JBH2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £315

<b>CONTROL BRACKETRY*</b>				
JBH1207		Standard Bracket RH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	○
JBH1208		Standard Bracket LH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	○
JBH1707		Swing Away RH Bracket	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	£135
JBH1708		Swing Away LH Bracket	(NC 2283, 2287, 2207, 2265, 2275, 1404)	£135
JBH1709		Matrx Resolve Swing in/out Bracket RH	(NC 2284, 2288, 2208, 2266, 2276, 1405)	£245
JBH1710		Matrx Resolve Swing in/out Bracket LH	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405)	£245
JBH2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405)	£495
JBH2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 2203, 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404)	£495
JBH2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2203, 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405)	○
JBH2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2128, 2108, 2219, 1221, 1312)	○
JBH2227		Bracket RH mounted for Attendant	(NC 2128, 2108, 2219, 1222, 1363)	£65
JBH2228		Bracket LH mounted for Attendant		£65
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228)	POA
<b>DUAL CONTROLS</b>				
JBH2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 2108, 1221, 1222, 2227, 2228)	£395
JBH2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2108, 2259, 2260)	£745
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2108, 2259, 2260, 0002, 2259, 2260, 1221, 1222)	£2,315
<b>DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)</b>				
JBH2259		Bracket RH, Attendant	(NC 2291, 2245, 1312)	○
JBH2260		Bracket LH, Attendant	(NC 2292, 2245, 1363)	○
JBH1221		Bracket RH, detachable	(NC 2291, 2219, 2227, 1312)	○
JBH1222		Bracket LH, detachable	(NC 2292, 2219, 2228, 1363)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
<b>COMPACT REMOTES</b>				
JBH2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£825
JBH2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289)	£3,190
JBH2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215)	£4,100
JBH2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215)	£2,285

<b>COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)</b>				
JBH2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBH2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBH2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBH2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
JBH2285		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JBH2286		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JBH2283		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBH2284		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBH2287		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBH2288		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2289		Individual mounting option per request (Special). Requires additional information/specification.		POA
<b>CHIN CONTROLS</b>				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2204)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX- IN500 module needed)	(NC 2201, 2203, 2204)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2204)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX- IN500 module needed)	(NC 2201, 2203, 2204)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204)	£2,015
<b>CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)</b>				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		£2,280
<b>HEAD CONTROLS</b>				
JBH2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£4,220
JBH2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£3,315
<b>FOOT CONTROL</b>				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2204)	£2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204)	£2,740
<b>FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)</b>				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification		POA
<b>TRAY CONTROL</b>				
JBH2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128)	£6,535

<b>TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)</b>			
JBH2275	<input type="checkbox"/>	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415) £495
JBH2276	<input type="checkbox"/>	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339) £495
JBZ2269	<input type="checkbox"/>	Individual mounting option per request (Special) Requires additional information/specification	POA
<b>WIRELESS ACCESSORIES FOR ASL CONTROLS</b>			
JBH2249	<input type="checkbox"/>	Wireless Triple Switch Receiver (ASL557-3)	£340
JBH2250	<input type="checkbox"/>	Wireless Mouse Emulator (ASL558)	£810
<b>LINX 10-WAY SWITCH</b>			
JBH2251	<input type="checkbox"/>	10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
<b>LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)</b>			
JBH2217	<input type="checkbox"/>	Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404) ○
JBH2218	<input type="checkbox"/>	Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405) ○
JBH2291	<input type="checkbox"/>	Bracket RH, mounted for attendant	(NC 2259, 1221, 2108) £65
JBH2292	<input type="checkbox"/>	Bracket LH, mounted for attendant	(NC 2260, 1222, 2108) £65
<b>LINX USB CHARGER PORT</b>			
JBH2253	<input type="checkbox"/>	USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £160
<b>LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)</b>			
JBH2281	<input type="checkbox"/>	RH Mount	○
JBH2282	<input type="checkbox"/>	LH Mount	○
<b>LINX INPUT MODULE</b>			
JBH2273	<input type="checkbox"/>	DLX-IN200 Module (Jack port for external power button)	£270
UNZ9249	<input type="checkbox"/>	DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	(NC 2201, 2203, 2203) £785
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
JBH2302	<input type="checkbox"/>	Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
JBH2300	<input type="checkbox"/>	To be enclosed only	○
JBZ2301	<input type="checkbox"/>	To be connected and programmed (Special)	POA
<b>BUDDY BUTTONS</b>			
JBH1685	<input type="checkbox"/>	Red Buddy Button	£90
JBH1686	<input type="checkbox"/>	Yellow Buddy Button	£90
JBH1687	<input type="checkbox"/>	Green Buddy Button	£90
JBH1690	<input type="checkbox"/>	Blue Buddy Button	£90
<b>EGG SWITCH</b>			
JBH2254	<input type="checkbox"/>	Egg Switch, black	£90
<b>PIKO BUTTONS</b>			
JBH1711	<input type="checkbox"/>	Red Piko Button	£160
JBH1712	<input type="checkbox"/>	Yellow Piko Button	£160
JBH1713	<input type="checkbox"/>	Green Piko Button	£160
JBH1714	<input type="checkbox"/>	Orange Piko Button	£160
JBH1715	<input type="checkbox"/>	Blue Piko Button	£160
JBH1716	<input type="checkbox"/>	Black Piko Button	£160

<b>BODYPOINT JOYSTICK KNOBS</b>			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
<b>5:1 SWITCH ADAPTOR</b>			
JBH1681		5:1 Switch Adapter ( 9-pin sub D to 5 jack ports)	£120
<b>SHROUD COLOUR*</b>			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300
<b>RIM COLOUR*</b>			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○ (NC 2872, 2873, 2874, 2875)
JBZ2887		Premium - Falcon White Rim Inserts	○ (NC 2872, 2873, 2874, 2875, 0299)
<b>FORKS*</b>			
JBH2631		Single sided fork	£250
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
<b>CASTORS*</b>			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
<b>CASTOR LOCKS</b>			
JBH1490		Castor Lock	£195
<b>DRIVE WHEELS*</b>			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
<b>C.T.C SUSPENSION*</b>			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
<b>SUSPENSION SETTING*</b>			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆

<b>MOTOR TYPE*</b>			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254) £515
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2201, 2203, 1437, 1254) £600
<b>BATTERY TYPE*</b>			
JBH1259		60 Ah Gel	◆
JBH1260		73.5 Ah Gel	(NC 1450) £100
<b>G TRAC</b>			
JBH1215		G-Trac Module	£400
<b>CHARGER*</b>			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	£155
<b>ANTI-TIPPERS*</b>			
JBH1205		Anti-tipper	◆
<b>LIGHTS</b>			
JBH1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors	(NC 2203, 1908, 1909) £450
JBH1254		Without Light system	(NC 1909) ○
<b>KERB RIDER</b>			
JBH1244		Kerb Rider with park position feature	(NC 0556, 0589, 1294) £215
<b>TRANSPORTATION KIT*</b>			
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 1091, 1092, 1093) ○
<b>EASY ACCESS FRONT TIE DOWNS</b>			
JBH2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
<b>FIXATION BELT*</b>			
JBH1329		Lap belt with Auto Lock Standard mount on backrest	○
JBH2113		Lap belt with Auto Lock Flexible mount in side rails	○
JBH2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBH2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBH2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBH2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBH2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBH2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90

<b>REAR VIEW MIRROR</b>			
JBH1415		Rear view mirror RH	(NC 2217, 2285, 2209, 2213, 2265, 2275, 2283, 2287) £45
JBH1339		Rear view mirror LH	(NC 2218, 2286, 2210, 2214, 2266, 2276, 2284, 2288) £45
<b>BAG HOOKS</b>			
JBH2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2172, 2173, 2227, 2228, 2233, 2240, 2257, 2291, 2292) £80
<b>CANE HOLDER</b>			
JBH1312		Cane holder RH	(NC 2227, 2259, 1221) £45
JBH1363		Cane holder LH	(NC 2228, 2260, 1222) £45
<b>WHEEL FENDERS</b>			
JBH1416		Fender for Castors	£80
JBH1417		Fender for Drive wheels	£160
<b>TRAY</b>			
JBH1404		Tray, Swing Away RH	(NC 1707, 2209, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217) £200
JBH1405		Tray, Swing Away LH	(NC 1708, 2210, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218) £200

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

JBH0002	AVIVA RX40 HD Modulite
JBH2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
JBH2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
JBH2123	480 mm Seat Plate
JBH2124	530 mm Seat Plate
JBH2106	Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)
JBH2132	410 mm
JBH2134	460 mm
JBH2142	Matrx Flo-tech Lite Cushion
JBH2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)
JBH2144	Matrx Vi Cushion (Mesh fabric)
JBH2148	Standard Cushion (90mm high) - For 550mm Seat Plate Only
JBH1980	Tapered Wedges (Pair)
JBH1981	Abductor
JBH1982	Lateral Wedge
JBH1983	Anterior Seat Wedges
JBH1980	Tapered Wedges (Pair)
JBH1981	Abductor
JBH1982	Lateral Wedge
JBH1983	Anterior Seat Wedges
JBH1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
JBH2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.
JBH2108	Backposts without Push Handles
JBH2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
JBH2179	Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).
JBH2172	Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
JBH2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
JBH2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBH2380	405 mm (16 inch) MAC BH
JBH2381	455 mm (18 inch) MAC BH
JBH2382	505 mm (20 inch) MAC BH

JBH2383	560 mm (22 inch) MAC BH
JBH2385	90 mm (3.5 inch) MAC Contour Depth LH
JBH2386	150 mm (6 inch) MAC Contour Depth LH
JBH2387	230 mm (9 inch) MAC Contour Depth LH
JBH2389	90 mm (3.5 inch) MAC Contour Depth RH
JBH2390	150 mm (6 inch) MAC Contour Depth RH
JBH2391	230 mm (9 inch) MAC Contour Depth RH
JBH2393	Backrest -25 mm (-1 inch) to seat width
JBH2394	Backrest matches seat width
JBH2395	Backrest +25 mm (+1 inch) to seat width
JBH2151	Black Nylon Fabric TR18
JBH2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
JBH2155	Black Spacer Fabric TR36
JBH2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
JBH2198	Padded Armrest with padded side guard (Water Fall Type)
JBH2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
JBH2192	Following Armrest. Height adjustable from 230 - 300 mm
JBH2194	Following Armrest. Height adjustable from 300 - 360 mm
JBH2195	Short Armpad (325 mm)
JBH2196	Long Armpad (375 mm)
JBH2197	Long Ergonomic Armpad (397 mm)
JBH2199	Soft Armpad (REA Azalea Wide Pad)
JBH2184	Hip Support LH, pad size 140 x 100 mm
JBH2186	Hip Support LH, pad size 230 x 100 mm
JBH2188	Hip Support LH, pad size 330 x 100 mm
JBH2185	Hip Support RH, pad size 140 x 100 mm
JBH2187	Hip Support RH, pad size 230 x 100 mm
JBH2189	Hip Support RH, pad size 330 x 100 mm
JBH1350	Neckrest depth adjustable
JBH1432	Headrest sideways adjustable
JBH1427	Headrest with cheek support
JBH2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
JBH2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
JBH2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
JBH1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBH1445	Lateral Support with extended arm, multi-adjustable pad
JBH1091	Manually settable Seat Tilt 0 - 15°
JBH1092	Powered "fixed pivot" Tilt 20°
JBH1093	Powered Tilt (0° - 30°) with CoG shift
JBH1094	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected
JBH0977	Powered Recliner 30°
JBH0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
JBH0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm

JBH0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
JBH0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBH0501	Without Leg Rests
JBH2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JBH0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
JBH0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0512	Fixed Footplate pair
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH0501	Without footplates for option 0501 without Leg Rests
JBH2201	DLX-REM110 (Drive Only)
JBH2203	DLX-REM211 (Drive & Actuators)
JBH2204	DLX-REM216 (Drive, Actuators & Lights)
JBH2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JBH2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JBH2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JBH1707	Swing Away RH Bracket
JBH1708	Swing Away LH Bracket
JBH1709	Matrx Resolve Swing in/out Bracket RH
JBH1710	Matrx Resolve Swing in/out Bracket LH
JBH2213	Centre Mount Swing Away to RH Nucleus Bracket
JBH2214	Centre Mount Swing Away to LH Nucleus Bracket
JBH2209	REM500 Swing Away RH Multi-Adjustable Bracket
JBH2210	REM500 Swing Away LH Multi-Adjustable Bracket
JBH2227	Bracket RH mounted for Attendant
JBH2228	Bracket LH mounted for Attendant
JBH2219	LiNX DLX-ACU200 (Compact Dual Control)
JBH2245	LiNX DLX-IDC (Intuitive Dual Control)
JBH2259	Bracket RH, Attendant
JBH2260	Bracket LH, Attendant
JBH1221	Bracket RH, detachable
JBH1222	Bracket LH, detachable
JBH2231	LiNX DLX-CR400 (Compact Remote)
JBH2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JBH2233	ASL MEC (Proportional MEC Joystick ASL130)
JBH2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount Swing Away to RH (Nucleus)

JBH2266	Centre Mount Swing Away to LH (Nucleus)
JBZ2215	Individual mounting option per request (Special). Requires additional information/specification.
JBH2285	Tray Mount RH (ASL603, 100mm x 150mm)
JBH2286	Tray Mount LH (ASL603, 100mm x 150mm)
JBH2283	Centre Mount Swing Away to RH (Nucleus)
JBH2284	Centre Mount Swing Away to LH (Nucleus)
JBH2287	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JBH2288	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JBZ2289	Individual mounting option per request (Special). Requires additional information/specification.
JBH2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
JBH2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
JBH2257	Four-Switch Proximity Array (ASL106)
JBH2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2253	USB Charger Port DLX-USB02-A
JBH2273	DLX-IN200 Module (Jack port for external power button)
JBH1686	Yellow Buddy Button
JBH1687	Green Buddy Button
JBH1690	Blue Buddy Button
JBH2254	Egg Switch, black
JBH1711	Red Piko Button
JBH1712	Yellow Piko Button
JBH1713	Green Piko Button
JBH1714	Orange Piko Button
JBH1715	Blue Piko Button
JBH1716	Black Piko Button
JBH2872	Mid range - Silver Shroud
JBH2873	Mid range - Bengal Red Shroud
JBH2874	Mid range - Mercury Blue Shroud
JBH2875	Mid range - Midnight Purple Shroud
JBZ0299	Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.
JBH1294	Indoor setting, reduced speed approx. 3 km/h
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBH1254	Without Light system
JBH1244	Kerb Rider with park position feature
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.

JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JBH2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBH2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBH2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBH2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBH2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBH2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):