



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

July 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Modulite

Standard Features

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths from 410mm - 510mm (16" - 20"), stepless adjustment
- Tension adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short amp pads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM110 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £7,250

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBG0001	AVIVA RX40 Modulite	£7,250
JBZ1300	Special AVIVA RX40 Modulite (Extended Lead Time). Please provide details in comments section.	POA
CHASSIS WIDTH*		
JBG2810	Narrow (560 mm/22")	◆
SEATING TYPE*		
JBG2100	Modulite Seat Unit	◆

SEAT FRAME*			
JBG2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	○
JBG2112		Telescopic Seat Frame + Sling Seat with Side Rails. Width 380 mm to 530 mm. +30 mm seat height. Additional 50 mm width by arm rest adjustment.	○
JBX2102		Telescopic Seat Frame and Plate with Side Rails width 530 to 550 mm	(NC 2122, 2123) £110
SEAT WIDTH*			
JBG2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	(NC 2174, 2175, 2176, 2148) ○
JBG2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	○
JBG2123		480 mm Seat Plate	(NC 2704, 2706, 2148) ○
JBG2124		530 mm Seat Plate	(NC 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2106, 2148) ○
JBG2104		550 mm Seat Plate	(NC 2132, 2134, 2142, 2143, 2144, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2183, 2192, 2193, 2194, 2198, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1496, 1445, 0556, 0589, 0526, 2701, 2702, 2703, 2704, 2705, 2706, 1404, 1405) ○
BODY FIT KIT			
JBG2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)	(NC 2192, 2194) ○
SEAT DEPTH*			
JBG2132		410 mm	(NC 2197, 2198) ○
JBG2134		460 mm	○
JBG2136		510 mm	○
SEAT CUSHION*			
JBG2142		Matrx Flo-tech Lite Cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) £120
JBG2143		Matrx Flo-tech Contour Visco; (2-way stretch material)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) £169
JBG2147		Matrx Libra Cushion, 90 mm thick (Extended lead time)	£410
JBG2144		Matrx Vi Cushion (Mesh fabric)	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983) £400
JBG2148		Standard Cushion (90mm high) - For 550mm Seat Plate Only	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) ○
JBG2140		Without seat cushion	(NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985) ○
1st POSITIONING ACCESORIES FOR MATRX LIBRA			
JBG1980		Tapered Wedges (Pair)	£30
JBG1981		Abductor	£23
JBG1982		Lateral Wedge	£20
JBG1983		Anterior Seat Wedges	£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
JBG1980		Tapered Wedges (Pair)	£30
JBG1981		Abductor	£23
JBG1982		Lateral Wedge	£20
JBG1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
JBG1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54

BACKREST SUPPORT*				
JBG2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.	(NC 2173, 2052, 2053)	○
JBG2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. Push Handles included.		○
PUSH HANDLES FOR BACKPOSTS				
JBG2108		Backposts without Push Handles	(NC 2219, 2220, 2227, 2228, 2245, 2291, 2292)	○
BACKREST UPHOLSTERY / CUSHION*				
JBG2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 2230, 2247, 2052, 2053)	◆
JBG2179		Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).	(NC 2192, 2193, 2194, 2199, 2155, 2230, 2247)	£155
JBG2172		Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2128)	£250
JBG2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2128)	£400
JBG2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen		£353
JBG2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen		£353
JBG2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen		£494
JBG2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen		£558
JBG2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen		£546
JBG2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445)	£475
JBG2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445)	£575
JBG2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2230, 2247, 2151, 2152, 2155)	○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
JBG2380		405 mm (16 inch) MAC BH	(NC 1445)	○
JBG2381		455 mm (18 inch) MAC BH	(NC 1445)	○
JBG2382		505 mm (20 inch) MAC BH	(NC 1445)	○
JBG2383		560 mm (22 inch) MAC BH	(NC 1445)	○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
JBG2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445)	£100
JBG2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445)	£100
JBG2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445)	£100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
JBG2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445)	£100
JBG2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445)	£100
JBG2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445)	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
JBG2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
JBG2394		Backrest matches seat width	(NC 1445)	○
JBG2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○

UPHOLSTERY COLOUR*			
JBG2151		Black Nylon Fabric TR18	◆
JBG2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	○
JBG2155		Black Spacer Fabric TR36 (NC 2171)	○
JBG2156		Black Nylon (For 550mm SW Only)	○
JBX1002		No Upholstery Colour required - No backrest upholstery required (NC 2179, 2173, 2174, 2175, 2176, 2177, 2178)	○
ARMRESTS*			
JBG2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	◆
JBG2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	£40
JBG2198		Padded Armrest with padded side guard (Water Fall Type) (NC 2195, 2196, 2197, 2199)	£120
JBG2183		Inrigged Armrests. Reduces SW by 25 mm each side	£30
JBG2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). (NC 1445, 2179)	£65
JBG2192		Following Armrest. Height adjustable from 230 - 300 mm (NC 2179)	£185
JBG2194		Following Armrest. Height adjustable from 300 - 360 mm (NC 2179)	£185
ARMPADS*			
JBG2195		Short Armpad (325 mm)	◆
JBG2196		Long Armpad (375 mm) (NC 2198)	○
JBG2197		Long Ergonomic Armpad (397 mm) (NC 2198, 2230, 2247)	○
JBG2199		Soft Armpad (REA Azalea Wide Pad)	£95
JBG3000		Please select with Option 2198 - Padded Armrest with padded side guard (Water Fall Type) (NC 2181, 2183, 2193, 2192)	○
HIP SUPPORTS LH			
JBG2184		Hip Support LH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBG2186		Hip Support LH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBG2188		Hip Support LH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
HIP SUPPORTS RH			
JBG2185		Hip Support RH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBG2187		Hip Support RH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
JBG2189		Hip Support RH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS			
JBG2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JBG2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○

HEADREST / NECKREST (Adequate headrest option must be selected)*				
JBG1350		Neckrest depth adjustable	(NC 2247)	£80
JBG1432		Headrest sideways adjustable	(NC 2247)	£225
JBG1427		Headrest with cheek support	(NC 2247)	£225
JBG2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	(NC 2247, 2170)	£325
JBG2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	(NC 2247, 2170)	£325
JBG2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware	(NC 2247, 2170)	£440
JBG1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171, 2174, 2175, 2176, 2177, 2178)	○
JBG1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		○
LATERAL SUPPORTS				
JBG1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£165
MANUAL / POWERED SEATING OPTIONS*				
JBG1091		Manually settable Seat Tilt 0 - 15°		○
JBG1092		Powered "fixed pivot" Tilt 20°	(NC 2201)	£400
JBG1093		Powered Tilt (0° - 30°) with CoG shift	(NC 2201)	£1,100
JBG2052		Powered Tilt (0° - 50°) with CoG shift	(NC 2201, 2171)	£1,450
JBG1094		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected	(NC 2201, 2220)	£1,900
JBG2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG. 6 point Tie downs must be selected	(NC 2201, 2220, 1910, 1912, 2171)	£2,500
SEAT TO GROUND HEIGHT*				
JBG2857		435 mm		○
JBG2858		460 mm		○
JBG2859		485 mm		○
BACKREST RECLINER*				
JBG0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°		◆
JBG0977		Powered Recliner 30°		£500
LEG REST TYPE*				
JBG0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 1244, 1416, 0501)	£370
JBG0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 1244, 1416, 0501)	£1,620
JBG0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 2365, 0559, 0516, 0517, 0501)	◆
JBG0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0559, 0516, 0517, 0554, 0555, 0501)	£120
JBG0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£235
JBG0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£660
JBG0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517)	○

FOOTPLATE TYPE*			
JBG2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£105
JBG0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBG0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)	£135
JBG0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBG0512		Fixed Footplate pair	(NC 0533, 0531) ◆
JBG0504		Angle Adjustable Footplate pair	(NC 0526) £75
JBG0513		Angle and depth adjustable Footplate pair	£120
JBG0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531) ○
AMPUTEE SUPPORT			
JBG0554		Amputee Support RH	£195
JBG0555		Amputee Support LH	£195
LINX CONTROLS*			
JBG2201		DLX-REM110 (Drive Only)	(NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908) ◆
JBG2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 1908, 1909) ○
JBG2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2270) £200
JBG2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2220) £425
JBG2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 2220, 2270, 1709, 1710) £765
REM400 TOGGLE OPTION			
JBG2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £315
CONTROL BRACKETRY*			
JBG1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBG1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBG1707		Swing Away RH Bracket	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) £135
JBG1708		Swing Away LH Bracket	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) £135
JBG1709		Matrx Resolve Swing in/out Bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 1404) £245
JBG1710		Matrx Resolve Swing in/out Bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 1405) £245
JBG2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) £495
JBG2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) £495
JBG2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2203, 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404) ○
JBG2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2203, 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405) ○
JBG2227		Bracket RH mounted for Attendant	(NC 2220, 2128, 2108, 2219, 1221, 1312) £65
JBG2228		Bracket LH mounted for Attendant	(NC 2220, 2128, 2108, 2219, 1222, 1363) £65
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
DUAL CONTROLS			
JBG2220		LiNX DLX-REM050 (Dual Control used as secondary remote)	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228, 9259) £395
JBG2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 2108, 1221, 1222, 2227, 2228, 9259) £395
JBG2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2108, 2259, 2260, 9259) £745
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2108, 2259, 2260, 0001, 2259, 2260, 1221, 1222) £2,315

DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
JBG2259		Dual Control Standard Bracket RH	(NC 2291, 2245, 1312)	○
JBG2260		Dual Control Standard Bracket LH	(NC 2292, 2245, 1363)	○
JBG1221		IDC Dual Control Detachable Bracket RH	(NC 2291, 2220, 2219, 2227, 1312)	○
JBG1222		IDC Dual Control Detachable Bracket LH	(NC 2292, 2220, 2219, 2228, 1363)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
JBG2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£825
JBG2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288, 2289, 2220)	£3,190
JBG2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215, 2220)	£4,100
JBG2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128, 2207, 2208, 2265, 2266, 2215, 2220)	£2,285
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
JBG2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBG2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBG2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBG2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
JBG2285		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JBG2286		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JBG2283		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBG2284		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBG2287		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JBG2288		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JBZ2289		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220)	£2,015

CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)			
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)	£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)	£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)	£2,280
HEAD CONTROLS			
JBG2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £4,220
JBG2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0001, 2220) £4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0001, 2220) £3,315
FOOT CONTROL			
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2201, 2203, 2204, 2220) £2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2201, 2203, 2204, 2220) £2,740
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)			
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.	POA
TRAY CONTROL			
JBG2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128) £6,535
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)			
JBG2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415) £495
JBG2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339) £495
JBZ2269		Individual mounting option per request (Special). Requires additional information/specification.	POA
WIRELESS ACCESSORIES FOR ASL CONTROLS			
JBG2249		Wireless Triple Switch Receiver (ASL557-3)	£340
JBG2250		Wireless Mouse Emulator (ASL558)	£810
LINX 10-WAY SWITCH			
JBG2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
JBG2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404) ○
JBG2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405) ○
JBG2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 2108) £65
JBG2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 2108) £65
LINX USB CHARGER PORT			
JBG2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £160
LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
JBG2281		RH Mount	○
JBG2282		LH Mount	○
LINX INPUT MODULE			
JBG2273		DLX-IN200 Module (Jack port for external power button)	£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	(NC 1094) £785

FURTHER LINX COMPONENTS / ACCESSORIES			
JBG2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
JBG2300		To be enclosed only	○
JBZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JBG1685		Red Buddy Button	£90
JBG1686		Yellow Buddy Button	£90
JBG1687		Green Buddy Button	£90
JBG1690		Blue Buddy Button	£90
EGG SWITCH			
JBG2254		Egg Switch, black	£90
PIKO BUTTONS			
JBG1711		Red Piko Button	£160
JBG1712		Yellow Piko Button	£160
JBG1713		Green Piko Button	£160
JBG1714		Orange Piko Button	£160
JBG1715		Blue Piko Button	£160
JBG1716		Black Piko Button	£160
BODYPOINT JOYSTICK KNOBS			
JBG2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JBG2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JBG2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JBG2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
5:1 SWITCH ADAPTOR			
JBG1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£120
SHROUD COLOUR*			
JBG2871		Mid range - Black/Anthracite Shroud	○
JBG2872		Mid range - Silver Shroud	○
JBG2873		Mid range - Bengal Red Shroud	○
JBG2874		Mid range - Mercury Blue Shroud	○
JBG2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300
RIM COLOUR*			
JBG2882		Mid range - Black Rim Inserts	○
JBG2883		Mid range - Silver Rim Inserts	○
JBG2884		Mid range - Bengal Red Rim Inserts	○
JBG2885		Mid range - Mercury Blue Rim Inserts	○
JBG2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
FORKS*			
JBG2631		Single sided fork	£250
JBG2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆

CASTORS*			
JBG0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	◆
JBG0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JBG0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JBG0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)	○
CASTOR LOCKS			
JBG1490		Castor Lock	£195
DRIVE WHEELS*			
JBG0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBG0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
JBG0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed	○
C.T.C SUSPENSION*			
JBG2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBG1910		Soft (NC 1094)	◆
JBG1912		Medium (NC 1094)	○
JBG1911		Firm	○
JBG1915		Very firm	○
MOTOR TYPE*			
JBG1294		Indoor setting, reduced speed approx. 3 km/h	○
JBG1282		6 km/h	◆
JBG1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 2203, 1437, 1254)	£515
JBG1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 0828, 2201, 2203, 1437, 1254)	£600
BATTERY TYPE*			
JBG1259		60 Ah Gel	◆
JBG1260		73.5 Ah Gel (NC 1450)	£100
G TRAC			
JBG1215		G-Trac Module	£400
CHARGER*			
JBG1343		Opti Charge Charger 8A (IP44)	◆
JBG1313		Opti Charge Charger 10A (IP44)	£155
ANTI-TIPPERS*			
JBG1205		Anti-tipper	◆
LIGHTS			
JBG1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors (NC 2203)	£450
JBG1254		Without Light system (NC 1909)	○
KERB RIDER			
JBG1244		Kerb Rider with park position feature (NC 0556, 0589, 1294)	£215

TRANSPORTATION KIT*			
JBG1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JBG1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBG1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 1091, 1092, 1093) ○
EASY ACCESS FRONT TIE DOWNS			
JBG2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
FIXATION BELT*			
JBG1329		Lap belt with Auto Lock Standard mount on backrest	○
JBG2113		Lap belt with Auto Lock Flexible mount in side rails	◆
JBG2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBG2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £95
JBG2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBG2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBG2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
JBG2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £90
REAR VIEW MIRROR			
JBG1415		Rear view mirror RH	(NC 2217, 2285, 2209, 2213, 2265, 2275, 2283, 2287) £45
JBG1339		Rear view mirror LH	(NC 2218, 2286, 2210, 2214, 2266, 2276, 2284, 2288) £45
BAG HOOKS			
JBG2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	(NC 2172, 2173, 2227, 2228, 2233, 2240, 2257, 2291, 2292) £80
CANE HOLDER			
JBG1312		Cane holder RH	(NC 2227, 2259, 1221) £45
JBG1363		Cane holder LH	(NC 2228, 2260, 1222) £45
WHEEL FENDERS			
JBG1416		Fender for Castors	£80
JBG1417		Fender for Drive wheels	£160
TRAY			
JBG1404		Tray, Swing Away RH	(NC 1707, 2209, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217) £200
JBG1405		Tray, Swing Away LH	(NC 1708, 2210, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218) £200

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBG0001	AVIVA RX40 Modulite
JBG2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
JBG2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
JBG2123	480 mm Seat Plate
JBG2124	530 mm Seat Plate
JBG2106	Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side)
JBG2132	410 mm
JBG2134	460 mm
JBG2142	Matrx Flo-tech Lite Cushion
JBG2143	Matrx Flo-tech Contour Visco; (2-way stretch material)
JBG2147	Matrx Libra Cushion, 90 mm thick (Extended lead time)
JBG2144	Matrx Vi Cushion (Mesh fabric)
JBG2148	Standard Cushion (90mm high) - For 550mm Seat Plate Only
JBG1980	Tapered Wedges (Pair)
JBG1981	Abductor
JBG1982	Lateral Wedge
JBG1983	Anterior Seat Wedges
JBG1980	Tapered Wedges (Pair)
JBG1981	Abductor
JBG1982	Lateral Wedge
JBG1983	Anterior Seat Wedges
JBG1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
JBG2108	Backposts without Push Handles
JBG2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
JBG2179	Flex3 Tension Adjustable Upholstery. Height 540 mm (seat plate to top of upholstery).
JBG2172	Flex3 adjustable Backrest Plates (without Backrest Cushion). Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBG2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion. Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBG2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
JBG2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBG2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBG2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBG2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBG2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
JBG2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]
JBG2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBG2380	405 mm (16 inch) MAC BH
JBG2381	455 mm (18 inch) MAC BH
JBG2382	505 mm (20 inch) MAC BH

JBG2383	560 mm (22 inch) MAC BH
JBG2385	90 mm (3.5 inch) MAC Contour Depth LH
JBG2386	150 mm (6 inch) MAC Contour Depth LH
JBG2387	230 mm (9 inch) MAC Contour Depth LH
JBG2389	90 mm (3.5 inch) MAC Contour Depth RH
JBG2390	150 mm (6 inch) MAC Contour Depth RH
JBG2391	230 mm (9 inch) MAC Contour Depth RH
JBG2393	Backrest -25 mm (-1 inch) to seat width
JBG2394	Backrest matches seat width
JBG2395	Backrest +25 mm (+1 inch) to seat width
JBG2151	Black Nylon Fabric TR18
JBG2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
JBG2155	Black Spacer Fabric TR36
JBG2181	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm
JBG2198	Padded Armrest with padded side guard (Water Fall Type)
JBG2183	Inrigged Armrests. Reduces SW by 25 mm each side
JBG2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm).
JBG2192	Following Armrest. Height adjustable from 230 - 300 mm
JBG2194	Following Armrest. Height adjustable from 300 - 360 mm
JBG2195	Short Armpad (325 mm)
JBG2196	Long Armpad (375 mm)
JBG2197	Long Ergonomic Armpad (397 mm)
JBG2199	Soft Armpad (REA Azalea Wide Pad)
JBG2184	Hip Support LH, pad size 140 x 100 mm
JBG2186	Hip Support LH, pad size 230 x 100 mm
JBG2188	Hip Support LH, pad size 330 x 100 mm
JBG2185	Hip Support RH, pad size 140 x 100 mm
JBG2187	Hip Support RH, pad size 230 x 100 mm
JBG2189	Hip Support RH, pad size 330 x 100 mm
JBG1350	Neckrest depth adjustable
JBG1432	Headrest sideways adjustable
JBG1427	Headrest with cheek support
JBG2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
JBG2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
JBG2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
JBG1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBG1445	Lateral Support with extended arm, multi-adjustable pad
JBG1091	Manually settable Seat Tilt 0 - 15°
JBG1092	Powered "fixed pivot" Tilt 20°
JBG1093	Powered Tilt (0° - 30°) with CoG shift
JBG2052	Powered Tilt (0° - 50°) with CoG shift
JBG1094	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift. 6 point Tie downs must be selected
JBG2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG. 6 point Tie downs must be selected

JBG0977	Powered Recliner 30°
JBG0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
JBG0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
JBG0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
JBG0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
JBG0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBG0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBG0501	Without Leg Rests
JBG2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JBG0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
JBG0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")
JBG0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBG0512	Fixed Footplate pair
JBG0504	Angle Adjustable Footplate pair
JBG0513	Angle and depth adjustable Footplate pair
JBG0501	Without footplates for option 0501 without Leg Rests
JBG0554	Amputee Support RH
JBG0555	Amputee Support LH
JBG2201	DLX-REM110 (Drive Only)
JBG2203	DLX-REM211 (Drive & Actuators)
JBG2204	DLX-REM216 (Drive, Actuators & Lights)
JBG2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JBG2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JBG2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JBG1707	Swing Away RH Bracket
JBG1708	Swing Away LH Bracket
JBG1709	Matrx Resolve Swing in/out Bracket RH
JBG1710	Matrx Resolve Swing in/out Bracket LH
JBG2213	Centre Mount Swing Away to RH Nucleus Bracket
JBG2214	Centre Mount Swing Away to LH Nucleus Bracket
JBG2209	REM500 Swing Away RH Multi-Adjustable Bracket
JBG2210	REM500 Swing Away LH Multi-Adjustable Bracket
JBG2227	Bracket RH mounted for Attendant
JBG2228	Bracket LH mounted for Attendant
JBG2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
JBG2219	LiNX DLX-ACU200 (Compact Dual Control)
JBG2245	LiNX DLX-IDC (Intuitive Dual Control)
UNZ9259	mo-vis Scoot Control (Only available via Special)
JBG2259	Dual Control Standard Bracket RH
JBG2260	Dual Control Standard Bracket LH
JBG1221	IDC Dual Control Detachable Bracket RH

JBG1222	IDC Dual Control Detachable Bracket LH
JBG2231	LiNX DLX-CR400 (Compact Remote)
JBG2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JBG2233	ASL MEC (Proportional MEC Joystick ASL130)
JBG2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
JBG2207	RH Mount
JBG2208	LH Mount
JBG2265	Centre Mount Swing Away to RH (Nucleus)
JBG2266	Centre Mount Swing Away to LH (Nucleus)
JBZ2215	Individual mounting option per request (Special). Requires additional information/specification.
JBG2285	Tray Mount RH (ASL603, 100mm x 150mm)
JBG2286	Tray Mount LH (ASL603, 100mm x 150mm)
JBG2283	Centre Mount Swing Away to RH (Nucleus)
JBG2284	Centre Mount Swing Away to LH (Nucleus)
JBG2287	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JBG2288	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JBZ2289	Individual mounting option per request (Special). Requires additional information/specification.
JBG2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
JBG2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
JBG2257	Four-Switch Proximity Array (ASL106)
JBG2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBG2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBG2217	Bracket RH, mounted for user
JBG2218	Bracket LH, mounted for user
JBG2291	Bracket RH, mounted for attendant
JBG2292	Bracket LH, mounted for attendant
JBG2253	USB Charger Port DLX-USB02-A
JBG2273	DLX-IN200 Module (Jack port for external power button)
JBG1686	Yellow Buddy Button
JBG1687	Green Buddy Button
JBG1690	Blue Buddy Button
JBG2254	Egg Switch, black
JBG1711	Red Piko Button
JBG1712	Yellow Piko Button
JBG1713	Green Piko Button
JBG1714	Orange Piko Button
JBG1715	Blue Piko Button
JBG1716	Black Piko Button
JBG0828	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
JBG1910	Soft
JBG1912	Medium
JBG1294	Indoor setting, reduced speed approx. 3 km/h
JBG1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights

JBG1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBG1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBG1254	Without Light system
JBG1244	Kerb Rider with park position feature
JBG1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.
JBG2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBG2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBG2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBG2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBG2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBG2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBG1415	Rear view mirror RH
JBG1339	Rear view mirror LH
JBG2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
JBG1312	Cane holder RH
JBG1363	Cane holder LH
JBG1416	Fender for Castors
JBG1404	Tray, Swing Away RH
JBG1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):