



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

July 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Standard

Standard Features

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Modulite Standard Seat
- One piece seat plate (380mm/15" or 430mm/17" SW), additional 50mm/2" arm rest width adjustment
- Seat depths from 410mm - 510mm (16" - 20")
- Tension adjustable backrest
- Standard foam cushion (50mm thick)
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 15°)
- Manual settable Backrest Recline (0°- 30° in 3.75° steps)
- 80° hangers with fixed footplates
- LiNX DLX-REM211 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension (Control - Traction – Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £4,650



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBF0001	AVIVA RX40 Standard	£4,650
JBZ1300	Special AVIVA RX40 Standard (Extended Lead Time). Please provide details in comments section.	POA
CHASSIS WIDTH*		
JBL2810	Chassis Narrow (560 mm/22")	◆
SEATING TYPE*		
JBF2105	Modulite Standard Seat (Standard seat plate and backrest frame in fixed width) (NC 1401)	◆

SEAT FRAME*			
JBF2110	<input type="checkbox"/>	One Piece Seat Plate	◆
SEAT WIDTH*			
JBF2121	<input type="checkbox"/>	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.	○
JBF2122	<input type="checkbox"/>	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.	◆
SEAT DEPTH*			
JBF2132	<input type="checkbox"/>	410 mm (NC 2197)	○
JBF2134	<input type="checkbox"/>	460 mm	○
JBF2136	<input type="checkbox"/>	510 mm	○
SEAT CUSHION*			
JBF2141	<input type="checkbox"/>	Standard Foam Cushion, 50 mm height	◆
JBF2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	£120
JBF2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco; (2-way stretch material)	£169
JBF2140	<input type="checkbox"/>	Without seat cushion	○
BACKREST SUPPORT*			
JBF2161	<input type="checkbox"/>	Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs	◆
BACKREST UPHOLSTERY / CUSHION*			
JBF2171	<input type="checkbox"/>	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	◆
JBF2170	<input type="checkbox"/>	Without Backrest Upholstery (Back Posts/Back Frame only) (NC 2151)	○
UPHOLSTERY COLOUR*			
JBF2151	<input type="checkbox"/>	Black Nylon Fabric TR18	◆
JBX1002	<input type="checkbox"/>	No Upholstery Colour required - No backrest upholstery required	○
ARMRESTS*			
JBF2191	<input type="checkbox"/>	T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm	◆
JBF2181	<input type="checkbox"/>	T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm	£40
JBF2183	<input type="checkbox"/>	Inrigged Armrests. Reduces SW by 25 mm each side	○
ARMPADS*			
JBF2195	<input type="checkbox"/>	Short Armpad (325 mm)	◆
JBF2196	<input type="checkbox"/>	Long Armpad (375 mm)	○
JBF2197	<input type="checkbox"/>	Long Ergonomic Armpad (397 mm)	○
JBF2199	<input type="checkbox"/>	Soft Armpad (REA Azalea Wide Pad)	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
JBF1350	<input type="checkbox"/>	Neckrest depth adjustable	£80
JBF1495	<input type="checkbox"/>	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
MANUAL / POWERED SEATING OPTIONS*			
JBF1091	<input type="checkbox"/>	Manually settable Seat Tilt 0 - 15°	○
JBF1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	◆

SEAT TO GROUND HEIGHT*			
JBF2857		435 mm	○
JBF2858		460 mm	○
JBF2859		485 mm	○
BACKREST RECLINER*			
JBF0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
LEG REST TYPE*			
JBF0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm (NC 0522)	◆
JBF0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	£120
JBF0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	£235
JBF0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad (NC 0522)	£660
FOOTPLATE TYPE*			
JBF0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBF0504		Angle Adjustable Footplate pair (NC 0526)	£75
JBF0513		Angle and depth adjustable Footplate pair	£120
JBF0522		Angle adjustable Footboard (NC 0531, 1244)	£165
LINX CONTROLS*			
JBF2203		DLX-REM211 (Drive & Actuators) (NC 1437, 1908, 1909)	◆
JBF2204		DLX-REM216 (Drive, Actuators & Lights)	£200
CONTROL BRACKETRY*			
JBF1207		Standard Bracket RH (NC 1404)	○
JBF1208		Standard Bracket LH (NC 1405)	○
JBF1707		Swing Away RH Bracket (NC 1404)	£135
JBF1708		Swing Away LH Bracket (NC 1405)	£135
JBF2227		Bracket RH mounted for Attendant (NC 2220, 2219, 1312)	£65
JBF2228		Bracket LH mounted for Attendant (NC 2220, 2219, 1363)	£65
JBZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
DUAL CONTROLS			
JBF2220		LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 2273, 2227, 2228)	£395
JBF2219		LiNX DLX-ACU200 (Compact Dual Control) (NC 2227, 2228)	£395
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)			
JBF2259		Dual Control Standard Bracket RH (NC 2291, 1312)	○
JBF2260		Dual Control Standard Bracket LH (NC 2292, 1363)	○
LINX 10-WAY SWITCH			
JBF2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
JBF2217		Bracket RH, mounted for user (NC 1404)	○
JBF2218		Bracket LH, mounted for user (NC 1405)	○
JBF2291		Bracket RH, mounted for attendant (NC 2259)	£65
JBF2292		Bracket LH, mounted for attendant (NC 2260)	£65
LINX IN200 MODULE			
JBF2273		DLX-IN200 Module (Jack port for external power button)	£270

FURTHER LINX COMPONENTS / ACCESSORIES			
JBF2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
JBF2300		To be enclosed only	○
JBZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JBF1685		Red Buddy Button	£90
JBF1686		Yellow Buddy Button	£90
JBF1687		Green Buddy Button	£90
JBF1690		Blue Buddy Button	£90
EGG SWITCH			
JBF2254		Egg Switch, black	£90
PIKO BUTTONS			
JBF1711		Red Piko Button	£160
JBF1712		Yellow Piko Button	£160
JBF1713		Green Piko Button	£160
JBF1714		Orange Piko Button	£160
JBF1715		Blue Piko Button	£160
JBF1716		Black Piko Button	£160
BODYPOINT JOYSTICK KNOBS			
JBF2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JBF2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JBF2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JBF2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
SHROUD COLOUR*			
JBF2871		Mid range - Black/Anthracite Shroud	○
JBF2872		Mid range - Silver Shroud	○
JBF2873		Mid range - Bengal Red Shroud	○
JBF2874		Mid range - Mercury Blue Shroud	○
JBF2875		Mid range - Midnight Purple Shroud	○
JBZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300
RIM COLOUR*			
JBF2882		Mid range - Black Rim Inserts	○
JBF2883		Mid range - Silver Rim Inserts	○
JBF2884		Mid range - Bengal Red Rim Inserts	○
JBF2885		Mid range - Mercury Blue Rim Inserts	○
JBF2886		Mid range - Midnight Purple Rim Inserts	○
JBZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
FORKS*			
JBF2631		Single sided fork	£250
JBF2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆

FRONT WHEEL*			
JBF0835	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	◆
JBF0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JBF0838	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JBF0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)	○
CASTOR LOCKS			
JBF1490	<input type="checkbox"/>	Castor Lock	£195
DRIVE WHEEL*			
JBF0704	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBF0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
JBF0706	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed	○
C.T.C SUSPENSION*			
JBF2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBF1910	<input type="checkbox"/>	Soft	◆
JBF1912	<input type="checkbox"/>	Medium	○
JBF1911	<input type="checkbox"/>	Firm	○
JBF1915	<input type="checkbox"/>	Very firm	○
MOTOR TYPE*			
JBF1294	<input type="checkbox"/>	Indoor setting, reduced speed approx. 3 km/h	○
JBF1282	<input type="checkbox"/>	6 km/h	◆
JBF1908	<input type="checkbox"/>	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 2203, 1254)	£515
JBF1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 0828, 2203, 1254)	£600
BATTERY TYPE*			
JBF1259	<input type="checkbox"/>	60 Ah Gel	◆
JBF1260	<input type="checkbox"/>	73.5 Ah Gel (NC 1450)	£100
G TRAC			
JBF1215	<input type="checkbox"/>	G-Trac Module	£400
CHARGER*			
JBF1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	◆
JBF1313	<input type="checkbox"/>	Opti Charge Charger 10A (IP44)	£155
ANTI-TIPPERS*			
JBF1205	<input type="checkbox"/>	Anti-tipper	◆
LIGHTS*			
JBF1437	<input type="checkbox"/>	LED Light system with indicators 30 Lux headlights, included with 10km/h motors (NC 2203)	£450
JBF1254	<input type="checkbox"/>	Without Light system (NC 1909)	○
KERB RIDER			
JBF1244	<input type="checkbox"/>	Kerb Rider with park position feature (NC 1294)	£215

TRANSPORTATION KIT*			
JBF1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JBF1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JBF1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 1091, 1092) ○
EASY ACCESS FRONT TIE DOWNS			
JBF2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
FIXATION BELT*			
JBF1329		Lap belt with Auto Lock Standard mount on backrest	◆
REAR VIEW MIRROR			
JBF1415		Rear view mirror RH	(NC 2217) £45
JBF1339		Rear view mirror LH	(NC 2218) £45
CANE HOLDER			
JBF1312		Cane holder RH	(NC 2227, 2259) £45
JBF1363		Cane holder LH	(NC 2228, 2260) £45
WHEEL FENDERS			
JBF1416		Fender for Castors	£80
JBF1417		Fender for Drive wheels	£160
TRAY			
JBF1404		Tray, Swing Away RH	(NC 1707, 2217) £200
JBF1405		Tray, Swing Away LH	(NC 1708, 2218) £200

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBF2151	Black Nylon Fabric TR18
JBF2197	Long Ergonomic Armpad (397 mm)
JBF1091	Manually settable Seat Tilt 0 - 15°
JBF1092	Powered "fixed pivot" Tilt 20°
JBF0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
JBF0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
JBF0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
JBF0522	Angle adjustable Footboard
JBF2203	DLX-REM211 (Drive & Actuators)
JBF1707	Swing Away RH Bracket
JBF1708	Swing Away LH Bracket
JBF2227	Bracket RH mounted for Attendant
JBF2228	Bracket LH mounted for Attendant
JBF2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
JBF2219	LiNX DLX-ACU200 (Compact Dual Control)
JBF2259	Dual Control Standard Bracket RH
JBF2260	Dual Control Standard Bracket LH
JBF2217	Bracket RH, mounted for user
JBF2218	Bracket LH, mounted for user
JBF2291	Bracket RH, mounted for attendant
JBF2292	Bracket LH, mounted for attendant
JBF2273	DLX-IN200 Module (Jack port for external power button)
JBF1686	Yellow Buddy Button
JBF1687	Green Buddy Button
JBF1690	Blue Buddy Button
JBF2254	Egg Switch, black
JBF1711	Red Piko Button
JBF1712	Yellow Piko Button
JBF1713	Green Piko Button
JBF1714	Orange Piko Button
JBF1715	Blue Piko Button
JBF1716	Black Piko Button
JBF0828	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
JBF1294	Indoor setting, reduced speed approx. 3 km/h
JBF1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBF1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBF1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JBF1254	Without Light system
JBF1244	Kerb Rider with park position feature
JBF1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.

JBF1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JBF1312	Cane holder RH
JBF1363	Cane holder LH
JBF1404	Tray, Swing Away RH
JBF1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):