



Yes, you can:

Power Invacare Ireland Prescription Form

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Ultra

Standard Features

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Ultra Low Maxx Seating System
- Seat widths of 405 mm - 505 mm / 16" - 20" & 480 mm - 555 mm / 19" - 22"
- Seat depths of 405 mm - 505 mm / 16" - 20" & 480 mm - 580 mm / 19" - 23"
- Backrest heights of 560mm - 660mm in 25mm increments
- Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5")
- Fixed STF height and seat angle setting
- Fixed Angle Backrest with Rehab backrest
- Manual Centre Mount legrest
- LiNX DLX-REM211 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic
- 14" (355 mm) x 3" (75 mm), Pneumatic drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: €11,937



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JAA0001	AVIVA RX40 Ultra	€11,937
JAZ1300	Special AVIVA RX40 Ultra (Extended Lead Time). Please provide details in comments section.	<input type="radio"/>
CHASSIS WIDTH*		
JAA2810	Narrow (560 mm/22")	<input type="radio"/>
JAA2830	Wide (635 mm/25")	<input type="radio"/>
SEATING TYPE*		
JAA2020	Ultra Low Maxx Powered Positioning System	<input checked="" type="checkbox"/>

SEAT ADJUSTMENTS*			
JAA2099		Fixed STF height and seat angle setting	(NC 2857, 2054, 2055, 2056, 1401) ◆
JAA2052		Ultra Low Maxx CG Powered Tilt (0 to 50°)	€1,735
JAA2053		Ultra Low Maxx CG Powered Tilt (0 to 45°) combined with 300 mm Seat Lifter. (Increased STF Height by 4 mm)	(NC 2313, 2314, 2315, 2263, 2264, 2261, 2262, 1494) €3,096
SEAT TO GROUND HEIGHT*			
JAA2856		420 mm	○
JAA2857		435 mm	(NC 2055, 2056, 2099) ○
JAA2858		460 mm	(NC 2055, 2056) ○
JAA2859		485 mm	(NC 2055, 2056) ○
ANTERIOR AND POSTERIOR TILT*			
JAA2057		Ultra Low Maxx CG Tilt Standard Mount. No anterior or posterior tilt	◆
JAZ2054		Ultra Low Maxx CG tilt 5° Fixed Anterior Assist mount. Max tilt angle reduced by 5°. With centre-mounted legrests only. Front STF height as defined above	(NC 2313, 2314, 2315, 2857, 2858, 2859, 2099) €458
JAZ2055		Ultra Low Maxx CG tilt 10° Fixed Anterior Assist mount. Max tilt angle reduced by 10°. With centre-mounted legrests only. Front STF height as defined above.	(NC 2312, 2313, 2314, 2315, 2857, 2858, 2859, 2099) €570
JAZ2056		Ultra Low Maxx CG tilt 5° Fixed Posterior-tilt Assist mount. Max tilt angle increased by 5°. Front STF height increased by 6 cm over chosen SH setting.	(NC 2313, 2314, 2315, 2857, 2858, 2859, 2099) €458
BACKREST RECLINER*			
JAA2069		Fixed Angle Setting (NO POWERED RECLINE)	(NC 2096, 2097, 2098, 2263, 2264, 2245) ◆
JAA2060		Powered Recline 90 to 170°, Standard Mount	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) €1,152
JAA2063		Powered Recline, 8° Precline mount, 82° to 162°	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) €1,201
JAA2066		Powered Recline, 30° Precline mount, 60° to 140°	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) €1,201
SEAT FRAME WIDTH RANGE*			
JAA2021		405 mm - 505 mm (increments of 25 mm) / 16" - 20"	(NC 2031) ○
JAA2003		480 mm - 555 mm (increments of 25 mm) / 19" - 22". (Additional Delivery Time)	(NC 2026, 2027) ○
SEAT FRAME DEPTH RANGE*			
JAA2008		405 mm - 505 mm (increments of 25 mm) / 16" - 20"	(NC 2041, 2042) ○
JAA2024		480 mm - 580 mm (increments of 25 mm) / 19" - 23". (Additional Delivery Time)	(NC 2036, 2037, 2810) ○
SEAT WIDTH*			
JAA2025		405 mm / 16"	(NC 2041, 2042, 0572, 0577, 0579, 0584, 2701, 2703, 2705) ○
JAA2026		430 mm / 17"	(NC 2041, 2042, 2701, 2703, 2705) ○
JAA2027		455 mm / 18"	(NC 2701, 2703, 2705) ○
JAA2028		480 mm / 19"	(NC 2704, 2706, 2550) ○
JAA2029		505 mm / 20"	(NC 2702, 2704, 2706) ○
JAA2030		530 mm / 21"	(NC 2035, 2036, 2145, 2702, 2704, 2706) ○
JAA2031		555 mm / 22"	(NC 2035, 2036, 2145, 2702, 2704, 2706) ○

SEAT DEPTH*			
JAA2035		405 mm / 16"	(NC 2030, 2031, 2701, 2702, 2703, 2704, 2705, 2706) <input type="radio"/>
JAA2036		430 mm / 17"	(NC 2030, 2031) <input type="radio"/>
JAA2037		455 mm / 18"	<input type="radio"/>
JAA2038		480 mm / 19"	<input type="radio"/>
JAA2039		505 mm / 20"	<input type="radio"/>
JAA2040		530 mm / 21"	<input type="radio"/>
JAA2041		555 mm / 22"	(NC 2025, 2026, 2145) <input type="radio"/>
JAA2042		580 mm / 23"	(NC 2025, 2026, 2145, 2146) <input type="radio"/>
SEAT CUSHION*			
JAA2144		Matrx Vi Cushion [Mesh fabric]	€473 <input type="radio"/>
JAA2145		Matrx Flo-Tech Lite Black smooth cover, 75 mm thick Up to 127 kg	(NC 2042, 2030, 2031) €154 <input type="radio"/>
JAA2146		Matrx Flo-tech Contour Black smooth cover, 100 mm thick Up to 191 kg	€217 <input type="radio"/>
JAA2147		Matrx Libra Cushion 90 mm thick Standard: up to 136 kg Extended lead time: + 4 days Additional Delivery Time + 4 days	€509 <input type="radio"/>
JAA2140		Without seat cushion	<input type="radio"/>
FIRST POSITIONING ACCESORIES FOR MATRX LIBRA AND FLO-TECH LITE			
JAA1980		Tapered Wedges (Pair)	€28 <input type="radio"/>
JAA1981		Abductor	€58 <input type="radio"/>
JAA1982		Lateral Wedge	€46 <input type="radio"/>
JAA1983		Anterior Seat Wedges	€43 <input type="radio"/>
POSITIONING ACCESSORIES FOR MATRX VI			
JAA1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	€71 <input type="radio"/>
BACK CANE HEIGHT - STRAIGHT OR ANGLED BACK CANES*			
JAA2510		Straight back cane height 455 mm / 18"	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563) <input type="radio"/>
JAA2511		Straight back cane height 505 mm / 20"	(NC 2543, 2051, 2320, 2552, 2553, 2555, 2557, 2177, 2563, 2559, 2178, 2175, 2561, 2562) <input type="radio"/>
JAA2512		Straight back cane height 555 mm / 22"	(NC 2320, 2553, 2555, 2177, 2563, 2178, 2175, 2561, 2562) <input type="radio"/>
JAA2513		Straight back cane height 610 mm / 24"	<input type="radio"/>
JAA2514		Angled (5°) back cane height 505 mm / 20"	(NC 2543, 2051, 2320, 2552, 2553, 2555, 2557, 2177, 2563, 2559, 2178, 2175, 2561, 2562) <input type="radio"/>
JAA2515		Angled (5°) back cane height 555 mm / 22"	(NC 2320, 2553, 2555, 2177, 2563, 2178, 2175, 2561, 2562) <input type="radio"/>
JAA2516		Angled (5°) back cane height 610 mm / 24"	<input type="radio"/>
JAA3007		Please select If powered recline is required	<input type="radio"/>
ANGLE SETTING FOR STRAIGHT BACK CANES*			
JAA2517		75°	(NC 2515, 2516) <input type="radio"/>
JAA2518		79°	(NC 2515, 2517) <input type="radio"/>
JAA2519		85°	(NC 2515, 2518) <input type="radio"/>
JAA2520		90°	(NC 2515, 2519) <input type="radio"/>
JAA2521		95°	(NC 2515, 2520) <input type="radio"/>
JAA2522		101°	(NC 2515, 2521) <input type="radio"/>
JAA2523		106°	(NC 2515, 2522) <input type="radio"/>
JAA2524		112°	(NC 2515, 2523) <input type="radio"/>
JAA2525		116°	(NC 2515, 2524) <input type="radio"/>
JAA3006		Please select if 5° Angle Back post or powered recline choosen	<input type="radio"/>

ANGLE SETTING FOR 5° ANGLED BACK CANES*				
JAA2526		81°	(NC 2511, 2512, 2513)	○
JAA2527		84°	(NC 2511, 2512, 2514)	○
JAA2528		91°	(NC 2511, 2512, 2515)	○
JAA2529		96°	(NC 2511, 2512, 2516)	○
JAA2530		101°	(NC 2511, 2512, 2517)	○
JAA2531		106°	(NC 2511, 2512, 2518)	○
JAA2532		111°	(NC 2511, 2512, 2519)	○
JAA2533		117°	(NC 2511, 2512, 2520)	○
JAA2534		121°	(NC 2511, 2512, 2521)	○
JAA3002		Straight Back post choosen or powered recline choosen		○
BACKREST SUPPORT*				
JAA2540		Standard Rehab Back Pan and Back Upholstery height 355 mm / 14"	(NC 2063, 2066)	○
JAA2541		Standard Rehab Back Pan and Back Upholstery height 405 mm / 16"	(NC 2063, 2066, 2510)	○
JAA2542		Standard Rehab Back Pan and Back Upholstery height 455 mm / 18"	(NC 2063, 2066, 2510, 2511, 2514)	○
JAA2543		Standard Rehab Back Pan and Back Upholstery height 505 mm / 20"	(NC 2063, 2066, 2510, 2511, 2512, 2514, 2515)	○
JAA2050		High Back, height 405 mm / 16" [Mesh fabric]		○
JAA2051		High Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	○
JAA2320		High Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	○
JAA2550		Matrx Elite E2 Back, height 355mm / 14" / 16 - 18" wide	(NC 2029)	€604
JAA2551		Matrx Elite E2 Back, height 405 mm / 16"		€604
JAA2552		Matrx Elite E2 Back, height 455 mm / 18"	(NC 2511, 2514)	€604
JAA2553		Matrx Elite E2 Back, height 505 mm / 20"	(NC 2511, 2512, 2514, 2515)	€604
JAA2554		Matrx Elite HD Back, height 405 mm /16" (Extended lead time: + 4 days)		€668
JAA2555		Matrx Elite HD Back, height 505 mm / 20" (Extended lead time: + 4 days)	(NC 2511, 2512, 2514, 2515)	€668
JAA2556		Matrx Elite TR Back, height 405 mm / 16" [Mesh fabric]		€651
JAA2557		Matrx Elite TR Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	€651
JAA2177		Matrx Elite TR Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	€651
JAA2563		Matrx Elite TR HD Back, height 505 mm / 20" [Mesh fabric] (Extended lead time: + 4 days)	(NC 2511, 2512, 2514, 2515)	€722
JAA2558		Matrx Elite Deep E2 Back, height 405 mm / 16" [Mesh fabric]		€722
JAA2559		Matrx Elite Deep E2 Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	€722
JAA2178		Matrx Elite Deep E2 Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	€722
JAA2174		Matrx PB Back, height 405 mm / 16" [Mesh fabric]		€651
JAA2175		Matrx PB Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	€651
JAA2560		Matrx PB Deep Back, height 405mm / 16"		€668
JAA2561		Matrx PB Deep Back, height 505 mm / 20"	(NC 2511, 2512, 2514, 2515)	€668
JAA2562		Matrx Elite Deep HD Back, height 505 mm / 20" (Extended lead time: + 4 days)	(NC 2511, 2512, 2514, 2515)	€805
JAA2168		Matrx MAC Back - Foam Cushion [Mesh fabric]		€650
JAA2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]		€750
JAA2570		Without Backrest Support, incl. Matrx interface	(NC 2322, 2323, 2324, 2325, 2326, 2327, 2230, 2247)	○

MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)			
JAA2380		405 mm (16 inch) MAC BH	<input type="radio"/>
JAA2381		455 mm (18 inch) MAC BH	<input type="radio"/>
JAA2382		505 mm (20 inch) MAC BH	<input type="radio"/>
JAA2383		560 mm (22 inch) MAC BH	<input type="radio"/>
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)			
JAA2385		90 mm (3.5 inch) MAC Contour Depth LH	€130
JAA2386		150 mm (6 inch) MAC Contour Depth LH	€130
JAA2387		230 mm (9 inch) MAC Contour Depth LH	€130
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)			
JAA2389		90 mm (3.5 inch) MAC Contour Depth RH	€130
JAA2390		150 mm (6 inch) MAC Contour Depth RH	€130
JAA2391		230 mm (9 inch) MAC Contour Depth RH	€130
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)			
JAA2393		Backrest -25 mm (-1 inch) to seat width	<input type="radio"/>
JAA2394		Backrest matches seat width	<input type="radio"/>
JAA2395		Backrest +25 mm (+1 inch) to seat width	<input type="radio"/>
BACKREST HEIGHT (NON-RECLINE)*			
JAA2504		455 mm / 18" ONLY WITH: Back height 355mm/14", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) <input type="radio"/>
JAA2505		485 mm / 19" ONLY WITH: Back height 355mm/14", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) <input type="radio"/>
JAA2506		505 mm / 20" ONLY WITH: Back height 405mm/16", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) <input type="radio"/>
JAA2501		535 mm / 21" ONLY WITH: Back height 405mm/16", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) <input type="radio"/>
JAA2070		555 mm / 22" ONLY WITH Back height 455mm/18", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2543, 2050, 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) <input type="radio"/>
JAA2071		585 mm / 23" ONLY WITH Back height 455mm/18", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2543, 2050, 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) <input type="radio"/>
JAA2072		610 mm / 24" ONLY WITH Back height 505mm/20", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2542, 2050, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) <input type="radio"/>
JAA2073		635 mm / 25" ONLY WITH Back height 505mm/20", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2542, 2050, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) <input type="radio"/>
JAA3020		Powered Recline or No backrest option selected	<input type="radio"/>

BACKREST HEIGHT (ESR RECLINE SYSTEMS)*			
JAA2506		505 mm / 20" ONLY WITH: Elite Back, height 355mm/14"(2550), MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2501		535 mm / 21" ONLY WITH: Elite Back, height 355mm/14"(2550), MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2070		555 mm / 22" ONLY WITH: Back height 405mm/16", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2071		585 mm / 23" ONLY WITH: Back height 405mm/16", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2072		610 mm / 24" ONLY WITH Back height 455mm/18", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563) ○
JAA2073		635 mm / 25" ONLY WITH Back height 455mm/18", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2074		660 mm / 26" ONLY WITH Back height 505mm/20", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA2502		685 mm / 27" ONLY WITH Back height 505mm/20", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA3020		NON -RECLINE BACKREST or No backrest SELECTED	(NC 2063, 2066, 2570) ○
ARMREST*			
JAA2307		Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5"). Not compatible with reclining backrest	(NC 2063, 2066, 2308, 2309) ◆
JAA2095		Cantilever Flip-Back Arm, width, height and angle adjustable, 240 - 330 mm (9.5" - 13")	(NC 2308, 2309) ○
JAA2096		Cantilever Flip-Back Arm, width, height and angle adjustable, 340 - 405 mm (13.5-16")	(NC 2308, 2309) ○
JAA2097		Two-post Recline Arm, Flip Back, adjustable height 230 - 330 mm (9-13")	○
JAA2098		Two-post Recline Arm, Flip Back, adjustable height 320 - 405 mm (12.5-16")	○
ARMREST ACCESSORIES			
JAA2308		1" Offset Armrest Spacer (pair), for two-post recline only	(NC 2096, 2307) €188
JAA2309		1/2" Offset Armrest Spacer (pair), for two-post recline only	(NC 2096, 2307) €188
ARMPAD LH*			
JAA2012		Desk arm pad with startex cover visco foam [255mm]; 40 mm thick	◆
JAA2013		Full length arm pad with startex cover visco foam [355mm], 40 mm thick	○
JAA2014		Ergonomic arm trough; 30 mm thick	(NC 1405, 2311, 2310, 2347, 2348, 2349, 2350) €154
JAA2015		without LH arm pad	○
ARMPAD RH*			
JAA2016		Desk arm pad with startex cover visco foam [255mm]; 40 mm thick	◆
JAA2017		Full length arm pad with startex cover visco foam [355mm]; 40 mm thick	○
JAA2018		Ergonomic arm trough, 30 mm thick	(NC 1405, 2311, 2310, 2351, 2352, 2353, 2354) €154
JAA2019		without RH arm pad	○
ARMPAD ACCESSORIES & MOUNTING OPTIONS LH			
JAA2080		Flat Hand Pad Left	(NC 2013, 2015) €154
JAA2081		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Left	(NC 2013, 2015) €68

ARMPAD ACCESSORIES & MOUNTING OPTIONS RH			
JAA2083		Flat Hand Pad Right (NC 2017, 2019)	€154
JAA2084		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Right (NC 2017, 2019)	€68
MULTI-AXIX UPPER EXTREMITY SUPPORTS LH			
JAA2082		Multi Axis Upper Extremity Support (MACES) Left, lowest armpad height will be increased by 43 mm	€387
MULTI-AXIX UPPER EXTREMITY SUPPORTS RH			
JAA2085		Multi Axis Upper Extremity Support (MACES) Right, lowest armpad height will be increased by 43 mm	€387
Height setting: ARMREST LH SEAT PAN TO TOP OF THE TUBE*			
O001		Factory setting	<input type="radio"/>
O002		Customized	<input type="radio"/>
O003		230 mm	<input type="radio"/>
O004		240 mm	<input type="radio"/>
O005		250 mm	<input type="radio"/>
O006		260 mm	<input type="radio"/>
O007		270 mm	<input type="radio"/>
O008		280 mm	<input type="radio"/>
O009		290 mm	<input type="radio"/>
O010		300 mm	<input type="radio"/>
O011		310 mm	<input type="radio"/>
O012		320 mm	<input type="radio"/>
O013		330 mm	<input type="radio"/>
O014		340 mm	<input type="radio"/>
O015		350 mm	<input type="radio"/>
O016		360 mm	<input type="radio"/>
O017		370 mm	<input type="radio"/>
O018		380 mm	<input type="radio"/>
O019		390 mm	<input type="radio"/>
O020		405 mm	<input type="radio"/>

Height setting: ARMREST RH SEAT PAN TO TOP OF THE TUBE*			
O021		Factory setting	<input type="radio"/>
O022		Customized	<input type="radio"/>
O023		230 mm	<input type="radio"/>
O024		240 mm	<input type="radio"/>
O025		250 mm	<input type="radio"/>
O026		260 mm	<input type="radio"/>
O027		270 mm	<input type="radio"/>
O028		280 mm	<input type="radio"/>
O029		290 mm	<input type="radio"/>
O030		300 mm	<input type="radio"/>
O031		310 mm	<input type="radio"/>
O032		320 mm	<input type="radio"/>
O033		330 mm	<input type="radio"/>
O034		340 mm	<input type="radio"/>
O035		350 mm	<input type="radio"/>
O036		360 mm	<input type="radio"/>
O037		370 mm	<input type="radio"/>
O038		380 mm	<input type="radio"/>
O039		390 mm	<input type="radio"/>
O040		405 mm	<input type="radio"/>

ELBOW BLOCK LH PADS			
JAA2347		Elbow Block XS (110 mm wide x 83 mm high)	€726
JAA2348		Elbow Block S (135 mm wide x 90 mm high)	€726
JAA2349		Elbow Block M (145 mm wide x 95 mm high)	€726
JAA2350		Elbow Block L (165 mm wide x 110 mm high)	€726

ELBOW BLOCK LH MOUNTING HARDWARE			
JAA2361		Vertical 3"/ Horizontal 5"	<input type="radio"/>
JAA2362		Vertical 3"/ Horizontal 7"	<input type="radio"/>
JAA2363		Vertical 5"/ Horizontal 5"	<input type="radio"/>
JAA2364		Vertical 5"/ Horizontal 7"	<input type="radio"/>

ELBOW BLOCKS RH PADS			
JAA2351		Elbow Block XS (110 mm wide x 83 mm high)	€726
JAA2352		Elbow Block S (135 mm wide x 90 mm high)	€726
JAA2353		Elbow Block M (145 mm wide x 95 mm high)	€726
JAA2354		Elbow Block L (165 mm wide x 110 mm high)	€726

ELBOW BLOCK RH MOUNTING HARDWARE			
JAA2371		Vertical 3"/ Horizontal 5"	<input type="radio"/>
JAA2372		Vertical 3"/ Horizontal 7"	<input type="radio"/>
JAA2373		Vertical 5"/ Horizontal 5"	<input type="radio"/>
JAA2374		Vertical 5"/ Horizontal 7"	<input type="radio"/>

HIP SUPPORTS LH				
JAA2184		Hip Support LH, Quick Release, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
JAA2114		Hip Support LH, Quick Release, pad size 180 x 150 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
JAA2186		Hip Support LH, Quick Release, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
JAA2188		Hip Support LH, Quick Release, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
HIP SUPPORTS RH				
JAA2185		Hip Support RH, Quick Release, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
JAA2115		Hip Support RH, Quick Release, pad size 180 x 150 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
JAA2187		Hip Support RH, Quick Release, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
JAA2189		Hip Support RH, Quick Release, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	€331
MOUNTING HARDWARE FOR HIP SUPPORTS				
JAA2164		Mounting hardware for hip supports (6,1" / 155 mm high)		○
JAA2165		Mounting hardware for hip supports (8,1" / 206 mm high)		○
HEADREST*				
JAA2092		Onyx Auto Style headrest with standard mounting hardware	(NC 2230, 2247, 2570)	€168
JAA2323		Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)	(NC 2230, 2247, 2570)	€386
JAA2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)	(NC 2230, 2247, 2570)	€386
JAA2322		Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)	(NC 2230, 2247, 2570)	€386
JAA2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23), complete with Elan Headrest Hardware (MEHW)	(NC 2230, 2247, 2570)	€386
JAA2326		Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover complete with Elan Headrest Hardware (MEHW)	(NC 2230, 2247, 2570)	€386
JAA2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover complete with Elan Headrest Hardware (MEHW)	(NC 2230, 2247, 2570)	€520
JAA1495		Head Control or equivalent Headrest, provided independently from this order. According to the requirements of ISO 7176-19, a headrest is compulsory for using the wheelchair as vehicle seat		○
LATERAL SUPPORTS LH				
JAA2343		Lateral Support Maxx Style Quick Release, pad size [110 mm wide x 83 mm high] LH	(NC 2559, 2178, 2560, 2561, 2562)	€560
JAA2345		Lateral Support Maxx Style Quick Release, pad size [135 mm wide x 90 mm high] LH	(NC 2559, 2178, 2560, 2561, 2562)	€560
JAA2043		Lateral Support Maxx Style Quick Release, pad size [145 mm wide x 95 mm high] LH	(NC 2559, 2178, 2560, 2561, 2562)	€560
JAA2045		Lateral Support Maxx Style Quick Release, pad size [165 mm wide x 110 mm high] LH	(NC 2559, 2178, 2560, 2561, 2562)	€560

LATERAL SUPPORTS RH			
JAA2344		Lateral Support Maxx Style Quick Release, pad size [110 mm wide x 83 mm high] RH	(NC 2559, 2178, 2560, 2561, 2562) €560
JAA2346		Lateral Support Maxx Style Quick Release, pad size [135 mm wide x 90 mm high] RH	(NC 2559, 2178, 2560, 2561, 2562) €560
JAA2044		Lateral Support Maxx Style Quick Release, pad size [145 mm wide x 95 mm high] RH	(NC 2559, 2178, 2560, 2561, 2562) €560
JAA2046		Lateral Support Maxx Style Quick Release, pad size [165 mm wide x 110 mm high] RH	(NC 2559, 2178, 2560, 2561, 2562) €560
LATERAL SUPPORT ACCESSORIES			
JAA2049		Additional telescopic link for Maxx Style lateral support, 25 mm, Pair	(NC 2559, 2178, 2560, 2561, 2562) €222
CENTRE MOUNTED LEGREST*			
JAA0586		Manual Centre Mounted Legrest, angle settable Length: 310 - 410 mm	(NC 0517, 1244, 0579, 0584, 0572, 0577) ◆
JAA0589		LNx Powered Centre Mounted Legrest, length compensating; Length: 340 - 410 mm	(NC 2262, 2263, 2264, 1244, 0579, 0584, 0572, 0577) €2,153
JAX3003		Without Centre Mount Legrests	○
FOOTPLATES / FOOT PLATFORM*			
JAA2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	(NC 3003) €135
JAA0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 3003) ○
JAA0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)	(NC 3003) €177
JAA2366		Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")	(NC 0586, 3003) €354
JAA0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	(NC 0586, 3003) €354
JAX3005		Please select If Centre mount footPlate/ platform is not required.	(NC 0589) ○
Height setting: CENTRE MOUNT LEGREST, OPTIONS 0586 & 0589 . SEAT PAN TO FOOTPLATE*			
0041		Factory setting	○
0042		Customized	○
0043		310 mm	○
0044		325 mm	○
0045		335 mm	○
0046		345 mm	○
0047		355 mm	○
0048		364 mm	○
0049		375 mm	○
0050		385 mm	○
0051		410 mm	○
0077		Please select if CENTRE MOUNT LEGREST OPTIONS 0586 & 0589 is not required	○
Height setting: CENTRE MOUNT LEGREST with TELESCOPING FOOT PLATFORM . SEAT PAN TO FOOTPLATE*			
0052		Factory setting	○
0053		Customized	○
0054		380 mm	○
0055		393 mm	○
0056		406 mm	○
0057		419 mm	○
0076		Please select if CENTRE MOUNT LEGREST with TELESCOPING FOOT PLATFORM is not required	○

SIDE MOUNTED HANGERS LH*			
JAA0539	<input type="checkbox"/>	Hanger ADE 0° - 80° LH (NC 2055, 2261, 2262, 2263, 2264)	€392
JAA0572	<input type="checkbox"/>	Pivot Plus Manual Elevating Legrest LH	€860
JAA0577	<input type="checkbox"/>	Pivot Plus Power Elevating Legrest LH (NC 2261, 2262, 2263, 2264)	€1,362
JAA0535	<input type="checkbox"/>	Without Hanger LH	<input type="radio"/>
SIDE MOUNTED HANGERS RH*			
JAA0549	<input type="checkbox"/>	Hanger ADE 0° - 80° RH (NC 2261, 2262, 2263, 2264)	€392
JAA0579	<input type="checkbox"/>	Pivot Plus Manual Elevating Legrest RH	€860
JAA0584	<input type="checkbox"/>	Pivot Plus Power Elevating Legrest RH (NC 2261, 2262, 2263, 2264)	€1,362
JAA0545	<input type="checkbox"/>	Without Hanger RH	<input type="radio"/>
INDIVIDUAL FOOTPLATE LH*			
JAA0509	<input type="checkbox"/>	Angle Adjustable Footplate LH	€21
JAA0511	<input type="checkbox"/>	Angle and depth adjustable Footplate LH	€79
JAX3012	<input type="checkbox"/>	Please select if Side mounted Legrest are not required	<input type="radio"/>
INDIVIDUAL FOOTPLATE RH*			
JAA0506	<input type="checkbox"/>	Angle Adjustable Footplate RH	€21
JAA0508	<input type="checkbox"/>	Angle and depth adjustable Footplate RH	€79
JAX3013	<input type="checkbox"/>	Please select if Side mounted Legrest are not required	<input type="radio"/>
Height setting: VARI LEGREST LH SEAT PAN TO FOOTPLATE*			
0080	<input type="checkbox"/>	Factory Setting	<input type="radio"/>
0081	<input type="checkbox"/>	Customized	<input type="radio"/>
0082	<input type="checkbox"/>	270 mm	<input type="radio"/>
0083	<input type="checkbox"/>	280 mm	<input type="radio"/>
0084	<input type="checkbox"/>	290 mm	<input type="radio"/>
0085	<input type="checkbox"/>	300 mm	<input type="radio"/>
0086	<input type="checkbox"/>	310 mm	<input type="radio"/>
0087	<input type="checkbox"/>	320 mm	<input type="radio"/>
0088	<input type="checkbox"/>	330 mm	<input type="radio"/>
0089	<input type="checkbox"/>	340 mm	<input type="radio"/>
0090	<input type="checkbox"/>	350 mm	<input type="radio"/>
0091	<input type="checkbox"/>	360 mm	<input type="radio"/>
0092	<input type="checkbox"/>	370 mm	<input type="radio"/>
0093	<input type="checkbox"/>	380 mm	<input type="radio"/>
0094	<input type="checkbox"/>	390 mm	<input type="radio"/>
0095	<input type="checkbox"/>	400 mm	<input type="radio"/>
0096	<input type="checkbox"/>	410 mm	<input type="radio"/>
0097	<input type="checkbox"/>	420 mm	<input type="radio"/>
0098	<input type="checkbox"/>	430 mm	<input type="radio"/>

Height setting: VARI LEGREST RH SEAT PAN TO FOOTPLATE*			
0099	<input type="checkbox"/>	Factory Setting	<input type="radio"/>
0100	<input type="checkbox"/>	Customized	<input type="radio"/>
0101	<input type="checkbox"/>	270 mm	<input type="radio"/>
0102	<input type="checkbox"/>	280 mm	<input type="radio"/>
0103	<input type="checkbox"/>	290 mm	<input type="radio"/>
0104	<input type="checkbox"/>	300 mm	<input type="radio"/>
0105	<input type="checkbox"/>	310 mm	<input type="radio"/>
0106	<input type="checkbox"/>	320 mm	<input type="radio"/>
0107	<input type="checkbox"/>	330 mm	<input type="radio"/>
0108	<input type="checkbox"/>	340 mm	<input type="radio"/>
0109	<input type="checkbox"/>	350 mm	<input type="radio"/>
0110	<input type="checkbox"/>	360 mm	<input type="radio"/>
0111	<input type="checkbox"/>	370 mm	<input type="radio"/>
0112	<input type="checkbox"/>	380 mm	<input type="radio"/>
0113	<input type="checkbox"/>	390 mm	<input type="radio"/>
0114	<input type="checkbox"/>	400 mm	<input type="radio"/>
0115	<input type="checkbox"/>	410 mm	<input type="radio"/>
0116	<input type="checkbox"/>	420 mm	<input type="radio"/>
0117	<input type="checkbox"/>	430 mm	<input type="radio"/>

Height setting: PIVOT PLUS LEGREST LH SEAT PAN TO FOOTPLATE*			
0058	<input type="checkbox"/>	Factory setting	<input type="radio"/>
0059	<input type="checkbox"/>	Customized	<input type="radio"/>
0060	<input type="checkbox"/>	402 mm	<input type="radio"/>
0061	<input type="checkbox"/>	414 mm	<input type="radio"/>
0062	<input type="checkbox"/>	426 mm	<input type="radio"/>
0063	<input type="checkbox"/>	438 mm	<input type="radio"/>
0064	<input type="checkbox"/>	450 mm	<input type="radio"/>
0065	<input type="checkbox"/>	462 mm	<input type="radio"/>
0066	<input type="checkbox"/>	474 mm	<input type="radio"/>
0080	<input type="checkbox"/>	Please select if PIVOT PLUS LEGREST is not required	<input type="radio"/>

Height setting: PIVOT PLUS LEGREST RH SEAT PAN TO FOOTPLATE*			
0067	<input type="checkbox"/>	Factory setting	<input type="radio"/>
0068	<input type="checkbox"/>	Customized	<input type="radio"/>
0069	<input type="checkbox"/>	402 mm	<input type="radio"/>
0070	<input type="checkbox"/>	414 mm	<input type="radio"/>
0071	<input type="checkbox"/>	426 mm	<input type="radio"/>
0072	<input type="checkbox"/>	438 mm	<input type="radio"/>
0073	<input type="checkbox"/>	450 mm	<input type="radio"/>
0074	<input type="checkbox"/>	462 mm	<input type="radio"/>
0075	<input type="checkbox"/>	474 mm	<input type="radio"/>
0079	<input type="checkbox"/>	Please select if PIVOT PLUS LEGREST is not required	<input type="radio"/>

CONTROLS*			
JAA2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2270, 2272, 1908, 1909, 2209, 2210) ◆
JAA2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2270, 2272, 2209, 2210) €198
JAA2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2210) €448
JAA2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1207, 1208, 1707, 1708, 1709, 1710) €1,019
REM400 TOGGLE OPTION			
JAA2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2204, 2206) €415
CONTROL BRACKETRY*			
JAA1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404, 2311) ○
JAA1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405, 2310) ○
JAA1707		Swing Away RH Bracket	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404, 2311) €175
JAA1708		Swing Away LH Bracket	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405, 2310) €175
JAA1709		Matrx Resolve Swing in/out Bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 2235, 2267, 2241, 1404) €281
JAA1710		Matrx Resolve Swing in/out Bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 2236, 2268, 2242, 1405) €281
JAA2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 2311) ○
JAA2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405, 2310) ○
JAA2075		Swing Away + height / width adjust. Mount on RH for user, Quad link design	€320
JAA2076		Swing Away + height / width adjust. Mount on LH for user, Quad link design	€320
JAA2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 2283, 2287, 2207, 2265, 2275, 2217, 1404, 1405, 2310, 2311, 1415) €654
JAA2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 2284, 2288, 2208, 2266, 2276, 2218, 1404, 1405, 2310, 2311, 1339) €654
JAA2227		Bracket RH mounted for Attendant	(NC 1221, 2206, 2219) €84
JAA2228		Bracket LH mounted for Attendant	(NC 1222, 2206, 2219) €84
JAZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
DUAL CONTROLS			
JAA2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1221, 1222, 2227, 2228, 9259) €518
JAA2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2259, 2260, 9259) €984
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2259, 2260, 0001, 2259, 2260, 1221, 1222) €3,020
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)			
JAA2259		Dual Control Standard Bracket RH	(NC 2245) ○
JAA2260		Dual Control Standard Bracket LH	(NC 2245) ○
JAA1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227) ○
JAA1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2228) ○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.	POA

COMPACT REMOTES				
JAA2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288)	€944
JAA2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288)	€1,260
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€2,630
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	€4,160
JAA2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2207, 2208, 2265, 2266, 2215)	€4,741
JAA2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2207, 2208, 2265, 2266, 2215)	€2,642
JAA2243		ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)	(NC 2203, 2204, 2255, 2256, 2207, 2208, 2265, 2266, 2215, 2285, 2286, 2283, 2284, 2287, 2288)	€3,394

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
JAA2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404, 2311)	○
JAA2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405, 2310)	○
JAA2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 2310, 2311, 1415)	€654
JAA2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 2310, 2311, 1339)	€654
JAA2285		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JAA2286		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JAA2283		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	€654
JAA2284		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	€654
JAA2287		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	€654
JAA2288		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	€654
JAA2235		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311)	○
JAA2236		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311)	○
JAA2267		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311)	€654
JAA2268		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311)	€654
JAA2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311)	€654
JAA2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311)	€654
JAZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA

CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	€2,630

CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)			
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)	€1,255
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)	€1,820
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)	€3,005
HEAD CONTROLS			
JAA2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2570) €4,881
JAA2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2570) €5,717
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 0001) €9,135
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 0001) €5,435
FOOT CONTROL			
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2203, 2204) €3,800
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204) €3,615
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)			
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.	POA
TRAY CONTROL			
JAA2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240) €7,560
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)			
JAA2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415) €654
JAA2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339) €654
JAZ2269		Individual mounting option per request (Special). Requires additional information/specification.	POA
SIP N PUFF HEAD CONTROL			
JAA2255		Sip'n Puff Head Control (ASL109), includes Input Module and Lip Switch (ASL 301 LIP Combines simple sip and puff controls with head movements	(NC 2203, 2204, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2243) €7,300
SIP N PUFF CONTROL & INPUT MODULE			
JAA2256		Sip'n Puff Control includes Sip'n Puff kit, Lip Switch ASL301LIP and Input module	(NC 2203, 2204, 2243) €1,221
WIRELESS ACCESSORIES FOR ASL CONTROLS			
JAA2249		Wireless Triple Switch Receiver (ASL557-3)	€448
JAA2250		Wireless Mouse Emulator (ASL558)	€1,072
LINX INTERFACE OPTIONS			
JAA2271		DLX-IN500 Module (sub D9, Sip'n Puff)	(NC 2255, 2256, 2203, 2204) €1,032
JAA2273		DLX-IN200 Module (Jack port for external power button + control inputs via DB9)	€316
JAA2272		DLX-OUT500 (allows environmental control access)	(NC 2204) €1,243
JAA1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	€159
SPECIAL ALTERNATIVE JOYSTICKS & INPUTS			
JAZ2216		Special alternative Joystick/Input and bracket as per request	POA
LINX 10-WAY SWITCH			
JAA2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	€663

LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)				
JAA2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404)	○
JAA2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405)	○
JAA2226		For attendant, mounted on back shroud	(NC 2218)	€84
TOGGLES / BUTTONS				
JAA2261		Stereo Toggle - direct control for tilt-only	(NC 0584, 0589, 2053, 2060, 2063, 2066)	€382
JAA2262		Stereo Button - direct control for tilt-only	(NC 0584, 0589, 2053, 2060, 2063, 2066)	€365
JAA2263		Four-Way Toggle - direct control for tilt & recline	(NC 2069, 0577, 0584, 0589)	€475
JAA2264		Four-Way Button - direct control for tilt & recline	(NC 2069, 0577, 0584, 0589)	€508
TOGGLES / BUTTONS MOUNT				
JAA2279		RH Mount (of Toggle/Button)		○
JAA2280		LH Mount (of Toggle/Button)		○
JAA2278		Without mounting		○
LiNX USB CHARGER PORT				
JAA2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706)	€217
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)				
JAA2281		RH Mount		○
JAA2282		LH Mount		○
REMOTE STOP SWITCH				
JAA2258		Remote Stop Switch (ASL504)		€1,378
FURTHER LiNX COMPONENTS / ACCESSORIES				
JAA2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
JAA2300		To be enclosed only		○
JAZ2301		To be connected and programmed (Special)		POA
BUDDY BUTTONS				
JAA1685		Red Buddy Button		€116
JAA1686		Yellow Buddy Button		€116
JAA1687		Green Buddy Button		€116
JAA1690		Blue Buddy Button		€116
EGG SWITCH				
JAA2254		Egg Switch, black		€116
PIKO BUTTONS				
JAA1711		Red Piko Button		€207
JAA1712		Yellow Piko Button		€207
JAA1713		Green Piko Button		€207
JAA1714		Orange Piko Button		€207
JAA1715		Blue Piko Button		€207
JAA1716		Black Piko Button		€207
BODYPOINT JOYSTICK KNOBS				
JAA2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A		€85
JAA2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A		€85
JAA2723		Domed Rubber, Joystick Handle, Bodypoint PC107A		€85
JAA2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A		€85

SHROUD COLOUR*			
JAA2871		Mid range - Black/Anthracite Shroud	○
JAA2872		Mid range - Silver Shroud	○
JAA2873		Mid range - Bengal Red Shroud	○
JAA2874		Mid range - Mercury Blue Shroud	○
JAA2875		Mid range - Midnight Purple Shroud	○
JAZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	€281
RIM COLOUR*			
JAA2882		Mid range - Black	○
JAA2883		Mid range - Silver	○
JAA2884		Mid range - Bengal Red	○
JAA2885		Mid range - Mercury Blue	○
JAA2886		Mid range - Midnight Purple	○
JAZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
FORKS*			
JAA2631		Single sided fork	€462
JAA2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
CASTORS*			
JAA0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	◆
JAA0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JAA0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JAA0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)	○
CASTOR LOCKS			
JAA1490		Castor Lock	€259
DRIVE WHEELS*			
JAA0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JAA0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
JAA0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed	○
C.T.C SUSPENSION*			
JAA2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JAA1910		Soft	○
JAA1912		Medium	○
JAA1911		Firm	○
JAA1915		Very firm	○

MOTOR TYPE*			
JAA1294		Indoor setting, reduced speed approx. 3 km/h	○
JAA1282		6 km/h	◆
JAA1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 1437, 1254, 2837) €713
JAA1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 1437, 1254, 2837, 0828) €857
BATTERY TYPE*			
JAA1259		60 Ah Gel	◆
JAA1260		73.5 Ah Gel	(NC 1450) €121
G TRAC			
JAA1215		G-Trac Module	€425
CHARGER*			
JAA1343		OptiCharge Charger 8A (IP44)	◆
JAA1313		OptiCharge Charger 10A (IP44)	€198
ANTI-TIPPERS*			
JAA1205		Anti-tipper	◆
LIGHTS			
JAA1437		LED Light system with indicators 30 Lux headlights, included with 10km/h motors	(NC 1908, 1909) €629
JAA2837		LED Light system extra bright with indicators 60 Lux headlights	(NC 1908, 1909) €188
JAA1254		Without Light system	(NC 1909) ○
KERB RIDER			
JAA1244		Kerb Rider with park position feature	(NC 0586, 0589, 1294) €278
TRANSPORTATION KIT*			
JAA1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	€769
JAA1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JAA1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 2052) ○
EASY ACCESS FRONT TIE DOWNS			
JAA2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	€64
FIXATION BELT*			
JAA2113		Lap belt with Auto Lock Flexible mount in side rails	◆
JAA2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) €101
JAA2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) €101
JAA2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) €104
JAA2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) €104
JAA2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) €104
JAA2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) €104

REAR VIEW MIRROR			
JAA1415		Rearview mirror RH	(NC 2265, 2267, 2275) €58
JAA1339		Rearview mirror LH	(NC 2266, 2268, 2276) €58
OTHER OPTIONS			
JAA1438		Luggage rack	€432
JAA1416		Fender for Castors	€116
JAA1417		Fender for Drive wheels	€222
JAA1404		Tray, Swing Away RH	(NC 1707, 2207, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) €219
JAA1405		Tray, Swing Away LH	(NC 1708, 2208, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) €219
JAA2310		Tray Flip and fold away 1/2 lap tray LH (14" W x 10" D)	(NC 1707, 2207, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) €969
JAA2311		Tray Flip and fold away 1/2 lap tray RH (14" W x 10" D)	(NC 1708, 2208, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) €969
JAZ2312		Vent Box, Low Profile Recline, Universal	(NC 2054, 2055, 2056) POA
JAZ2314		Vent Box, Low Profile Recline, Trilogy	(NC 2054, 2055, 2056) POA
JAZ2315		Vent Box, Low Profile Recline, LTV – Easy Fit	(NC 2054, 2055, 2056) POA
JAZ2313		Oxygen Holder, Low Profile Recline	(NC 2054, 2055, 2056) POA

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JAA0001	AVIVA RX40 Ultra
JAA2810	Narrow (560 mm/22")
JAA2099	Fixed STF height and seat angle setting
JAA2052	Ultra Low Maxx CG Powered Tilt (0 to 50°)
JAA2053	Ultra Low Maxx CG Powered Tilt (0 to 45°) combined with 300 mm Seat Lifter. (Increased STF Height by 4 mm)
JAA2857	435 mm
JAA2858	460 mm
JAA2859	485 mm
JAZ2054	Ultra Low Maxx CG tilt 5° Fixed Anterior Assist mount. Max tilt angle reduced by 5°. With centre-mounted legrests only. Front STF height as defined above
JAZ2055	Ultra Low Maxx CG tilt 10° Fixed Anterior Assist mount. Max tilt angle reduced by 10°. With centre-mounted legrests only. Front STF height as defined above.
JAZ2056	Ultra Low Maxx CG tilt 5° Fixed Posterior-tilt Assist mount. Max tilt angle increased by 5°. Front STF height increased by 6 cm over chosen SH setting.
JAA2069	Fixed Angle Setting (NO POWERED RECLINE)
JAA2060	Powered Recline 90 to 170° , Standard Mount
JAA2063	Powered Recline, 8° Precline mount, 82° to 162°
JAA2066	Powered Recline, 30° Precline mount, 60° to 140°
JAA2025	405 mm / 16"
JAA2026	430 mm / 17"
JAA2027	455 mm / 18"
JAA2028	480 mm / 19"
JAA2029	505 mm / 20"
JAA2030	530 mm / 21"
JAA2031	555 mm / 22"
JAA2035	405 mm / 16"
JAA2036	430 mm / 17"
JAA2037	455 mm / 18"
JAA2041	555 mm / 22"
JAA2042	580 mm / 23"
JAA2145	Matrix Flo-Tech Lite Black smooth cover, 75 mm thick Up to 127 kg
JAA2146	Matrix Flo-tech Contour Black smooth cover, 100 mm thick Up to 191 kg
JAA2510	Straight back cane height 455 mm / 18"
JAA2511	Straight back cane height 505 mm / 20"
JAA2512	Straight back cane height 555 mm / 22"
JAA2513	Straight back cane height 610 mm / 24"
JAA2514	Angled (5°) back cane height 505 mm / 20"
JAA2515	Angled (5°) back cane height 555 mm / 22"
JAA2516	Angled (5°) back cane height 610 mm / 24"
JAA2517	75°
JAA2518	79°
JAA2519	85°
JAA2520	90°

JAA2521	95°
JAA2522	101°
JAA2523	106°
JAA2524	112°
JAA2540	Standard Rehab Back Pan and Back Upholstery height 355 mm / 14"
JAA2541	Standard Rehab Back Pan and Back Upholstery height 405 mm / 16"
JAA2542	Standard Rehab Back Pan and Back Upholstery height 455 mm / 18"
JAA2543	Standard Rehab Back Pan and Back Upholstery height 505 mm / 20"
JAA2050	High Back, height 405 mm / 16" [Mesh fabric]
JAA2051	High Back, height 455 mm / 18" [Mesh fabric]
JAA2320	High Back, height 505 mm / 20" [Mesh fabric]
JAA2550	Matrx Elite E2 Back, height 355mm / 14" / 16 - 18" wide
JAA2551	Matrx Elite E2 Back, height 405 mm / 16"
JAA2552	Matrx Elite E2 Back, height 455 mm / 18"
JAA2553	Matrx Elite E2 Back, height 505 mm / 20"
JAA2554	Matrx Elite HD Back, height 405 mm / 16" (Extended lead time: + 4 days)
JAA2555	Matrx Elite HD Back, height 505 mm / 20" (Extended lead time: + 4 days)
JAA2556	Matrx Elite TR Back, height 405 mm / 16" [Mesh fabric]
JAA2557	Matrx Elite TR Back, height 455 mm / 18" [Mesh fabric]
JAA2177	Matrx Elite TR Back, height 505 mm / 20" [Mesh fabric]
JAA2563	Matrx Elite TR HD Back, height 505 mm / 20" [Mesh fabric] (Extended lead time: + 4 days)
JAA2558	Matrx Elite Deep E2 Back, height 405 mm / 16" [Mesh fabric]
JAA2559	Matrx Elite Deep E2 Back, height 455 mm / 18" [Mesh fabric]
JAA2178	Matrx Elite Deep E2 Back, height 505 mm / 20" [Mesh fabric]
JAA2174	Matrx PB Back, height 405 mm / 16" [Mesh fabric]
JAA2175	Matrx PB Back, height 505 mm / 20" [Mesh fabric]
JAA2560	Matrx PB Deep Back, height 405mm / 16"
JAA2561	Matrx PB Deep Back, height 505 mm / 20"
JAA2562	Matrx Elite Deep HD Back, height 505 mm / 20" (Extended lead time: + 4 days)
JAA2570	Without Backrest Support, incl. Matrx interface
JAA2307	Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5"). Not compatible with reclining backrest
JAA2096	Cantilever Flip-Back Arm, width, height and angle adjustable, 340 - 405 mm (13.5-16")
JAA2097	Two-post Recline Arm, Flip Back, adjustable height 230 - 330 mm (9-13")
JAA2098	Two-post Recline Arm, Flip Back, adjustable height 320 - 405 mm (12.5-16")
JAA2308	1" Offset Armrest Spacer (pair), for two-post recline only
JAA2309	1/2" Offset Armrest Spacer (pair), for two-post recline only
JAA2013	Full length arm pad with startex cover visco foam [355mm], 40 mm thick
JAA2015	without LH arm pad
JAA2017	Full length arm pad with startex cover visco foam [355mm]; 40 mm thick
JAA2019	without RH arm pad

JAA2347	Elbow Block XS (110 mm wide x 83 mm high)
JAA2348	Elbow Block S (135 mm wide x 90 mm high)
JAA2349	Elbow Block M (145 mm wide x 95 mm high)
JAA2350	Elbow Block L (165 mm wide x 110 mm high)
JAA2351	Elbow Block XS (110 mm wide x 83 mm high)
JAA2352	Elbow Block S (135 mm wide x 90 mm high)
JAA2353	Elbow Block M (145 mm wide x 95 mm high)
JAA2354	Elbow Block L (165 mm wide x 110 mm high)
JAA2114	Hip Support LH, Quick Release, pad size 180 x 150 mm
JAA2186	Hip Support LH, Quick Release, pad size 230 x 100 mm
JAA2188	Hip Support LH, Quick Release, pad size 330 x 100 mm
JAA2185	Hip Support RH, Quick Release, pad size 140 x 100 mm
JAA2115	Hip Support RH, Quick Release, pad size 180 x 150 mm
JAA2187	Hip Support RH, Quick Release, pad size 230 x 100 mm
JAA2189	Hip Support RH, Quick Release, pad size 330 x 100 mm
JAA2092	Onyx Auto Style headrest with standard mounting hardware
JAA2323	Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2322	Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23), complete with Elan Headrest Hardware (MEHW)
JAA2326	Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover complete with Elan Headrest Hardware (MEHW)
JAA2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover complete with Elan Headrest Hardware (MEHW)
JAA0586	Manual Centre Mounted Legrest, angle settable Length: 310 - 410 mm
JAA0589	LNx Powered Centre Mounted Legrest, length compensating; Length: 340 - 410 mm
JAX3003	Without Centre Mount Legrests
JAA0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JAA0572	Pivot Plus Manual Elevating Legrest LH
JAA0577	Pivot Plus Power Elevating Legrest LH
JAA0579	Pivot Plus Manual Elevating Legrest RH
JAA0584	Pivot Plus Power Elevating Legrest RH
JAA2203	DLX-REM211 (Drive & Actuators)
JAA2204	DLX-REM216 (Drive, Actuators & Lights)
JAA2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JAA2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JAA2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JAA1207	Standard Bracket RH
JAA1208	Standard Bracket LH
JAA1707	Swing Away RH Bracket
JAA1708	Swing Away LH Bracket
JAA1709	Matrx Resolve Swing in/out Bracket RH
JAA1710	Matrx Resolve Swing in/out Bracket LH
JAA2209	REM500 Swing Away RH Multi-Adjustable Bracket

JAA2210	REM500 Swing Away LH Multi-Adjustable Bracket
JAA2213	Centre Mount Swing Away to RH Nucleus Bracket
JAA2214	Centre Mount Swing Away to LH Nucleus Bracket
JAA2227	Bracket RH mounted for Attendant
JAA2228	Bracket LH mounted for Attendant
JAA2219	LiNX DLX-ACU200 (Compact Dual Control)
JAA2245	LiNX DLX-IDC (Intuitive Dual Control)
UNZ9259	mo-vis Scoot Control (Only available via Special)
JAA2259	Dual Control Standard Bracket RH
JAA2260	Dual Control Standard Bracket LH
JAA1221	IDC Dual Control Detachable Bracket RH
JAA1222	IDC Dual Control Detachable Bracket LH
JAA2231	LiNX DLX-CR400 (Compact Remote)
JAA2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JAA2233	ASL MEC (Proportional MEC Joystick ASL130)
JAA2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
JAA2243	ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)
JAA2207	RH Mount
JAA2208	LH Mount
JAA2265	Centre Mount Swing Away to RH (Nucleus)
JAA2266	Centre Mount Swing Away to LH (Nucleus)
JAA2285	Tray Mount RH (ASL603, 100mm x 150mm)
JAA2286	Tray Mount LH (ASL603, 100mm x 150mm)
JAA2283	Centre Mount Swing Away to RH (Nucleus)
JAA2284	Centre Mount Swing Away to LH (Nucleus)
JAA2287	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAA2288	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAA2235	Tray Mount RH (ASL603, 100mm x 150mm)
JAA2236	Tray Mount LH (ASL603, 100mm x 150mm)
JAA2267	Centre Mount Swing Away to RH (Nucleus)
JAA2268	Centre Mount Swing Away to LH (Nucleus)
JAA2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAA2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAZ2215	Individual mounting option per request (Special). Requires additional information/specification.
JAA2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
JAA2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
JAA2257	Four-Switch Proximity Array (ASL106)
JAA2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JAA2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JAA2255	Sip'n Puff Head Control (ASL109), includes Input Module and Lip Switch (ASL 301 LIP Combines simple sip and puff controls with head movements)
JAA2256	Sip'n Puff Control includes Sip'n Puff kit, Lip Switch ASL301LIP and Input module
JAA2271	DLX-IN500 Module (sub D9, Sip'n Puff)
JAA2272	DLX-OUT500 (allows environmental control access)

JAA2217	Bracket RH, mounted for user
JAA2218	Bracket LH, mounted for user
JAA2261	Stereo Toggle - direct control for tilt-only
JAA2262	Stereo Button - direct control for tilt-only
JAA2263	Four-Way Toggle - direct control for tilt & recline
JAA2264	Four-Way Button - direct control for tilt & recline
JAA2253	USB Charger Port DLX-USB02-A
JAA1686	Yellow Buddy Button
JAA1687	Green Buddy Button
JAA1690	Blue Buddy Button
JAA2254	Egg Switch, black
JAA1711	Red Piko Button
JAA1712	Yellow Piko Button
JAA1713	Green Piko Button
JAA1714	Orange Piko Button
JAA1715	Blue Piko Button
JAA1716	Black Piko Button
JAA0828	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
JAA1294	Indoor setting, reduced speed approx. 3 km/h
JAA1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JAA1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JAA1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JAA2837	LED Light system extra bright with indicators 60 Lux headlights
JAA1254	Without Light system
JAA1244	Kerb Rider with park position feature
JAA1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.
JAA1494	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.
JAA1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JAA2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JAA2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JAA2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JAA2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JAA2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JAA2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JAA1415	Rearview mirror RH
JAA1339	Rearview mirror LH
JAA1404	Tray, Swing Away RH
JAA1405	Tray, Swing Away LH
JAA2310	Tray Flip and fold away 1/2 lap tray LH (14" W x 10" D)
JAA2311	Tray Flip and fold away 1/2 lap tray RH (14" W x 10" D)

JAZ2312	Vent Box, Low Profile Recline, Universal
JAZ2314	Vent Box, Low Profile Recline, Trilogy
JAZ2315	Vent Box, Low Profile Recline, LTV – Easy Fit
JAZ2313	Oxygen Holder, Low Profile Recline



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):