



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

July 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Ultra

Standard Features

- Max User Weight 136kg (21 Stone)
- Crash tested ISO7176 - 19
- Ultra Low Maxx Seating System
- Seat widths of 405 mm - 505 mm / 16" - 20" & 480 mm - 555 mm / 19" - 22"
- Seat depths of 405 mm - 505 mm / 16" - 20" & 480 mm - 580 mm / 19" - 23"
- Backrest heights of 560mm - 660mm in 25mm increments
- Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5")
- Fixed STF height and seat angle setting
- Fixed Angle Backrest with Rehab backrest
- Manual Centre Mount legrest
- LiNX DLX-REM211 with standard bracket
- Choice of 5 shroud colours and 5 rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic
- 14" (355 mm) x 3" (75 mm), Pneumatic drive wheels
- CTC Suspension (Control - Traction - Comfort)
- 6 km/h (4mph) 4 pole motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £9,500



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JAA0001	AVIVA RX40 Ultra	£9,500
JAZ1300	Special AVIVA RX40 Ultra (Extended Lead Time). Please provide details in comments section.	<input type="radio"/>
CHASSIS WIDTH*		
JAA2810	Narrow (560 mm/22")	<input type="radio"/>
JAA2830	Wide (635 mm/25")	<input type="radio"/>
SEATING TYPE*		
JAA2020	Ultra Low Maxx Powered Positioning System	<input checked="" type="checkbox"/>

SEAT ADJUSTMENTS*			
JAA2099		Fixed STF height and seat angle setting	(NC 2857, 2054, 2055, 2056, 1401) ◆
JAA2052		Ultra Low Maxx CG Powered Tilt (0 to 50°)	£1,450
JAA2053		Ultra Low Maxx CG Powered Tilt (0 to 45°) combined with 300 mm Seat Lifter. (Increased STF Height by 4 mm)	(NC 2313, 2314, 2315, 2263, 2264, 2261, 2262, 1494) £2,500
SEAT TO GROUND HEIGHT*			
JAA2856		420 mm	○
JAA2857		435 mm	(NC 2055, 2056, 2099) ○
JAA2858		460 mm	(NC 2055, 2056) ○
JAA2859		485 mm	(NC 2055, 2056) ○
ANTERIOR AND POSTERIOR TILT*			
JAA2057		Ultra Low Maxx CG Tilt Standard Mount. No anterior or posterior tilt	◆
JAZ2054		Ultra Low Maxx CG tilt 5° Fixed Anterior Assist mount. Max tilt angle reduced by 5°. With centre-mounted legrests only. Front STF height as defined above	(NC 2313, 2314, 2315, 2857, 2858, 2859, 2099) £320
JAZ2055		Ultra Low Maxx CG tilt 10° Fixed Anterior Assist mount. Max tilt angle reduced by 10°. With centre-mounted legrests only. Front STF height as defined above.	(NC 2312, 2313, 2314, 2315, 2857, 2858, 2859, 2099) £400
JAZ2056		Ultra Low Maxx CG tilt 5° Fixed Posterior-tilt Assist mount. Max tilt angle increased by 5°. Front STF height increased by 6 cm over chosen SH setting.	(NC 2313, 2314, 2315, 0592, 2857, 2858, 2859, 2099) £320
BACKREST RECLINER*			
JAA2069		Fixed Angle Setting (NO POWERED RECLINE)	(NC 2096, 2097, 2098, 2263, 2264, 2245) ◆
JAA2060		Powered Recline 90 to 170° , Standard Mount	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) £875
JAA2063		Powered Recline, 8° Precline mount, 82° to 162°	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) £915
JAA2066		Powered Recline, 30° Precline mount, 60° to 140°	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) £915
SEAT FRAME WIDTH RANGE*			
JAA2021		405 mm - 505 mm (increments of 25 mm) / 16" - 20"	(NC 2031) ○
JAA2003		480 mm - 555 mm (increments of 25 mm) / 19" - 22". (Additional Delivery Time)	(NC 2026, 2027) ○
SEAT FRAME DEPTH RANGE*			
JAA2008		405 mm - 505 mm (increments of 25 mm) / 16" - 20"	(NC 2041, 2042) ○
JAA2024		480 mm - 580 mm (increments of 25 mm) / 19" - 23". (Additional Delivery Time)	(NC 2036, 2037, 2810) ○
SEAT WIDTH*			
JAA2025		405 mm / 16"	(NC 2041, 2042, 0572, 0577, 0579, 0584, 2701, 2703, 2705) ○
JAA2026		430 mm / 17"	(NC 2041, 2042, 2701, 2703, 2705) ○
JAA2027		455 mm / 18"	(NC 2701, 2703, 2705) ○
JAA2028		480 mm / 19"	(NC 2704, 2706, 2550) ○
JAA2029		505 mm / 20"	(NC 2702, 2704, 2706) ○
JAA2030		530 mm / 21"	(NC 2035, 2036, 2145, 2702, 2704, 2706) ○
JAA2031		555 mm / 22"	(NC 2035, 2036, 2145, 2702, 2704, 2706) ○

SEAT DEPTH*			
JAA2035		405 mm / 16"	(NC 2030, 2031, 2701, 2702, 2703, 2704, 2705, 2706) ○
JAA2036		430 mm / 17"	(NC 2030, 2031) ○
JAA2037		455 mm / 18"	○
JAA2038		480 mm / 19"	○
JAA2039		505 mm / 20"	○
JAA2040		530 mm / 21"	○
JAA2041		555 mm / 22"	(NC 2025, 2026, 2145) ○
JAA2042		580 mm / 23"	(NC 2025, 2026, 2145, 2146) ○
SEAT CUSHION*			
JAA2144		Matrx Vi Cushion [Mesh fabric]	£400
JAA2145		Matrx Flo-Tech Lite Black smooth cover, 75 mm thick Up to 127 kg	(NC 2042, 2030, 2031) £120
JAA2146		Matrx Flo-tech Contour Black smooth cover, 100 mm thick Up to 191 kg	£169
JAA2147		Matrx Libra Cushion 90 mm thick Standard: up to 136 kg Extended lead time: + 4 days Additional Delivery Time + 4 days	£410
JAA2140		Without seat cushion	○
1st FIRST POSITIONING ACCESORIES FOR MATRX LIBRA AND FLO-TECH LITE			
JAA1980		Tapered Wedges (Pair)	£30
JAA1981		Abductor	£23
JAA1982		Lateral Wedge	£20
JAA1983		Anterior Seat Wedges	£27
2nd FIRST POSITIONING ACCESORIES FOR MATRX LIBRA AND FLO-TECH LITE			
JAA1980		Tapered Wedges (Pair)	£30
JAA1981		Abductor	£23
JAA1982		Lateral Wedge	£20
JAA1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
JAA1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACK CANE HEIGHT - STRAIGHT OR ANGLED BACK CANES*			
JAA2510		Straight back cane height 455 mm / 18"	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563) ○
JAA2511		Straight back cane height 505 mm / 20"	(NC 2543, 2051, 2320, 2552, 2553, 2555, 2557, 2177, 2563, 2559, 2178, 2175, 2561, 2562) ○
JAA2512		Straight back cane height 555 mm / 22"	(NC 2320, 2553, 2555, 2177, 2563, 2178, 2175, 2561, 2562) ○
JAA2513		Straight back cane height 610 mm / 24"	○
JAA2514		Angled (5°) back cane height 505 mm / 20"	(NC 2543, 2051, 2320, 2552, 2553, 2555, 2557, 2177, 2563, 2559, 2178, 2175, 2561, 2562) ○
JAA2515		Angled (5°) back cane height 555 mm / 22"	(NC 2320, 2553, 2555, 2177, 2563, 2178, 2175, 2561, 2562) ○
JAA2516		Angled (5°) back cane height 610 mm / 24"	○
JAA3007		Please select If powered recline is required	○

ANGLE SETTING FOR STRAIGHT BACK CANES*				
JAA2517	<input type="checkbox"/>	75°	(NC 2515, 2516)	<input type="radio"/>
JAA2518	<input type="checkbox"/>	79°	(NC 2515, 2517)	<input type="radio"/>
JAA2519	<input type="checkbox"/>	85°	(NC 2515, 2518)	<input type="radio"/>
JAA2520	<input type="checkbox"/>	90°	(NC 2515, 2519)	<input type="radio"/>
JAA2521	<input type="checkbox"/>	95°	(NC 2515, 2520)	<input type="radio"/>
JAA2522	<input type="checkbox"/>	101°	(NC 2515, 2521)	<input type="radio"/>
JAA2523	<input type="checkbox"/>	106°	(NC 2515, 2522)	<input type="radio"/>
JAA2524	<input type="checkbox"/>	112°	(NC 2515, 2523)	<input type="radio"/>
JAA2525	<input type="checkbox"/>	116°	(NC 2515, 2524)	<input type="radio"/>
JAA3006	<input type="checkbox"/>	Please select if 5° Angle Back post or powered recline choosen		<input type="radio"/>
ANGLE SETTING FOR 5° ANGLED BACK CANES*				
JAA2526	<input type="checkbox"/>	81°	(NC 2511, 2512, 2513)	<input type="radio"/>
JAA2527	<input type="checkbox"/>	84°	(NC 2511, 2512, 2514)	<input type="radio"/>
JAA2528	<input type="checkbox"/>	91°	(NC 2511, 2512, 2515)	<input type="radio"/>
JAA2529	<input type="checkbox"/>	96°	(NC 2511, 2512, 2516)	<input type="radio"/>
JAA2530	<input type="checkbox"/>	101°	(NC 2511, 2512, 2517)	<input type="radio"/>
JAA2531	<input type="checkbox"/>	106°	(NC 2511, 2512, 2518)	<input type="radio"/>
JAA2532	<input type="checkbox"/>	111°	(NC 2511, 2512, 2519)	<input type="radio"/>
JAA2533	<input type="checkbox"/>	117°	(NC 2511, 2512, 2520)	<input type="radio"/>
JAA2534	<input type="checkbox"/>	121°	(NC 2511, 2512, 2521)	<input type="radio"/>
JAA3002	<input type="checkbox"/>	Straight Back post choosen or powered recline choosen		<input type="radio"/>

BACKREST SUPPORT*				
JAA2540		Standard Rehab Back Pan and Back Upholstery height 355 mm / 14"	(NC 2063, 2066)	○
JAA2541		Standard Rehab Back Pan and Back Upholstery height 405 mm / 16"	(NC 2063, 2066, 2510)	○
JAA2542		Standard Rehab Back Pan and Back Upholstery height 455 mm / 18"	(NC 2063, 2066, 2510, 2511, 2514)	○
JAA2543		Standard Rehab Back Pan and Back Upholstery height 505 mm / 20"	(NC 2063, 2066, 2510, 2511, 2512, 2514, 2515)	○
JAA2050		High Back, height 405 mm / 16" [Mesh fabric]		○
JAA2051		High Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	○
JAA2320		High Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	○
JAA2550		Matrx Elite E2 Back, height 355mm / 14" / 16 - 18" wide	(NC 2029)	£460
JAA2551		Matrx Elite E2 Back, height 405 mm / 16"		£460
JAA2552		Matrx Elite E2 Back, height 455 mm / 18"	(NC 2511, 2514)	£460
JAA2553		Matrx Elite E2 Back, height 505 mm / 20"	(NC 2511, 2512, 2514, 2515)	£460
JAA2554		Matrx Elite HD Back, height 405 mm /16" (Extended lead time: + 4 days)		£509
JAA2555		Matrx Elite HD Back, height 505 mm / 20" (Extended lead time: + 4 days)	(NC 2511, 2512, 2514, 2515)	£509
JAA2556		Matrx Elite TR Back, height 405 mm / 16" [Mesh fabric]		£496
JAA2557		Matrx Elite TR Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	£496
JAA2177		Matrx Elite TR Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	£496
JAA2563		Matrx Elite TR HD Back, height 505 mm / 20" [Mesh fabric] (Extended lead time: + 4 days)	(NC 2511, 2512, 2514, 2515)	£550
JAA2558		Matrx Elite Deep E2 Back, height 405 mm / 16" [Mesh fabric]		£550
JAA2559		Matrx Elite Deep E2 Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	£550
JAA2178		Matrx Elite Deep E2 Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	£550
JAA2174		Matrx PB Back, height 405 mm / 16" [Mesh fabric]		£496
JAA2175		Matrx PB Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	£496
JAA2560		Matrx PB Deep Back, height 405mm / 16"		£509
JAA2561		Matrx PB Deep Back, height 505 mm / 20"	(NC 2511, 2512, 2514, 2515)	£509
JAA2562		Matrx Elite Deep HD Back, height 505 mm / 20" (Extended lead time: + 4 days)	(NC 2511, 2512, 2514, 2515)	£613
JAA2168		Matrx MAC Back - Foam Cushion [Mesh fabric]		£475
JAA2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]		£575
JAA2570		Without Backrest Support, incl. Matrx interface	(NC 2322, 2323, 2324, 2325, 2326, 2327, 2230, 2247)	○
MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
JAA2380		405 mm (16 inch) MAC BH		○
JAA2381		455 mm (18 inch) MAC BH		○
JAA2382		505 mm (20 inch) MAC BH		○
JAA2383		560 mm (22 inch) MAC BH		○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
JAA2385		90 mm (3.5 inch) MAC Contour Depth LH		£100
JAA2386		150 mm (6 inch) MAC Contour Depth LH		£100
JAA2387		230 mm (9 inch) MAC Contour Depth LH		£100

MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)			
JAA2389		90 mm (3.5 inch) MAC Contour Depth RH	£100
JAA2390		150 mm (6 inch) MAC Contour Depth RH	£100
JAA2391		230 mm (9 inch) MAC Contour Depth RH	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)			
JAA2393		Backrest -25 mm (-1 inch) to seat width	○
JAA2394		Backrest matches seat width	○
JAA2395		Backrest +25 mm (+1 inch) to seat width	○
BACKREST HEIGHT (NON-RECLINE)*			
JAA2504		455 mm / 18" ONLY WITH: Back height 355mm/14", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2505		485 mm / 19" ONLY WITH: Back height 355mm/14", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2506		505 mm / 20" ONLY WITH: Back height 405mm/16", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2501		535 mm / 21" ONLY WITH: Back height 405mm/16", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2070		555 mm / 22" ONLY WITH Back height 455mm/18", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2543, 2050, 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2071		585 mm / 23" ONLY WITH Back height 455mm/18", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2543, 2050, 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2072		610 mm / 24" ONLY WITH Back height 505mm/20", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2542, 2050, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA2073		635 mm / 25" ONLY WITH Back height 505mm/20", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2542, 2050, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA3020		Powered Recline or No backrest option selected	○
BACKREST HEIGHT (ESR RECLINE SYSTEMS)*			
JAA2506		505 mm / 20" ONLY WITH: Elite Back, height 355mm/14"(2550), MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2501		535 mm / 21" ONLY WITH: Elite Back, height 355mm/14"(2550), MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2070		555 mm / 22" ONLY WITH: Back height 405mm/16", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2071		585 mm / 23" ONLY WITH: Back height 405mm/16", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2072		610 mm / 24" ONLY WITH Back height 455mm/18", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563) ○
JAA2073		635 mm / 25" ONLY WITH Back height 455mm/18", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2074		660 mm / 26" ONLY WITH Back height 505mm/20", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA2502		685 mm / 27" ONLY WITH Back height 505mm/20", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA3020		NON -RECLINE BACKREST or No backrest SELECTED	(NC 2063, 2066, 2570) ○

ARMREST*			
JAA2307		Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5"). Not compatible with reclining backrest	(NC 2063, 2066, 2308, 2309) ◆
JAA2095		Cantilever Flip-Back Arm, width, height and angle adjustable, 240 - 330 mm (9.5" - 13")	(NC 2308, 2309) ○
JAA2096		Cantilever Flip-Back Arm, width, height and angle adjustable, 340 - 405 mm (13.5-16")	(NC 2308, 2309) ○
JAA2097		Two-post Recline Arm, Flip Back, adjustable height 230 - 330 mm (9-13")	○
JAA2098		Two-post Recline Arm, Flip Back, adjustable height 320 - 405 mm (12.5-16")	○
ARMREST ACCESSORIES			
JAA2308		1" Offset Armrest Spacer (pair), for two-post recline only	(NC 2096, 2307) £145
JAA2309		1/2" Offset Armrest Spacer (pair), for two-post recline only	(NC 2096, 2307) £145
ARMPAD LH*			
JAA2012		Desk arm pad with startex cover visco foam [255mm]; 40 mm thick	◆
JAA2013		Full length arm pad with startex cover visco foam [355mm], 40 mm thick	○
JAA2014		Ergonomic arm trough; 30 mm thick	(NC 1405, 2311, 2310, 2347, 2348, 2349, 2350) £110
JAA2015		without LH arm pad	○
ARMPAD RH*			
JAA2016		Desk arm pad with startex cover visco foam [255mm]; 40 mm thick	◆
JAA2017		Full length arm pad with startex cover visco foam [355mm]; 40 mm thick	○
JAA2018		Ergonomic arm trough, 30 mm thick	(NC 1405, 2311, 2310, 2351, 2352, 2353, 2354) £110
JAA2019		without RH arm pad	○
ARMPAD ACCESSORIES & MOUNTING OPTIONS LH			
JAA2080		Flat Hand Pad Left	(NC 2013, 2015) £110
JAA2081		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Left	(NC 2013, 2015) £55
ARMPAD ACCESSORIES & MOUNTING OPTIONS RH			
JAA2083		Flat Hand Pad Right	(NC 2017, 2019) £110
JAA2084		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Right	(NC 2017, 2019) £55
MULTI-AXIX UPPER EXTREMITY SUPPORTS LH			
JAA2082		Multi Axis Upper Extremity Support (MACES) Left, lowest armpad height will be increased by 43 mm	£270
MULTI-AXIX UPPER EXTREMITY SUPPORTS RH			
JAA2085		Multi Axis Upper Extremity Support (MACES) Right, lowest armpad height will be increased by 43 mm	£270

Height setting: ARMREST LH SEAT PAN TO TOP OF THE TUBE*			
O001		Factory setting	<input type="radio"/>
O002		Customized	<input type="radio"/>
O003		230 mm	<input type="radio"/>
O004		240 mm	<input type="radio"/>
O005		250 mm	<input type="radio"/>
O006		260 mm	<input type="radio"/>
O007		270 mm	<input type="radio"/>
O008		280 mm	<input type="radio"/>
O009		290 mm	<input type="radio"/>
O010		300 mm	<input type="radio"/>
O011		310 mm	<input type="radio"/>
O012		320 mm	<input type="radio"/>
O013		330 mm	<input type="radio"/>
O014		340 mm	<input type="radio"/>
O015		350 mm	<input type="radio"/>
O016		360 mm	<input type="radio"/>
O017		370 mm	<input type="radio"/>
O018		380 mm	<input type="radio"/>
O019		390 mm	<input type="radio"/>
O020		405 mm	<input type="radio"/>

Height setting: ARMREST RH SEAT PAN TO TOP OF THE TUBE*			
O021		Factory setting	<input type="radio"/>
O022		Customized	<input type="radio"/>
O023		230 mm	<input type="radio"/>
O024		240 mm	<input type="radio"/>
O025		250 mm	<input type="radio"/>
O026		260 mm	<input type="radio"/>
O027		270 mm	<input type="radio"/>
O028		280 mm	<input type="radio"/>
O029		290 mm	<input type="radio"/>
O030		300 mm	<input type="radio"/>
O031		310 mm	<input type="radio"/>
O032		320 mm	<input type="radio"/>
O033		330 mm	<input type="radio"/>
O034		340 mm	<input type="radio"/>
O035		350 mm	<input type="radio"/>
O036		360 mm	<input type="radio"/>
O037		370 mm	<input type="radio"/>
O038		380 mm	<input type="radio"/>
O039		390 mm	<input type="radio"/>
O040		405 mm	<input type="radio"/>

ELBOW BLOCK LH PADS			
JAA2347		Elbow Block XS (110 mm wide x 83 mm high)	£475
JAA2348		Elbow Block S (135 mm wide x 90 mm high)	£475
JAA2349		Elbow Block M (145 mm wide x 95 mm high)	£475
JAA2350		Elbow Block L (165 mm wide x 110 mm high)	£475
ELBOW BLOCK LH MOUNTING HARDWARE			
JAA2361		Vertical 3"/ Horizontal 5"	○
JAA2362		Vertical 3"/ Horizontal 7"	○
JAA2363		Vertical 5"/ Horizontal 5"	○
JAA2364		Vertical 5"/ Horizontal 7"	○
ELBOW BLOCKS RH PADS			
JAA2351		Elbow Block XS (110 mm wide x 83 mm high)	£475
JAA2352		Elbow Block S (135 mm wide x 90 mm high)	£475
JAA2353		Elbow Block M (145 mm wide x 95 mm high)	£475
JAA2354		Elbow Block L (165 mm wide x 110 mm high)	£475
ELBOW BLOCK RH MOUNTING HARDWARE			
JAA2371		Vertical 3"/ Horizontal 5"	○
JAA2372		Vertical 3"/ Horizontal 7"	○
JAA2373		Vertical 5"/ Horizontal 5"	○
JAA2374		Vertical 5"/ Horizontal 7"	○
HIP SUPPORTS LH			
JAA2184		Hip Support LH, Quick Release, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
JAA2114		Hip Support LH, Quick Release, pad size 180 x 150 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
JAA2186		Hip Support LH, Quick Release, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
JAA2188		Hip Support LH, Quick Release, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
HIP SUPPORTS RH			
JAA2185		Hip Support RH, Quick Release, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
JAA2115		Hip Support RH, Quick Release, pad size 180 x 150 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
JAA2187		Hip Support RH, Quick Release, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
JAA2189		Hip Support RH, Quick Release, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS			
JAA2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JAA2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○

HEADREST*			
JAA2092		Onyx Auto Style headrest with standard mounting hardware (NC 2230, 2247, 2570)	£120
JAA2323		Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£325
JAA2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£325
JAA2322		Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£325
JAA2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£325
JAA2326		Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£325
JAA2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£440
JAA1495		Head Control or equivalent Headrest, provided independently from this order. According to the requirements of ISO 7176-19, a headrest is compulsory for using the wheelchair as vehicle seat	○
LATERAL SUPPORTS LH			
JAA2343		Lateral Support Maxx Style Quick Release, pad size [110 mm wide x 83 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£395
JAA2345		Lateral Support Maxx Style Quick Release, pad size [135 mm wide x 90 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£395
JAA2043		Lateral Support Maxx Style Quick Release, pad size [145 mm wide x 95 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£395
JAA2045		Lateral Support Maxx Style Quick Release, pad size [165 mm wide x 110 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£395
LATERAL SUPPORTS RH			
JAA2344		Lateral Support Maxx Style Quick Release, pad size [110 mm wide x 83 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£395
JAA2346		Lateral Support Maxx Style Quick Release, pad size [135 mm wide x 90 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£395
JAA2044		Lateral Support Maxx Style Quick Release, pad size [145 mm wide x 95 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£395
JAA2046		Lateral Support Maxx Style Quick Release, pad size [165 mm wide x 110 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£395
LATERAL SUPPORT ACCESSORIES			
JAA2049		Additional telescopic link for Maxx Style lateral support, 25 mm, Pair (NC 2559, 2178, 2560, 2561, 2562)	£160
CENTRE MOUNTED LEGREST*			
JAA0586		Manual Centre Mounted Legrest, angle settable Length: 310 - 410 mm (NC 0517, 1244, 0579, 0584, 0572, 0577)	◆
JAA0589		LNx Powered Centre Mounted Legrest, length compensating; Length: 340 - 410 mm (NC 2262, 2263, 2264, 1244, 0579, 0584, 0572, 0577)	£1,620
JAZ0592		LNx Powered Centre Mounted Legrest with telescoping foot platform (lowers to the ground; Length: 380 - 419 mm). Selectable length depends on STF height. Custom-made, requires customized setting! Custom-made, requires customized setting! (NC 0516, 0517, 2056, 2261, 2262, 2263, 2264, 2365, 2366, 1244, 0579, 0584, 0572, 0577)	£2,450
JAX3003		Without Centre Mount Legrests	○

FOOTPLATES / FOOT PLATFORM*			
JAA2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5") (NC 3003)	£105
JAA0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5") (NC 3003)	○
JAA0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5") (NC 3003)	£135
JAA2366		Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5") (NC 0586, 3003)	£270
JAA0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75") (NC 0586, 3003)	£270
JAX3005		Please select If Centre mount footPlate/ platform is not required. (NC 0589)	○
Height setting: CENTRE MOUNT LEGREST, OPTIONS 0586 & 0589 . SEAT PAN TO FOOTPLATE*			
0041		Factory setting	○
0042		Customized	○
0043		310 mm	○
0044		325 mm	○
0045		335 mm	○
0046		345 mm	○
0047		355 mm	○
0048		364 mm	○
0049		375 mm	○
0050		385 mm	○
0051		410 mm	○
0077		Please select if CENTRE MOUNT LEGREST OPTIONS 0586 & 0589 is not required	○
Height setting: CENTRE MOUNT LEGREST with TELESCOPING FOOT PLATFORM . SEAT PAN TO FOOTPLATE*			
0052		Factory setting	○
0053		Customized	○
0054		380 mm	○
0055		393 mm	○
0056		406 mm	○
0057		419 mm	○
0076		Please select if CENTRE MOUNT LEGREST with TELESCOPING FOOT PLATFORM is not required	○
SIDE MOUNTED HANGERS LH*			
JAA0539		Hanger ADE 0° - 80° LH (NC 2055, 2261, 2262, 2263, 2264)	£330
JAA0572		Pivot Plus Manual Elevating Legrest LH	£675
JAA0577		Pivot Plus Power Elevating Legrest LH (NC 2261, 2262, 2263, 2264)	£1,050
JAA0535		Without Hanger LH	○
SIDE MOUNTED HANGERS RH*			
JAA0549		Hanger ADE 0° - 80° RH (NC 2261, 2262, 2263, 2264)	£330
JAA0579		Pivot Plus Manual Elevating Legrest RH	£675
JAA0584		Pivot Plus Power Elevating Legrest RH (NC 2261, 2262, 2263, 2264)	£1,050
JAA0545		Without Hanger RH	○
INDIVIDUAL FOOTPLATE LH*			
JAA0509		Angle Adjustable Footplate LH	£75
JAA0511		Angle and depth adjustable Footplate LH	£60
JAX3012		Please select if Side mounted Legrest are not required	○

INDIVIDUAL FOOTPLATE RH*		
JAA0506	<input type="checkbox"/>	Angle Adjustable Footplate RH £75
JAA0508	<input type="checkbox"/>	Angle and depth adjustable Footplate RH £60
JAX3013	<input type="checkbox"/>	Please select if Side mounted Legrest are not required <input type="radio"/>
Height setting: VARI LEGREST LH SEAT PAN TO FOOTPLATE*		
0080	<input type="checkbox"/>	Factory Setting <input type="radio"/>
0081	<input type="checkbox"/>	Customized <input type="radio"/>
0082	<input type="checkbox"/>	270 mm <input type="radio"/>
0083	<input type="checkbox"/>	280 mm <input type="radio"/>
0084	<input type="checkbox"/>	290 mm <input type="radio"/>
0085	<input type="checkbox"/>	300 mm <input type="radio"/>
0086	<input type="checkbox"/>	310 mm <input type="radio"/>
0087	<input type="checkbox"/>	320 mm <input type="radio"/>
0088	<input type="checkbox"/>	330 mm <input type="radio"/>
0089	<input type="checkbox"/>	340 mm <input type="radio"/>
0090	<input type="checkbox"/>	350 mm <input type="radio"/>
0091	<input type="checkbox"/>	360 mm <input type="radio"/>
0092	<input type="checkbox"/>	370 mm <input type="radio"/>
0093	<input type="checkbox"/>	380 mm <input type="radio"/>
0094	<input type="checkbox"/>	390 mm <input type="radio"/>
0095	<input type="checkbox"/>	400 mm <input type="radio"/>
0096	<input type="checkbox"/>	410 mm <input type="radio"/>
0097	<input type="checkbox"/>	420 mm <input type="radio"/>
0098	<input type="checkbox"/>	430 mm <input type="radio"/>
Height setting: VARI LEGREST RH SEAT PAN TO FOOTPLATE*		
0099	<input type="checkbox"/>	Factory Setting <input type="radio"/>
0100	<input type="checkbox"/>	Customized <input type="radio"/>
0101	<input type="checkbox"/>	270 mm <input type="radio"/>
0102	<input type="checkbox"/>	280 mm <input type="radio"/>
0103	<input type="checkbox"/>	290 mm <input type="radio"/>
0104	<input type="checkbox"/>	300 mm <input type="radio"/>
0105	<input type="checkbox"/>	310 mm <input type="radio"/>
0106	<input type="checkbox"/>	320 mm <input type="radio"/>
0107	<input type="checkbox"/>	330 mm <input type="radio"/>
0108	<input type="checkbox"/>	340 mm <input type="radio"/>
0109	<input type="checkbox"/>	350 mm <input type="radio"/>
0110	<input type="checkbox"/>	360 mm <input type="radio"/>
0111	<input type="checkbox"/>	370 mm <input type="radio"/>
0112	<input type="checkbox"/>	380 mm <input type="radio"/>
0113	<input type="checkbox"/>	390 mm <input type="radio"/>
0114	<input type="checkbox"/>	400 mm <input type="radio"/>
0115	<input type="checkbox"/>	410 mm <input type="radio"/>
0116	<input type="checkbox"/>	420 mm <input type="radio"/>
0117	<input type="checkbox"/>	430 mm <input type="radio"/>

Height setting: PIVOT PLUS LEGREST LH SEAT PAN TO FOOTPLATE*			
0058		Factory setting	○
0059		Customized	○
0060		402 mm	○
0061		414 mm	○
0062		426 mm	○
0063		438 mm	○
0064		450 mm	○
0065		462 mm	○
0066		474 mm	○
0080		Please select if PIVOT PLUS LEGREST is not required	○
Height setting: PIVOT PLUS LEGREST RH SEAT PAN TO FOOTPLATE*			
0067		Factory setting	○
0068		Customized	○
0069		402 mm	○
0070		414 mm	○
0071		426 mm	○
0072		438 mm	○
0073		450 mm	○
0074		462 mm	○
0075		474 mm	○
0079		Please select if PIVOT PLUS LEGREST is not required	○
CONTROLS*			
JAA2203		DLX-REM211 (Drive & Actuators)	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2270, 2272, 1908, 1909, 2209, 2210) ◆
JAA2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2270, 2272, 2209, 2210) £200
JAA2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2210) £425
JAA2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1207, 1208, 1707, 1708, 1709, 1710) £765
REM400 TOGGLE OPTION			
JAA2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2204, 2206) £315

CONTROL BRACKETRY*				
JAA1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404, 2311)	○
JAA1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405, 2310)	○
JAA1707		Swing Away RH Bracket	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404, 2311)	£135
JAA1708		Swing Away LH Bracket	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405, 2310)	£135
JAA1709		Matrx Resolve Swing in/out Bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 2235, 2267, 2241, 1404)	£245
JAA1710		Matrx Resolve Swing in/out Bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 2236, 2268, 2242, 1405)	£245
JAA2209		REM500 Swing Away RH Multi-Adjustable Bracket	(NC 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 2311)	○
JAA2210		REM500 Swing Away LH Multi-Adjustable Bracket	(NC 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405, 2310)	○
JAA2075		Swing Away + height / width adjust. Mount on RH for user, Quad link design		£225
JAA2076		Swing Away + height / width adjust. Mount on LH for user, Quad link design		£225
JAA2213		Centre Mount Swing Away to RH Nucleus Bracket	(NC 2283, 2287, 2207, 2265, 2275, 2217, 1404, 1405, 2310, 2311, 1415)	£495
JAA2214		Centre Mount Swing Away to LH Nucleus Bracket	(NC 2284, 2288, 2208, 2266, 2276, 2218, 1404, 1405, 2310, 2311, 1339)	£495
JAA2227		Bracket RH mounted for Attendant	(NC 1221, 2206, 2219)	£65
JAA2228		Bracket LH mounted for Attendant	(NC 1222, 2206, 2219)	£65
JAZ2239		Individual mounting option per request (Special). Requires additional information/specification.		POA
DUAL CONTROLS				
JAA2219		LiNX DLX-ACU200 (Compact Dual Control)	(NC 1221, 1222, 2227, 2228, 9259)	£395
JAA2245		LiNX DLX-IDC (Intuitive Dual Control)	(NC 2259, 2260, 9259)	£745
UNZ9259		mo-vis Scoot Control (Only available via Special)	(NC 2259, 2260, 0001, 2259, 2260, 1221, 1222)	£2,315
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)				
JAA2259		Dual Control Standard Bracket RH	(NC 2245)	○
JAA2260		Dual Control Standard Bracket LH	(NC 2245)	○
JAA1221		IDC Dual Control Detachable Bracket RH	(NC 2219, 2227)	○
JAA1222		IDC Dual Control Detachable Bracket LH	(NC 2219, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification.		POA
COMPACT REMOTES				
JAA2231		LiNX DLX-CR400 (Compact Remote)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288)	£825
JAA2232		LiNX DLX-CR400LF (Compact Remote Low Force)	(NC 2203, 2204, 2285, 2286, 2283, 2284, 2287, 2288)	£1,075
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£2,015
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 0001, 2207, 2208, 2265, 2266, 2285, 2286, 2283, 2284, 2287, 2288)	£3,190
JAA2233		ASL MEC (Proportional MEC Joystick ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2207, 2208, 2265, 2266, 2215)	£4,100
JAA2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2207, 2208, 2265, 2266, 2215)	£2,285
JAA2243		ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)	(NC 2203, 2204, 2255, 2256, 2207, 2208, 2265, 2266, 2215, 2285, 2286, 2283, 2284, 2287, 2288)	£2,365

COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)				
JAA2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404, 2311)	○
JAA2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405, 2310)	○
JAA2265		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 2310, 2311, 1415)	£495
JAA2266		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 2310, 2311, 1339)	£495
JAA2285		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JAA2286		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JAA2283		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JAA2284		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JAA2287		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£495
JAA2288		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£495
JAA2235		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311)	○
JAA2236		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311)	○
JAA2267		Centre Mount Swing Away to RH (Nucleus)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311)	£495
JAA2268		Centre Mount Swing Away to LH (Nucleus)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311)	£495
JAA2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311)	£495
JAA2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311)	£495
JAZ2215		Individual mounting option per request (Special). Requires additional information/specification.		POA
CHIN CONTROLS				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 2203, 2204)	£2,600
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,730
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 2203, 2204)	£2,320
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,520
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£1,395
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,015
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)		£960
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)		£1,380
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)		£2,280
HEAD CONTROLS				
JAA2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2570)	£4,220
JAA2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2570)	£4,945
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 0001)	£4,775
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2175, 0001)	£3,315
FOOT CONTROL				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 2203, 2204)	£2,880
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204)	£2,740
FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.		POA

TRAY CONTROL			
JAA2257	<input type="checkbox"/>	Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240) £6,535
TRAY CONTROL BRACKETRY (Only applicable if Tray Control is selected)			
JAA2275	<input type="checkbox"/>	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415) £495
JAA2276	<input type="checkbox"/>	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339) £495
JAZ2269	<input type="checkbox"/>	Individual mounting option per request (Special). Requires additional information/specification.	POA
SIP N PUFF HEAD CONTROL			
JAA2255	<input type="checkbox"/>	Sip'n Puff Head Control (ASL109), includes Input Module and Lip Switch (ASL 301 LIP Combines simple sip and puff controls with head movements	(NC 2203, 2204, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2243) £6,120
SIP N PUFF CONTROL & INPUT MODULE			
JAA2256	<input type="checkbox"/>	Sip'n Puff Control includes Sip'n Puff kit, Lip Switch ASL301LIP and Input module	(NC 2203, 2204, 2243) £855
WIRELESS ACCESSORIES FOR ASL CONTROLS			
JAA2249	<input type="checkbox"/>	Wireless Triple Switch Receiver (ASL557-3)	£340
JAA2250	<input type="checkbox"/>	Wireless Mouse Emulator (ASL558)	£810
LINX INTERFACE OPTIONS			
JAA2271	<input type="checkbox"/>	DLX-IN500 Module (sub D9, Sip'n Puff)	(NC 2255, 2256, 2203, 2204) £785
JAA2273	<input type="checkbox"/>	DLX-IN200 Module (Jack port for external power button + control inputs via DB9)	£270
JAA2272	<input type="checkbox"/>	DLX-OUT500 (allows environmental control access)	(NC 2204) £785
JAA1681	<input type="checkbox"/>	5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£120
SPECIAL ALTERNATIVE JOYSTICKS & INPUTS			
JAZ2216	<input type="checkbox"/>	Special alternative Joystick/Input and bracket as per request	POA
LINX 10-WAY SWITCH			
JAA2251	<input type="checkbox"/>	10-Way Switch DLX-SW10-A for direct Powered Seating Access	£500
LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)			
JAA2217	<input type="checkbox"/>	Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404) ○
JAA2218	<input type="checkbox"/>	Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405) ○
JAA2226	<input type="checkbox"/>	For attendant, mounted on back shroud	(NC 2218) £65
TOGGLES / BUTTONS			
JAA2261	<input type="checkbox"/>	Stereo Toggle - direct control for tilt-only	(NC 0584, 0589, 0592, 2053, 2060, 2063, 2066) £270
JAA2262	<input type="checkbox"/>	Stereo Button - direct control for tilt-only	(NC 0584, 0589, 0592, 2053, 2060, 2063, 2066) £255
JAA2263	<input type="checkbox"/>	Four-Way Toggle - direct control for tilt & recline	(NC 2069, 0577, 0584, 0589, 0592) £330
JAA2264	<input type="checkbox"/>	Four-Way Button - direct control for tilt & recline	(NC 2069, 0577, 0584, 0589, 0592) £355
TOGGLES / BUTTONS MOUNT			
JAA2279	<input type="checkbox"/>	RH Mount (of Toggle/Button)	○
JAA2280	<input type="checkbox"/>	LH Mount (of Toggle/Button)	○
JAA2278	<input type="checkbox"/>	Without mounting	○
LINX USB CHARGER PORT			
JAA2253	<input type="checkbox"/>	USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £160
LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)			
JAA2281	<input type="checkbox"/>	RH Mount	○
JAA2282	<input type="checkbox"/>	LH Mount	○

FURTHER LINX COMPONENTS / ACCESSORIES			
JAA2302		Not Necessary	(NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716) ◆
JAA2300		To be enclosed only	○
JAZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JAA1685		Red Buddy Button	£90
JAA1686		Yellow Buddy Button	£90
JAA1687		Green Buddy Button	£90
JAA1690		Blue Buddy Button	£90
EGG SWITCH			
JAA2254		Egg Switch, black	£90
PIKO BUTTONS			
JAA1711		Red Piko Button	£160
JAA1712		Yellow Piko Button	£160
JAA1713		Green Piko Button	£160
JAA1714		Orange Piko Button	£160
JAA1715		Blue Piko Button	£160
JAA1716		Black Piko Button	£160
BODYPOINT JOYSTICK KNOBS			
JAA2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£75
JAA2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£75
JAA2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
JAA2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75
SHROUD COLOUR*			
JAA2871		Mid range - Black/Anthracite Shroud	○
JAA2872		Mid range - Silver Shroud	○
JAA2873		Mid range - Bengal Red Shroud	○
JAA2874		Mid range - Mercury Blue Shroud	○
JAA2875		Mid range - Midnight Purple Shroud	○
JAZ0299		Premium - Individual Shroud Colour. RAL number required Extended lead lime: + 15 days.	£300
RIM COLOUR*			
JAA2882		Mid range - Black	○
JAA2883		Mid range - Silver	○
JAA2884		Mid range - Bengal Red	○
JAA2885		Mid range - Mercury Blue	○
JAA2886		Mid range - Midnight Purple	○
JAZ2887		Premium - Individual Rim Insert Colour. RAL number required Extended lead lime: + 15 days. Price Included with Premium Shroud Kit	○
FORKS*			
JAA2631		Single sided fork	£250
JAA2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆

CASTORS*			
JAA0835	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	◆
JAA0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JAA0838	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JAA0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)	○
CASTOR LOCKS			
JAA1490	<input type="checkbox"/>	Castor Lock	£195
DRIVE WHEELS*			
JAA0704	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JAA0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
JAA0706	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed	○
C.T.C SUSPENSION*			
JAA2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JAA1910	<input type="checkbox"/>	Soft	○
JAA1912	<input type="checkbox"/>	Medium	○
JAA1911	<input type="checkbox"/>	Firm	○
JAA1915	<input type="checkbox"/>	Very firm	○
MOTOR TYPE*			
JAA1294	<input type="checkbox"/>	Indoor setting, reduced speed approx. 3 km/h	○
JAA1282	<input type="checkbox"/>	6 km/h	◆
JAA1908	<input type="checkbox"/>	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 1437, 1254)	£515
JAA1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 1437, 1254, 0828)	£600
BATTERY TYPE*			
JAA1259	<input type="checkbox"/>	60 Ah Gel	◆
JAA1260	<input type="checkbox"/>	73.5 Ah Gel (NC 1450)	£100
G TRAC			
JAA1215	<input type="checkbox"/>	G-Trac Module	£400
CHARGER*			
JAA1343	<input type="checkbox"/>	OptiCharge Charger 8A (IP44)	◆
JAA1313	<input type="checkbox"/>	OptiCharge Charger 10A (IP44)	£155
ANTI-TIPPERS*			
JAA1205	<input type="checkbox"/>	Anti-tipper	◆
LIGHTS			
JAA1437	<input type="checkbox"/>	LED Light system with indicators 30 Lux headlights, included with 10km/h motors (NC 1908, 1909)	£450
JAA1254	<input type="checkbox"/>	Without Light system (NC 1909)	○
KERB RIDER			
JAA1244	<input type="checkbox"/>	Kerb Rider with park position feature (NC 0586, 0589, 0592, 1294)	£215

TRANSPORTATION KIT*			
JAA1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
JAA1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.	◆
JAA1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.	(NC 2052) ○
EASY ACCESS FRONT TIE DOWNS			
JAA2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down. For transport worthy versions (1494/1401) only	£55
FIXATION BELT*			
JAA2113		Lap belt with Auto Lock Flexible mount in side rails	◆
JAA2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) £95
JAA2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) £95
JAA2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) £90
JAA2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) £90
JAA2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) £90
JAA2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) £90
REAR VIEW MIRROR			
JAA1415		Rearview mirror RH	(NC 2265, 2267, 2275) £45
JAA1339		Rearview mirror LH	(NC 2266, 2268, 2276) £45
OTHER OPTIONS			
JAA1438		Luggage rack	£350
JAA1416		Fender for Castors	£80
JAA1417		Fender for Drive wheels	£160
JAA1404		Tray, Swing Away RH	(NC 1707, 2207, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £200
JAA1405		Tray, Swing Away LH	(NC 1708, 2208, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £200
JAA2310		Tray Flip and fold away 1/2 lap tray LH (14" W x 10" D)	(NC 1707, 2207, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £735
JAA2311		Tray Flip and fold away 1/2 lap tray RH (14" W x 10" D)	(NC 1708, 2208, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £735
JAZ2312		Vent Box, Low Profile Recline, Universal	(NC 2054, 2055, 2056) POA
JAZ2314		Vent Box, Low Profile Recline, Trilogy	(NC 2054, 2055, 2056) POA
JAZ2315		Vent Box, Low Profile Recline, LTV – Easy Fit	(NC 2054, 2055, 2056) POA
JAZ2313		Oxygen Holder, Low Profile Recline	(NC 2054, 2055, 2056) POA

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JAA0001	AVIVA RX40 Ultra
JAA2810	Narrow (560 mm/22")
JAA2099	Fixed STF height and seat angle setting
JAA2052	Ultra Low Maxx CG Powered Tilt (0 to 50°)
JAA2053	Ultra Low Maxx CG Powered Tilt (0 to 45°) combined with 300 mm Seat Lifter. (Increased STF Height by 4 mm)
JAA2857	435 mm
JAA2858	460 mm
JAA2859	485 mm
JAZ2054	Ultra Low Maxx CG tilt 5° Fixed Anterior Assist mount. Max tilt angle reduced by 5°. With centre-mounted legrests only. Front STF height as defined above
JAZ2055	Ultra Low Maxx CG tilt 10° Fixed Anterior Assist mount. Max tilt angle reduced by 10°. With centre-mounted legrests only. Front STF height as defined above.
JAZ2056	Ultra Low Maxx CG tilt 5° Fixed Posterior-tilt Assist mount. Max tilt angle increased by 5°. Front STF height increased by 6 cm over chosen SH setting.
JAA2069	Fixed Angle Setting (NO POWERED RECLINE)
JAA2060	Powered Recline 90 to 170° , Standard Mount
JAA2063	Powered Recline, 8° Precline mount, 82° to 162°
JAA2066	Powered Recline, 30° Precline mount, 60° to 140°
JAA2025	405 mm / 16"
JAA2026	430 mm / 17"
JAA2027	455 mm / 18"
JAA2028	480 mm / 19"
JAA2029	505 mm / 20"
JAA2030	530 mm / 21"
JAA2031	555 mm / 22"
JAA2035	405 mm / 16"
JAA2036	430 mm / 17"
JAA2037	455 mm / 18"
JAA2041	555 mm / 22"
JAA2042	580 mm / 23"
JAA2145	Matrix Flo-Tech Lite Black smooth cover, 75 mm thick Up to 127 kg
JAA2146	Matrix Flo-tech Contour Black smooth cover, 100 mm thick Up to 191 kg
JAA2510	Straight back cane height 455 mm / 18"
JAA2511	Straight back cane height 505 mm / 20"
JAA2512	Straight back cane height 555 mm / 22"
JAA2513	Straight back cane height 610 mm / 24"
JAA2514	Angled (5°) back cane height 505 mm / 20"
JAA2515	Angled (5°) back cane height 555 mm / 22"
JAA2516	Angled (5°) back cane height 610 mm / 24"
JAA2517	75°
JAA2518	79°
JAA2519	85°
JAA2520	90°

JAA2521	95°
JAA2522	101°
JAA2523	106°
JAA2524	112°
JAA2540	Standard Rehab Back Pan and Back Upholstery height 355 mm / 14"
JAA2541	Standard Rehab Back Pan and Back Upholstery height 405 mm / 16"
JAA2542	Standard Rehab Back Pan and Back Upholstery height 455 mm / 18"
JAA2543	Standard Rehab Back Pan and Back Upholstery height 505 mm / 20"
JAA2050	High Back, height 405 mm / 16" [Mesh fabric]
JAA2051	High Back, height 455 mm / 18" [Mesh fabric]
JAA2320	High Back, height 505 mm / 20" [Mesh fabric]
JAA2550	Matrx Elite E2 Back, height 355mm / 14" / 16 - 18" wide
JAA2551	Matrx Elite E2 Back, height 405 mm / 16"
JAA2552	Matrx Elite E2 Back, height 455 mm / 18"
JAA2553	Matrx Elite E2 Back, height 505 mm / 20"
JAA2554	Matrx Elite HD Back, height 405 mm / 16" (Extended lead time: + 4 days)
JAA2555	Matrx Elite HD Back, height 505 mm / 20" (Extended lead time: + 4 days)
JAA2556	Matrx Elite TR Back, height 405 mm / 16" [Mesh fabric]
JAA2557	Matrx Elite TR Back, height 455 mm / 18" [Mesh fabric]
JAA2177	Matrx Elite TR Back, height 505 mm / 20" [Mesh fabric]
JAA2563	Matrx Elite TR HD Back, height 505 mm / 20" [Mesh fabric] (Extended lead time: + 4 days)
JAA2558	Matrx Elite Deep E2 Back, height 405 mm / 16" [Mesh fabric]
JAA2559	Matrx Elite Deep E2 Back, height 455 mm / 18" [Mesh fabric]
JAA2178	Matrx Elite Deep E2 Back, height 505 mm / 20" [Mesh fabric]
JAA2174	Matrx PB Back, height 405 mm / 16" [Mesh fabric]
JAA2175	Matrx PB Back, height 505 mm / 20" [Mesh fabric]
JAA2560	Matrx PB Deep Back, height 405mm / 16"
JAA2561	Matrx PB Deep Back, height 505 mm / 20"
JAA2562	Matrx Elite Deep HD Back, height 505 mm / 20" (Extended lead time: + 4 days)
JAA2570	Without Backrest Support, incl. Matrx interface
JAA2307	Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5"). Not compatible with reclining backrest
JAA2096	Cantilever Flip-Back Arm, width, height and angle adjustable, 340 - 405 mm (13.5-16")
JAA2097	Two-post Recline Arm, Flip Back, adjustable height 230 - 330 mm (9-13")
JAA2098	Two-post Recline Arm, Flip Back, adjustable height 320 - 405 mm (12.5-16")
JAA2308	1" Offset Armrest Spacer (pair), for two-post recline only
JAA2309	1/2" Offset Armrest Spacer (pair), for two-post recline only
JAA2013	Full length arm pad with startex cover visco foam [355mm], 40 mm thick
JAA2015	without LH arm pad
JAA2017	Full length arm pad with startex cover visco foam [355mm]; 40 mm thick
JAA2019	without RH arm pad

JAA2347	Elbow Block XS (110 mm wide x 83 mm high)
JAA2348	Elbow Block S (135 mm wide x 90 mm high)
JAA2349	Elbow Block M (145 mm wide x 95 mm high)
JAA2350	Elbow Block L (165 mm wide x 110 mm high)
JAA2351	Elbow Block XS (110 mm wide x 83 mm high)
JAA2352	Elbow Block S (135 mm wide x 90 mm high)
JAA2353	Elbow Block M (145 mm wide x 95 mm high)
JAA2354	Elbow Block L (165 mm wide x 110 mm high)
JAA2114	Hip Support LH, Quick Release, pad size 180 x 150 mm
JAA2186	Hip Support LH, Quick Release, pad size 230 x 100 mm
JAA2188	Hip Support LH, Quick Release, pad size 330 x 100 mm
JAA2185	Hip Support RH, Quick Release, pad size 140 x 100 mm
JAA2115	Hip Support RH, Quick Release, pad size 180 x 150 mm
JAA2187	Hip Support RH, Quick Release, pad size 230 x 100 mm
JAA2189	Hip Support RH, Quick Release, pad size 330 x 100 mm
JAA2092	Onyx Auto Style headrest with standard mounting hardware
JAA2323	Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2322	Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23), complete with Elan Headrest Hardware (MEHW)
JAA2326	Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover complete with Elan Headrest Hardware (MEHW)
JAA2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover complete with Elan Headrest Hardware (MEHW)
JAA0586	Manual Centre Mounted Legrest, angle settable Length: 310 - 410 mm
JAA0589	LNx Powered Centre Mounted Legrest, length compensating; Length: 340 - 410 mm
JAZ0592	LNx Powered Centre Mounted Legrest with telescoping foot platform (lowers to the ground; Length: 380 - 419 mm). Selectable length depends on STF height. Custom-made, requires customized setting! Custom-made, requires customized setting!
JAX3003	Without Centre Mount Legrests
JAA2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JAA0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
JAA2366	Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")
JAA0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JAA0572	Pivot Plus Manual Elevating Legrest LH
JAA0577	Pivot Plus Power Elevating Legrest LH
JAA0579	Pivot Plus Manual Elevating Legrest RH
JAA0584	Pivot Plus Power Elevating Legrest RH
JAA2203	DLX-REM211 (Drive & Actuators)
JAA2204	DLX-REM216 (Drive, Actuators & Lights)
JAA2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
JAA2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
JAA2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
JAA1207	Standard Bracket RH
JAA1208	Standard Bracket LH

JAA1707	Swing Away RH Bracket
JAA1708	Swing Away LH Bracket
JAA1709	Matrx Resolve Swing in/out Bracket RH
JAA1710	Matrx Resolve Swing in/out Bracket LH
JAA2209	REM500 Swing Away RH Multi-Adjustable Bracket
JAA2210	REM500 Swing Away LH Multi-Adjustable Bracket
JAA2213	Centre Mount Swing Away to RH Nucleus Bracket
JAA2214	Centre Mount Swing Away to LH Nucleus Bracket
JAA2227	Bracket RH mounted for Attendant
JAA2228	Bracket LH mounted for Attendant
JAA2219	LiNX DLX-ACU200 (Compact Dual Control)
JAA2245	LiNX DLX-IDC (Intuitive Dual Control)
UNZ9259	mo-vis Scoot Control (Only available via Special)
JAA2259	Dual Control Standard Bracket RH
JAA2260	Dual Control Standard Bracket LH
JAA1221	IDC Dual Control Detachable Bracket RH
JAA1222	IDC Dual Control Detachable Bracket LH
JAA2231	LiNX DLX-CR400 (Compact Remote)
JAA2232	LiNX DLX-CR400LF (Compact Remote Low Force)
JAA2233	ASL MEC (Proportional MEC Joystick ASL130)
JAA2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
JAA2243	ASL Paediatric (Proportional Joystick ASL132, incl. DLX-IN500)
JAA2207	RH Mount
JAA2208	LH Mount
JAA2265	Centre Mount Swing Away to RH (Nucleus)
JAA2266	Centre Mount Swing Away to LH (Nucleus)
JAA2285	Tray Mount RH (ASL603, 100mm x 150mm)
JAA2286	Tray Mount LH (ASL603, 100mm x 150mm)
JAA2283	Centre Mount Swing Away to RH (Nucleus)
JAA2284	Centre Mount Swing Away to LH (Nucleus)
JAA2287	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAA2288	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAA2235	Tray Mount RH (ASL603, 100mm x 150mm)
JAA2236	Tray Mount LH (ASL603, 100mm x 150mm)
JAA2267	Centre Mount Swing Away to RH (Nucleus)
JAA2268	Centre Mount Swing Away to LH (Nucleus)
JAA2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAA2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
JAZ2215	Individual mounting option per request (Special). Requires additional information/specification.
JAA2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
JAA2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
JAA2257	Four-Switch Proximity Array (ASL106)
JAA2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JAA2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)

JAA2255	Sip'n Puff Head Control (ASL109), includes Input Module and Lip Switch (ASL 301 LIP Combines simple sip and puff controls with head movements)
JAA2256	Sip'n Puff Control includes Sip'n Puff kit, Lip Switch ASL301LIP and Input module
JAA2271	DLX-IN500 Module (sub D9, Sip'n Puff)
JAA2272	DLX-OUT500 (allows environmental control access)
JAA2217	Bracket RH, mounted for user
JAA2218	Bracket LH, mounted for user
JAA2261	Stereo Toggle - direct control for tilt-only
JAA2262	Stereo Button - direct control for tilt-only
JAA2263	Four-Way Toggle - direct control for tilt & recline
JAA2264	Four-Way Button - direct control for tilt & recline
JAA2253	USB Charger Port DLX-USB02-A
JAA1686	Yellow Buddy Button
JAA1687	Green Buddy Button
JAA1690	Blue Buddy Button
JAA2254	Egg Switch, black
JAA1711	Red Piko Button
JAA1712	Yellow Piko Button
JAA1713	Green Piko Button
JAA1714	Orange Piko Button
JAA1715	Blue Piko Button
JAA1716	Black Piko Button
JAA0828	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
JAA1294	Indoor setting, reduced speed approx. 3 km/h
JAA1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JAA1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JAA1437	LED Light system with indicators 30 Lux headlights, included with 10km/h motors
JAA1254	Without Light system
JAA1244	Kerb Rider with park position feature
JAA1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.
JAA1494	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory. Not compatible with Seat lifter.
JAA1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Required for variants with combined Lift & Tilt module. Adequate headrest is compulsory.
JAA2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JAA2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JAA2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JAA2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JAA2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JAA2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JAA1415	Rearview mirror RH
JAA1339	Rearview mirror LH

JAA1404	Tray, Swing Away RH
JAA1405	Tray, Swing Away LH
JAA2310	Tray Flip and fold away 1/2 lap tray LH (14" W x 10" D)
JAA2311	Tray Flip and fold away 1/2 lap tray RH (14" W x 10" D)
JAZ2312	Vent Box, Low Profile Recline, Universal
JAZ2314	Vent Box, Low Profile Recline, Trilogy
JAZ2315	Vent Box, Low Profile Recline, LTV – Easy Fit
JAZ2313	Oxygen Holder, Low Profile Recline



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):