



AVIVA RX10

POWER SERIES



Access all areas

Ultra-compact rear wheel drive powerchair



| AVIVA RX10

Control, confidence and independence every day.

Introducing the AVIVA RX10 ultra-compact rear-wheel drive power chair designed to support effortless control in any environment.

Narrow, intuitive and reassuring to drive, it's built to deliver comfort and reliable performance.



| AVIVA RX10

Design. Technology. Performance.

The new exciting member of the **AVIVA RX** range. The RX10's compact base supports stable, controlled performance, helping you drive with confidence and in comfort everyday.



539 mm*

Driving unit width



Ultra-compact design

With the shortest and narrowest wheelbase in the series, the **AVIVA RX10** negotiates indoor spaces with ease. Moving from room to room, around corners or turning 360°, the ultra-compact design makes manoeuvrability effortless in the tightest of spaces.



Comfort and control

14" drive wheels, a choice of **8 or 9"** castors and a 75A power module ensure you drive in comfort and control at all times. From indoor surfaces to outdoor terrain, the **RX10** offers a simple mobility solution without compromise.



Minimalist design

With simplicity at its core, the minimalist design of the **AVIVA RX10** has specifically been developed, to give you discreet and dependable mobility wherever you go.





| AVIVA RX10

Easy set up with *Action*® seating.

The **AVIVA RX10** has been designed in combination with the market leading Action seating options. Making individual positioning and adjustments quick and easy for optimal comfort.



0-20°

Seat tilt

Powered available as an option

0-30°

Manual backrest recline

Powered available as an option



Choice of seat widths

430 or 480 mm plus each size is adjustable +25mm per side on each armrest and legrest.



Tension adjustable backrest

contours individual shape for optimal comfort and support.



60 mm foam seat cushion

for added comfort.



Swing-away 80° legrest for ease of transfer

Adjusts in length 245 - 370 mm.



Height and depth adjustable armrest

240 - 315 mm height up to 90 mm depth.



Adjustable seat depth

410 - 510 mm.



| AVIVA RX10

Precise control with **LiNX**® technology.

Designed for ease-of-use and precise control, the **AVIVA RX10** features LiNX technology to help make everyday navigation easy.



LED Battery gauge

gives reassurance of charge wherever you go.



Speed dial

for independent and tailored control.



Swing-away control arm

making joystick positioning and access up to and under desks easy.

Depending on configuration the RX10 is available with LiNX REM060, REM150, REM210 or REM215.



AVIVA RX10

Technical data.

For further information about this product, including the user manual, please visit our website www.invacare.eu.com.

Colour



Black

804 - 809 mm*

Turning radius
*depending on the configuration



Seat width

430 - 480 mm
+ 25 mm
adjustment
per side in the
armrests and
legrests.



Seat depth

410 - 510 mm
In 25 mm
increments.



Seat height

451 - 471 mm



Backrest height

510 mm



Armrest height

240 - 315 mm



Footrest length

245 - 370 mm



Seat angle

0° - 20°
(manual)
0° - 20°
(powered)



Backrest angle

90° - 120°



Total height

1350 mm



Total length inc.
footrests

1035 mm



Total length
without footrests

824 mm



Driving unit
width

539 mm



Max. user
weight

136 kg



Battery
capacity

51Ah (Gel)



Max. climbable
obstacle height

50 mm



Max. safe slope

6° (10.5%)



Speed

6 or 10 km/h



Driving range

29.2 km (6 km/h)
25.6 km (10 km/h)



Total weight

113.6 kg



Drive wheels

14"



Front castors

8 or 9"

AVIVA RX

AVIVA FX

TDX SP2

Esprit Action

Kompas

Storm 4

P O W E R S E R I E S

Tel: +41 61 487 70 70

hqeuropa@invacare.com

www.invacare.eu.com

© 2026 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. AVIVA RX10_EU-22052026

