



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Manual Retail Prescription Form

10 DECEMBER 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Action 1R Stock RRP

### Standard Features

- Max user weight 125kg (19.5 stone)
- Crash tested ISO 7176-Part 19
- Transit (attendant) propulsion or self (occupant) propulsion
- Seat width range 15" - 20" (380mm - 505mm)
- Seat depth 16" (405mm) or 17" (430mm) hook and loop adjustability
- Seat Height 20" (510mm)
- Half folding 7° backrest 16.5" (420mm)
- Black nylon upholstery
- Swing away and detachable legrest hangers with flip up footplates & heel loops
- Full length, fixed height, swing away armrest
- 24" x 1 3/8" quick release rear wheel, black solid tyre or
- 12" x 2" fixed axle, black solid tyre
- 8" x 1 1/4" (200mm x 25mm) solid front castor
- Lapbelt (adjustable both sides)
- Push to lock brakes for Self propel wheels
- Footbrakes for Transit wheel



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>Action 1 Transit</b>		
DK1681332	380mm (15")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681333	405mm (16")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681334	430mm (17")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681335	455mm (18")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681336	480mm (19")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681337	510mm (20")SW x 405mm to 430mm (16"-17") ADJ SD	£290

Action 1 Self Propel			
DK1681326		380mm (15")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681327		405mm (16")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681328		430mm (17")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681329		455mm (18")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681330		480mm (19")SW x 405mm to 430mm (16"-17") ADJ SD	£290
DK1681331		510mm (20")SW x 405mm to 430mm (16"-17") ADJ SD	£290

			<b>TOTAL PRICE</b>	
--	--	--	--------------------	--

**Non Conforming Options:**



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):