



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

APRIL 2024
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 2NG Stock (New - Transit Lite)

Standard Features

- A2 Transit Lite
- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Specific Transit Frame
- Anthracite Black
- Choice of 4 seat widths 16"-19" (406mm-457mm)
- Depth adjustable 16"-17" (406mm-432mm)
- Black padded nylon seat upholstery
- Half folding height adjustable backrest 17"-20" (430mm-510mm)
- Lap belt
- Swing in/out 80° detachable hanger with adjustable heel strap
- Black 7" x 1.1/4"(178mm x 32mm) welded castor housing with aluminum fork
- Attendant brake with parking function
- Push to lock brakes
-
- Net Price: £195.00
- Price Break - Transit Version only
- 1-2 units - £195.00
- 3-4 units - £185.00
- 5+ units - £175.00



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
A2NG TRANSIT LITE		
BQ1673109	A2NG 15" W X 16" to 17" D	£195
BQ1673110	A2NG 16" W X 16" to 17" D	£195
BQ1673111	A2NG 17" W X 16" to 17" D	£195
BQ1673112	A2NG 18" W X 16" to 17" D	£195
BQ1673113	A2NG 19" W X 16" to 17" D	£195

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):