



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 2NG Stock spec 1

Standard Features

- Action 2 Self Propel and Transit Version
- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 20
- Lap belt (adjustable both sides)
- Anti-tippers (24" wheels only)
- Anthracite Black
- Choice of 5 seat widths 15"-19" (381mm-457mm)
- Depth adjustable 16"-17" (406mm-432mm)
- Black padded nylon seat upholstery
- Half folding 7° backrest 17" (430mm)
- Swing in/out 80° detachable hanger with adjustable heel strap
- Removable, lightweight and height adjustable full arm (200mm > 260mm)
- Black 24" x 1.3/8" quick release rear wheel solid tyre or
- 12" X 2" fixed axle, black solid tyre wheel forward position
- Black 8" x 1.1/4" (200mm x 25mm) solid front castor aluminium fork
- Push to lock brakes



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
A2NG SELF PROPEL		
BQ1529465-AF	A2NG 15" W X 16" to 17" D	€348
BQ1529466-AF	A2NG 16" W X 16" to 17" D	€348
BQ1529469-AF	A2NG 17" W X 16" to 17" D	€348
BQ1529470-AF	A2NG 18" W X 16" to 17" D	€348
BQ1529485-AF	A2NG 19" W X 16" to 17" D	€348

A2NG TRANSIT MK 1

BQ1529486- AF		A2NG 15" W X 16" to 17" D	€348
BQ1529487- AF		A2NG 16" W X 16" to 17" D	€348
BQ1529490- AF		A2NG 17" W X 16" to 17" D	€348
BQ1529491- AF		A2NG 18" W X 16" to 17" D	€348
BQ1529492- AF		A2NG 19" W X 16" to 17" D	€348

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):