



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

DECEMBER 2025  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Action 3NG Comfort

### Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion or Attendant Propulsion
- Choice of 4 seat widths (16" - 20") 405mm to 505mm
- Choice of 2 seat depths (16" or 18") 400mm or 450mm
- Choice of 3 seat height (16" - 18") 410mm to 460mm
- Black nylon seat upholstery
- Reclining (Ratchet) backrest 0° - 30° with Laguna shell 20° (510mm)
- Comfort Black multi stretch PU coated fabric Passad 2 backrest
- Elevating legrest hangers 0° to 75° including calf pad
- Flip up, angle and depth adjustable footplates
- Removable, height and depth adjustable comfort armrests
- Black multi stretch PU coated fabric trunk supports and headrest
- Choice of 8 colours (with 7 special colour options)
- Push or Pull to lock brakes (configuration dependant)
- Comfort transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

**Standard Base price: €1,498**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description		
<b>BASE PRODUCT*</b>			
CEG0005	ACTION 3NG COMFORT	(NC 0852)	<input type="radio"/>
CEGSPEC	SPECIAL		POA
<b>PROPULSION*</b>			
CEG0100	SELF PROPULSION	(NC 1062, 1074)	€1,498
CEG0110	TRANSIT PROPULSION	(NC 0290, 0300, 1068, 1150, 1180, 1260, 1261, 1275, 1390)	€1,498

<b>SEAT WIDTH*</b>			
CEG0190		16" / 405 MM SEAT WIDTH	○
CEG0200		17" / 430 MM SEAT WIDTH	◆
CEG0210		18" / 455 MM SEAT WIDTH	○
CEG0230		20" / 505 MM SEAT WIDTH	○
<b>SEAT DEPTH*</b>			
CEG0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEG0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
<b>SEAT HEIGHT*</b>			
CEG0280		18" / 460 MM SEAT HEIGHT	◆
CEG0290		17" / 435 MM SEAT HEIGHT (NC 1010, 1062)	○
CEG0300		16" / 410 MM SEAT HEIGHT (NC 1010, 1062, 1370, 1440, 1441)	○
<b>SEAT UPHOLSTERY*</b>			
CEG0340		NYLON WITH HOOK AND LOOP STRIPS	◆
CEG0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM	○
CEG0372		CUSHION COMFORT MISTRAL 2 LAYER (INCLUDING SEAT UPHOLSTERY 3 VELCROS) (NC 0852, 0869, 0871, 0873, 0874, 0876)	○
<b>BACKREST*</b>			
CEG0441		RECLINING (RATCHET) 0/30° LAGUNA SHELL 20"	◆
<b>BACKREST UPHOLSTERY*</b>			
CEG0490		COMFORT BACKREST LAGUNA BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT	○
CEG0491		COMFORT BACKREST PASSAD 2 BLACK MULTI STRETCH PU COATED 20" / 510 MM HEIGHT	◆
<b>LEGREST HANGER RH*</b>			
CEG0520		ELEVATING LEGREST RH SWING IN/SWING OUT (NC 66666)	◆
CEG0545		AMPUTEE LEGREST RH (NC 0610, 1111)	€150
<b>LEGREST HANGER LH*</b>			
CEG0570		ELEVATING LEGREST LH SWING IN/SWING OUT (NC 55555)	◆
CEG0595		AMPUTEE LEGREST LH (NC 0630, 2222)	€150
<b>FOOTREST RH*</b>			
CEG0610		FLIP UP DEPTH AND ANGLE ADJ RH	◆
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) (NC 0545)	◆
CEE66666		NO FOOTREST	○
<b>FOOTREST LH*</b>			
CEG0630		FLIP UP DEPTH AND ANGLE ADJ LH	◆
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY) (NC 0595)	◆
CEE55555		NO FOOTREST	○
<b>ARMREST RH*</b>			
CEG0691		REMOVABLE RH (N°4) COMFORT ARM	◆
CEG0732		REMOVABLE RH (N°4) HEMI ARM (NC 1481)	€69
<b>ARMREST LH*</b>			
CEG0791		REMOVABLE LH (N°4) COMFORT ARM	◆
CEG0832		REMOVABLE LH N°4 HEMI ARM (NC 1481)	€69

SEAT CUSHION			
CEG0869		MATRX FLO-TECH LITE CUSHION - Low Risk	€168
CEG0871		MATRX FLO-TECH CONTOUR VISCO CUSHION NG - HIGH RISK	€236
CEG0873		MATRX FLO-TECH IMAGE CUSHION - High Risk	€382
CEG0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk	€566
CEG0876		MATRX LIBRA CUSHION - Very High Risk	€645
SEAT CUSHION DEPTH			
CEG0878		FOR 400 MM SEAT DEPTH	○
CEG0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEG0852		SEAT RIGIDIZER	€79
CEG0881		TRUNK SUPPORT COMFORT BLACK MULTI STRETCH PU COATED FABRIC (PAIR) ONLY FOR LAGUNA SHELL	◆
CEG0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC	◆
CEG0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	€189
CEG0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	€189
POSITIONING ACCESSORIES FOR MATRX LIBRA & FLOTECH			
CEG0891		TAPERED WEDGES (PAIR)	€28
CEG0892		ABDUCTOR	€58
CEG0893		LATERAL WEDGE	€46
CEG0894		ANTERIOR SEAT WEDGE	€43
COLOUR*			
CEG0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEG0964		PEARL WHITE	○
CEG0965		ELECTRIC GREEN	○
CEG0966		MATT BLACK	○
CEG0969		RED MATT	○
CEG0972		WATER BLUE	○
CEG0974		BLUE	○
CEG0976		DARK PURPLE	○
CEG0980		ANTHRACITE BLACK	○
CEG0991		SPECIAL COLOUR ON DEMAND ONLY	€238
SWO3333		DARK RED - SPECIAL	€238
SWO4444		PEARL GREY - SPECIAL	€238
SWO5555		INTENSE YELLOW	€238
SWO6666		VIBRANT ORANGE	€238
SWO7777		POLAR BLUE	€238
SWO8888		CANDY PINK	€238
CASTOR HOUSING*			
CEG0995		STANDARD	◆
CEG0996		LIGHT WEIGHT	€35

FRONT CASTORS*			
CEG1010		8" X 1.1/4" BLACK SOLID TYRE	◆
CEG1050		6" X 1.1/4" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CEG1062		12" (NC 1071, 1072, 1084, 1088, 1140, 1150, 1180, 1260, 1261, 1390)	○
CEG1068		24"	◆
REAR WHEELS*			
CEG1071		STANDARD (BASIC) 4 POINTS FIXING HANDRIM (BLACK HUB, COLOURLESS RIM AND ZINC PLATTED SPOKES ) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1170, 1280)	◆
CEG1072		REINFORCED SELF PROPEL 6 POINTS FIXING HANDRIM (COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1280)	€24
CEG1073		DRUM BRAKE SELF-PROPEL 22 AND 24" 4 PARTS FIXING HANDRIM COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES. TRANSIT 12" - BLACK PLASTIC HUB AND RIM DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1201, 1270, 1275)	○
CEG1074		COMPOSITE (BLACK PLASTIC HUB AND RIM) 7 POINTS FIXING HANDRIM TRANSIT BLACK PLASTIC HUB & RIM DISTANCE WHEEL-HANDRIM 14MM (NC 1088, 1260, 1261, 1280)	○
TYRES*			
CEG1082		SOLID BLACK 1.3/8" / 35MM	◆
CEG1084		PNEUMATIC PROFILE BLACK 1.3/8" / 35MM	€24
CEG1088		PNEUMATIC LOW PROFILE BLACK 1" / 25MM	○
HANDRIM			
CEG1140		ALUMINIUM PRE-ANODIZED HANDRIM	◆
CEG1150		HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC	€40
REAR WHEEL AXLE*			
CEG1170		FIXED AXLE	○
CEG1180		Q/R AXLE	◆
REAR WHEEL POSITION*			
CEG1200		PASSIVE REAR WHEELS POSITION	◆
CEG1201		DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)	€84
PUMP			
CEG1230		PUMP AND SUPPORT	€16
CHASSIS OPTION			
CEG1241		FRONT WHEELS EXTENSION	€30
CEG1260		SPOKE GUARDS PAIR 24" BLK (NC 1261)	€93
CEG1261		SPOKE GUARDS PAIR CLEAR	€118
BRAKES*			
CEG1270		PUSH TO LOCK BRAKES	◆
CEG1275		PULL TO LOCK BRAKES	○
CEG1280		ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)	€243
COMFORT TRANSPORTATION KIT*			
CEG1346		COMFORT TRANSPORTATION KIT (NC 1363, 1364)	◆

SAFETY OPTIONS			
CEG1361		BODYPOINT 2-PT POSTURE BELT M (NC 1362, 1363, 1364)	€171
CEG1362		BODYPOINT 2-PT POSTURE BELT L (NC 1363, 1364)	€171
CEG1363		BODYPOINT 4-PT POSTURE BELT MEDIUM (NC 1364)	€161
CEG1364		BODYPOINT 2-PT POSTURE BELT LARGE	€161
CEG1370		REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441)	○
CEG1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	€85
CEG1390		BRAKE EXTENSION PAIR	€16
OTHER OPTIONS			
CEG1430		CANE HOLDER	€28
CEG1440		TIPPING AID RH	€16
CEG1441		TIPPING AID LH	€16
CEG1481		COMFORT ACRYLIC TRAY	€125
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	€105

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

CEG0290	17" / 435 MM SEAT HEIGHT
CEG0300	16" / 410 MM SEAT HEIGHT
CEG0545	AMPUTEE LEGREST RH
CEG0595	AMPUTEE LEGREST LH
CEG0610	FLIP UP DEPTH AND ANGLE ADJ RH
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE66666	NO FOOTREST
CEG0630	FLIP UP DEPTH AND ANGLE ADJ LH
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE55555	NO FOOTREST
CEG0869	MATRX FLO-TECH LITE CUSHION - Low Risk
CEG0871	MATRX FLO-TECH CONTOUR VISCO CUSHION NG - HIGH RISK
CEG0873	MATRX FLO-TECH IMAGE CUSHION - High Risk
CEG0874	MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk
CEG0876	MATRX LIBRA CUSHION - Very High Risk
CEG0878	FOR 400 MM SEAT DEPTH
CEG0879	FOR 450 MM SEAT DEPTH
CEG0852	SEAT RIGIDIZER
CEG1010	8" X 1.1/4" BLACK SOLID TYRE
CEG1062	12"
CEG1068	24"
CEG1071	STANDARD (BASIC) 4 POINTS FIXING HANDRIM (BLACK HUB, COLOURLESS RIM AND ZINC PLATTED SPOKES ) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)
CEG1072	REINFORCED SELF PROPEL 6 POINTS FIXING HANDRIM (COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)
CEG1074	COMPOSITE (BLACK PLASTIC HUB AND RIM) 7 POINTS FIXING HANDRIM TRANSIT BLACK PLASTIC HUB & RIM DISTANCE WHEEL-HANDRIM 14MM
CEG1084	PNEUMATIC PROFILE BLACK 1.3/8" / 35MM
CEG1088	PNEUMATIC LOW PROFILE BLACK 1" / 25MM
CEG1140	ALUMINIUM PRE-ANODIZED HANDRIM
CEG1150	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC
CEG1170	FIXED AXLE
CEG1180	Q/R AXLE
CEG1201	DUAL AMPUTEE REAR WHEEL POSITION
CEG1260	SPOKE GUARDS PAIR 24" BLK
CEG1261	SPOKE GUARDS PAIR CLEAR
CEG1270	PUSH TO LOCK BRAKES
CEG1275	PULL TO LOCK BRAKES
CEG1280	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)
CEG1362	BODYPOINT 2-PT POSTURE BELT L
CEG1363	BODYPOINT 4-PT POSTURE BELT MEDIUM
CEG1364	BODYPOINT 2-PT POSTURE BELT LARGE
CEG1370	REAR ANTI-TIPPERS PAIR
CEG1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID

---

CEG1390	BRAKE EXTENSION PAIR
CEG1440	TIPPING AID RH
CEG1441	TIPPING AID LH
CEG1481	COMFORT ACRYLIC TRAY

---



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):