



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

APRIL 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG DHP

Standard Features

- Max user weight 125kg
- Crash tested ISO 7176-Part 19
- Dual Hand Rim Self Propulsion (DHP)
- Choice of 6 seat widths (15" - 20") 380mm to 505mm
- Choice of 2 seat depths (16" or 18") 400mm or 450mm
- Choice of 3 seat heights (17" - 19") 435mm to 485mm
- Black Nylon Upholstery
- Fixed 7° (17") 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in/Swing out legrest hangers
- Flip up footpates
- Swing away and detachable desk armrests
- Quick release axle
- Choice of 8 standard colours (with 7 special colour options)
- Push or Pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: £985



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description		
BASE PRODUCT*			
CEH0006	ACTION 3NG DHR	(NC 1084)	<input type="radio"/>
CEHSPEC	SPECIAL		POA
PROPULSION*			
CEH0140	DUAL HANDRIM RH PROPULSION	(NC 1330)	£985
CEH0150	DUAL HANDRIM LH PROPULSION	(NC 1320)	£985

SEAT WIDTH*			
CEH0180		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm) (NC 0852)	○
CEH0190		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEH0200		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEH0210		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEH0220		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEH0230		20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CEH0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEH0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
CEH0252		16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH (NC 0330, 0350)	£21
SEAT HEIGHT*			
CEH0270		19" / 485 MM SEAT HEIGHT	○
CEH0280		18" / 460 MM SEAT HEIGHT	◆
CEH0290		17" / 435 MM SEAT HEIGHT (NC 1010, 1030)	○
SEAT UPHOLSTERY*			
CEH0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEH0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEH0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM	○
CEH0350		PADDED NYLON	£30
BACKREST*			
CEH0380		FIXED 0° BACKREST 2 HEIGHTS 16" / 400 MM or 17" / 430 MM (NC 0471, 0477, 0945)	○
CEH0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT (NC 0471, 0477, 0945)	○
CEH0395		FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£67
CEH0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES (NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	£80
CEH0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES (NC 0471, 0477, 0945)	◆
CEH0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES (NC 0471, 0477, 0945)	○
CEH0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE (NC 0471, 0477, 0945)	£43
CEH0432		ANGLE ADJUSTABLE AND FOLD DOWN (-15° to + 15°) 17" (NC 0465, 0471, 0474, 0478, 0480, 0481, 0882)	£86
CEH0455		FOLDING 7° BACKREST 17" (NC 0471, 0477, 0945)	£62
CEH0460		FOLDING 7° BACKREST 20" (NC 0471, 0477, 0945)	£62
BACKREST UPHOLSTERY*			
CEH0465		WITHOUT BACK UPHOLSTERY (NC 0882)	○
CEH0471		STANDARD NYLON SLING (JEMIMA) (NC 0882, 0905, 0906, 0907)	○
CEH0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	○
CEH0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER	£55
CEH0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)	£62
CEH0480		CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE (NC 0940, 0945)	£120
CEH0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE (NC 0940, 0945)	£120

LEGREST HANGER RH*				
CEH0510		80° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CEH0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 66666)	£116
CEH0545		AMPUTEES LEGREST RH	(NC 0600, 0610, 0611, 0930, 0936, 1111)	£133
LEGREST HANGER LH*				
CEH0560		80° LH SWING IN/SWING OUT	(NC 2222)	◆
CEH0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 55555)	£116
CEH0595		AMPUTEES LEGREST LH	(NC 0611, 0620, 0630, 0930, 0936, 2222)	£133
FOOTREST RH*				
CEH0600		FLIP UP RH FOOTPLATE	(NC 0937, 0938)	◆
CEH0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		£43
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		£43
CEH0611		FLIP UP ONE PIECE	(NC 0620, 0630, 0935, 0936, 0937, 0938)	£78
CEE66666		NO FOOTREST		○
FOOTREST LH*				
CEF0620		FLIP UP LH FOOTPLATE	(NC 0937, 0938)	◆
CEH0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		£43
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		£43
CEE55555		NO FOOTREST		○
ARMREST RH*				
CEH0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796, 1480, 1483)	○
CEH0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796)	○
CEH0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 1480, 1483)	◆
CEH0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796)	○
CEH0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)	(NC 0740, 0750, 0760, 0770, 0795, 0796, 1480, 1483)	£56
CEH0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH	(NC 0740, 0750, 0760, 0770, 0795, 0796, 1483)	£56
CEH0695		REMOVABLE RH (N°4) A DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 1480, 1483)	£60
CEH0696		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 1483)	£60
CEH0732		REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED	(NC 1480, 1483)	£59

ARMREST LH*			
CEH0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)	○
CEH0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEH0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT (NC 1480, 1484)	◆
CEH0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT	○
CEH0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM) (NC 1480, 1484)	£56
CEH0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH (NC 1484)	£56
CEH0795		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.) AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)	£60
CEH0796		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED (NC 1484)	£60
CEH0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED (NC 1480, 1484)	£59
SEAT CUSHION			
CEH0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER	£24
CEH0869		MATRX FLO-TECH LITE CUSHION - Low Risk	£121
CEH0871		MATRX FLO-TECH CONTOUR VISCO CUSHION - High Risk	£169
CEH0873		MATRX FLO-TECH IMAGE CUSHION - High Risk	£288
CEH0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk	£398
CEH0876		MATRX LIBRA CUSHION - Very High Risk	£410
SEAT CUSHION DEPTH			
CEH0878		FOR 400 MM SEAT DEPTH	○
CEH0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEH0852		SEAT RIGIDIZER	£39
CEH0882		LATERAL POSITIONING PADS (PAIR)	£61
CEH0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	£131
CEH0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	£159
CEH0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	£183
CEH0930		CALF STRAP	◆
CEH0935		LONGER FOOTREST TUBE (NC 0936)	£30
CEH0936		SHORTER FOOTREST SUPPORT	£42
CEH0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	£114
CEH0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	£114
CEH0940		BACKBRACE (NC 0945)	£34
CEH0945		FOLDABLE BACKBRACE	£38

COLOUR*			
CEH0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEH0964	<input type="checkbox"/>	PEARL WHITE	○
CEH0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEH0966	<input type="checkbox"/>	MATT BLACK	○
CEH0969	<input type="checkbox"/>	RED MATT	○
CEH0972	<input type="checkbox"/>	WATER BLUE	○
CEG0974	<input type="checkbox"/>	BLUE	○
CEG0976	<input type="checkbox"/>	DARK PURPLE	○
CEH0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
CEH0991	<input type="checkbox"/>	SPECIAL COLOUR ON DEMAND ONLY	£169
SWO3333	<input type="checkbox"/>	DARK RED - SPECIAL	£169
SWO4444	<input type="checkbox"/>	PEARL GREY - SPECIAL	£169
SWO5555	<input type="checkbox"/>	INTENSE YELLOW	£169
SWO6666	<input type="checkbox"/>	VIBRANT ORANGE	£169
SWO7777	<input type="checkbox"/>	POLAR BLUE	£169
SWO8888	<input type="checkbox"/>	CANDY PINK	£169
CASTOR HOUSING*			
CEH0995	<input type="checkbox"/>	STANDARD	◆
CEH0996	<input type="checkbox"/>	LIGHT WEIGHT	£28
FRONT CASTORS*			
CEH1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆
CEH1030	<input type="checkbox"/>	8" X 2" BLACK SOLID TYRE	£28
CEH1035	<input type="checkbox"/>	6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD	£28
CEH1050	<input type="checkbox"/>	6" X 1.1/4" BLACK SOLID TYRE	£28
REAR WHEEL SIZE*			
CEH1068	<input type="checkbox"/>	24"	◆
REAR WHEELS*			
CEH1075	<input type="checkbox"/>	ONE ARM DRIVE (NC 1084)	◆
TYRES*			
CEH1082	<input type="checkbox"/>	SOLID BLACK 1.3/8" / 35MM	£21
CEH1084	<input type="checkbox"/>	PNEUMATIC PROFILE BLACK 1.3/8" / 35MM	◆
REAR WHEEL AXLE*			
CEH1180	<input type="checkbox"/>	Q/R AXLE	◆
REAR WHEEL POSITION*			
CEH1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	○
CEH1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	◆
PUMP			
CEH1230	<input type="checkbox"/>	PUMP AND SUPPORT	£14
CHASSIS OPTION			
CEH1241	<input type="checkbox"/>	FRONT WHEELS EXTENSION	£22
BRAKES*			
CEH1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
CEH1275	<input type="checkbox"/>	PULL TO LOCK BRAKES	○
CEH1320	<input type="checkbox"/>	OAD WIRE CONTROLLED BRAKE RH (NC 1390)	£45
CEH1330	<input type="checkbox"/>	OAD WIRE CONTROLLED BRAKE LH (NC 1390)	£45

TRANSPORTATION KIT*			
CEH1345		TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CEH1360		LAP BELT	(NC 1361, 1362, 1363, 1364) ◆
CEH1361		BODYPOINT 2-PT POSTURE BELT M	(NC 1362, 1363, 1364) £83
CEH1362		BODYPOINT 2-PT POSTURE BELT L	(NC 1363, 1364) £83
CEH1363		BODYPOINT 4-PT POSTURE BELT M	(NC 1364) £98
CEH1364		BODYPOINT 4-PT POSTURE BELT L	£98
CEH1370		REAR ANTI-TIPPERS PAIR	(NC 1375, 1440, 1441) £36
CEH1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID	(NC 1440, 1441) £72
CEH1380		REFLECTORS KIT	£8
CEH1390		BRAKE EXTENSION PAIR	£14
OTHER OPTIONS			
CEH1430		CANE HOLDER	£23
CEH1440		TIPPING AID RH	£14
CEH1441		TIPPING AID LH	£14
CEH1480		ACRYLIC TRAY	(NC 1483, 1484, 1485) £116
CEH1483		HALF TABLE TRAY RH	(NC 1484) £89
CEH1484		HALF TABLE TRAY LH	£89
CEH1485		COMFORT PADDING FOR HALF TABLE TRAY	£44
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£107
1595002		ViaGo Push aid. Includes as standard: Drive Unit, Charger, Control unit & Brackets (Note: 10 minute fitting time)	£1,995

TOTAL PRICE	
--------------------	--

Non Conforming Options:

CEH0330	NYLON W/O HOOK AND LOOP STRIPS
CEH0350	PADDED NYLON
CEH0465	WITHOUT BACK UPHOLSTERY
CEH0471	STANDARD NYLON SLING (JEMIMA)
CEH0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEH0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEH0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEH0480	CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE
CEH0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEH0600	FLIP UP RH FOOTPLATE
CEH0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEH0611	FLIP UP ONE PIECE
CEE66666	NO FOOTREST
CEF0620	FLIP UP LH FOOTPLATE
CEH0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE55555	NO FOOTREST
CEH0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEH0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEH0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEH0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEH0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CEH0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CEH0795	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.) AUTOLOCK SAFETY DEVICE INCLUDED
CEH0796	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEH0878	FOR 400 MM SEAT DEPTH
CEH0879	FOR 450 MM SEAT DEPTH
CEH0852	SEAT RIGIDIZER
CEH0882	LATERAL POSITIONING PADS (PAIR)
CEH0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CEH0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEH0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEH0930	CALF STRAP
CEH0935	LONGER FOOTREST TUBE
CEH0936	SHORTER FOOTREST SUPPORT
CEH0937	BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED
CEH0938	BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED
CEH0940	BACKBRACE
CEH0945	FOLDABLE BACKBRACE

CEH1010	8" X 1.1/4" BLACK SOLID TYRE
CEH1030	8" X 2" BLACK SOLID TYRE
CEH1084	PNEUMATIC PROFILE BLACK 1.3/8" / 35MM
CEH1320	OAD WIRE CONTROLLED BRAKE RH
CEH1330	OAD WIRE CONTROLLED BRAKE LH
CEH1361	BODYPOINT 2-PT POSTURE BELT M
CEH1362	BODYPOINT 2-PT POSTURE BELT L
CEH1363	BODYPOINT 4-PT POSTURE BELT M
CEH1364	BODYPOINT 4-PT POSTURE BELT L
CEH1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEH1390	BRAKE EXTENSION PAIR
CEH1440	TIPPING AID RH
CEH1441	TIPPING AID LH
CEH1480	ACRYLIC TRAY
CEH1483	HALF TABLE TRAY RH
CEH1484	HALF TABLE TRAY LH
CEH1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):