



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Light

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion with 24" Quick release rear wheels
- 8" X 1.1/4" castors, black non marking
- Choice of 6 seat widths 380mm to 505mm (15" - 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 6 seat heights 385mm to 510mm (15" - 20")
- Black nylon upholstery
- Folding 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / Swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: €1,038



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code **Description**

BASE PRODUCT*

CES0017	☐	ACTION 3NG LIGHT	☐
CESSPEC	☐	SPECIAL	POA

PROPULSION*				
CES0100		SELF PROPULSION		€1,038
SEAT WIDTH*				
CES44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	(NC 0660, 0670, 0680, 0690, 0695, 0696, 0732, 0738, 0760, 0770, 0780, 0790, 0795, 0796, 0832, 0838, 0849, 0869, 0871, 0873, 0874, 0876, 0878, 0879, 0891, 0892, 0893, 0894)	○
CES0180		15" / 380 MM SEAT WIDTH	(NC 0852)	○
CES0190		16" / 405 MM SEAT WIDTH		○
CES0200		17" / 430 MM SEAT WIDTH		◆
CES0210		18" / 455 MM SEAT WIDTH		○
CES0220		19" / 480 MM SEAT WIDTH		○
CES0230		20" / 505 MM SEAT WIDTH		○
SEAT DEPTH*				
CES0240		16" / 400 MM SEAT DEPTH	(NC 0879)	○
CES0250		18" / 450 MM SEAT DEPTH	(NC 0878)	◆
CES0252		16" TO 18" / 406 TO 457 MM ADJUSTABLE SEAT DEPTH	(NC 0330, 0350)	€24
SEAT HEIGHT*				
CES0260		20" / 510 MM SEAT HEIGHT	(NC 0511, 0561, 0738, 1035, 1050, 1054, 1056, 1066)	○
CES0270		19" / 485 MM SEAT HEIGHT	(NC 0511, 0561, 1054)	○
CES0280		18" / 460 MM SEAT HEIGHT		◆
CES0290		17" / 435 MM SEAT HEIGHT	(NC 1010, 1030)	○
CES0300		16" / 410 MM SEAT HEIGHT	(NC 0935, 1010, 1030, 1056, 1370, 1440, 1441)	○
CES0310		15" / 385 MM SEAT HEIGHT ,	(NC 0511, 0561, 0660, 0670, 0738, 0760, 0770, 0838, 0935, 1010, 1030, 1035, 1050, 1056, 1068, 1370, 1440, 1441)	○
SEAT UPHOLSTERY*				
CES0330		NYLON W/O HOOK AND LOOP STRIPS		◆
CES0340		NYLON WITH HOOK AND LOOP STRIPS		○
CES0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM		○
CES0350		PADDED NYLON		€36
BACKREST*				
CES0380		FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM	(NC 0471, 0477, 0945)	○
CES0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT	(NC 0471, 0477, 0945)	○
CES0395		FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	€80
CES0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	€105
CES0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	(NC 0471, 0477, 0945)	◆
CES0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES	(NC 0471, 0477, 0945)	○
CES0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE	(NC 0471, 0477, 0945)	€50
CES0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HT ADJ TO 16" / 400 MM FIX PUSH HANDLES	(NC 0471, 0477, 0480, 0481, 0905, 0906, 0907, 0940)	○
CES0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17"	(NC 0465, 0471, 0474, 0478, 0480, 0481, 0882, 0952)	€102
CES0455		FOLDING 7° BACKREST 17"	(NC 0471, 0477, 0945)	€74
CES0460		FOLDING 7° BACKREST 20"	(NC 0471, 0477, 0945)	€74

BACKREST UPHOLSTERY*				
CES0465		WITHOUT BACK UPHOLSTERY	(NC 0882)	○
CES0471		STANDARD NYLON SLING (JEMIMA)	(NC 0882, 0905, 0906, 0907)	○
CES0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER		○
CES0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER		€63
CES0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	€71
CES0480		CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE	(NC 0940, 0945)	€142
CES0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	€142
LEGREST HANGER RH*				
CES0510		80° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CES0511		90° RH SWING IN/SWING OUT	(NC 0600, 0936, 1010, 1030, 1035, 1050, 1056, 1111, 66666)	○
CES0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 66666)	€150
CES0545		AMPUTEE LEGREST RH	(NC 0600, 0610, 0611, 0930, 0936, 1111)	€159
LEGREST HANGER LH*				
CES0560		80° LH SWING IN/SWING OUT	(NC 2222)	◆
CES0561		90° LH SWING IN/SWING OUT	(NC 0620, 0936, 1010, 1030, 1035, 1050, 1056, 2222)	○
CES0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 55555)	€150
CES0595		AMPUTEE LEGREST LH	(NC 0611, 0620, 0630, 0930, 0936, 2222)	€159
FOOTREST RH*				
CES0600		FLIP UP RH FOOTPLATE	(NC 0937, 0938)	◆
CES0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		€50
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		€50
CES0611		FLIP UP ONE PIECE	(NC 0620, 0630, 0935, 0936, 0937, 0938)	€93
CES66666		NO FOOTREST		○
FOOTREST LH*				
CES0620		FLIP UP LH FOOTPLATE	(NC 0937, 0938)	◆
CES0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		€50
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		€50
CES55555		NO FOOTREST		○

ARMREST RH*				
CES0640		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 1480, 1483)	○
CES0650		REMOVABLE RH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838)	○
CES0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0952, 0955, 1480, 1483)	◆
CES0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0952, 0955)	○
CES0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0952, 0955, 1480, 1483)	€67
CES0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0952, 0955, 1483)	€67
CES0695		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0952, 0955, 1480, 1483)	€71
CES0696		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0952, 0955, 1483)	€71
CES0732		REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION	(NC 0838, 0952, 1480, 1483)	€69
CES0738		REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0795, 0796, 0832, 0952, 0955, 1066, 1200, 1480, 1483)	€41
ARMREST LH*				
CES0740		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT	(NC 1480, 1484)	○
CES0750		REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT		○
CES0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	(NC 0952, 0955, 1480, 1484)	◆
CES0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0952, 0955)	○
CES0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)	(NC 0952, 0955, 1480, 1484)	€67
CES0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH	(NC 0952, 0955, 1484)	€67
CES0795		REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.	(NC 0952, 0955, 1480, 1484)	€71
CES0796		REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.	(NC 0952, 0955, 1484)	€71
CES0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION	(NC 0952, 1480, 1484)	€69
CES0838		REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0952, 0955, 1066, 1200, 1484)	€41
SEAT CUSHION				
CES0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER	(NC 0952)	€27
CES0869		MATRX FLO-TECH LITE CUSHION - Low Risk	(NC 0952)	€168
CES0871		MATRX FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK	(NC 0952)	€236
CES0873		MATRX FLO-TECH IMAGE CUSHION - High Risk	(NC 0952)	€382
CES0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk	(NC 0952)	€566
CES0876		MATRX LIBRA CUSHION - Very High Risk	(NC 0952)	€645

SEAT CUSHION DEPTH			
CES0878		FOR 400 MM SEAT DEPTH	○
CES0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CES0852		SEAT RIGIDIZER (NC 0952)	€79
CES0882		LATERAL POSITIONING PADS (PAIR)	€72
CES0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	€158
CES0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	€190
CES0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	€217
SWO3333		MULT ADJUSTABLE LATERALS	€340
CES0930		CALF STRAP	◆
CES0935		LONGER FOOTREST TUBE	€36
CES0936		SHORTER FOOTREST SUPPORT	€48
CES0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	€189
CES0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	€189
CES0940		BACKBRACE (NC 0945)	€41
CES0945		FOLDABLE BACKBRACE	€44
CES0952		IN-RIGGED ARMREST (NC 1480, 1483, 1484)	€48
CES0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○
POSITIONING ACCESSORIES FOR MATRX LIBRA & FLOTECH			
CES0891		TAPERED WEDGES (PAIR)	€28
CES0892		ABDUCTOR	€58
CES0893		LATERAL WEDGE	€46
CES0894		ANTERIOR SEAT WEDGE	€43
COLOUR*			
CES0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CES0964		PEARL WHITE	○
CES0965		ELECTRIC GREEN	○
CES0966		MATT BLACK	○
CES0969		RED MATT	○
CES09722		WATER BLUE	○
CEA0974		BLUE	○
CEA0976		DARK PURPLE	○
CES0980		ANTHRACITE BLACK	○
CES0991		SPECIAL COLOUR ON DEMAND ONLY	€238
SWO4444		DARK RED - SPECIAL	€238
SWO5555		PEARL GREY - SPECIAL	€238
SWO6666		INTENSE YELLOW	€238
SWO7777		VIBRANT ORANGE	€238
SWO8888		POLAR BLUE	€238
SWO9999		CANDY PINK	€238
CASTOR HOUSING*			
CES0996		LIGHT WEIGHT	◆

FRONT CASTORS*				
CES1010		8" X 1.1/4" BLACK SOLID TYRE		◆
CES1030		8" X 2" BLACK SOLID TYRE		€32
CES1035		6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD		€32
CES1050		6" X 1.1/4" BLACK SOLID TYRE		€32
CES1054		5" X 1" BLACK SOLID TYRE		€32
CES1056		5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD		€32
REAR WHEEL SIZE*				
CES1066		22"	(NC 1092)	○
CES1068		24"		◆
REAR WHEELS*				
CES1076		LIGHTWEIGHT 6 PARTS FIXING HANDRIM (BLACK HUB AND RIM, 24 BLACK RADIAL SPOKES) DISTANCE WHEEL -HANDRIM 11 OR 16MM (SET AT 11MM)	(NC 1230)	◆
TYRES*				
CES1082		SOLID BLACK		◆
CES1084		PNEUMATIC PROFILE BLACK		○
CES1088		PNEUMATIC LOW PROFILE BLACK		○
CES1092		PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM		€41
HANDRIM*				
CES1140		ALUMINIUM PRE-ANODIZED HANDRIM		◆
CES1150		HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC		€40
REAR WHEEL AXLE*				
CES1180		Q/R AXLE		◆
REAR WHEEL POSITION*				
CES1190		STANDARD REAR WHEELS POSITION		◆
CES1200		PASSIVE REAR WHEELS POSITION		○
PUMP				
CES1230		PUMP AND SUPPORT		€16
CHASSIS OPTION				
CES1241		FRONT WHEELS EXTENSION		€30
BRAKES*				
CES1270		PUSH TO LOCK BRAKES		◆
CES1275		PULL TO LOCK BRAKES		○
TRANSPORTATION KIT*				
CES1345		TRANSPORTATION BRACKET		◆
SAFETY OPTIONS*				
CES1360		LAP BELT	(NC 1361, 1362, 1363, 1364)	◆
CES1361		BODYPOINT 2-PT POSTURE BELT M	(NC 1362, 1363, 1364)	€171
CES1362		BODYPOINT 2-PT POSTURE BELT L	(NC 1363, 1364)	€171
CES1363		BODYPOINT 4-PT POSTURE BELT M	(NC 1364)	€161
CES1364		BODYPOINT 4-PT POSTURE BELT L		€161
CES1370		REAR ANTI-TIPPERS PAIR	(NC 1375, 1440, 1441)	€43
CES1375		REAR ANTI-TIPPERS PAIR WITH TIPPING AID	(NC 1440, 1441)	€85
CES1380		REFLECTORS KIT		€10
CES1390		BRAKE EXTENSION PAIR		€16

OTHER OPTIONS			
CES1430		CANE HOLDER	€28
CES1440		TIPPING AID RH	€16
CES1441		TIPPING AID LH	€16
CES1480		ACRYLIC TRAY (NC 1483, 1484, 1485)	€150
CES1483		HALF TABLE TRAY RH (NC 1484)	€62
CES1484		HALF TABLE TRAY LH	€62
CES1485		COMFORT PADDING FOR HALF TABLE TRAY	€30
SWO10002		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	€105

			TOTAL PRICE	
--	--	--	-------------	--

Non Conforming Options:	
CES0330	NYLON W/O HOOK AND LOOP STRIPS
CES0350	PADDED NYLON
CES0465	WITHOUT BACK UPHOLSTERY
CES0471	STANDARD NYLON SLING (JEMIMA)
CES0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CES0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CES0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CES0480	CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE
CES0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CES0511	90° RH SWING IN/SWING OUT
CES0561	90° LH SWING IN/SWING OUT
CES0600	FLIP UP RH FOOTPLATE
CES0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CES0611	FLIP UP ONE PIECE
CES66666	NO FOOTREST
CES0620	FLIP UP LH FOOTPLATE
CES0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CES55555	NO FOOTREST
CES0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CES0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CES0680	SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CES0690	SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CES0695	REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.
CES0696	REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.
CES0732	REMOVABLE RH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION
CES0738	REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)
CES0740	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. 200 TO 280 MM HEIGHT
CES0750	REMOVABLE LH (N°1) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. 200 TO 280 MM HEIGHT
CES0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CES0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CES0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170-295 MM)
CES0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170-295 MM) AND DEPTH
CES0795	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ.
CES0796	REMOVABLE LH (N°4) AUTOLOCK SAFETY DEVICE INCLUDED FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ.
CES0832	REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION

CES0838	REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)
CES0849	BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER
CES0869	MATRX FLO-TECH LITE CUSHION - Low Risk
CES0871	MATRX FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK
CES0873	MATRX FLO-TECH IMAGE CUSHION - High Risk
CES0874	MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk
CES0876	MATRX LIBRA CUSHION - Very High Risk
CES0878	FOR 400 MM SEAT DEPTH
CES0879	FOR 450 MM SEAT DEPTH
CES0852	SEAT RIGIDIZER
CES0882	LATERAL POSITIONING PADS (PAIR)
CES0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CES0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CES0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CES0930	CALF STRAP
CES0935	LONGER FOOTREST TUBE
CES0936	SHORTER FOOTREST SUPPORT
CES0937	BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED
CES0938	BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED
CES0940	BACKBRACE
CES0945	FOLDABLE BACKBRACE
CES0952	IN-RIGGED ARMREST
CES0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CES0891	TAPERED WEDGES (PAIR)
CES0892	ABDUCTOR
CES0893	LATERAL WEDGE
CES0894	ANTERIOR SEAT WEDGE
CES1010	8" X 1.1/4" BLACK SOLID TYRE
CES1030	8" X 2" BLACK SOLID TYRE
CES1035	6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD
CES1050	6" X 1.1/4" BLACK SOLID TYRE
CES1054	5" X 1" BLACK SOLID TYRE
CES1056	5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD
CES1066	22"
CES1068	24"
CES1092	PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM
CES1200	PASSIVE REAR WHEELS POSITION
CES1230	PUMP AND SUPPORT
CES1361	BODYPOINT 2-PT POSTURE BELT M
CES1362	BODYPOINT 2-PT POSTURE BELT L
CES1363	BODYPOINT 4-PT POSTURE BELT M
CES1364	BODYPOINT 4-PT POSTURE BELT L
CES1370	REAR ANTI-TIPPERS PAIR
CES1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID

CES1440	TIPPING AID RH
CES1441	TIPPING AID LH
CES1480	ACRYLIC TRAY
CES1483	HALF TABLE TRAY RH
CES1484	HALF TABLE TRAY LH
CES1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):