



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Self Propel Compact

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion with 20" Quick release rear wheels
- Choice of 3 seat widths 380mm to 430mm (15" - 17")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 5 seat heights 360mm to 460mm (14" - 18")
- Black nylon upholstery
- Folding 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 90° Swing in / Swing out legrest hangers
- Flip up, Angle and Depth Adjustable Footplates
- Swing away and detachable desk armrests
- Choice of 8 standard colours (with 7 special colour options)
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: €966



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEN0012	ACTION 3NG COMPACT	☐
CENSPEC	SPECIAL	POA
PROPULSION*		
CEN0100	SELF PROPULSION	€966

SEAT WIDTH*				
CEN44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 (lapbelt adjustment 420mm - 1520mm)	(NC 0660, 0670, 0695, 0696, 0760, 0770, 0795, 0796, 0849, 0869, 0871, 0873, 0874, 0876, 0878, 0879, 0891, 0892, 0893, 0894)	○
CEN0180		15" / 380 MM SEAT WIDTH	(NC 0852)	○
CEN0190		16" / 405 MM SEAT WIDTH		◆
CEN0200		17" / 430 MM SEAT WIDTH		○
SEAT DEPTH*				
CEN0240		16" / 400 MM SEAT DEPTH	(NC 0879)	◆
CEN0250		18" / 450 MM SEAT DEPTH	(NC 0878)	○
SEAT HEIGHT*				
CEN0280		18" / 460 MM SEAT HEIGHT		○
CEN0290		17" / 435 MM SEAT HEIGHT		◆
CEN0300		16" / 410 MM SEAT HEIGHT	(NC 1370, 1440, 1441)	○
CEN0310		15" / 385 MM SEAT HEIGHT ,	(NC 0511, 0561, 0660, 0670, 0760, 0770, 1370, 1440, 1441)	○
CEN0315		14" / 360 MM SEAT HEIGHT, 0° SEAT ANGLE (AVAILABLE WITH LOW FORK, 5" FRONT WHEEL AND 20" REAR WHEELS)	(NC 0511, 0561, 0611, 0660, 0670, 0760, 0770, 1066, 1375, 1440, 1441)	○
SEAT UPHOLSTERY*				
CEN0330		NYLON W/O HOOK AND LOOP STRIPS		◆
CEN0340		NYLON WITH HOOK AND LOOP STRIPS		○
CEN0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM		○

BACKREST*				
CEN0380		FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM (DELIVERED 17" / 430 MM) FIXED PUSH HANDLES	(NC 0471, 0477, 0945)	○
CEN0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT FIXED PUSH HANDLES	(NC 0471, 0477, 0945)	○
CEN0395		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 MM TO 18" / 460 MM (DELIVERED 15" / 380 MM) FIXED PUSH HANDLES (follow upholstery height)	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	€80
CEN0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM (DELIVERED 15" / 380 MM) HEIGHT ADJUSTABLE PUSH HANDLES (independantly from upholstery height)	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	€105
CEN0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM (DELIVERED 17" / 430 MM) FIXED PUSH HANDLES	(NC 0471, 0477, 0945)	◆
CEN0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIXED PUSH HANDLES	(NC 0471, 0477, 0945)	○
CEN0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HEIGHT ADJUSTABLE PUSHHANDLE (6" / 150 MM)	(NC 0471, 0477, 0945)	€50
CEN0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HEIGHT (ADJUSTABLE TO 16" / 406 MM) FIXED PUSH HANDLES (SHORT)	(NC 0471, 0477, 0480, 0481, 0905, 0906, 0907, 0940)	○
CEN0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) FIXED 7°, HEIGHT 17" / 425MM WITH FIXED LONG PUSHHANDLES HEIGHT 17" / 425 MM 18" / 450 MM SET AT 425MM / 17" (2 positions, 1" / 25mm steps)	(NC 0465, 0471, 0474, 0478, 0480, 0481, 0882, 0952)	€102
CEN0445		RECLINING (GAS SPRING) 0/30° BACKREST INCL. BACKRACE AND PASSIVE REAR WHEEL POS. 3 HEIGHTS 15" to 17" / 390/410/430 MM (DELIVERED 17"/ 430 MM) ANTITIPPERS REQUIRED	(NC 0471, 0477, 0940, 0945, 1190)	€226
CEN0450		RECLINING (GAS SPRING) 0/30° BACKREST INCL. BACKRACE AND PASSIVE REAR WHEEL POS. 3 HEIGHTS 18" to 20" / 470/490/510 MM (DELIVERED 20" / 510 MM) ANTITIPPERS REQUIRED	(NC 0471, 0477, 0940, 0945, 1190)	€226
CEN0455		FOLDING 7° BACKREST 17" / 430 MM HEIGHT FIXED PUSH HANDLES	(NC 0471, 0477, 0945)	€74
CEN0460		FOLDING 7° BACKREST 20" / 510 MM HEIGHT FIXED PUSH HANDLES	(NC 0471, 0477, 0945)	€74
BACKREST UPHOLSTERY*				
CEN0465		WITHOUT BACK UPHOLSTERY	(NC 0882)	○
CEN0471		STANDARD NYLON SLING (JEMIMA)	(NC 0882, 0905, 0906, 0907)	○
CEN0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER		○
CEN0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER		€63
CEN0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	€71
CEN0480		CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE	(NC 0940, 0945)	€142
CEN0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	€142
LEGREST HANGER RH*				
CEN0511		90° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CEN0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 66666)	€150
CEN0545		AMPUTEE LEGREST RH	(NC 0610, 0611, 0930, 1111)	€159
LEGREST HANGER LH*				
CEN0561		90° LH SWING IN/SWING OUT	(NC 2222)	◆
CEN0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 0611, 0930, 55555)	€150
CEN0595		AMPUTEE LEGREST LH	(NC 0611, 0630, 0930, 2222)	€159

FOOTREST RH*				
CEN0610		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		○
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		○
CEN0611		FLIP UP ONE PIECE	(NC 0630, 0938)	€93
CEN66666		NO FOOTREST		○
FOOTREST LH*				
CEN0630		FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.		○
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		○
CEN55555		NO FOOTREST		○
ARMREST RH*				
CEN0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0795, 0796, 1480, 1483)	○
CEN0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0795, 0796)	○
CEN0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0795, 0796, 0952, 0955, 1480, 1483)	◆
CEN0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0795, 0796, 0952, 0955)	○
CEN0695		REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0952, 0955, 1480, 1483)	€71
CEN0696		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0952, 0955, 1483)	€71
ARMREST LH*				
CEN0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 1480, 1484)	○
CEN0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED		○
CEN0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	(NC 0952, 0955, 1480, 1484)	◆
CEN0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0952, 0955)	○
CEN0795		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0952, 0955, 1480, 1484)	€71
CEN0796		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0952, 0955, 1484)	€71

SEAT CUSHION			
CEN0849		BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER (NC 0952)	€27
CEN0869		MATRX FLO-TECH LITE CUSHION - Low Risk (NC 0952)	€168
CEN0871		MATRX FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK (NC 0952)	€236
CEN0873		MATRX FLO-TECH IMAGE CUSHION - High Risk (NC 0952)	€382
CEN0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk (NC 0952)	€566
CEN0876		MATRX LIBRA CUSHION - Very High Risk (NC 0952)	€645
SEAT CUSHION DEPTH			
CEN0878		FOR 400 MM SEAT DEPTH	○
CEN0879		FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEN0852		SEAT RIGIDIZER (NC 0952)	€79
CEN0882		LATERAL POSITIONING PADS (PAIR)	€72
CEN0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	€158
CEN0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	€190
CEN0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	€217
CEN0930		CALF STRAP	○
CEN0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	€189
CEN0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	€189
CEN0940		BACKBRACE (NC 0945)	€41
CEN0945		FOLDABLE BACKBRACE	€44
CEN0952		IN-RIGGED ARMREST (NC 1480, 1483, 1484)	€48
CEN0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○
POSITIONING ACCESSORIES FOR MATRX LIBRA & FLOTECH			
CEN0891		TAPERED WEDGES (PAIR)	€28
CEN0892		ABDUCTOR	€58
CEN0893		LATERAL WEDGE	€46
CEN0894		ANTERIOR SEAT WEDGE	€43

COLOUR*			
CEN0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEN0964		PEARL WHITE	○
CEN0965		ELECTRIC GREEN	○
CEN0966		MATT BLACK	○
CEN0969		RED MATT	○
CEN0972		WATER BLUE	○
CEG0974		BLUE	○
CEG0976		DARK PURPLE	○
CEN0980		ANTHRACITE BLACK	○
CEN0991		SPECIAL COLOUR ON DEMAND ONLY	€238
SWO4444		DARK RED - SPECIAL	€238
SWO5555		PEARL GREY - SPECIAL	€238
SWO6666		INTENSE YELLOW	€238
SWO7777		VIBRANT ORANGE	€238
SWO8888		POLAR BLUE	€238
SWO9999		CANDY PINK	€238
CASTOR HOUSING*			
CEN0995		STANDARD	◆
CEN0996		LIGHT WEIGHT	€35
FRONT CASTORS*			
CEN1054		5" X 1" BLACK SOLID TYRE	◆
REAR WHEEL SIZE*			
CEN1064		20" (NC 1071, 1076, 1088)	◆
CEN1066		22"	○
REAR WHEELS*			
CEN1071		STANDARD 4 PARTS FIXING HANDRIM (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)	◆
CEN1072		REINFORCED 6 PARTS FIXING HANDRIM (COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) DISTANCE WHEEL- HANDRIM 11 OR 16MM (SET AT 11MM)	€24
CEN1076		LIGHTWEIGHT 6 PARTS FIXING HANDRIM (BLACK HUB AND RIM, 24 BLACK RADIAL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)	(NC 1230) €70
TYRES*			
CEN1082		SOLID BLACK	◆
CEN1084		PNEUMATIC PROFILE BLACK	€24
CEN1088		PNEUMATIC LOW PROFILE BLACK	○
HANDRIM*			
CEN1140		ALUMINIUM PRE-ANODIZED HANDRIM	◆
REAR WHEEL AXLE*			
CEN1180		Q/R AXLE	◆
REAR WHEEL POSITION*			
CEN1190		STANDARD REAR WHEELS POSITION	◆
CEN1200		PASSIVE REAR WHEELS POSITION	○

PUMP			
CEN1230	☐	PUMP AND SUPPORT	€16
CHASSIS OPTION			
CEN1241	☐	FRONT WHEELS EXTENSION	€30
BRAKES*			
CEN1270	☐	PUSH TO LOCK BRAKES	◆
CEN1275	☐	PULL TO LOCK BRAKES	○
TRANSPORTATION KIT*			
CEN1345	☐	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CEN1360	☐	LAP BELT (NC 1361, 1362, 1363, 1364)	◆
CEN1361	☐	BODYPOINT 2-PT POSTURE BELT M (NC 1362, 1363, 1364)	€171
CEN1362	☐	BODYPOINT 2-PT POSTURE BELT L (NC 1363, 1364)	€171
CEN1363	☐	BODYPOINT 4-PT POSTURE BELT M (NC 1364)	€161
CEN1364	☐	BODYPOINT 4-PT POSTURE BELT L	€161
CEN1370	☐	REAR ANTI-TIPPERS PAIR (NC 1375, 1440, 1441)	€43
CEN1375	☐	REAR ANTI-TIPPERS PAIR WITH TIPPING AID (NC 1440, 1441)	€85
CEN1380	☐	REFLECTORS KIT	€10
CEN1390	☐	BRAKE EXTENSION PAIR	€16
OTHER OPTIONS			
CEN1440	☐	TIPPING AID RH	€15
CEN1441	☐	TIPPING AID LH	€15
CEN1480	☐	ACRYLIC TRAY (NC 1483, 1484, 1485)	€150
CEN1483	☐	HALF TABLE TRAY RH (NC 1484)	€62
CEN1484	☐	HALF TABLE TRAY LH	€62
CEN1485	☐	COMFORT PADDING FOR HALF TABLE TRAY	€30
SWO10002	☐	ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	€105

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:	
CEN0465	WITHOUT BACK UPHOLSTERY
CEN0471	STANDARD NYLON SLING (JEMIMA)
CEN0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEN0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEN0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEN0480	CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE
CEN0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEN0511	90° RH SWING IN/SWING OUT
CEN0561	90° LH SWING IN/SWING OUT
CEN0610	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEN0611	FLIP UP ONE PIECE
CEN66666	NO FOOTREST
CEN0630	FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEN55555	NO FOOTREST
CEN0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CEN0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CEN0695	REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEN0696	REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEN0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEN0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEN0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEN0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEN0795	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210-345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEN0796	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210-345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEN0849	BASIC 50 MM FOAM CUSHION - NYLON AND ANTI-SLIP COVER
CEN0869	MATRX FLO-TECH LITE CUSHION - Low Risk
CEN0871	MATRX FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK
CEN0873	MATRX FLO-TECH IMAGE CUSHION - High Risk
CEN0874	MATRX FLO-TECH SOLUTION XTRA CUSHION - High Risk
CEN0876	MATRX LIBRA CUSHION - Very High Risk
CEN0878	FOR 400 MM SEAT DEPTH
CEN0879	FOR 450 MM SEAT DEPTH
CEN0852	SEAT RIGIDIZER
CEN0882	LATERAL POSITIONING PADS (PAIR)
CEN0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CEN0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC

CEN0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEN0930	CALF STRAP
CEN0938	BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED
CEN0940	BACKBRACE
CEN0945	FOLDABLE BACKBRACE
CEN0952	IN-RIGGED ARMREST
CEN0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEN0891	TAPERED WEDGES (PAIR)
CEN0892	ABDUCTOR
CEN0893	LATERAL WEDGE
CEN0894	ANTERIOR SEAT WEDGE
CEN1066	22"
CEN1071	STANDARD 4 PARTS FIXING HANDRIM (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)
CEN1076	LIGHTWEIGHT 6 PARTS FIXING HANDRIM (BLACK HUB AND RIM, 24 BLACK RADIAL SPOKES) DISTANCE WHEEL- HANDRIM 11 OR 16MM (SET AT 11MM)
CEN1088	PNEUMATIC LOW PROFILE BLACK
CEN1190	STANDARD REAR WHEELS POSITION
CEN1230	PUMP AND SUPPORT
CEN1361	BODYPOINT 2-PT POSTURE BELT M
CEN1362	BODYPOINT 2-PT POSTURE BELT L
CEN1363	BODYPOINT 4-PT POSTURE BELT M
CEN1364	BODYPOINT 4-PT POSTURE BELT L
CEN1370	REAR ANTI-TIPPERS PAIR
CEN1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID
CEN1440	TIPPING AID RH
CEN1441	TIPPING AID LH
CEN1480	ACRYLIC TRAY
CEN1483	HALF TABLE TRAY RH
CEN1484	HALF TABLE TRAY LH
CEN1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):