



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 3NG Self Propel

Standard Features

- Max user weight 125kg (19.7 stone)
- Crash tested ISO 7176-Part 19
- Self Propulsion with 22" or 24" Quick release rear wheels
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 6 seat heights 385mm to 510mm (15" to 20")
- Black nylon upholstery
- Folding 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / Swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and lap belt (adjustable both sides)
- Adjustable heel straps
- User tool kit

Standard Base price: €866



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEE0003	ACTION 3NG SELF PROPEL	☐
CEESPEC	SPECIAL	POA
PROPULSION*		
CEE0100	SELF PROPULSION	€866

SEAT WIDTH*				
CEE44444		13" / 330 MM SEAT WIDTH - SELECT SEAT WIDTH OPTION CEE0180 AND IN-RIGGED ARMREST OPTION CEE0952 NEEDS TO BE SELECTED (lapbelt adjustment 420mm - 1520mm)	(NC 0660, 0670, 0680, 0690, 0695, 0696, 0732, 0738, 0760, 0770, 0780, 0790, 0795, 0796, 0832, 0838, 0849, 0869, 0871, 0873, 0874, 0876, 0878, 0879, 0891, 0892, 0893, 0894)	○
CEE0180		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	(NC 0852)	○
CEE0190		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)		○
CEE0200		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)		◆
CEE0210		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)		○
CEE0220		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)		○
CEE0230		20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)		○
SEAT DEPTH*				
CEE0240		16" / 400 MM SEAT DEPTH	(NC 0879)	○
CEE0250		18" / 450 MM SEAT DEPTH	(NC 0878)	◆
CEE0252		16" TO 18" / 400 TO 450 MM ADJUSTABLE SEAT DEPTH	(NC 0330, 0350)	€24
SEAT HEIGHT*				
CEE0260		20" / 510 MM SEAT HEIGHT	(NC 0511, 0561, 0738, 1035, 1050, 1054, 1056, 1066)	○
CEE0270		19" / 485 MM SEAT HEIGHT	(NC 0511, 0561, 1054)	○
CEE0280		18" / 460 MM SEAT HEIGHT		◆
CEE0290		17" / 435 MM SEAT HEIGHT	(NC 1010, 1030)	○
CEE0300		16" / 410 MM SEAT HEIGHT	(NC 0935, 1010, 1030, 1056, 1370, 1440, 1441)	○
CEE0310		15" / 385 MM SEAT HEIGHT ,	(NC 0511, 0561, 0660, 0670, 0738, 0760, 0770, 0838, 0935, 1010, 1030, 1035, 1050, 1056, 1068, 1073, 1280, 1370, 1440, 1441)	○
SEAT UPHOLSTERY*				
CEE0330		NYLON W/O HOOK AND LOOP STRIPS		◆
CEE0340		NYLON WITH HOOK AND LOOP STRIPS		○
CEE0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM		○
CEE0350		PADDED NYLON		€36

BACKREST - Anti-tippers should be selected with reclining backrests 0435, 0435, 0440, 0445 and 0450*				
CEE0380		FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM	(NC 0471, 0477, 0945)	○
CEE0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT	(NC 0471, 0477, 0945)	○
CEE0395		FIXED 7° BACKREST HT ADJ 15" / 380 MM TO 18" / 460 MM FIX PUSH HANDLES	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	€80
CEE0396		FIXED 7° BACKREST HEIGHT ADJUST 15" / 380 TO 18" / 460 MM HT ADJ PUSH HANDLES	(NC 0474, 0477, 0478, 0481, 0905, 0906, 0907, 0940)	€105
CEE0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	(NC 0471, 0477, 0945)	◆
CEE0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES	(NC 0471, 0477, 0945)	○
CEE0425		FIXED 7° BACKREST 17" / 430 MM HEIGHT HT ADJ PUSHHANDLE	(NC 0471, 0477, 0945)	€50
CEE0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HT ADJ TO 16" / 400 MM FIX PUSH HANDLES	(NC 0471, 0477, 0480, 0481, 0905, 0906, 0907, 0940, 1280)	○
CEE0432		ANGLE ADJUSTABLE AND FOLD DOWN (15° to + 15°) 17"	(NC 0465, 0471, 0474, 0478, 0480, 0481, 0882, 0952)	€102
CEE0435		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 15" 16" 17"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190)	€152
CEE0445		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 15" 16" 17"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190)	€226
CEE0450		RECLINING (GAS SPRING) 0/30° 3 HEIGHTS 18" 19" 20"	(NC 0471, 0477, 0738, 0838, 0940, 0945, 1190)	€226
CEE0455		FOLDING 7° BACKREST 17"	(NC 0471, 0477, 0945)	€74
CEE0460		FOLDING 7° BACKREST 20"	(NC 0471, 0477, 0945)	€74
BACKREST UPHOLSTERY*				
CEE0465		WITHOUT BACK UPHOLSTERY	(NC 0882)	○
CEE0471		STANDARD NYLON SLING (JEMIMA)	(NC 0882, 0905, 0906, 0907)	○
CEE0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER		○
CEE0477		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER		€63
CEE0478		CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	€71
CEE0480		CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE	(NC 0940, 0945)	€142
CEE0481		CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE	(NC 0940, 0945)	€142
LEGREST HANGER RH*				
CEE0510		80° RH SWING IN/SWING OUT	(NC 1111, 66666)	◆
CEE0511		90° STRAIGHT RH SWING IN/ OUT HANGER BLACK COLOUR ONLY	(NC 0600, 0936, 1010, 1030, 1035, 1050, 1056, 1111, 66666)	○
CEE0520		ELEVATING LEGREST RH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 66666)	€150
CEE0545		AMPUTEE LEGREST RH	(NC 0600, 0610, 0611, 0930, 0936, 1111)	€159
LEGREST HANGER LH*				
CEE0560		80° LH SWING IN/SWING OUT	(NC 2222)	◆
CEE0561		90° STRAIGHT LH SWING IN/ OUT HANGER BLACK COLOUR ONLY	(NC 0620, 0936, 1010, 1030, 1035, 1050, 1056, 2222)	○
CEE0570		ELEVATING LEGREST LH SWING IN/SWING OUT	(NC 0611, 0930, 0935, 0936, 55555)	€150
CEE0595		AMPUTEE LEGREST LH	(NC 0611, 0620, 0630, 0930, 0936, 2222)	€159

FOOTREST RH*				
CEE0600		FLIP UP RH FOOTPLATE	(NC 0937, 0938)	◆
CEE0610		FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.		€50
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		€50
CEE0611		FLIP UP ONE PIECE	(NC 0620, 0630, 0935, 0936, 0937, 0938)	€93
CEE66666		NO FOOTREST		○
FOOTREST LH*				
CEE0620		FLIP UP LH FOOTPLATE	(NC 0937, 0938)	◆
CEE0630		FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.		€50
SWO2222		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)		€50
CEE55555		NO FOOTREST		○
ARMREST RH*				
CEE0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838, 1480, 1483)	○
CEE0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0760, 0770, 0780, 0790, 0795, 0796, 0838)	○
CEE0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0952, 0955, 1480, 1483)	◆
CEE0670		SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0740, 0750, 0780, 0790, 0795, 0796, 0838, 0952, 0955)	○
CEE0680		SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0952, 0955, 1480, 1483)	€67
CEE0690		SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH	(NC 0740, 0750, 0760, 0770, 0795, 0796, 0838, 0952, 0955, 1483)	€67
CEE0695		REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0952, 0955, 1480, 1483)	€71
CEE0696		REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0838, 0952, 0955, 1483)	€71
CEE0732		REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0838, 0952, 1480, 1483)	€69
CEE0738		REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0740, 0750, 0760, 0770, 0780, 0790, 0795, 0796, 0832, 0952, 0955, 1066, 1200, 1201, 1480, 1483)	€41

ARMREST LH*				
CEE0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	(NC 1480, 1484)	○
CEE0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED		○
CEE0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	(NC 0952, 0955, 1480, 1484)	◆
CEE0770		SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT	(NC 0952, 0955)	○
CEE0780		SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)	(NC 0952, 0955, 1480, 1484)	€67
CEE0790		SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH	(NC 0952, 0955, 1484)	€67
CEE0795		REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0952, 0955, 1480, 1484)	€71
CEE0796		REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED	(NC 0952, 0955, 1484)	€71
CEE0832		REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION	(NC 0952, 1480, 1484)	€69
CEE0838		REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)	(NC 0952, 0955, 1066, 1200, 1201, 1484)	€41
SEAT CUSHION				
CEE0849		BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER	(NC 0952)	€27
CEE0869		MATRX FLO-TECH LITE CUSHION - Low Risk	(NC 0952)	€168
CEE0871		MATRX FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK	(NC 0952)	€236
CEE0873		MATRX FLO-TECH IMAGE CUSHION - High Risk	(NC 0952)	€382
CEE0874		MATRX FLO-TECH SOLUTION XTRA CUSHION - Very High Risk	(NC 0952)	€566
CEE0876		MATRX LIBRA CUSHION - Very High Risk	(NC 0952)	€645
SEAT CUSHION DEPTH				
CEE0878		FOR 400 MM SEAT DEPTH		○
CEE0879		FOR 450 MM SEAT DEPTH		○
POSITIONING ACCESSORIES FOR MATRIX LIBRA & FLOTECH				
CEE0891		TAPERED WEDGES (PAIR)		€28
CEE0892		ABDUCTOR		€58
CEE0893		LATERAL WEDGE		€46
CEE0894		ANTERIOR SEAT WEDGE		€43

OTHER SEATING OPTIONS			
CEE0852		SEAT RIGIDIZER (NC 0952)	€79
CEE0882		LATERAL POSITIONING PADS (PAIR)	€72
CEE0905		NECKREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0906, 0907, 0940, 0945)	€158
CEE0906		HEADREST BLACK MULTI STRETCH PU COATED FABRIC (NC 0907, 0940, 0945)	€190
CEE0907		HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC (NC 0940, 0945)	€217
SWO3333		MULT ADJUSTABLE LATERALS	€340
CEE0930		CALF STRAP	◆
CEE0935		LONGER FOOTREST TUBE	€36
CEE0936		SHORTER FOOTREST SUPPORT	€48
CEE0937		BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED	€189
CEE0938		BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED	€189
CEE0940		BACKBRACE (NC 0945)	€41
CEE0945		FOLDABLE BACKBRACE	€44
CEE0952		INRIGGED ARMREST (NC 1480, 1483, 1484)	€48
CEE0955		W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○
COLOUR*			
CEE0961		SLATE GREY (RAL7024 FINE TEXTURE)	○
CEE0964		PEARL WHITE	○
CEE0965		ELECTRIC GREEN	○
CEE0966		MATT BLACK	○
CEE0969		RED MATT	○
CEE0972		WATER BLUE	○
CEA0974		BLUE	○
CEA0976		DARK PURPLE	○
CEE0980		ANTHRACITE BLACK	○
CEE0991		SPECIAL COLOUR ON DEMAND ONLY	€238
SWO4444		DARK RED SPECIAL	€238
SWO5555		PEARL GREY SPECIAL	€238
SWO6666		INTENSE YELLOW	€238
SWO7777		VIBRANT ORANGE	€238
SWO8888		POLAR BLUE	€238
SWO9999		CANDY PINK	€238
CASTOR HOUSING*			
CEE0995		STANDARD	◆
CEE0996		LIGHT WEIGHT	€35
FRONT CASTORS*			
CEE1010		8" X 1.1/4" BLACK SOLID TYRE	◆
CEE1030		8" X 2" BLACK SOLID TYRE	€32
CEE1035		6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD	€32
CEE1050		6" X 1.1/4" BLACK SOLID TYRE	€32
CEE1054		5" X 1" BLACK SOLID TYRE	€32
CEE1056		5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD	€32

REAR WHEEL SIZE*			
CEE1066		22" (NC 1092, 1260)	○
CEE1068		24"	◆
REAR WHEELS*			
CEE1071		STANDARD 4 PARTS FIXING HANDRIM (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1092, 1170, 1280)	◆
CEE1072		REINFORCED 6 PARTS FIXING HANDRIM (COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1092, 1280)	€24
CEE1073		DRUM BRAKE SELFPROPEL 22 AND 24" - 4 PARTS FIXING HANDRIM - COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES. TRANSIT 12" BLACK PLASTIC HUB AND RIM - DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1092, 1201, 1270, 1275)	○
CEE1076		LIGHTWEIGHT 6 PARTS FIXING HANDRIM (BLACK HUB AND RIM, 24 BLACK RADIAL SPOKES) DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM) (NC 1170, 1201, 1230, 1260, 1261, 1280)	€70
TYRES*			
CEE1082		SOLID BLACK 1.3/8" / 35MM	◆
CEE1084		PNEUMATIC PROFILE BLACK 1.3/8" / 35MM	€24
CEE1088		PNEUMATIC LOW PROFILE BLACK 1" / 25MM	○
CEE1092		PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM	€41
HANDRIM*			
CEE1140		ALUMINIUM PREANODIZED HANDRIM	◆
CEE1150		HANDRIM COATED IN BLACK ANTISLIP PLASTIC	€40
REAR WHEEL AXLE*			
CEE1170		FIXED AXLE	○
CEE1180		Q/R AXLE	◆
REAR WHEEL POSITION*			
CEE1190		STANDARD REAR WHEELS POSITION	◆
CEE1200		PASSIVE REAR WHEELS POSITION	○
CEE1201		DUAL AMPUTEE REAR WHEEL POSITION (NC 1275, 1280)	€84
PUMP			
CEE1230		PUMP AND SUPPORT	€16
CHASSIS OPTION			
CEE1241		FRONT WHEELS EXTENSION	€30
CEE1260		SPOKE GUARDS PAIR 24" BLK (NC 1261)	€93
CEE1261		SPOKE GUARDS PAIR CLEAR	€118
BRAKES*			
CEE1270		PUSH TO LOCK BRAKES	◆
CEE1275		PULL TO LOCK BRAKES	○
CEE1280		ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)	€243
TRANSPORTATION KIT*			
CEE1345		TRANSPORTATION BRACKETS	◆

SAFETY OPTIONS*				
CEE1360		LAP BELT	(NC 1361, 1362, 1363, 1364)	◆
CEE1361		BODYPOINT 2-PT POSTURE BELT M	(NC 1362, 1363, 1364)	€171
CEE1362		BODYPOINT 2-PT POSTURE BELT L	(NC 1363, 1364)	€171
CEE1363		BODYPOINT 4-PT POSTURE BELT M	(NC 1364)	€161
CEE1364		BODYPOINT 4-PT POSTURE BELT L		€161
CEE1370		REAR ANTITIPPERS PAIR	(NC 1375, 1440, 1441)	€43
CEE1375		REAR ANTITIPPERS PAIR WITH TIPPING AID	(NC 1440, 1441)	€85
CEE1380		REFLECTORS KIT		€10
CEE1390		BRAKE EXTENSION PAIR		€16
OTHER OPTIONS				
CEE1430		CANE HOLDER		€28
CEE1440		TIPPING AID RH		€16
CEE1441		TIPPING AID LH		€16
CEE1480		ACRYLIC TRAY	(NC 1483, 1484, 1485)	€150
CEE1483		HALF TABLE TRAY RH	(NC 1484)	€62
CEE1484		HALF TABLE TRAY LH		€62
CEE1485		COMFORT PADDING FOR HALF TABLE TRAY		€30
SWO10002		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618		€105

TOTAL PRICE				
-------------	--	--	--	--

Non Conforming Options:	
CEE0330	NYLON W/O HOOK AND LOOP STRIPS
CEE0350	PADDED NYLON
CEE0465	WITHOUT BACK UPHOLSTERY
CEE0471	STANDARD NYLON SLING (JEMIMA)
CEE0474	LUMBAR SUPPORT BI MATERIAL PADDED COVER
CEE0477	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER
CEE0478	CONTOUR ADJUSTABLE (3+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEE0480	CONTOUR ADJUSTABLE (7+1 STRAPS) PADDED SHAPED COVER INCL BACKBRACE
CEE0481	CONTOUR ADJUSTABLE (7+1 STRAPS) BI MATERIAL PADDED COVER INCL BACKBRACE
CEE0511	90° STRAIGHT RH SWING IN/ OUT HANGER BLACK COLOUR ONLY
CEE0561	90° STRAIGHT LH SWING IN/ OUT HANGER BLACK COLOUR ONLY
CEE0600	FLIP UP RH FOOTPLATE
CEE0610	FLIP UP DEPTH AND ANGLE AJUSTABLE RH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE0611	FLIP UP ONE PIECE
CEE66666	NO FOOTREST
CEE0620	FLIP UP LH FOOTPLATE
CEE0630	FLIP UP DEPTH AND ANGLE AJUSTABLE LH FOOTPLATE NB: IF SELECTED WITH A FIXED FOOTREST FOOTPLATE WILL BE ANGLE AND DEPTH ADJ AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJ NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO2222	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
CEE55555	NO FOOTREST
CEE0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CEE0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0680	SWING AWAY DETACHABLE RH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)
CEE0690	SWING AWAY DETACHABLE RH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH
CEE0695	REMOVABLE RH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEE0696	REMOVABLE RH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEE0732	REMOVABLE RH (N°4) ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION AUTOLOCK SAFETY DEVICE INCLUDED
CEE0738	REMOVABLE MUDGUARD RH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)
CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0780	SWING AWAY DETACHABLE LH (N°3) DESK ADJ. IN HEIGHT (170 to 295 MM)
CEE0790	SWING AWAY DETACHABLE LH (N°3) FULL LENGTH ADJ. IN HEIGHT (170 to 295 MM) AND DEPTH
CEE0795	REMOVABLE LH (N°4) DESK ADJ. IN HEIGHT (210345 MM) + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEE0796	REMOVABLE LH (N°4) FULL LENGTH ADJ. IN HEIGHT (210345 MM) AND DEPTH + SKIRTGUARD HEIGHT ADJ. AUTOLOCK SAFETY DEVICE INCLUDED
CEE0832	REMOVABLE LH N°4 ARMPAD HEMI BLACK MULTI STRETCH PU COATED FABRIC MULTI POSITION ARMPAD HEMI MULTI POSITION

CEE0838	REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) HEIGHT ADJUSTABLE (155 TO 220 MM)
CEE0849	BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER
CEE0869	MATRX FLO-TECH LITE CUSHION - Low Risk
CEE0871	MATRX FLO-TECH CONTOUR VISCO NG CUSHION - HIGH RISK
CEE0873	MATRX FLO-TECH IMAGE CUSHION - High Risk
CEE0874	MATRX FLO-TECH SOLUTION XTRA CUSHION - Very High Risk
CEE0876	MATRX LIBRA CUSHION - Very High Risk
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE0891	TAPERED WEDGES (PAIR)
CEE0892	ABDUCTOR
CEE0893	LATERAL WEDGE
CEE0894	ANTERIOR SEAT WEDGE
CEE0852	SEAT RIGIDIZER
CEE0882	LATERAL POSITIONING PADS (PAIR)
CEE0905	NECKREST BLACK MULTI STRETCH PU COATED FABRIC
CEE0906	HEADREST BLACK MULTI STRETCH PU COATED FABRIC
CEE0907	HEADREST WITH CHEEK SUPPORT BLACK MULTI STRETCH PU COATED FABRIC
CEE0930	CALF STRAP
CEE0935	LONGER FOOTREST TUBE
CEE0936	SHORTER FOOTREST SUPPORT
CEE0937	BODYPOINT ANKLE HUGGERS (PAIR) SIZE M - UNMOUNTED
CEE0938	BODYPOINT ANKLE HUGGERS (PAIR) SIZE L - UNMOUNTED
CEE0940	BACKBRACE
CEE0945	FOLDABLE BACKBRACE
CEE0952	INRIGGED ARMREST
CEE0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEE1010	8" X 1.1/4" BLACK SOLID TYRE
CEE1030	8" X 2" BLACK SOLID TYRE
CEE1035	6" X 1.1/4" BLACK SOLID TYRE MEDIUM HARD
CEE1050	6" X 1.1/4" BLACK SOLID TYRE
CEE1054	5" X 1" BLACK SOLID TYRE
CEE1056	5.1/2" X 1.3/4" BLACK SOLID TYRE MEDIUM HARD
CEE1066	22"
CEE1068	24"
CEE1073	DRUM BRAKE SELFPROPEL 22 AND 24" - 4 PARTS FIXING HANDRIM - COLOURLESS HUB, CLEAR ANODIZED RIM AND STAINLESS STEEL SPOKES. TRANSIT 12" BLACK PLASTIC HUB AND RIM - DISTANCE WHEEL-HANDRIM 11 OR 16MM (SET AT 11MM)
CEE1092	PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM
CEE1170	FIXED AXLE
CEE1190	STANDARD REAR WHEELS POSITION
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1201	DUAL AMPUTEE REAR WHEEL POSITION
CEE1230	PUMP AND SUPPORT
CEE1260	SPOKE GUARDS PAIR 24" BLK

CEE1261	SPOKE GUARDS PAIR CLEAR
CEE1270	PUSH TO LOCK BRAKES
CEE1275	PULL TO LOCK BRAKES
CEE1280	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)
CEE1361	BODYPOINT 2-PT POSTURE BELT M
CEE1362	BODYPOINT 2-PT POSTURE BELT L
CEE1363	BODYPOINT 4-PT POSTURE BELT M
CEE1364	BODYPOINT 4-PT POSTURE BELT L
CEE1370	REAR ANTITIPPERS PAIR
CEE1375	REAR ANTITIPPERS PAIR WITH TIPPING AID
CEE1440	TIPPING AID RH
CEE1441	TIPPING AID LH
CEE1480	ACRYLIC TRAY
CEE1483	HALF TABLE TRAY RH
CEE1484	HALF TABLE TRAY LH
CEE1485	COMFORT PADDING FOR HALF TABLE TRAY



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):