



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

DECEMBER 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Action 5 Rigid

Standard Features

- Max user weight 130kg
- Crash tested ISO 7176 - Part 9
- Self Propulsion with 22" or 24" quick release rear wheels
- Choice of 7 seat widths 15" - 21" (380mm to 530mm)
- Choice of 7 seat depths 14" - 20" (350mm to 500mm)
- Choice of 8 front seat to floor heights 14" - 21" (360mm to 535mm)
- Choice of 7 rear seat to floor heights 14" - 20" (360mm to 510mm)
- Tension adjustable nylon seat upholstery
- Fixed backrest (90° to the seat rail)
- Fixed backrest, nylon and mesh padded cover
- Fixed front end or 80° swing in/swing out legrest hangers
- Flip up footplates
- sideguards / skirtguards
- Choice of colours and 2 component accents
- 2° Camber
- Push to lock brakes
- Transportation kit with lap belt
- Adjustable heel straps
- User tool kit



Standard Base price: £1,288



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SJE0016	ACTION 5 RIGID	(NC 1510) ☐
SJESPEC	SPECIAL	POA
FRAME*		
SJE0200	SWING AWAY FRAME	(NC 3230, 3440, 3450) £1,288
SJE0210	FIXED FRONT END 80 DEGREE ANGLE FRAME	(NC 1111, 3010, 3020, 3030, 3040, 3050, 3110, 3120, 3130, 3140, 3150, 3333, 3400, 3410, 3420, 3430) £1,288
PROPULSION*		
SJE0300	SELF PROPULSION	☒

SEAT WIDTH*			
SJE1030		15" SW	○
SJE1040		16" SW	○
SJE1050		17" SW	○
SJE1060		18" SW	○
SJE1070		19" SW	○
SJE1075		20" SW (NC 1100, 1110, 1170)	○
SJE1078		21" SW (NC 1100, 1110, 1170)	○
SEAT DEPTH*			
SJE1100		14" SD (NC 1270, 3040, 3140, 7080, 7100, 7110, 7600, 7610)	○
SJE1110		15" SD (NC 1270, 7600)	○
SJE1120		16" SD (NC 1270)	○
SJE1130		17" SD	○
SJE1140		18" SD	○
SJE1150		19" SD	○
SJE1160		20" SD	○
SJE1170		14" - 16" ADJ SD (NC 1510, 7600)	○
SJE1180		16" - 18" ADJ SD (NC 1510)	○
SJE1190		18" - 20" ADJ SD (NC 1510)	○
FRONT SEAT HEIGHT*			
SJE1200		14" FSH (NC 1310, 1320, 1330, 1340, 1350, 1360, 3420, 3430, 3450, 7040, 7055, 7070, 7080, 7100, 7110)	○
SJE1210		15" FSH (NC 1320, 1330, 1340, 1350, 1360, 3430, 7000, 7040, 7070, 7080, 7100, 7110)	○
SJE1220		16" FSH (NC 1330, 1340, 1350, 1360, 7080, 7100, 7110)	○
SJE1230		17" FSH (NC 1300, 1340, 1350, 1360, 7110)	○
SJE1240		18" FSH (NC 1300, 1310, 1350, 1360)	○
SJE1250		19" FSH (NC 1300, 1310, 1320, 1360)	○
SJE1260		20" FSH (NC 1300, 1310, 1320, 1330, 7000, 7040)	○
SJE1270		21" FSH (NC 1300, 1310, 1320, 1330, 1340, 3030, 3130, 7000, 7040, 7055, 7070, 7080, 7100)	○
REAR SEAT HEIGHT*			
SJE1300		14" RSH (NC 7250, 7260)	○
SJE1310		15" RSH	○
SJE1320		16" RSH	○
SJE1330		17" RSH	○
SJE1340		18" RSH	○
SJE1350		19" RSH (NC 7200)	○
SJE1360		20" RSH (NC 7200, 7230, 7240)	○
SEAT UPHOLSTERY*			
SJE1510		NYLON FIXED	○
SJE1520		TENSION AND DEPTH ADJUSTABLE	○
BACKREST ATTACHMENT*			
SJE2000		FIXED 90°	○
SJE2020		ANGLE ADJUSTABLE (-15° TO +15°) AND FOLDING	£100

BACKREST FIXED HEIGHT - INCLUDING FIXED PUSH HANDLES*				
SJE2143		14" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 2505, 8160, 8200)	○
SJE2150		15" BACK HT	(NC 2300, 2310, 2400, 2410, 2420, 8200)	○
SJE2153		16" BACK HT	(NC 2300, 2310, 2400, 2410, 2420)	○
SJE2160		17" BACK HT	(NC 2300, 2310, 2400, 2410, 2420)	○
SJE2163		18" BACK HT	(NC 2300, 2310, 2400, 2410, 2420)	○
SJE2170		19" BACK HT	(NC 2300, 2310, 2400, 2410, 2420)	○
SJE2173		20" BACK HT	(NC 2300, 2310, 2400, 2410, 2420)	○
SJE7777		No Fixed backrest, use only when selecting height adjustable backrest, option 2300 or 2310	(NC 8888)	○
BACKREST HEIGHT ADJUSTABLE - PUSH HANDLES MUST BE SELECTED IN THE SECTION BELOW*				
SJE2300		HT ADJ 14.5" - 16"	(NC 2505)	£25
SJE2310		HT ADJ 16" - 19.5"		£25
SJE8888		No height adjustable backrest, use only when selecting Fixed backrest, options 2143-2173		○
PUSH HANDLES - MUST BE SELECTED IN COMBINATION WITH HEIGHT ADJUSTABLE BACKRESTS				
SJE2400		W/O PUSH HANDLES	(NC 8160, 9100, 9200)	○
SJE2410		REAR SET	(NC 2505)	£49
SJE2420		INTEGRATED		£81
BACKREST UPHOLSTERY*				
SJE2500		W/O BACK UPH - DELIVERED WITH SIMPLE SLEEVE		○
SJE2505		FIXED NYLON AND MESH	(NC 5600)	◆
SJE2510		CONTOUR ADJUSTABLE		£74
LEGREST HANGER RH*				
SJE3000		W/O HANGER RH	(NC 5400)	○
SJE3010		70° RH	(NC 1111, 3440, 3450, 22222)	○
SJE3020		80° RH	(NC 1111, 3440, 3450, 22222)	○
SJE3030		90° RH	(NC 1111, 3200, 3440, 3450, 7080, 7100, 7110, 22222)	○
SJE3040		ELEVATING LEGREST RH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM	(NC 3230, 3400, 3440, 3450, 5400, 22222, 8120)	£116
SJE3050		AMPUTEE LEGREST RH	(NC 1111, 3200, 3210, 3230, 5400, 8100, 8120)	£133
LEGREST HANGER LH*				
SJE3100		W/O HANGER LH	(NC 5400)	○
SJE3110		70° LH	(NC 3333, 3440, 3450, 44444)	○
SJE3120		80° LH	(NC 3333, 3440, 3450, 44444)	○
SJE3130		90° LH	(NC 3300, 3333, 3440, 3450, 44444, 7080, 7100, 7110)	○
SJE3140		ELEVATING LEGREST LH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM	(NC 3400, 3440, 3450, 44444, 5400, 8120)	£116
SJE3150		AMPUTEE LEGREST LH	(NC 3300, 3310, 3333, 5400, 8100, 8120)	£133

FOOTREST RH*			
SJE3200		FLIP UP RH (NC 3410)	◆
SJE3210		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.	£43
SWO1111		FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	£43
SJE3230		ONE PIECE (NC 3400, 3410, 3420, 3430)	£186
SJE22222		NO FOOTPLATE REQUIRED	○
FOOTREST LH*			
SJE3300		FLIP UP LH (NC 3410)	◆
SJE3310		FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.	£43
SWO3333		FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)	£43
SJE44444		NO FOOTPLATE REQUIRED	○
FOOTREST TUBE*			
SJE3400		SWING AWAY FRAME HIGH MOUNTED TUBE	○
SJE3410		SWING AWAY FRAME SHORT TUBE	○
SJE3420		SWING AWAY FRAME STANDARD TUBE	○
SJE3430		SWING AWAY FRAME LONG TUBE	○
SJE3440		FIXED FRAME STANDARD TUBE	○
SJE3450		FIXED FRAME LONG TUBE	○
MUDGUARDS/SIDEGUARDS*			
SJE3700		W/O SIDEGUARD/MUDGUARD	○
SJE3750		SIDEGUARDS (NC 4010, 4020, 4030, 4040, 4110, 4120, 4130, 4140)	◆
SJE3800		MUDGUARDS (NC 4010, 4020, 4030, 4040, 4110, 4120, 4130, 4140, 8160)	£74
ARMREST RH*			
SJE4000		W/O ARMREST	○
SJE4010		T-ARM DESK ARMPAD RH (NC 7600, 7610)	£74
SJE4020		T-ARM FULL ARMPAD RH	£74
SJE4030		SWING AWAY DESK ARMPAD RH	£99
SJE4040		SWING AWAY FULL ARMPAD RH	£99
ARMREST LH*			
SJE4100		W/O ARMREST	○
SJE4110		T-ARM DESK ARMPAD LH (NC 7600, 7610)	£74
SJE4120		T-ARM FULL ARMPAD LH	£74
SJE4130		SWING AWAY DESK ARMPAD LH	£99
SJE4140		SWING AWAY FULL ARMPAD LH	£99

SEAT CUSHION*			
SJE55555		WITHOUT SEAT CUSHION	○
SJE5000		BASIC FOAM CUSHION	£21
SEATING OPTIONS			
SJE5400		CALF STRAP	○
SJE5600		FIXED BACKBRACE	£21
FRAME COLOUR*			
SJE6000		PEARL WHITE	○
SJE6010		ANTHRACITE BLACK	○
SJE6015		BLUE	○
SJE6040		ELECTRIC GREEN	○
SJE6045		DARK PURPLE	○
SJE6060		PEARL GREY	○
SJE6070		MATT BLACK	○
SJE6099		SPECIAL COLOUR	£169
SWO10002		SLATE GREY - SPECIAL COLOUR 3 WEEK LEAD TIME	£169
SWO10003		RED MATT - SPECIAL COLOUR 3 WEEK LEAD TIME	£169
SWO10004		POLAR BLUE - SPECIAL COLOUR 3 WEEK LEAD TIME	£169
SWO10005		CANDY PINK - SPECIAL COLOUR 3 WEEK LEAD TIME	£169
ACCENT COLOUR*			
SJE6100		WITHOUT ACCENT COLOUR SILVER FORK AND REARWHEEL SUPPORT	○
SJE6110		BLUE ACCENT	○
SJE6200		PURPLE ACCENT	○
FRONT CASTORS*			
SJE7000		4" X 1.1/4 SOFTROLL	£37
SJE7040		4" X 1.3/8" ALU SOFTROLL	£37
SJE7055		5" X 1" SOLID	◆
SJE7070		5" X 1.3/8" ALU SOFTROLL	£37
SJE7080		6" X 1" SOLID	○
SJE7100		6" X 1.3/8" ALU SOFTROLL	£37
SJE7110		7" X 1.1/4" SOLID	○
REAR WHEELS*			
SJE7200		22" STANDARD, CROSS SPOKED WHEEL (NC 7320, 7410)	○
SJE7230		24" STANDARD, CROSS SPOKED WHEEL (NC 7320)	◆
SJE7240		24" LIGHTWEIGHT WHEEL (NC 7410, 8000, 8160)	£59
SJE7250		25" LIGHTWEIGHT WHEEL (NC 7302, 7304, 7315, 7410, 8000, 8160)	£59
SJE7260		25" SPINERGY WHEELS (NC 7302, 7304, 7315, 7410, 7420, 8000, 8160)	£561
REAR WHEEL TYRE*			
SJE7302		SOLID BLACK 1.3/8"	◆
SJE7304		PNEUMATIC TYRE BLACK 1.3/8"	○
SJE7315		PNEUMATIC LOW PROFILE BLACK 1"	○
SJE7320		PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM	£70

HANDRIM*			
SJE7400		ALUMINIUM HANDRIM	◆
SJE7410		STAINLESS STEEL	○
SJE7420		HANDRIM COATED IN BLACK ANTI-SLIP	£49
REAR WHEEL AXLE*			
SJE7500		Q/R AXLE	◆
REAR WHEEL POSITION*			
SJE7600		VERY ACTIVE	○
SJE7610		ACTIVE	○
SJE7620		MEDIUM ACTIVE	○
SJE7630		STANDARD	◆
CAMBER ANGLE*			
SJE7700		0° CAMBER	○
SJE7710		1° CAMBER	○
SJE7720		2° CAMBER	◆
SJE7730		3° CAMBER	○
SJE7740		4° CAMBER	○
CHASSIS OPTION			
SJE8000		SPOKE PROTECTORS	£81
PUMP			
SJE8030		PUMP - HIGH PRESSURE	£41
BRAKES*			
SJE8100		PUSH TO LOCK BRAKES	◆
SJE8110		PULL TO LOCK BRAKES	○
SJE8120		PUSH TO LOCK BRAKES WITH EXTENTION	£17
SJE8130		PULL TO LOCK BRAKES WITH EXTENTION	£17
SJE8160		ATTENDANT HUB BRAKES INCLUDES STANDARD USER PUSH TO LOCK BRAKE (OR PULL VERSION IF IT APPEARS AS THE ONLY CHOICE ABOVE)	£231
TRANSPORTATION KIT			
SJE8200		TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
SJE8300		LAP BELT	◆
SJE8500		REAR ANTI-TIPPER ACTIVE RH (NC 9200)	£64
SJE8550		REAR ANTI-TIPPER ACTIVE LH (NC 9100)	£64
SJE8600		REFLECTORS KIT (WHEELS)	£17
ATTENDANT OPTION			
SJE9000		CANE HOLDER	£23
SJE9100		TIPPING AID LH	£49
SJE9200		TIPPING AID RH	£49
OTHER OPTIONS			
SWO10001		ABDUCTED LEGRESTS (PR) 0° TO 30° - INTERNAL REFERENCE DE259618	£107
1595002		ViaGo Push aid. Includes as standard: Drive Unit, Charger, Control unit & Brackets (Note: 10 minute fitting time)	£1,995

TOTAL PRICE	
--------------------	--

Non Conforming Options:	
SJE1100	14" SD
SJE1110	15" SD
SJE1170	14" - 16" ADJ SD
SJE1270	21" FSH
SJE1300	14" RSH
SJE1310	15" RSH
SJE1320	16" RSH
SJE1330	17" RSH
SJE1340	18" RSH
SJE1350	19" RSH
SJE1360	20" RSH
SJE1510	NYLON FIXED
SJE2300	HT ADJ 14.5" - 16"
SJE2310	HT ADJ 16" - 19.5"
SJE8888	No height adjustable backrest, use only when selecting Fixed backrest, options 2143-2173
SJE2400	W/O PUSH HANDLES
SJE2410	REAR SET
SJE2420	INTEGRATED
SJE2505	FIXED NYLON AND MESH
SJE3010	70° RH
SJE3020	80° RH
SJE3030	90° RH
SJE3040	ELEVATING LEGREST RH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM
SJE3050	AMPUTEE LEGREST RH
SJE3110	70° LH
SJE3120	80° LH
SJE3130	90° LH
SJE3140	ELEVATING LEGREST LH, WITH CALF PAD INCLUDED - WITH SEAT HEIGHT 360 AND 385 MM, LEGREST MAY BE IN CONTACT WITH THE GROUND IF UL IS MORE THAN 400MM
SJE3150	AMPUTEE LEGREST LH
SJE3200	FLIP UP RH
SJE3210	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO1111	FLIP UP DEPTH AND ANGLE ADJUSTABLE RH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
SJE3230	ONE PIECE
SJE22222	NO FOOTPLATE REQUIRED
SJE3300	FLIP UP LH
SJE3310	FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: IF SELECTED WITH A FIXED FOOTREST - FOOTPLATE WILL BE ANGLE AND DEPTH ADJUSTABLE AND COMPATIBLE WITH ANKLE HUGGERS. IF SELECTED WITH AN ELEVATING LEGREST - FOOTPLATE WILL BE ANGLE, DEPTH AND WIDTH ADJUSTABLE - NOT COMPATIBLE WITH ANKLE HUGGERS.
SWO3333	FLIP UP DEPTH AND ANGLE ADJUSTABLE LH FOOTPLATE - NB: COMPATIBLE WITH ANKLE HUGGERS (FOR ELEVATING LEGRESTS ONLY)
SJE44444	NO FOOTPLATE REQUIRED

SJE3400	SWING AWAY FRAME HIGH MOUNTED TUBE
SJE3410	SWING AWAY FRAME SHORT TUBE
SJE3420	SWING AWAY FRAME STANDARD TUBE
SJE3430	SWING AWAY FRAME LONG TUBE
SJE3440	FIXED FRAME STANDARD TUBE
SJE3450	FIXED FRAME LONG TUBE
SJE4010	T-ARM DESK ARMPAD RH
SJE4020	T-ARM FULL ARMPAD RH
SJE4030	SWING AWAY DESK ARMPAD RH
SJE4040	SWING AWAY FULL ARMPAD RH
SJE4110	T-ARM DESK ARMPAD LH
SJE4120	T-ARM FULL ARMPAD LH
SJE4130	SWING AWAY DESK ARMPAD LH
SJE4140	SWING AWAY FULL ARMPAD LH
SJE5400	CALF STRAP
SJE5600	FIXED BACKBRACE
SJE7000	4" X 1.1/4 SOFTROLL
SJE7040	4" X 1.3/8" ALU SOFTROLL
SJE7055	5" X 1" SOLID
SJE7070	5" X 1.3/8" ALU SOFTROLL
SJE7080	6" X 1" SOLID
SJE7100	6" X 1.3/8" ALU SOFTROLL
SJE7110	7" X 1.1/4" SOLID
SJE7200	22" STANDARD, CROSS SPOKED WHEEL
SJE7230	24" STANDARD, CROSS SPOKED WHEEL
SJE7240	24" LIGHTWEIGHT WHEEL
SJE7250	25" LIGHTWEIGHT WHEEL
SJE7260	25" SPINERGY WHEELS
SJE7302	SOLID BLACK 1.3/8"
SJE7304	PNEUMATIC TYRE BLACK 1.3/8"
SJE7315	PNEUMATIC LOW PROFILE BLACK 1"
SJE7320	PNEUMATIC AND PUNCTURE RESISTANT MAXXIS ALPHA PROTECT TYRE 1" / 25 MM
SJE7410	STAINLESS STEEL
SJE7420	HANDRIM COATED IN BLACK ANTI-SLIP
SJE7600	VERY ACTIVE
SJE7610	ACTIVE
SJE8000	SPOKE PROTECTORS
SJE8100	PUSH TO LOCK BRAKES
SJE8120	PUSH TO LOCK BRAKES WITH EXTENTION
SJE8160	ATTENDANT HUB BRAKES INCLUDES STANDARD USER PUSH TO LOCK BRAKE (OR PULL VERSION IF IT APPEARS AS THE ONLY CHOICE ABOVE)
SJE8200	TRANSPORTATION BRACKETS

SJE9100	TIPPING AID LH
SJE9200	TIPPING AID RH



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):