

E-MOTION® M25



Please refer to www.invacare.ie for local distributor.

Please fill in all specifications with each order.

Account code

Customer order number.....

Company name

Address



- e-motion M25** €6,447
- Included as standard:
 - 2 e-motion drive wheels 22"
 - Speed: 6km/h 24"
 - Weight: Only 7.8 kg per wheel 26"
 - Brushed stainless steel push rims
 - Charger
 - Lithium-ion batteries in the wheel hubs
 - Pneumatic tyres with puncture protection

- Ergonomic Control System (ECS)** 1592486 €523
- Switches e-motion on and off
 - Two remote assistance levels
 - Roll back delay
 - Battery monitor
 - Self diagnosis



- Anti-tippers with jack-up function (pair)** 1489214 €382
- NB Always use anti-tippers (Alber or existing)

Mobility Plus Package

Download from App Store or Playstore as an In-App purchase

- Remote control
- Increase speed to 8.5km/h
- Cruise control function,
- EasyNavi wheelchair-friendly navigation
- ECS functions
- Push counter - evaluates your driving behaviour during journeys



Brackets required for installation onto wheelchair

- 1467178 STD
- Omit brackets €96

Fitting on wheelchair at Invacare factory & delivery* €86

Wheelchair order number.....

Collection & fitting to third party wheelchair & delivery** €86 + €110


Fitting by trained dealer** €0

*Only when ordered with an Invacare wheelchair.
**Mounting kits can only be sent to selected trained dealers.

Other

- Equipped with vehicle retention system
- Rear wheel extension available

Available wheel position:



- 1=rear
- 2=centre
- 3=front

Wheelchair details

Wheelchair manufacturer (e.g. Kuschall).....

Serial no.....

Year.....

Model (e.g. Compact).....

Seat width

Seat height (e.g.Hemi)

Degree of wheel camber

(Drum brakes must be removed)

E-MOTION® M25



Accessories

Push rims

- Stainless steel (standard) 1592526 €0
- 22" or 24" push rim cover, Tetra (pair) 1592527 €170
- CarboLife Curve L Black Powder coated 1592528 €369
- CarboLife Quadro Tetra-Gripp coated 1592520 €446



Tyres (pairs)

- Schwalbe Rightrun (Standard) Black 1592530 €0
- Schwalbe Marathon Plus Evolution with SmartGuard puncture protection, Black 1565894 €47
- Pro Speed puncture-proof solid PU 1" Black (+ 300g per wheel, not available on 22" or 26") 1592532 €125
- Puncture proof with inlay 1 3/8" (+ 500g per wheel, not available on 26") 1592533 €125

Push rim covers (retrofit)


- Push rim cover 22" 1592425 €183
- Push rim cover 24" 1592421 €183
- Push rim cover 26" 1592426 €183

Spoke Covers (pairs)

- 22" spoke cover Design 'Comic' 1592266 €183
- 24" spoke cover Design 'Sphere' 1592268 €190
- 26" spoke cover, Design 'Sphere' 1592269 €198
- 24" spoke cover, Plain 1592553 €190



Extra options

- Bracket for remote control  1592232 €190
 - LH RH
- Brake extension* 1467780 €34
- Transport axle 1591362 €93
- Ergonomic stubaxle 1591500 €178



Do you want to use both your standard and e-motion wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for e-motion * / **

Retrofit, stub axle to 1/2" diameter €134

*If the manual wheels are used in combination with e-motion, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.

**If drum brakes are installed, these must be removed as they are not compatible with e-motion. Standard lever brakes are used instead.

**maybe required depending on wheelchair type*

Comments (mounting etc.).....

.....

.....

.....