



Invacare Ltd  
 Unit 4 Pencoed Technology Park  
 Pencoed, Bridgend  
 CF35 5AQ  
 T: 01656 776222  
 E: ordersuk@invacare.com  
 www.invacare.co.uk



Please fill in all specifications with each order.

Account code

Customer order number

Company name

Address



**e-fix E35** **£4800**

- Included as standard: 22"  24"
- 2 e-fix drive wheels
  - LED display unit
  - Mobility Plus package
  - Range 16km (10 miles)
  - Max. user weight 120kg
  - 6.0 Ah battery (2kg) 7.5 Ah battery (2kg) £235
  - Wheel weight 7.8kg
  - Total weight 18.9kg

**e-fix E36** **£5565**

- Included as standard: 24"
- 2 e-fix drive wheels
  - LED display unit
  - Mobility Plus package
  - Range 20km (12.4 miles)
  - Max. user weight 160kg
  - 7.5 Ah battery (2kg) 6Ah battery (2kg) -£235
  - Wheel weight 8kg
  - Total weight 19.3kg

**Wheelchair details**

Wheelchair manufacturer (e.g. Kuschall)

Serial no Year

Model (e.g. Compact)

Seat width  cm

Seat height (e.g. Hemi)  cm

Armrest & Brake type

Degree of wheel camber °

*(Drum brakes must be removed)*

**Brackets required for installation onto wheelchair**

1457841  **£94**

Omit brackets

Fitting on wheelchair at Invacare factory & delivery\*  **£70**

Wheelchair order number

Collection & fitting to third party wheelchair & delivery\*\*  **£70 + £90**

Fitting by trained dealer\*\*\*  **£0**

**Control Unit Mounting**

Right side

Left side  **£63**

Right & Left (pluggable)

Bracket on: Frame mount  Armrest mount

*\*Only when ordered with an Invacare wheelchair.*

*\*\*Pending compatibility*

*\*\*\*Mounting kits can only be sent to selected trained dealers.*

**Do you want to use both your standard and e-fix wheels in the future? If so, please order the retrofit stub axle.**

**Manual wheel axle retrofit for e-fix \* / \*\***

Retrofit, stub axle to 1/2" diameter  **£96**

*\*If the manual wheels are used in combination with e-fix, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.*

*\*\*If drum brakes are installed, these must be removed as they are not compatible with e-fix. Standard lever brakes are used instead.*



## Accessories

Anti-tippers with jack-up function (pair)\* 1489214 £263

### Control Unit Accessories

Swivel arm 1467901 £334

Protection Bumper 1565836 £48

### Joystick Options

Option special 40mm/1.57" 1565669 £25

Option soft 70mm/2.76" 1566196 £25

Option T - handle 1565670 £82

*\*NB Always use anti-tippers (Alber or existing)*

### Tray table incl. control unit bracket

Control unit centered SW 38-43cm 1566207 £235

Control unit centered SW 43-48cm 1566208 £235

Installation kit for control unit with existing tray table 1566206 £138

### Extension for speed pre-selection

Set for speed pre-selection 1566197 £38

### Extension for coupling ring

Set for one drive wheel (single) 1566179 £96

Set for two drive wheels (pair) 1566178 £192

### Spoke Covers (pair)

22" wheel, decor 'Alber' 1566201 £60

22" wheel, transparent 1566203 £60

24" wheel, decor 'Alber' 1566200 £60

24" wheel, transparent 1566202 £60

### Attendant Control Bracket

Right 1457754 £117

Left 1457755 £117

### Intuitive Attendant Control

Control unit intuitive attendant control 1566165 £428

Adapter cable intuitive attendant control 1566177 £88

Additional handle 1566166 £42

Opti-box 1566184 £329

Speed/Stop limit switch 1467698 £70

### Battery pack

Spare battery 6.0 Ah 1565639 £517

Spare battery 7.5 Ah 1566153 £679

Transport bag for battery 1566168 £68

Attachment set for transport bag 1566169 £20

### Battery holder fixation

Fixed assembly set for cross brace 1566164 £163

Fixed assembly set for cross tube 1566163 £140

### Mobility Plus Package\*

Allows the individual to connect to their smartphone via Bluetooth and access functions which include battery status, speed, range, drive mode. It also allows speed up to 8km/h, access to the easyNavi navigation system, remote control function and on Android devices, the ability to monitor incoming phone calls and SMS via the e-fix display screen.

*\*iPhone and Android compatible*


## Other

External charger socket incl. clip 1566177 £91

Brake extension 1467780 £29

Rear wheel extension available

Available wheel position:



1=rear

2=centre

3=front

*\*NB Always use anti-tippers (Alber or existing)*

## Comments (mounting etc.)