



Please refer to www.invacare.ie for local distributor.

Please fill in all specifications with each order.

Account code

Customer order number.....

Company name

Address

.....



e-fix E35 €6,045

Included as standard: 22" 24"

- 2 e-fix drive wheels
- LED display unit
- Mobility Plus package
- Range 9.9 miles
- Max. user weight 120kg
- 6.0 Ah battery (2kg) 7.5 Ah battery (2kg) €343
- Wheel weight 7.8kg
- Total weight 18.9kg

e-fix E36 €7,079

Included as standard: 24"

- 2 e-fix drive wheels
- LED display unit
- Mobility Plus package
- Range 16 miles
- Max. user weight 160kg
- 7.5 Ah battery (2kg) 6Ah battery (2kg) €343
- Wheel weight 8kg
- Total weight 19.3kg

Wheelchair details

Wheelchair manufacturer (e.g. Kuschall).....

Serial no..... Year.....

Model (e.g. Compact).....

Seat width 

Seat height (e.g. Hemi) 

Armrest & Brake type

Degree of wheel camber

(Drum brakes must be removed)

Brackets required for installation onto wheelchair

1457841 €133
Omit brackets

Fitting on wheelchair at Invacare factory & delivery* €86
Wheelchair order number

Collection & fitting to third party wheelchair & delivery** €86 + €110

Fitting by trained dealer** €0

Control Unit Mounting

Right side
Left side
Right & Left (pluggable) €90
Bracket on: Frame tube Side part

*Only when ordered with an Invacare wheelchair.
**Mounting kits can only be sent to selected trained dealers.

Do you want to use both your standard and e-fix wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for e-fix * / **
Retrofit, stub axle to 1/2" diameter €131

*If the manual wheels are used in combination with e-fix, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.
**If drum brakes are installed, these must be removed as they are not compatible with e-fix. Standard lever brakes are used instead.



Accessories

Anti-tippers with jack-up function (pair)* 1489214 €382

Control Unit Accessories

Swivel arm 1467901 €486

Protection Bumper 1565836 €71

Joystick Options

Option special 40mm/1.57" 1565669 €35

Option soft 70mm/2.76" 1566196 €35

Option T - handle 1565670 €119

**NB Always use anti-tippers (Alber or existing)*

Tray table incl. control unit bracket

Control unit centered SW 38-43cm 1566207 €341

Control unit centered SW 43-48cm 1566208 €341

Installation kit for control unit with existing tray table 1566206 €200

Extension for speed pre-selection

Set for speed pre-selection 1566197 €54

Extension for coupling ring

Set for one drive wheel (single) 1566179 €140

Set for two drive wheels (pair) 1566178 €278

Spoke Covers (pair)

22" wheel, decor 'Alber' 1566201 €85

22" wheel, transparent 1566203 €85

24" wheel, decor 'Alber' 1566200 €85

24" wheel, transparent 1566202 €85


Other

External charger socket incl. clip 1566177 €127

Brake extension 1467780 €44

Rear wheel extension available

Available wheel position:



1=rear

2=center

3=front

**NB Always use anti-tippers (Alber or existing)*

Attendant Control Bracket

Right 1457754 €169

Left 1457755 €169

Intuitive Attendant Control

Control unit intuitive attendant control 1566165 €62

Adapter cable intuitive attendant control 1566177 €127

Additional handle 1566166 €61

Opti-box 1566184 €464

Speed/Stop limit switch 1467698 €96

Battery pack

Spare battery 6.0 Ah 1565639 €1110

Spare battery 7.5 Ah 1566153 €1,463

Transport bag for battery 1566168 €100

Attachment set for transport bag 1566169 €29

Battery holder fixation

Fixed assembly set for cross brace 1566164 €236

Fixed assembly set for cross tube 1566163 €202

Mobility Plus Package*

Allows the individual to connect to their smartphone via Bluetooth and access functions which include battery status, speed, range, drive mode. It also allows speed up to 8km/h, remote control function and on Android devices, the ability to monitor incoming phone calls and SMS via the e-fix display screen.

**iPhone and Android compatible*

Comments (mounting etc.)

.....

.....

.....

.....

.....

.....

.....