



Invacare Ireland Ltd
 Unit 5 Seatown Business Park,
 Seatown Road, Swords, Co. Dublin
 K67 K271
 T: 01 810 7084 F: 01 810 7085
 E: ordersireland@invacare.com
 www.invacare.ie



Please fill in all specifications with each order.

Account code

Customer order number.....

Company name

Address

.....

twion M24 €6,282

Included as standard:

- 2 twion push rim 24" drive wheels
- Speed: 6km/h to 10km/h
- Weight: 6.3kg per wheel
- Built-in lithium-ion batteries in the wheel hubs
- Charger with magnetic connectors
- Mobility Plus Package

twion T24 €6

Standard options same as above

- | | | |
|---------|--|--------------------------|
| 1591624 | White cover, Tyre Marathon, Curve L push rim | <input type="checkbox"/> |
| 1591623 | White cover, Tyre Schwalbe - One, Curve L push rim | <input type="checkbox"/> |
| 1591622 | Black cover, Tyre Marathon, Curve L push rim | <input type="checkbox"/> |
| 1591620 | Black cover, Tyre Schwalbe - One, Curve L push rim | <input type="checkbox"/> |
| 1592494 | Black cover, Puncture-proof, Curve L push rim | <input type="checkbox"/> |
| 1592495 | White cover, Puncture-proof, Curve L push rim | <input type="checkbox"/> |

Wheelchair details

Wheelchair manufacturer (e.g. Kuschall)

Wheelchair serial number

Model (e.g. Compact)

Seat width 

Seat height (e.g. Hemi) 

Degree of wheel camber °

(Drum brakes must be removed)



Brackets required for installation onto wheelchair

- | | | |
|---------------------|--------------------------|------|
| Bracket as standard | <input type="checkbox"/> | STD |
| Spare parts set | <input type="checkbox"/> | €290 |
| Omit brackets | <input type="checkbox"/> | -€95 |

Fitting twion on wheelchair at Invacare factory & delivery* €82

Wheelchair order number

Collection & fitting to third party wheelchair & delivery €82 + €105

Fitting by trained dealer** €0

**Only when ordered with an Invacare wheelchair.*

***Mounting kits can only be sent to selected trained dealers.*

Do you want to use both your standard and twion wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for twion * / **

Retrofit, stub axle to 1/2" diameter €125



**If the manual wheels are used in combination with twion, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.*

***If drum brakes are installed, these must be removed as they are not compatible with e-fix. Standard lever brakes are used instead.*



Accessories

Anti-tippers with jack-up function (pair)* 1489214 €357

Push rims - M24 ONLY

Matt aluminium, coated black STD
 Stainless steel, brushed (+400g each) €0

**NB Always use anti-tippers (Alber or existing)*



Covers (pairs) - M24 ONLY

White covers STD
 Black covers €37



Tyres (pairs) - M24 ONLY

Schwalbe Marathon Plus Evolution STD
 Schwalbe One (-300g each) €59
 Solid infill black tyre 1/38 €59



Extra options

Transportation axle 1591362 €87
 Loop axle 1591500 €166
 Brake extension* 1467780 €32



**maybe required depending on wheelchair type*

Mobility Plus Package*


Comprises of two wheel hub outer rings which allows the individual to connect via Bluetooth and the licence code information. Functions now include: Cruise control function, battery status, speed and distance, ECO / SPORT mode, adjust speed up to 10km/h, BlueDrive remote control mode, trip profile record.

**iPhone and Android compatible*

Other

External charger socket incl. clip 1566177 €119
 Rear wheel extension available

Available wheel position:



1=rear
 2=centre
 3=front

**NB Always use anti-tippers (Alber or existing)*

Comments (mounting etc.).....

.....

.....

.....

.....

.....

.....

.....

.....