



Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus,
 Seatown Road, Swords, Co Dublin
 K67 K271
 T: 01 810 7084 F: 01 810 7085
 E: ordersireland@invacare.com
 www.invacare.ie



Please fill in all specifications with each order.

Account code

Customer order number.....

Company name

Address

.....



e-fix E35 €5,067

Included as standard:

- 2 e-fix drive wheels 22" 24"
- LED display unit
- Mobility Plus package
- Range 9.9 miles
- Max. user weight 120kg
- 6.0 Ah battery (2kg) 7.5 Ah battery (2kg) €288
- Wheel weight 7.8kg
- Total weight 18.9kg

e-fix E36 €5,937

Included as standard:

- 2 e-fix drive wheels 24"
- LED display unit
- Mobility Plus package
- Range 16 miles
- Max. user weight 160kg
- 7.5 Ah battery (2kg) 6Ah battery (2kg) -€288
- Wheel weight 8kg
- Total weight 19.3kg

Wheelchair details

Wheelchair manufacturer (e.g. Kuschall).....

Serial no..... Year.....

Model (e.g. Compact).....

Seat width 

Seat height (e.g. Hemi) 

Armrest & Brake type

Degree of wheel camber

(Drum brakes must be removed)

Brackets required for installation onto wheelchair

1457841 €110
 Omit brackets

Fitting e-fix on wheelchair at Invacare factory* €56
 Wheelchair order number

Fitting to third party wheelchair €56
 Fitting by trained dealer** €0

Control Unit Mounting

Right side
 Left side
 Right & Left (pluggable) €74
 Bracket on: Frame tube Side part

*Only when ordered with an Invacare wheelchair.
 **Mounting kits can only be sent to selected trained dealers.

Do you want to use both your standard and e-fix wheels in the future? If so, please order the retrofit stub axle.

Manual wheel axle retrofit for e-fix * / **

Retrofit, stub axle to 1/2" diameter €109

*If the manual wheels are used in combination with e-fix, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.

**If drum brakes are installed, these must be removed as they are not compatible with e-fix. Standard lever brakes are used instead.



Accessories

Anti-tippers with jack-up function (pair)* 1489214 €309

Control Unit Accessories

Swivel arm 1467901 €392

Protection Bumper 1565836 €57

Joystick Options

Option special 40mm/1.57" 1565669 €29

Option soft 70mm/2.76" 1566196 €29

Option T - handle 1565670 €96

**NB Always use anti-tippers (Alber or existing)*

Tray table incl. control unit bracket

Control unit centered SW 38-43cm 1566207 €275

Control unit centered SW 43-48cm 1566208 €275

Installation kit for control unit with existing tray table 1566206 €162

Extension for speed pre-selection

Set for speed pre-selection 1566197 €44

Extension for coupling ring

Set for one drive wheel (single) 1566179 €112

Set for two drive wheels (pair) 1566178 €225

Spoke Covers (pair)

22" wheel, decor 'Alber' 1566201 €69

22" wheel, transparent 1566203 €69

24" wheel, decor 'Alber' 1566200 €69

24" wheel, transparent 1566202 €69


Other

External charger socket incl. clip 1566177 €106

Brake extension 1467780 €35

Rear wheel extension available

Available wheel position:



1=rear

2=center

3=front

**NB Always use anti-tippers (Alber or existing)*

Attendant Control Bracket

Right 1457754 €137

Left 1457755 €137

Intuitive Attendant Control

Control unit intuitive attendant control 1566165 €502

Adapter cable intuitive attendant control 1566177 €103

Additional handle 1566166 €49

Opti-box 1566184 €375

Speed/Stop limit switch 1467698 €78

Battery pack

Spare battery 6.0 Ah 1565639 €853

Spare battery 7.5 Ah 1566153 €1,126

Transport bag for battery 1566168 €80

Attachment set for transport bag 1566169 €23

Battery holder fixation

Fixed assembly set for cross brace 1566164 €191

Fixed assembly set for cross tube 1566163 €164

Mobility Plus Package*

Allows the individual to connect to their smartphone via Bluetooth and access functions which include battery status, speed, range, drive mode. It also allows speed up to 8km/h, access to the easyNavi navigation system, remote control function and on Android devices, the ability to monitor incoming phone calls and SMS via the e-fix display screen.

**iPhone and Android compatible*

Comments (mounting etc.)

.....

.....

.....

.....

.....

.....

.....