



Invacare® Top End® Pro™-2 Wheelchair

Product Description:

The Pro™-2 Basketball Wheelchair is designed for aspiring athletes who want durability and quick, easy adjustments, all at a great price. Adjustments that can take up to four hours in other chairs, take only 15 minutes in the Pro 2 Basketball wheelchair. Plus, it is so easy to order that you only need to measure seat width. Now that is easy!

About New Smart Forms:

Top End is excited to introduce a new way to order Top End products. Our new Smart Forms provide an easy and electronic way to fill out order forms, save copies, and submit to them via e-mail. Please see below for the features and useful information.

How to use:

- Use any Adobe Acrobat application to electronically fill out the order form
- Use the default selections to order a product with the most common configuration
- Use the "**Reset**" button located in the forms footer to reset selections to their defaults
- Use the "**Clear All**" button located in the forms footer to clear selections, including defaults. This feature allows for a blank order form to be printed and completed by hand
- Use the "**E-mail to Top End order box**" button located on the last page of the form to e-mail the completed form directly to the topendorders@invacare.com e-mail inbox
- Forms can be completed electronically, printed, and e-mailed or faxed to topendorders@invacare.com or 727-522-1007 if access to Acrobat Reader for Windows is unavailable

(Acrobat Reader for Windows with java-script enabled is required to save and attach to e-mail)

Features:

- Tabulated options for an organized and standardized work-flow
- New more specific option codes which will appear on confirmations to improve customer / manufacturing communication
- Common special requests have been given option codes
- Built in exclusions to help eliminate errors
- Ability to save the form after it is filled out
- Ability to e-mail the form inside Acrobat Reader rather than saving a copy and attaching to an e-mail
- Ability to clear the entire form so it can be printed and completed by hand in the field when a computer is not available
- Ability to reset form to default selections



Invacare® Top End® Pro™-2 Wheelchair



Order Request Type:			
<input type="radio"/> Quote	<input type="radio"/> CAD/Quote	<input type="radio"/> Order	PO # <input style="width: 100px;" type="text"/>
Provider / Customer Information:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight:	
Contact:	Address Line 2:	Leg Length:	
Phone #:	City:	Hip Width:	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	
		Email: <input style="width: 100px;" type="text"/>	

Complete chair starts at: \$1,995.00

Base Model:		
Item #:	Description:	Price:
<input type="radio"/> 10014	Top End Pro-2 All Sport	\$1995
<input type="radio"/> 10026	Factory Assembly	\$100

Optional Wing: (See Definitions)		
Item #:	Description:	Price:
<input type="radio"/> 25074	Clamp-on Basketball Wing Set	\$300

Sizing:		
Item #:	Description:	Price:
Seat Width:		
<input type="radio"/> 15005	10" Seat Width - Gecko Green	\$0
<input type="radio"/> 15009	12" Seat Width - Safety Yellow	\$0
<input type="radio"/> 15013	14" Seat Width - White Snow	\$0
<input type="radio"/> 15017	16" Seat Width - Simple Red	\$0
<input type="radio"/> 15021	18" Seat Width - Simple Blue	\$0
Seat Depth: (See Definitions)		
<input type="radio"/> 15357	Adjustable Seat Depth	\$0
Seat Height:		
<input type="radio"/> 15096	Adjustable front seat height: 18" to 20" Adjustable rear seat height 15" to 20"	\$0
Back Height:		
<input type="radio"/> 15136	Adjustable back height 11" to 14"	\$0

Wheelbase: (Front to Rear Axle)		
Item #:	Description:	Price:
<input type="radio"/> 75068	17" Wheelbase	\$0

Wheel Camber Angle:		
Item #:	Description:	Price:
<input type="radio"/> 45008	15 Degree Camber	\$0

Rear Wheel Options:		
Item #:	Description:	Price:
Size and Style:		
<input type="radio"/> 50002	25" Sun High Performance wheel	\$0
Handrim Size/Material:		
<input type="radio"/> 50301	25" Aluminum Anodized handrim	\$0
Handrim Style/Mounting:		
<input type="radio"/> 50308	Standard Mount handrim	\$0
Tires:		
<input type="radio"/> 50201	25" Primo, Grey tire	\$0

Center of Gravity: (front of backrest to center of axle)		
Item #:	Description:	Price:
<input type="radio"/> 45073	Adjustable 2.5" to 5" center of gravity	\$0

Back Type:		
Item #:	Description:	Price:
<input type="radio"/> 15159	Angle Adjustable, Non-Folding backrest	\$0

Side Guards:		
Item #:	Description:	Price:
<input type="radio"/> 20014	Screw-on, Aluminum sideguards	\$0

Anti-Tip:		
Item #:	Description:	Price:
<input type="radio"/> 20030	Single telescoping, removable, height adjustable anti-tip	\$0

Footrest:		
Item #:	Description:	Price:
<input type="radio"/> 60003	Clamp-on tubular, height adjustable footrest	\$0
<input type="radio"/> 60030	Click Strap, Single Pad for footrest	\$75



Invacare® Top End® Pro™-2 Wheelchair

Axles:		
Item #:	Description:	Price:
<input type="radio"/> 50400	Quick Release steel 1/2" axles	\$0

Upholstery:		
Item #:	Description:	Price:
Material/Color:		
<input type="radio"/> 30001	Nylon, Black upholstery	\$0
Seat Style:		
<input type="radio"/> 30019	Screw-On seat upholstery	\$0
Back Style:		
<input type="radio"/> 30012	Adjustable Tension back upholstery	\$0
Options:		
<input type="radio"/> 25063	2" Soft Foam Seat Cushion	\$0
<input type="radio"/> 25060	Hook and Loop Positioning Strap	\$0
<input type="radio"/> 30005	2" Soft Cushion	\$75
<input type="radio"/> 30006	3" Soft Cushion	\$75
<input type="radio"/> 30008	2" Hard Cushion	\$125
<input type="radio"/> 30009	3" Hard Cushion	\$125
<input type="radio"/> 30010	4" Hard Cushion	\$125

Options and Accessories:		
Item #:	Description:	Price:
Safety:		
<input type="radio"/> 25066	Calf Strap	\$0
<input type="radio"/> 25065	Rear Frame Protector	\$0
Positioning:		
<input type="radio"/> 25001	Click Strap set, 8" Small	\$75
<input type="radio"/> 25002	Click Strap set, 10" Medium	\$80
<input type="radio"/> 25003	Click Strap set, 12" Large	\$86
<input type="radio"/> 25012	Metal O-Ring Positioning Strap	\$60
<input type="radio"/> 25013	Hook and Loop Positioning Strap	\$40
Accessories:		
<input type="radio"/> 25011	Basketball Carrying Bag	\$75
Apparel:		
<input type="radio"/> 25009	T-Shirt, Large	\$10
<input type="radio"/> 25010	T-Shirt, Extra Large	\$10

Product Illustration:
Pro-2:



Definitions:	
Factory Assembly:	Attachment of seat frame to base frame, installation of anti-tip and footrests. Only assembly required upon arrival is attachment of rear wheels.
Adjustable Seat Depth:	All Invacare Top End Pro-2 All Sport Wheelchairs have adjustable seat depths. For seat widths of 10" and 12" the seat depth will be adjustable to depths of 10" and 12". For seat widths of 14", 16", and 18" the seat depth will be adjustable to depths of 13", 15" and 17".
Clamp-on Basketball Wing:	The Clamp-on Basketball Wing set components extend the protective bumper around the front of the chair. This is typically desired for use in Wheelchair Basketball. This option is not assembled on the chair and requires minimal assembly upon arrival.

Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	