


# Aquatec® | **Bathing Transfer aids**



Aquatec **Trans**

- Compatible with the *Aquatec* range of bath lifters
- Large comfortable seating surface
- Safe: turntable is only unlocked by the weight of the user
- Seat surface turns and slides easily
- Soft, skin-friendly cover: easy-care, wipe clean
- Also available in Grey

**3 year warranty**



Aquatec **Disk**

- Large comfortable seating surface
- Turns very easily
- Soft, skin-friendly cover: easy-care, wipes clean
- The *Aquatec Disk* can be used with the *Aquatec Orca* bath lifter range.
- Also available in Grey

**3 year warranty**

**AQUATEC®**



**Yes, you can.®**



Aquatec **Step**

- Stable frame
- Feet with non-slip rubber caps
- Non-slip, removable cover, machine-washable

**3 year warranty**

**150kg/24 stone weight limit**

## Technical data

Order numbers	Aquatec <b>Trans</b>	Aquatec <b>Disk</b>	Aquatec <b>Step</b>
Turning and transfer aid, blue	<b>4.03.002</b>		
Turning and transfer aid, grey	<b>4.03.008</b>		
Rotating disk, blue		<b>4.02.002</b>	
Rotating disk, grey		<b>4.02.008</b>	
Step, blue			<b>3.00.002</b>
Dimensions			
Diameter	365 mm	365 mm	–
Length	530 mm	–	400 mm
Width	–	–	300 mm
Height	25 mm	25 mm	165 mm
More info			
Weight	2.2 kg	1.6 kg	1.5 kg
Load capacity	130 kg (21 stone)	130 kg (21 stone)	150 kg (24 stone)
Colours	Blue, grey	Blue, grey	Blue, white
Warranty	3 Years	3 Years	3 Years