



aquatec®

 **Pure Cala**

RAISED TOILET SEAT

Discreet design for
everyday routines



Introducing the Aquatec **Pure Cala**®

The Aquatec **Pure Cala** is a new raised toilet seat, designed to support with personal care. Simple to use and built for lasting stability, it makes daily care routines feel natural, discreet, and reassuring.



Modern design

Pure Cala features a simple, modern design that is perfect for any home.



Quality assured

All products from Aquatec are designed, made and tested to the highest standards.



Simple and discreet

Designed to support everyday independence with ease and dignity.



Trusted by experts

Aquatec is a trusted brand used for decades by professionals.

Blends
effortlessly into
everyday life.





Everyday independence

The first in Aquatec's **Pure** range, **Pure Cala** is a secure, simple-to-fit 10cm raised toilet seat designed for comfort and everyday support.

 **Simple to set-up**
Tool-free set up makes installation quick and easy.

 **Easy-grip fixing**
The contoured, textured locking wheel makes mounting easier for those with limited dexterity.



[Click here to watch the video >>](#)



 **Optional lid**
Choose Pure Cala with or without lid.

 **Optional soft seat**
For increased comfort.

 **Optimal height**
Increased 10cm height for easier sit-to-stand movement.

 **Secure and supportive**
Built-in locking system with anti-slip rubber inserts provides stability.

 **Shaped for comfort**
Gently contoured seat gives optimal support and comfort.

 **Universal fit**
Compatible with a wide range of toilet bowls and bathrooms.



Front and back aperture



Stain resistant and easy to clean



Discreet unobtrusive design



Designed to fit into any home

The Aquatec **Pure Cala** easily and effortlessly fits into home or healthcare settings. It's purposely designed to break associated stigma and promotes everyday independence.

aquatec®

Aquatec continues to be a trusted name in personal care solutions. For decades, safety, simplicity and enduring quality have become its hallmark.



Toileting



Showering



Bathing



Safety

Technical specifications



Width of opening

20.5 cm



Total width

36.5 cm



Depth of opening

25.5 cm



Total depth

40.5 cm



Raised seat height from toilet bowl

10 cm



Total product height

15.5 cm
includes the
assembled
locking-system



Weight

1.1kg (with lid)
0.7kg (without lid)



Load capacity

250kg



White

Colour



5 yr

Service Life

Made in Europe



We're committed to using more eco-conscious materials, reducing plastic where possible, and exploring smarter packaging solutions.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

© 2025 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. **Aquatec Pure Cala_ EU-10/25**

