



aquatec[®]

Pure Range

Discreet design for
everyday routines





Designed to fit into any home

The Aquatec **Pure Range** brings together thoughtfully designed products that support personal care in a way that feels natural, discreet, and dependable. Created for home and healthcare settings and designed to break associated stigma and promote everyday independence.

[Click here to discover the range >>](#)

aquatec®

Aquatec continues to be a trusted name in personal care solutions. For decades, safety, simplicity and enduring quality have become its hallmark.



Toileting



Showering



Bathing



Safety







Pure Sol, Bay & Siro Shower stools

The Aquatec **Pure Sol, Bay & Siro** shower stools are designed to increase comfort, safety, and independence.

Stylish and compact, the stools are designed to fit neatly into corners or small spaces.



Pure Sol
Round shower stool



Pure Bay
Rectangular shower stool

Key features

-  **Lightweight & portable**
Allows for easy lifting and movement.
-  **Adjusts in seat height**
Personalised positioning helps relieve, back, hip and knee strain.
-  **Enhanced comfort**
Smooth seat edges helps reduce pressure and improve comfort.
-  **Anodized aluminium legs**
Offer stability, durability, and corrosion resistance.
-  **Stable on uneven surfaces**
Usable with legs set at different height* and fitted with non-slip rubber tips for secure use.



Pure Siro
Triangular shower stool

Easy assembly and care



Tool-free assembly and height adjustment
Simplifies setup and use.



Easy-to-handle locking system
For secure, quick assembly and disassembly.



Maintenance-free and easy to clean
Suitable for daily use in home or healthcare environments.

*For a horizontal seating position only.



Pure Isla Shower chairs

The Aquatec **Pure Isla** shower chairs are designed to provide a safe and dignified seating solution for personal care.

Design inspired by balance and movement, combining a sense of calm with strength and stability.



Pure Isla
Seat only



Pure Isla
Seat with
armrests



Pure Isla
Seat with
backrest

Easy assembly and care



Tool-free assembly and height adjustment
Simplifies setup and use.



Easy-to-handle locking system
For secure, quick assembly and disassembly.



Maintenance-free and easy to clean
Suitable for daily use in home or healthcare environments.

Key features



-  **Lightweight & portable**
Allows for easy lifting and movement.
-  **Adjusts in seat height and anterior tilt**
Personalised positioning helps relieve, back, hip and knee strain.
-  **Enhanced comfort**
Smooth seat edges helps reduce pressure and improve comfort.
-  **U-shaped seat**
Enhances accessibility. Can be positioned front or rear, depending on personal needs.
-  **Anodized aluminium legs**
Offer stability, durability, and corrosion resistance.
-  **Stable on uneven surfaces**
Usable with legs set at different height and fitted with non-slip rubber tips for secure use.

Pure Isla
Seat with armrests
and backrest



Pure Azure Commode

The Aquatec **Pure Azure** commode has a 3-in-1 design, for use as a shower chair, toilet frame or commode. Ensuring comfort and security as personal needs change over time.



Pure Azure
Configured as shower chair




Pure Azure
Configured as toilet frame

Key features




-  **Lightweight & portable**
Allows for easy lifting and movement.
-  **Adjusts in seat height and anterior tilt**
Personalised positioning helps relieve, back, hip and knee strain.
-  **Enhanced comfort**
Smooth seat edges helps reduce pressure and improve comfort.
-  **Wide hygiene recess**
For convenient use as an alternative to the toilet bowl.
-  **Versatile toilet pan holder**
Allows the pan to be removed from front and back.
-  **Anodized aluminium legs**
Offer stability, durability, and corrosion resistance.
-  **Stable on uneven surfaces**
Usable with legs set at different height and fitted with non-slip rubber tips for secure use.

Easy assembly and care

 **Tool-free assembly and height adjustment**
Simplifies setup and use.

 **Easy-to-handle locking system**
For secure, quick assembly and disassembly.

 **Maintenance-free and easy to clean**
Suitable for daily use in home or healthcare environments.

Pure Azure
Configured as commode



🔧 Pure Cala Raised toilet seat

The Aquatec **Pure Cala** raised toilet seat is designed to support with personal care. Simple to use and built for lasting stability, it makes daily care routines feel natural, discreet, and reassuring.



Simple to set-up
Tool-free set up makes installation quick and easy.



Easy-grip fixing
The contoured, textured locking wheel makes mounting easier for those with limited dexterity.



[Click here to watch the video >>](#)



Key features



Optional lid
Choose Pure Cala with or without lid.



Optional soft seat
For increased comfort.



Optimal height
Increased 10cm height for easier sit-to-stand movement.



Secure and supportive
Built-in locking system with anti-slip rubber inserts provides stability.



Shaped for comfort
Gently contoured seat gives optimal support and comfort.



Universal fit
Compatible with a wide range of toilet bowls and bathrooms.



Front and back aperture



Stain resistant and easy to clean



Discreet unobtrusive design



Technical specifications

Pure Sol, Bay & Siro Shower Stools



Seat width



Seat depth



Seat height

Pure Sol	390 mm	390 mm	420 – 570 mm
Pure Bay	435 mm	385 mm	420 – 570 mm
Pure Siro	455mm (at widest point)	440 mm	420 – 570 mm



Weight



Load capacity

Pure Sol	2.1 kg	150 kg
Pure Bay	2.3 kg	180 kg
Pure Siro	2.2 kg	130 kg



Colour



Service Life

Made in Europe

Pure Isla Shower Chair



Seat width



Seat depth



Seat height



Total height

500 mm	420 mm	420 – 570 mm	780 – 930 mm
--------	--------	--------------	--------------



Backrest height



Backrest width



Armrest height



Armrest width

360 mm	435 mm	210 mm	250 mm
--------	--------	--------	--------



Total width



Weight



Load capacity

615 mm	2.8 kg (seat) 3.5 kg (seat+armrest) 3.7 kg (seat+backrest) 4.4 kg (complete chair)	160 kg
--------	---	--------



Colour



Service Life

Made in Europe



Technical specifications

Pure Azure Commode



Seat width

500 mm



Seat depth

420 mm



Seat height

420 – 570 mm



Total height

780 – 930 mm



Backrest height

360 mm



Backrest width

435 mm



Armrest height

210 mm



Armrest width

250 mm



Total width

615 mm



Hygiene recess width

195 mm



Hygiene recess depth

235 mm



Weight

5.1 kg (with pan)



White

Colour



Service Life

Made in Europe



Load capacity

160 kg

Pure Cala Raised Toilet Seat



Width of opening

205 mm



Total width

365 mm




Depth of opening

255 mm



Total depth

405 mm



Raised seat height from toilet bowl

100 mm



Total product height

155 mm
includes the assembled locking-system



Weight

1.1kg (with lid)
0.7kg (without lid)



Load capacity

250kg



White

Colour



Service Life

Made in Europe



We're committed to using more eco-conscious materials, reducing plastic where possible, and exploring smarter packaging solutions.

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland