



aquatec®

Pure Range

Discreet design for everyday routines





aquatec



Designed to fit into any home

The Aquatec **Pure Range** brings together thoughtfully designed products that support personal care in a way that feels natural, discreet, and dependable. Created for home and healthcare settings and designed to break associated stigma and promote everyday independence.



Scan here to
discover the
range

aquatec®

Aquatec continues to be a trusted name in personal care solutions. For decades, safety, simplicity and enduring quality have become its hallmark.



Toileting



Showering



Bathing



Safety

Pure Sol, Bay & Siro Shower stools

The Aquatec **Pure Sol, Bay & Siro** shower stools are designed to increase comfort, safety, and independence.

Stylish and compact, the stools are designed to fit neatly into corners or small spaces.



Pure Sol
Round
shower stool



Pure Bay
Rectangular
shower stool

Easy assembly and care



**Tool-free assembly
and height adjustment**
Simplifies setup
and use.



**Easy-to-handle
locking system**
For secure, quick
assembly and
disassembly.



**Maintenance-free
and easy to clean**
Suitable for daily use
in home or healthcare
environments.




Key features



Lightweight & portable

Allows for easy lifting and movement.




Adjusts in seat height

Personalised positioning helps relieve, back, hip and knee strain.



Enhanced comfort

Smooth seat edges helps reduce pressure and improve comfort.



Anodized aluminium legs

Offer stability, durability, and corrosion resistance.



Stable on uneven surfaces

Usable with legs set at different height and fitted with non-slip rubber tips for secure use.

Pure Siro
Triangular
shower stool

Pure Isla Shower chairs

The Aquatec **Pure Isla** shower chairs are designed to provide a safe and dignified seating solution for personal care.

Design inspired by balance and movement, combining a sense of calm with strength and stability.



Pure Isla
Seat only



Pure Isla
Seat with
armrests



Pure Isla
Seat with
backrest

Easy assembly and care



**Tool-free assembly
and height adjustment**
Simplifies setup
and use.



**Easy-to-handle
locking system**
For secure, quick
assembly and
disassembly.



**Maintenance-free
and easy to clean**
Suitable for daily use
in home or healthcare
environments.



Key features



Lightweight & portable

Allows for easy lifting and movement.



Adjusts in seat height and anterior tilt

Personalised positioning helps relieve, back, hip and knee strain.



Enhanced comfort

Smooth seat edges helps reduce pressure and improve comfort.



U-shaped seat

Enhances accessibility. Can be positioned front or rear, depending on personal needs.



Anodized aluminium legs

Offer stability, durability, and corrosion resistance.



Stable on uneven surfaces

Usable with legs set at different height and fitted with non-slip rubber tips for secure use.

Pure Isla
Seat with armrests and backrest

Pure Azure Commode

The Aquatec **Pure Azure** commode has a 3-in-1 design, for use as a shower chair, toilet frame or commode.

Ensuring comfort and security as personal needs change over time.



Pure Azure
Configured as
shower chair



Pure Azure
Configured as
toilet frame

Easy assembly and care



**Tool-free assembly
and height adjustment**
Simplifies setup
and use.



**Easy-to-handle
locking system**
For secure, quick
assembly and
disassembly.



**Maintenance-free
and easy to clean**
Suitable for daily use
in home or healthcare
environments.



Key features



Lightweight & portable

Allows for easy lifting and movement.



Adjusts in seat height and anterior tilt

Personalised positioning helps relieve, back, hip and knee strain.



Enhanced comfort

Smooth seat edges helps reduce pressure and improve comfort.



Wide hygiene recess

For convenient use as an alternative to the toilet bowl.



Versatile toilet pan holder

Allows the pan to be removed from front and back.



Anodized aluminium legs

Offer stability, durability, and corrosion resistance.



Stable on uneven surfaces

Usable with legs set at different height and fitted with non-slip rubber tips for secure use.

Pure Azure
Configured
as commode

Pure Cala Raised toilet seat

The Aquatec **Pure Cala** raised toilet seat is designed to support with personal care. Simple to use and built for lasting stability, it makes daily care routines feel natural, discreet, and reassuring.



Simple to set-up

Tool-free set up makes installation quick and easy.



Easy-grip fixing

The contoured, textured locking wheel makes mounting easier for those with limited dexterity.



[Click here to watch the video >>](#)



Key features



Optional lid

Choose Pure Cala with or without lid.



Optional soft seat

For increased comfort.



Optimal height

Increased 10cm height for easier sit-to-stand movement.



Secure and supportive

Built-in locking system with anti-slip rubber inserts provides stability.



Shaped for comfort

Gently contoured seat gives optimal support and comfort.



Universal fit

Compatible with a wide range of toilet bowls and bathrooms.



Front and back aperture



Stain resistant and easy to clean



Discreet unobtrusive design

Technical specifications

Pure Sol, Bay & Siro Shower Stools



Seat width

Pure Sol	390 mm
Pure Bay	435 mm
Pure Siro	455mm (at widest point)



Seat depth

Pure Sol	390 mm
Pure Bay	385 mm
Pure Siro	440 mm



Seat height

Pure Sol	420 – 570 mm
Pure Bay	420 – 570 mm
Pure Siro	420 – 570 mm



Weight

Pure Sol	2.1 kg
Pure Bay	2.3 kg
Pure Siro	2.2 kg



Load capacity

Pure Sol	150 kg
Pure Bay	180 kg
Pure Siro	130 kg

White

Colour

5 yr

Service Life

Made in Europe



Pure Isla Shower Chair



Seat width

500 mm



Seat depth

420 mm



Seat height

420 – 570 mm



Total height

780 – 930 mm



Backrest height

360 mm



Backrest width

435 mm



Armrest height

210 mm



Armrest width

250 mm



Total width

615 mm



Weight

2.8 kg (seat)
3.5 kg (seat+armrest)
3.7 kg (seat+backrest)
4.4 kg (complete chair)



Load capacity

160 kg

White

Colour

5 yr

Service Life

**Made in
Europe**

Technical specifications

Pure Azure Commode



Seat width

500 mm



Seat depth

420 mm



Seat height

420 – 570 mm



Total height

780 – 930 mm



Backrest height

360 mm



Backrest width

435 mm



Armrest height

210 mm



Armrest width

250 mm



Total width

615 mm



Hygiene recess width

195 mm




Hygiene recess depth

235 mm



Weight

5.1 kg (with pan)



White

Colour



Service Life

Made in Europe



Load capacity

160 kg



Pure Cala Raised Toilet Seat



Width of opening

205 mm



Total width

365 mm



Depth of opening

255 mm



Total depth

405 mm



Raised seat height from toilet bowl

100 mm



Total product height

155 mm
includes the assembled locking-system



Weight

1.1kg (with lid)
0.7kg (without lid)



Load capacity

250kg

White

Colour

5 yr

Service Life

Made in Europe



We're committed to using more eco-conscious materials, reducing plastic where possible, and exploring smarter packaging solutions.

T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare UK & Ireland
Unit 4 Pencoed
Technology Park
Bridgend CF35 5AQ