



Certificate in terms of transportation

We hereby confirm, that the non-spillable valve-regulated lead-acid gel batteries with the capacities:

52 Ah 12 V (TNE12-52) TN Power
60 Ah 12 V MK Powered
74 Ah 12 V MK Powered

used in the product

Power Wheelchair Invacare® Aviva

are classified in terms of transportation as follows:

IATA/ICAO: NOT RESTRICTED per Special Provision A67
DOT: NOT RESTRICTED per 49 CFR 173.159
IMDG: NOT RESTRICTED per Special Provision 238

Isny, 14. Jul. 2020

Invacare GmbH

A handwritten signature in blue ink, appearing to read "M. Schreiber", is written over a light blue horizontal line.

i. V. Michael Schreiber
Quality Manager