



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Power Retail Prescription Form

SEPTEMBER 2024  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX20 Modulite

### Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 15°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,321**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
JBC0000	AVIVA RX 20 Base Model Max user weight up to 136 kg	£5,321
JBX1300	Special Extended lead time: + 10 days	POA
<b>CHASSIS WIDTH*</b>		
JBL2810	Chassis Narrow (560 mm/22")	◆
<b>SEATING TYPE*</b>		
JBC2100	Invacare® Modulite Seat Unit	◆

<b>SEAT FRAME*</b>			
JBC2111		Telescopic Seat Frame and Plate with Side Rails	◆
JBC2112		Telescopic Seat Frame + Sling Seat with Side Rails	○
<b>SEAT WIDTH*</b>			
JBC2121		380 mm Seat Plate (NC 2174, 2175, 2176)	○
JBC2122		430 mm Seat Plate	○
JBC2123		480 mm Seat Plate (NC 2704, 2706)	○
JBC2124		530 mm Seat Plate (NC 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2106)	○
<b>Body Fit Kit</b>			
JBC2106		Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side) (NC 2192, 2194)	○
<b>SEAT DEPTH*</b>			
JBC2132		Seat Depth 410 mm (NC 2197, 2198, 2701, 2702, 2703, 2704, 2705, 2706)	○
JBC2134		Seat Depth 460 mm	○
JBC2136		Seat Depth 510 mm	○
<b>SEAT CUSHION*</b>			
JBC2142		Matrx Flo-tech Cushion Flo-tech Lite	£118
JBC2143		Matrx Flo-tech Contour Visco	£165
JBC2147		Matrx Libra Cushion	£389
JBC2144		Matrx Vi Cushion	£360
JBC2140		Without seat cushion	○
<b>BACKREST SUPPORT*</b>			
JBC2163		Backposts short (480 mm) (NC 2173)	○
JBC2162		Backposts tall (540 mm)	○
<b>PUSH HANDLES FOR BACKPOSTS</b>			
JBC2108		Backposts without push handles (NC 2219, 2220, 2227, 2228, 2245, 2291, 2292)	○
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBC2171		Tension Adjustable Upholstery (NC 2230, 2247)	◆
JBC2179		Flex3 Tension Adjustable (NC 2192, 2193, 2194, 2199, 2155, 2230, 2247)	£153
JBC2172		Flex3 adjustable Backrest Plates without Backrest Cushion (NC 2128)	£250
JBC2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion (NC 2128)	£406
JBC2174		Matrx Posture Back, height 400 mm (up to 136 kg max user weight)	£558
JBC2175		Matrx Posture Back, height 500 mm (up to 136 kg max user weight)	£558
JBC2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£558
JBC2177		Matrx Elite TR Back, height 500 mm	£558
JBC2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£648
JBC2170		Without Backrest Upholstery (Back Posts/Back Frame only) (NC 2324, 2325, 2327, 2230, 2247)	○
<b>UPHOLSTERY TYPE*</b>			
JBC2151		Black Fabric Upholstery	◆
JBC2152		Black multi stretch PU coated fabric Upholstery	○
JBC2155		Black Spacer Fabric Upholstery (NC 2171)	○
JBX3002		No Upholstry Colour required - No backrest upholstery (NC 2179, 2173, 2174, 2175, 2176, 2177, 2178)	○

<b>ARMREST (Pair)*</b>			
JBC2191		T-Arm detachable, Standard	◆
JBC2181		T-Arm detachable, Extended height	£41
JBC2198		Padded Armrest with padded side guard (WATER FALL TYPE) (NC 2195, 2196, 2197, 2199)	£109
JBC2183		Inrigged Armrests Reduces SW by 25 mm each side	○
JBC2193		Flip-up Armrest Up to 150 kg max user weight (NC 1445, 2179)	£65
JBC2192		Following Armrest Height adjustable from 230 - 300 mm (NC 2179, 2106)	£185
JBC2194		Following Armrest Height adjustable from 300 - 360 mm (NC 2179, 2106)	£185
<b>ARMPAD (Pair)*</b>			
JBC2195		Short Armpad (325 mm)	◆
JBC2196		Long Armpad (375 mm) (NC 2198)	○
JBC2197		Long Ergonomic Armpad (397 mm) (NC 2198, 2230, 2247)	○
JBC2199		Soft Armpad	£96
JBX3000		Please select with Option 2198 - Padded Armrest with padded side guard (WATER FALL TYPE) (NC 2181, 2183, 2193, 2192)	○
<b>HIP SUPPORTS LH</b>			
JBC2184		Hip Support LH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
JBC2186		Hip Support LH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
JBC2188		Hip Support LH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
<b>HIP SUPPORTS RH</b>			
JBC2185		Hip Support RH, pad size 140 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
JBC2187		Hip Support RH, pad size 230 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
JBC2189		Hip Support RH, pad size 330 x 100 mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
JBC2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JBC2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
JBC1350		Neckrest depth adjustable (NC 2247)	£81
JBC1432		Headrest sideways adjustable (NC 2247)	£226
JBC1427		Headrest with cheek support (NC 2247)	£226
JBC2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (NC 2247, 2170)	£361
JBC2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (NC 2247, 2170)	£361
JBC2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard (NC 2247, 2170)	£489
JBC1496		Non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2171, 2174, 2175, 2176, 2177, 2178, 1288, 1289)	○
JBC1495		Non Invacare Headrest	○
<b>LATERAL SUPPORTS</b>			
JBC1445		Lateral Support with extended arm, multi-adjustable pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£168

<b>SEAT LIFTER / SEAT TILT*</b>			
JBC1091		Manually settable Seat Tilt 0 - 15°	○
JBC1092		Powered "fixed pivot" Tilt 20° (NC 2201)	£404
JBC1093		Powered Tilt (0° - 30°) with CoG shift (NC 2201)	£1,073
JBC2052		Electric Tilt (0°-50°) with COG shift (NC 2201, 2171)	£1,322
JBC1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs must be selected (NC 2201, 2220)	£1,893
JBC2053		Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift, 6 point Tie downs must be selected (NC 2201, 2220, 1910, 1912, 2171)	£2,359
<b>SEAT TO GROUND HEIGHT*</b>			
JBC2857		Seat-to-Floor Height 435 mm	○
JBC2858		Seat-to-Floor Height 460 mm	○
JBC2859		Seat-to-Floor Height 485 mm	○
<b>BACKREST RECLINER*</b>			
JBC0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBC0977		Powered Recliner 30°	£505
<b>CENTRE MOUNTED LEGREST*</b>			
JBC0556		Centre mount legrest, angle settable Length: 280 - 385 mm (NC 0517, 1244, 1416)	£373
JBC0589		LNx Power centre mount legrest, length compensating Length: 340 - 410 mm (NC 2201, 1244, 1416)	£1,640
JBC2855		Without Legrests	○
<b>CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*</b>			
JBC2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£103
JBC0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBC0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	£135
JBC0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBX3001		Please select if Centre mount footplate is not required. (NC 0589)	○
<b>SIDE MOUNTED HANGERS LH</b>			
JBC0555		Amputee Support LH	£194
<b>SIDE MOUNTED HANGERS RH</b>			
JBC0554		Amputee Support RH	£194
<b>FOOTREST HANGER TYPE*</b>			
JBC0526		Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
JBC0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools	£120
JBC0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad	£234
JBC0531		Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad (NC 2201)	£669
JBC0501		Without Hanger pair (NC 0504, 0513)	○
<b>FOOTPLATE TYPE*</b>			
JBC0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBC0504		Angle Adjustable Footplate pair (NC 0526)	£32
JBC0513		Angle and depth adjustable Footplate pair	£119
JBX3003		No Side Mounted Footplate required (NC 0503, 0533, 0531)	○

<b>CONTROLS*</b>			
JBC2201		Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module	(NC 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908) ◆
JBC2203		Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908) ○
JBC2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258) £151
JBC2205		Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets	(NC 2220) £341
JBC2206		Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233	(NC 2220, 2270, 1709, 1710) £776
<b>REM400 TOGGLE OPTION</b>			
JBC2270		Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)	(NC 2203, 2204, 2206) £316
<b>CONTROL BRACKETRY*</b>			
JBC1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBC1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBC1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) £132
JBC1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) £132
JBC1709		Matrx Resolve Swing in/out bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 1404) £219
JBC1710		Matrx Resolve Swing in/out bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 1405) £219
JBC2209		Swing away RH, multi-adjustable	(NC 2203, 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404) ○
JBC2210		Swing away LH, multi-adjustable	(NC 2203, 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405) ○
JBC2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) £498
JBC2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) £498
JBC2227		Bracket RH, mounted for Attendant	(NC 1287, 2220, 2128, 2108, 2219, 1221, 1312) £64
JBC2228		Bracket LH, mounted for Attendant	(NC 1286, 2220, 2128, 2108, 2219, 1222, 1363) £64
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA
<b>DUAL CONTROL</b>			
JBC2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228) £395
JBC2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 2108, 1221, 1222, 2227, 2228) £395
<b>INVACARE LINX IDC DUAL CONTROL</b>			
JBC2245		Invacare LiNX IDC Dual control	(NC 2108, 2259, 2260) £751
<b>DUAL CONTROL BRACKETRY</b>			
JBC2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBC2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
<b>INVACARE LINX IDC DUAL CONTROL BRACKETRY</b>			
JBC1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBC1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○
<b>LINX COMPACT REMOTE</b>			
JBC2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204) £718
JBC2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204) £960

LINX COMPACT REMOTE BRACKET				
JBC2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404)	○
JBC2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405)	○
JBC2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£498
JBC2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£498
JBX2215		Individual mounting option per request (Special) Requires additional information/specification		POA
CHIN CONTROL				
UNI9250		Mo-vis Micro joystick		POA
UNI9251		Mo-vis Multi joystick		POA
CHIN CONTROL BRACKETRY MANUAL				
JBC1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 1496)	○
JBC1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 1496)	○
CHIN CONTROL BRACKETRY POWERED				
JBC1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 2291, 1312, 1415)	£2,183
JBC1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 2292, 1363, 1339)	£2,183
LINX HEAD CONTROL				
JBC2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£3,719
JBC2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£4,355
LINX FOUR -SWITCH PROXIMITY ARRAY				
JBC2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128)	£5,760
JBC2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£498
JBC2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£498
JBX2269		Individual mounting option per request (Special) Requires additional information/specification		POA
LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL				
JBC2249		Wireless Triple Switch Receiver (ASL557-3)		£341
JBC2250		Wireless Mouse Emulator (ASL558)		£817
LINX MICRO EXTREMITY CONTROL BRACKET				
JBC2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128)	£3,612
LINX HEAVY DUTY REMOTE				
JBC2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128)	£2,013
LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET				
JBC2285		Tray Mount RH (100x150 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	○
JBC2286		Tray Mount LH (100x150 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	○
JBC2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£498
JBC2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£498
JBC2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2209, 2213, 2217, 1404, 1405, 1415)	£498
JBC2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2210, 2214, 2218, 1404, 1405, 1339)	£498
JBX2289		Individual mounting option per request (Special) Requires additional information/specification		POA



<b>LINX 10-WAY SWITCH</b>			
JBC2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£505
<b>LINX 10-WAY SWITCH BRACKET</b>			
JBC2217		Bracket RH, mounted for user (NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404)	○
JBC2218		Bracket LH, mounted for user (NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405)	○
JBC2291		Bracket RH, mounted for attendant (NC 2259, 1221, 1287, 2108)	£64
JBC2292		Bracket LH, mounted for attendant (NC 2260, 1222, 1288, 2108)	£64
<b>LINX USB CHARGER PORT</b>			
JBC2253		USB Charger Port DLX-USB02-A (NC 2702, 2703, 2704, 2705, 2706)	£165
<b>LINX USB CHARGER PORT BRACKET</b>			
JBC2281		RH Mount	○
JBC2282		LH Mount	○
<b>REMOTE STOP SWITCH</b>			
JBC2258		Remote Stop Switch (ASL504) (NC 2203, 2204)	£1,050
<b>LINX INPUT MODULE</b>			
JBC2273		IN200 Module (Jack port for external power button )	£242
UNI9252		IN500 Module (for DB9 Drive input) (NC 2201, 2203, 2203)	£785
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
JBC2300		To be enclosed only	○
JBX2301		To be connected and programmed (Special)	POA
<b>BUDDY BUTTONS</b>			
JBC1685		Red Buddy Button	£88
JBC1686		Yellow Buddy Button	£88
JBC1687		Green Buddy Button	£88
JBC1690		Blue Buddy Button	£88
<b>EGG SWITCH</b>			
JBC2254		Egg Switch, black	£88
<b>PIKO BUTTONS</b>			
JBC1711		Red Piko Button	£158
JBC1712		Yellow Piko Button	£158
JBC1713		Green Piko Button	£158
JBC1714		Orange Piko Button	£158
JBC1715		Blue Piko Button	£158
JBC1716		Black Piko Button	£158
<b>BODYPOINT JOYSTICK KNOBS</b>			
JBC2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£63
JBC2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£63
JBC2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£63
JBC2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£63
<b>5:1 SWITCH ADAPTOR</b>			
JBC1681		5:1 Switch Adapter ( 9-pin sub D to 5 jack ports)	£121

<b>SHROUD COLOUR*</b>			
JBC2871		Mid range - Black/Anthracite Shroud	○
JBC2872		Mid range - Silver Shroud	○
JBC2873		Mid range - Bengal Red Shroud	○
JBC2874		Mid range - Mercury Blue Shroud	○
JBC2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days	(NC 2883, 2884, 2885, 2886, 2888, 2889, 2890, 2891) £214
JBG2865		Premium - Racing Green Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2889, 2890, 2891) £104
JBG2866		Premium - Deep Ocean Blue Satin Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2890, 2891) £104
JBG2867		Premium - Frosty Blue Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2891) £104
JBG2868		Premium - Gold Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2890) £104
<b>RIM COLOUR*</b>			
JBC2882		Mid range - Black Rim Inserts	○
JBC2883		Mid range - Silver Rim Inserts	○
JBC2884		Mid range - Bengal Red Rim Inserts	○
JBC2885		Mid range - Mercury Blue Rim Inserts	○
JBC2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	(NC 2872, 2873, 2874, 2875, 2865, 2866, 2867, 2868) ○
JBG2888		Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2866, 2867, 2868) ○
JBG2889		Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2867, 2868) ○
JBG2890		Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2868) ○
JBG2891		Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2867) ○
<b>FORKS*</b>			
JBC2631		Single sided fork	£352
JBC2632		Double sided fork Without castor rim insert	◆
<b>FRONT WHEEL*</b>			
JBC0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBC0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JBC0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JBC0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	◆
<b>CASTOR LOCKS</b>			
JBC1490		Castor Lock	£197
<b>DRIVE WHEEL*</b>			
JBC0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBC0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
<b>C.T.C SUSPENSION*</b>			
JBC2870		CTC Suspension Control - Traction - Comfort	◆



<b>SUSPENSION SETTING*</b>				
JBC1910	<input type="checkbox"/>	Soft	(NC 1094)	◆
JBC1912	<input type="checkbox"/>	Medium	(NC 1094)	○
JBC1911	<input type="checkbox"/>	Firm		○
JBC1915	<input type="checkbox"/>	Very firm		○
<b>MOTOR TYPE*</b>				
JBC1294	<input type="checkbox"/>	Indoor setting, reduced speed approx. 3 km/h		○
JBC1282	<input type="checkbox"/>	6 km/h		◆
JBC1908	<input type="checkbox"/>	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254)	£543
<b>BATTERY TYPE*</b>				
JBC1259	<input type="checkbox"/>	60 Ah Gel		◆
<b>G TRAC</b>				
JBC1215	<input type="checkbox"/>	G-Trac Module		£323
<b>CHARGER*</b>				
JBC1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)		◆
JBC1313	<input type="checkbox"/>	Opti Charge Charger 10A (IP44)		£151
<b>ANTI-TIPPERS*</b>				
JBC1205	<input type="checkbox"/>	Anti-tipper		◆
<b>LIGHTS</b>				
JBC1437	<input type="checkbox"/>	LED Light system with Indicators 30 Lux headlights	(NC 2203, 1908)	£480
JBC1254	<input type="checkbox"/>	Without Light system		○
<b>KERB RIDER</b>				
JBC1244	<input type="checkbox"/>	Kerb Rider With park position feature	(NC 0556, 0589, 1294)	£212
<b>TRANSPORTATION KIT*</b>				
JBC1450	<input type="checkbox"/>	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering		£586
JBC1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory. Not compatible with Seat lifter		◆
JBC1401	<input type="checkbox"/>	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) required for variants with combined Lift & Tilt module Adequate headrest is compulsory	(NC 1091, 1092, 1093)	○
<b>EASY ACCESS FRONT TIE DOWNS</b>				
JBC2894	<input type="checkbox"/>	Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transport worthy versions (1494/1401) only		£49
<b>LAP BELT*</b>				
JBC1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest		○
JBC2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails		◆

POSITIONING BELT				
JBC2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£77
JBC2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£77
JBC2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
JBC2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
JBC2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
JBC2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
REAR VIEW MIRROR				
JBC1415		Rear view mirror RH	(NC 2217, 2285, 2209, 2213, 2265, 2275, 2283, 2287)	£44
JBC1339		Rear view mirror LH	(NC 2218, 2286, 2210, 2214, 2266, 2276, 2284, 2288)	£44
OTHER OPTIONS				
JBC2128		Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2172, 2173, 2227, 2228, 2233, 2240, 2257, 2291, 2292)	£77
JBC1312		Cane holder RH	(NC 2227, 2259, 1221)	£42
JBC1363		Cane holder LH	(NC 2228, 2260, 1222)	£42
JBC1416		Fender for Castors		£88
JBC1417		Fender for Drive wheels		£169
JBC1404		Tray, Swing Away RH	(NC 1707, 2209, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217)	£167
JBC1405		Tray, Swing Away LH	(NC 1708, 2210, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218)	£167

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

JBC2121	380 mm Seat Plate
JBC2122	430 mm Seat Plate
JBC2123	480 mm Seat Plate
JBC2124	530 mm Seat Plate
JBC2106	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)
JBC2132	Seat Depth 410 mm
JBC2143	Matrx Flo-tech Contour Visco
JBC2108	Backposts without push handles
JBC2171	Tension Adjustable Upholstery
JBC2179	Flex3 Tension Adjustable
JBC2172	Flex3 adjustable Backrest Plates without Backrest Cushion
JBC2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
JBC2174	Matrx Posture Back, height 400 mm (up to 136 kg max user weight)
JBC2175	Matrx Posture Back, height 500 mm (up to 136 kg max user weight)
JBC2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBC2177	Matrx Elite TR Back, height 500 mm
JBC2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBC2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBC2155	Black Spacer Fabric Upholstery
JBC2181	T-Arm detachable, Extended height
JBC2198	Padded Armrest with padded side guard (WATER FALL TYPE)
JBC2183	Inrigged Armrests Reduces SW by 25 mm each side
JBC2193	Flip-up Armrest Up to 150 kg max user weight
JBC2192	Following Armrest Height adjustable from 230 - 300 mm
JBC2194	Following Armrest Height adjustable from 300 - 360 mm
JBC2195	Short Armpad (325 mm)
JBC2196	Long Armpad (375 mm)
JBC2197	Long Ergonomic Armpad (397 mm)
JBC2199	Soft Armpad
JBC2184	Hip Support LH, pad size 140 x 100 mm
JBC2186	Hip Support LH, pad size 230 x 100 mm
JBC2188	Hip Support LH, pad size 330 x 100 mm
JBC2185	Hip Support RH, pad size 140 x 100 mm
JBC2187	Hip Support RH, pad size 230 x 100 mm
JBC2189	Hip Support RH, pad size 330 x 100 mm
JBC1350	Neckrest depth adjustable
JBC1432	Headrest sideways adjustable
JBC1427	Headrest with cheek support
JBC2324	Elan Standard Pad 255 mm / 10" wide w/standard cover

JBC2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover
JBC2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard
JBC1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBC1445	Lateral Support with extended arm, multi-adjustable pad
JBC1091	Manually settable Seat Tilt 0 - 15°
JBC1092	Powered "fixed pivot" Tilt 20°
JBC1093	Powered Tilt (0° - 30°) with CoG shift
JBC1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift, 6 point Tie downs must be selected
JBC0977	Powered Recliner 30°
JBC0556	Centre mount legrest, angle settable Length: 280 - 385 mm
JBC0589	LNX Power centre mount legrest, length compensating Length: 340 - 410 mm
JBC0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBC0526	Standard hanger 80° pair Lower leg length: 350 - 450 mm
JBC0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBC0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBC0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad
JBC0504	Angle Adjustable Footplate pair
JBC0513	Angle and depth adjustable Footplate pair
JBC2201	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBC2203	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBC2204	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBC2205	Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets
JBC2206	Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233
JBC2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JBC1707	Swing Away Bracket RH
JBC1708	Swing Away Bracket LH
JBC1709	Matrx Resolve Swing in/out bracket RH
JBC1710	Matrx Resolve Swing in/out bracket LH
JBC2209	Swing away RH, multi-adjustable
JBC2210	Swing away LH, multi-adjustable
JBC2213	Centre Mount, swing-away to RH, Nucleus
JBC2214	Centre Mount, swing-away to LH, Nucleus
JBC2227	Bracket RH, mounted for Attendant
JBC2228	Bracket LH, mounted for Attendant
JBC2220	Dual Control LiNX (DLX-REM050 as secondary remote)
JBC2219	Compact Dual Control LiNX (DLX-ACU200)
JBC2245	Invacare LiNX IDC Dual control
JBC2259	Bracket RH, Attendant
JBC2260	Bracket LH, Attendant
JBC1221	Bracket RH, detachable

JBC1222	Bracket LH, detachable
JBC2231	Proportional Compact Remote DLX-CR400
JBC2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBC2207	RH Mount
JBC2208	LH Mount
JBC2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBC2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBC1289	Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest
JBC1288	Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest
JBC1287	Swing Away Bracket RH, powered
JBC1286	Swing Away Bracket LH, powered
JBC2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBC2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBC2257	Four-Switch Proximity Array (ASL106)
JBC2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBC2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBC2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JBC2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBC2285	Tray Mount RH (100x150 mm)
JBC2286	Tray Mount LH (100x150 mm)
JBC2283	Centre Mount, swing-away to RH, Nucleus
JBC2284	Centre Mount, swing-away to LH, Nucleus
JBC2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBC2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBC2217	Bracket RH, mounted for user
JBC2218	Bracket LH, mounted for user
JBC2291	Bracket RH, mounted for attendant
JBC2292	Bracket LH, mounted for attendant
JBC2253	USB Charger Port DLX-USB02-A
JBC2258	Remote Stop Switch (ASL504)
JBC2273	IN200 Module (Jack port for external power button )
JBC2872	Mid range - Silver Shroud
JBC2873	Mid range - Bengal Red Shroud
JBC2874	Mid range - Mercury Blue Shroud
JBC2875	Mid range - Midnight Purple Shroud
JBX0299	Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days
JBG2865	Premium - Racing Green Shroud
JBG2866	Premium - Deep Ocean Blue Satin Shroud
JBG2867	Premium - Frosty Blue Shroud
JBG2868	Premium - Gold Shroud
JBC2883	Mid range - Silver Rim Inserts
JBC2884	Mid range - Bengal Red Rim Inserts

JBC2885	Mid range - Mercury Blue Rim Inserts
JBC2886	Mid range - Midnight Purple Rim Inserts
JBX2887	Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days
JBG2888	Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit
JBG2889	Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit
JBG2890	Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit
JBG2891	Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit
JBC1910	Soft
JBC1912	Medium
JBC1294	Indoor setting, reduced speed approx. 3 km/h
JBC1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBC1437	LED Light system with Indicators 30 Lux headlights
JBC1254	Without Light system
JBC1244	Kerb Rider With park position feature
JBC2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBC2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBC2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBC2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBC2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBC2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBC1415	Rear view mirror RH
JBC1339	Rear view mirror LH
JBC2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBC1312	Cane holder RH
JBC1363	Cane holder LH
JBC1416	Fender for Castors
JBC1404	Tray, Swing Away RH
JBC1405	Tray, Swing Away LH





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):