



**Yes, you can:**

# Power Retail Prescription Form

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2025  
 Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## AVIVA RX40 Standard Seat

### Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- One piece seat plate width 380mm (15") or 430mm (17") with standard foam cushion (50mm, 2")
- seat depths 410mm - 510mm (16" - 20")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manually settable Seat Tilt 0 - 15°
- Manual settable backrest Recline (0° - 30°, in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX Control REM211 with standard Bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic
- 14" (355 mm) x 3" (75 mm), Pneumatic drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole Durawatt motors and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £4,521**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

| Code                  | Description   |        |
|-----------------------|---|--------|
| <b>BASE MODEL*</b>    |   |        |
| JBF0001               | AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries | £4,521 |
| JBX1300               | Special Extended lead time: + 10 days   | POA    |
| <b>CHASSIS WIDTH*</b> |   |        |
| JBL2810               | Chassis Narrow (560 mm/22")   | ◆      |
| <b>SEATING TYPE*</b>  |   |        |
| JBF2105               | Modulite Standard Seat Standard Seat Plate and Backrest Frame in fixed Width (NC 1401)  | ◆      |

| <b>SEAT FRAME*</b>  |                          |   |      |
|---|--------------------------|---|------|
| JBF2110   | <input type="checkbox"/> | One Piece Seat Plate  | ◆    |
| <b>SEAT WIDTH*</b>  |                          |   |      |
| JBF2121   | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests  | ○    |
| JBF2122   | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests  | ◆    |
| <b>SEAT DEPTH*</b>  |                          |   |      |
| JBF2132   | <input type="checkbox"/> | Seat Depth 410 mm (NC 2197)   | ○    |
| JBF2134   | <input type="checkbox"/> | Seat Depth 460 mm   | ○    |
| JBF2136   | <input type="checkbox"/> | Seat Depth 510 mm   | ○    |
| <b>SEAT CUSHION*</b>  |                          |   |      |
| JBF2141   | <input type="checkbox"/> | Standard Foam Cushion 50 mm height  | ◆    |
| JBF2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite  | £117 |
| JBF2143   | <input type="checkbox"/> | Matrx Flo-tech Contour Visco; (2-way stretch material)  | £165 |
| JBF2140   | <input type="checkbox"/> | Without seat cushion  | ○    |
| <b>BACKREST SUPPORT*</b>  |                          |   |      |
| JBF2161   | <input type="checkbox"/> | Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs   | ◆    |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |   |      |
| JBF2171   | <input type="checkbox"/> | Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm        | ◆    |
| JBF2170   | <input type="checkbox"/> | Without Backrest Upholstery (Back Posts/Back Frame only)  | ○    |
| <b>UPHOLSTERY TYPE*</b>   |                          |   |      |
| JBF2151   | <input type="checkbox"/> | Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)   | ◆    |
| <b>ARMREST (Pair)*</b>  |                          |   |      |
| JBF2191   | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm   | ◆    |
| JBF2181   | <input type="checkbox"/> | T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm  | £44  |
| JBF2183   | <input type="checkbox"/> | Inrigged Armrests Reduces SW by 25 mm each side   | ○    |
| <b>ARMPAD (Pair)*</b>   |                          |   |      |
| JBF2195   | <input type="checkbox"/> | Short Armpad (325 mm)   | ◆    |
| JBF2196   | <input type="checkbox"/> | Long Armpad (375 mm)  | ○    |
| JBF2197   | <input type="checkbox"/> | Long Ergonomic Armpad (397 mm)  | ○    |
| JBF2199   | <input type="checkbox"/> | Soft Armpad<br>In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote            | ○    |
| <b>HEADREST / NECKREST - Adequate headrest option must be selected*</b> |                          |   |      |
| JBF1350   | <input type="checkbox"/> | Neckrest depth adjustable   | £81  |
| JBF1495   | <input type="checkbox"/> | Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts | ○    |
| <b>SEAT LIFTER / SEAT TILT*</b>   |                          |   |      |
| JBF1091   | <input type="checkbox"/> | Manually settable Seat Tilt 0 - 15°   | ○    |
| JBF1092   | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20°  | ◆    |

| <b>SEAT T GROUND HEIGHT*</b>      |  |  |      |
|-----------------------------------|--|--|------|
| JBF2857                           |  | Seat-to-Floor Height 435 mm  | ○    |
| JBF2858                           |  | Seat-to-Floor Height 460 mm  | ○    |
| JBF2859                           |  | Seat-to-Floor Height 485 mm  | ○    |
| <b>BACKREST RECLINER*</b>         |  |  |      |
| JBF0975                           |  | Manual Settable Recliner (in 3.75° steps) 90° - 120°   | ◆    |
| <b>FOOTREST HANGER TYPE*</b>      |  |  |      |
| JBF0526                           |  | Standard hanger 80° pair Lower leg length: 350 - 450 mm (NC 0522)  | ◆    |
| JBF0503                           |  | Hanger Vari F pair Lower leg length: 290 - 460 mm<br>Angle adjustment: 70 - 0° with tools                                | £121 |
| JBF0533                           |  | Hanger Vari A pair Lower leg length: 290 - 460 mm<br>Angle adjustment: 70 - 0° with handwheel and calf pad               | £234 |
| JBF0531                           |  | Hanger ADE pair Lower leg length: 290 - 460 mm (NC 0522)<br>Powered Angle adjustment: 80 - 0° with calf pad              | £669 |
| <b>FOOTPLATE TYPE*</b>            |  |  |      |
| JBF0512                           |  | Fixed Footplate pair (NC 0533, 0531)   | ◆    |
| JBF0504                           |  | Angle Adjustable Footplate pair (NC 0526)  | £32  |
| JBF0513                           |  | Angle and depth adjustable Footplate pair  | £118 |
| JBF0522                           |  | Angle adjustable Footboard (NC 0531, 1244)   | £154 |
| <b>CONTROLS*</b>                  |  |  |      |
| JBF2203                           |  | Remote LiNX (NC 1437, 1908, 1909)<br>DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators | ◆    |
| JBF2204                           |  | Remote LiNX<br>DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights      | £151 |
| <b>CONTROL BRACKETRY*</b>         |  |  |      |
| JBF1207                           |  | Standard Bracket RH (NC 1404)  | ○    |
| JBF1208                           |  | Standard Bracket LH (NC 1405)  | ○    |
| JBF1707                           |  | Swing Away Bracket RH (NC 1404)  | £133 |
| JBF1708                           |  | Swing Away Bracket LH (NC 1405)  | £133 |
| JBF2227                           |  | Bracket RH, mounted for Attendant (NC 2220, 2219, 1312)  | £64  |
| JBF2228                           |  | Bracket LH, mounted for Attendant (NC 2220, 2219, 1363)  | £64  |
| JBX2239                           |  | Individual mounting option per request (Special)<br>Requires additional information/specification                        | POA  |
| <b>DUAL CONTROL</b>               |  |  |      |
| JBF2220                           |  | Dual Control LiNX (DLX-REM050 as secondary remote) (NC 2273, 2227, 2228)   | £395 |
| JBF2219                           |  | Compact Dual Control LiNX (DLX-ACU200) (NC 2227, 2228)   | £395 |
| <b>DUAL CONTROL BRACKETRY</b>     |  |  |      |
| JBF2259                           |  | Bracket RH, Attendant (NC 2291, 1312)  | ○    |
| JBF2260                           |  | Bracket LH, Attendant (NC 2292, 1363)  | ○    |
| <b>LINX 10-WAY SWITCH</b>         |  |  |      |
| JBF2251                           |  | 10-Way Switch DLX-SW10-A for direct Powered Seating Access   | £505 |
| <b>LINX 10-WAY SWITCH BRACKET</b> |  |  |      |
| JBF2217                           |  | Bracket RH, mounted for user (NC 1404)   | ○    |
| JBF2218                           |  | Bracket LH, mounted for user (NC 1405)   | ○    |
| JBF2291                           |  | Bracket RH, mounted for attendant (NC 2259)  | £64  |
| JBF2292                           |  | Bracket LH, mounted for attendant (NC 2260)  | £64  |
| <b>LINX IN200 MODULE</b>          |  |  |      |
| JBF2273                           |  | IN200 Module (Jack port for external power button + control inputs via DB9)  | £241 |

| <b>FURTHER LINX COMPONENTS / ACCESSORIES</b> |  |  |      |
|--|--|--|------|
| JBF2300                                      |  | To be enclosed only  | ○    |
| JBX2301                                      |  | To be connected and programmed (Special)   | POA  |
| <b>BUDDY BUTTONS</b>                         |  |  |      |
| JBF1685                                      |  | Red Buddy Button   | £88  |
| JBF1686                                      |  | Yellow Buddy Button  | £88  |
| JBF1687                                      |  | Green Buddy Button   | £88  |
| JBF1690                                      |  | Blue Buddy Button  | £88  |
| <b>EGG SWITCH</b>                            |  |  |      |
| JBF2254                                      |  | Egg Switch, black  | £88  |
| <b>PIKO BUTTONS</b>                          |  |  |      |
| JBF1711                                      |  | Red Piko Button  | £158 |
| JBF1712                                      |  | Yellow Piko Button   | £158 |
| JBF1713                                      |  | Green Piko Button  | £158 |
| JBF1714                                      |  | Orange Piko Button   | £158 |
| JBF1715                                      |  | Blue Piko Button   | £158 |
| JBF1716                                      |  | Black Piko Button  | £158 |
| <b>BODYPOINT JOYSTICK KNOBS</b>              |  |  |      |
| JBF2721                                      |  | 3 inch U-Shaped, Joystick Handle, Bodypoint PC101A   | £63  |
| JBF2722                                      |  | 4 inch U-Shaped, Joystick Handle, Bodypoint PC102A   | £63  |
| JBF2723                                      |  | Domed Rubber, Joystick Handle, Bodypoint PC107A  | £63  |
| JBF2724                                      |  | Mushroom Grooved, Joystick Handle, Bodypoint PC110A  | £63  |
| <b>SHROUD COLOUR*</b>                        |  |  |      |
| JBF2871                                      |  | Mid range - Black/Anthracite Shroud  | ○    |
| JBF2872                                      |  | Mid range - Silver Shroud  | ○    |
| JBF2873                                      |  | Mid range - Bengal Red Shroud  | ○    |
| JBF2874                                      |  | Mid range - Mercury Blue Shroud  | ○    |
| JBF2875                                      |  | Mid range - Midnight Purple Shroud   | ○    |
| JBX0299                                      |  | Premium - Individual Shroud Colour RAL number required<br>Extended lead lime: + 15 days  | £240 |
| <b>RIM COLOUR*</b>                           |  |  |      |
| JBF2882                                      |  | Mid range - Black Rim Inserts  | ○    |
| JBF2883                                      |  | Mid range - Silver Rim Inserts   | ○    |
| JBF2884                                      |  | Mid range - Bengal Red Rim Inserts   | ○    |
| JBF2885                                      |  | Mid range - Mercury Blue Rim Inserts   | ○    |
| JBF2886                                      |  | Mid range - Midnight Purple Rim Inserts  | ○    |
| JBX2887                                      |  | Premium - Individual Rim Insert Colour RAL number required<br>Extended lead lime: + 15 days Price Included with Premium Shroud Kit | ○    |
| <b>FORKS*</b>                                |  |  |      |
| JBF2631                                      |  | Single sided fork  | £352 |
| JBF2632                                      |  | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)  | ◆    |

| FRONT WHEEL*        |                          |  |      |
|---------------------|--------------------------|--|------|
| JBF0835             | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic   | ◆    |
| JBF0836             | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre   | ○    |
| JBF0838             | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed   | ○    |
| JBF0828             | <input type="checkbox"/> | 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)   | ○    |
| CASTOR LOCKS        |                          |  |      |
| JBF1490             | <input type="checkbox"/> | Castor Lock  | £197 |
| DRIVE WHEEL*        |                          |  |      |
| JBF0704             | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Pneumatic  | ◆    |
| JBF0705             | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert   | ○    |
| JBF0706             | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed  | ○    |
| C.T.C SUSPENSION*   |                          |  |      |
| JBF2870             | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort  | ◆    |
| SUSPENSION SETTING* |                          |  |      |
| JBF1910             | <input type="checkbox"/> | Soft (suggested for up to 70 kg User Weight)   | ◆    |
| JBF1912             | <input type="checkbox"/> | Medium (suggested for up to 90 kg User Weight)   | ○    |
| JBF1911             | <input type="checkbox"/> | Firm (suggested for up to 110 kg User Weight)  | ○    |
| JBF1915             | <input type="checkbox"/> | Very firm (suggested for up to 136 kg User Weight)   | ○    |
| MOTOR TYPE*         |                          |  |      |
| JBF1294             | <input type="checkbox"/> | Indoor setting, reduced speed approx. 3 km/h   | ○    |
| JBF1282             | <input type="checkbox"/> | 6 km/h   | ◆    |
| JBF1908             | <input type="checkbox"/> | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights (NC 2203, 1254, 2837)                        | £543 |
| JBF1909             | <input type="checkbox"/> | 12 km/h Ultra High Speed Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights (NC 0828, 2203, 1254, 2837) | £653 |
| BATTERY TYPE*       |                          |  |      |
| JBF1259             | <input type="checkbox"/> | 60 Ah Gel  | ◆    |
| JBF1260             | <input type="checkbox"/> | 73,5 Ah Gel (NC 1450)  | £92  |
| G TRAC              |                          |  |      |
| JBF1215             | <input type="checkbox"/> | G-Trac Module  | £324 |
| CHARGER*            |                          |  |      |
| JBF1343             | <input type="checkbox"/> | Opti Charge Charger 8A (IP44)  | ◆    |
| JBF1313             | <input type="checkbox"/> | Opti Charge Charger 10A (IP44)   | £151 |
| ANTI-TIPPERS*       |                          |  |      |
| JBF1205             | <input type="checkbox"/> | Anti-tipper  | ◆    |
| LIGHTS*             |                          |  |      |
| JBF1437             | <input type="checkbox"/> | LED Light system with Indicators 30 Lux headlights (NC 2203)   | £479 |
| JBF2837             | <input type="checkbox"/> | LED Light system extra bright with indicators 60 Lux headlights (NC 2203)  | £143 |
| JBF1254             | <input type="checkbox"/> | Without Light system (NC 1909)   | ○    |
| KERB RIDER          |                          |  |      |
| JBF1244             | <input type="checkbox"/> | Kerb Rider With park position feature (NC 1294)  | £212 |

**TRANSPORTATION KIT\***

|         |  |  |                   |
|---------|--|--|-------------------|
| JBF1450 |  | Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box<br>Docking Station for car implementation has to be ordered separately with DAHL Engineering | £586              |
| JBF1494 |  | 4-Point tie-down<br>(Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory   | ◆                 |
| JBF1401 |  | 6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory                          | (NC 1091, 1092) ○ |

**EASY ACCESS FRONT TIE DOWNS**

|         |  |  |     |
|---------|--|--|-----|
| JBF2894 |  | Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transport worthy versions (1494/1401) only | £49 |
|---------|--|--|-----|

**LAP BELT\***

|         |  |  |   |
|---------|--|--|---|
| JBF1329 |  | Lap belt with Auto Lock Standard mount on backrest | ◆ |
|---------|--|--|---|

**REAR VIEW MIRROR**

|         |  |                     |               |
|---------|--|---------------------|---------------|
| JBF1415 |  | Rear view mirror RH | (NC 2217) £44 |
| JBF1339 |  | Rear view mirror LH | (NC 2218) £44 |

**OTHER OPTIONS**

|         |  |                         |                      |
|---------|--|-------------------------|----------------------|
| JBF1312 |  | Cane holder RH          | (NC 2227, 2259) £42  |
| JBF1363 |  | Cane holder LH          | (NC 2228, 2260) £42  |
| JBF1416 |  | Fender for Castors      | £88                  |
| JBF1417 |  | Fender for Drive wheels | £169                 |
| JBF1404 |  | Tray, Swing Away RH     | (NC 1707, 2217) £167 |
| JBF1405 |  | Tray, Swing Away LH     | (NC 1708, 2218) £167 |

**TOTAL PRICE**

**Non Conforming Options:**

|         |   |
|---------|---|
| JBF2197 | Long Ergonomic Armpad (397 mm)  |
| JBF1091 | Manually settable Seat Tilt 0 - 15°   |
| JBF1092 | Powered "fixed pivot" Tilt 20°  |
| JBF0526 | Standard hanger 80° pair Lower leg length: 350 - 450 mm   |
| JBF0533 | Hanger Vari A pair Lower leg length: 290 - 460 mm<br>Angle adjustment: 70 - 0° with handwheel and calf pad  |
| JBF0531 | Hanger ADE pair Lower leg length: 290 - 460 mm<br>Powered Angle adjustment: 80 - 0° with calf pad   |
| JBF0522 | Angle adjustable Footboard  |
| JBF2203 | Remote LiNX<br>DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators  |
| JBF1707 | Swing Away Bracket RH   |
| JBF1708 | Swing Away Bracket LH   |
| JBF2227 | Bracket RH, mounted for Attendant   |
| JBF2228 | Bracket LH, mounted for Attendant   |
| JBF2220 | Dual Control LiNX (DLX-REM050 as secondary remote)  |
| JBF2219 | Compact Dual Control LiNX (DLX-ACU200)  |
| JBF2259 | Bracket RH, Attendant   |
| JBF2260 | Bracket LH, Attendant   |
| JBF2217 | Bracket RH, mounted for user  |
| JBF2218 | Bracket LH, mounted for user  |
| JBF2291 | Bracket RH, mounted for attendant   |
| JBF2292 | Bracket LH, mounted for attendant   |
| JBF2273 | IN200 Module (Jack port for external power button + control inputs via DB9)   |
| JBF0828 | 8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre  |
| JBF1294 | Indoor setting, reduced speed approx. 3 km/h  |
| JBF1908 | 10 km/h Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights   |
| JBF1909 | 12 km/h Ultra High Speed Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights  |
| JBF1437 | LED Light system with Indicators 30 Lux headlights  |
| JBF2837 | LED Light system extra bright with indicators 60 Lux headlights   |
| JBF1254 | Without Light system  |
| JBF1244 | Kerb Rider With park position feature   |
| JBF1450 | Preparation for DAHL-Docking station<br>(Adequate headrest is compulsory) Includes a mounting bracket in the battery box<br>Docking Station for car implementation has to be ordered separately with DAHL Engineering |
| JBF1401 | 6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with<br>combined Lift & Tilt module Adequate headrest is compulsory                          |
| JBF1312 | Cane holder RH  |
| JBF1363 | Cane holder LH  |
| JBF1404 | Tray, Swing Away RH   |
| JBF1405 | Tray, Swing Away LH   |



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_